

SALVIN[®]
2025 COMPREHENSIVE CATALOG

**Everything For Your
Implant Practice
But The Implants[®]**



**Dr. Michael A. Pikos**

Oral Surgeon
Palm Harbor, FL

"Salvin Dental is an important resource for our practice and our course participants. They are great people to work with and won't let you make a mistake."

**Dr. Scott Froum**

Periodontist
New York, NY

"Not only does Salvin have a great product line, but their support of academics and advancement of dental technology is unparalleled."

**Dr. Hom-Lay Wang**

Periodontist
Ann Arbor, MI

"Salvin Dental has an excellent product line. Their commitment to quality and service has been a valuable resource for both my practice and my students."

**Dr. Shaun Young**

Oral Surgeon
Tampa, FL

"Salvin continues to be a significant portion of our regenerative portfolio due to its consistency. They have always proven to be at the top when it comes to grafting materials and our practice continues to see great results."

**Dr. Michael Sonick**

Periodontist
Fairfield, CT

"There's a rare, consistent quality to Salvin products and service that's hard to find anywhere else. I use the Salvin catalog confidently for the implant courses and seminars I teach."

**Dr. Alvaro Ordonez**

General Dentist
Coral Gables, FL

"The team at Salvin Dental is always looking out for my best interest. Relationships like that are hard to find and impossible to replace."

**Dr. Dwayne Karateew**

Periodontist/Prosthodontist
Chicago, IL

"Every time I place an order with Salvin Dental my expectations for quality and service are exceeded. I consider them to be an invaluable partner in my practice and an asset to the profession."

**Dr. Randy Resnik**

Prosthodontist
Pittsburgh, PA

"The entire Salvin team is always willing to go above and beyond. I know I can count on them to get me exactly what I need, when I need it."

**Dr. James Wilson**

Periodontist
Tampa Bay, FL

"For over 20 years Salvin has been my most valuable corporate partner. They are customer satisfaction oriented and have a comprehensive portfolio. Their culture frowns on high-pressure salesmanship making them seem more like friends than sales reps."

**Dr. Scott Ganz**

Prosthodontist
Fort Lee, NJ

"As implant dentistry has evolved, so has Salvin Dental as a valuable resource to my practice. They are always well informed, and provide a personal touch that is very much appreciated by me and my staff."

**Dr. Dmitriy Ivanov**

Periodontist,
Redondo Beach, CA

"We have been using Salvin Dental products in our periodontal practices for over 10 years with the highest success. The quality of their products is exceptional, and they consistently deliver outstanding results for our patients."

**Dr. Mia Geisinger**

Periodontist
Birmingham, AL

"As an educator and a researcher, I am consistently impressed by Salvin Dental's support for residents and new practitioners. They have been an amazing resource for our practice and our students!"

**Dr. Brian Jackson**

General Dentist
Utica, NY

"I take pride in being able to say that I have been with Salvin Dental Specialties since its very beginning. Today, Salvin is a critical player on my team. I could not imagine performing implant surgery without Salvin."

**Dr. Jaime Lozada**

Prosthodontist
Loma Linda, CA

"The Salvin catalog is a snapshot of the very best products for the up-to-date implant practice. They always have the newest technology, and going the extra mile is the standard procedure at Salvin."

**Dr. Pat Allen**

Periodontist
Dallas, TX

"It's the quality people behind the quality products that sets Salvin Dental apart from the rest. Their knowledgeable, helpful, and friendly staff makes working with them a very pleasant experience."

**Dr. Christopher Richardson**

Periodontist
Richmond, VA

"Salvin Dental Specialties is so much more than a dental supply company. I consider Salvin and their team members to be an essential and collaborative partner in the care I provide my patients."

**Dr. Russ Baer**

General Dentist
Chicago, IL

"There are three things you need in Implant Dentistry: a well trained dentist, a good implant system, and your Salvin catalog. Salvin's products and customer support are an irreplaceable part of my implant practice."

**Dr. Fadi Hasan**

Periodontist,
McLean, VA

"What truly sets Salvin apart is their support team and management, who provide a personalized, one-of-a-kind support experience at any time of day or night. Their commitment to excellence makes them an invaluable partner in our daily practice."

**Dr. Tim Kosinski**

General Dentist
Bingham Farms, MI

"Working with the Salvin team for the past 30 years, I know that I can count on them for top of the line products and world class service. They are a valuable resource for my practice."

**Dr. Gordon Christensen**

Prosthodontist
Provo, UT

"I consider Salvin Dental to be a premium company in providing implant accessory supplies. The Salvin catalog is not only a resource for my practice, but also for the doctors who attend our courses."

**Dr. Jennifer Hirsch Doobrow**

Periodontist
Cullman, AL

"Loyalty, integrity, and excellence are what Salvin Dental embodies. Their commitment to impeccable service and outstanding quality is why they are a leader in the dental implant community."



Dr. George Kotsakis
Periodontist
Hoboken, NJ

"It's always a pleasure to work with Salvin Dental because it's always a partnership rather than a business transaction. They genuinely care and always deliver. I would highly recommend Salvin."



Dr. Natalie Wong
Prosthodontist
Toronto, Ontario, Canada

"Delivering excellent quality instrumentation with great service is what you can expect from Salvin. Professional, knowledgeable, and courteous, I consider them an integral part of my dental implant team."



Dr. Robert Levine
Periodontist
Philadelphia, PA

"Working with Salvin is like working with family. You get service, great products and ability to talk directly to staff who are incredibly helpful, knowledgeable & personable. I highly recommend Salvin and their quality products."



Dr. Burt Melton
Prosthodontist
Albuquerque, NM

"Salvin Dental has all the things you need for a day-to-day implant practice as well as all the special things you need from time to time. Salvin's products are always available for you when you need them in a hurry."



Dr. Bach Le
Oral Surgeon
Whittier, CA

"Salvin Dental has one of the most professional sales teams in the dental industry. They know that my time is valuable and always deliver more than they promise."



Dr. Lee Silverstein
Periodontist
Atlanta, GA

"When Salvin says everything for the implant practice but the implants, they mean it! Salvin is our prime resource for expertise and incredible service. They are a doctor's true friend in our industry."



Dr. Joseph Cristoforo
Periodontist
Madison, WI

"Knowledge. Service. Commitment. Salvin Consistently provides vast product knowledge, unparalleled service, and commitment to enhancing the clinical experience for both the provider and patient."



Dr. Todd Engel
General Dentist
Charlotte, NC

"Salvin Dental consistently delivers great products and excellent service. They have become an important resource for my practice."



Dr. Peter Nordland
Periodontist
La Jolla, CA

"The best aspect of Salvin Dental is the personal attention and knowledge that all of the staff members have whenever I call. I am always greeted with a kind, courteous, knowledgeable representative."



Dr. George Romanos
Periodontist
Stony Brook, NY

"The Salvin Group is the best partner of my practice due to the immediate customer service and professionalism. Their products are of high quality and a joy to work with them."



Dr. Carlos Rubio
General Dentist
Algodones, Mexico

"Salvin has been a major part in my professional growth as a dentist and a teacher. I have always been able to count on them to have the best quality material, as well as the service. Thanks Salvin for being there along the way."



Dr. Richard Rasmussen
Periodontist
Tampa, FL

"Salvin has a great team with one thing in mind...To make my practice and me their first priority, every time I call. I know that when I call Salvin, I've called the right place."



Dr. Sayward E. Duggan
Periodontist
Yorktown, VA

"Salvin consistently continues to set the bar high in the dental industry with excellent products, competitive prices, and above all superior client services. Thank you Salvin for your loyal professional support."



Dr. Tyler Borg
Periodontist
Northglenn, CO

"The team at Salvin Dental offers consistent support, continued product innovation, and top notch customer service. All in the pursuit of the ideal doctor and patient experience alike."



Dr. Robert Heller
Prosthodontist
Lewis Center, OH

"I confidently recommend Salvin Dental to our Midwest Institute students. Their knowledge and helpfulness make them a pleasure to do business with."



Dr. George Duello
Periodontist
St. Louis, MO

"These guys are more than just suppliers to my practice. Salvin Dental's team are proven consultants who help us every time we call. These professionals really know implant dentistry."



Dr. Craig M. Misch
Oral Surgeon/Prosthodontist
Sarasota, FL

"Salvin has a thorough product selection with great personalized customer service. I keep going back to Salvin for whatever I need for my implant practice and highly recommend them to my colleagues."



Dr. Rick Ferguson
General Dentist
Pembroke Pines, FL

"Dealing with Salvin is a first-class experience. They have a reputation for always putting their customers first. They provide great service, value, and expertise to my practice and the students at our course."



Dr. Jack Krauser
Periodontist
Boca Raton, FL

"The Salvin catalog is perhaps the most valuable tool in implant dentistry. Their unique products and excellent service are a tremendous resource for my private practice and for our students at the school."



Dr. Jason Stoner
Periodontist
Dublin, OH

"Every time I call Salvin, I can count on speaking with a knowledgeable person with the answers I need. They know that my time is valuable and always deliver what they promise."



Dr. Adam S. Folek
General Dentist
Norfolk, VA

"Salvin has provided expertise and experience necessary to take my implant practice to the next level. They do everything they can to ensure my success. Salvin helps me equip my offices to offer the highest level of treatment to my patients."



Dr. John Russo
Periodontist
Sarasota, FL

"Salvin Dental is a great resource for my practice and especially for the doctors who take my courses. Salvin has what you need whether you're getting started or doing advanced procedures."



Dr. Michael Pruett
General Dentist
Augusta, GA

"I love the way that Salvin does business. Their willingness to go the extra mile and their commitment to serving makes it a pleasure to deal with them."



Dr. Nick Shumaker
Periodontist
Fort Collins, CO

"When it comes to quality products and great customer service, nobody does it better than Salvin!"



Dr. Fred Sakamoto
Periodontist
Blacklick, OH

"Salvin Dental runs a first class operation with excellent customer service and products. We consider Salvin as part of our team with periodontal, implant and related surgical cases."



Dr. Canyon Keeney
Oral Surgeon
Mechanicsville, VA

"Salvin Dental has one of the most professional sales teams in the dental industry. They know that my time is valuable and always deliver more than they promise."



Dr. Lauralee Nygaard
Periodontist
Spokane Valley, WA

"When Salvin says everything for the implant practice but the implants, they mean it! Salvin is our prime resource for expertise and incredible service. They are a doctor's true friend in our industry."



Dr. Etienne Lacrampe
Periodontist
Walnut Creek, CA

"I use Salvin for all my biologics. There is no one I trust more to provide the best results for both my patients and my referring dentists."



Dr. Richard Martin
Oral Surgeon
Lewisville, TX

"When you work with Salvin Dental, you do not feel like you are working with a company, you feel like you are working with your family. Their attention to detail and reliability is unmatched in the dental field."



Dr. Eric Hamrick
Periodontist
Greenville, SC

"For over twenty-five years, Salvin has been our practice's go-to company for surgical supplies and instruments. Their service, vast inventory, and prices are comparable to none."



Dr. Avi Schetritt
Periodontist
Miami, FL

"Great products, great service, and great people! Their professionalism is what sets Salvin apart. They deliver every time!"



Dr. Bradley Ross
Periodontist
Miami, FL

"Salvin provides top quality surgical armamentarium and biomaterials, but their people truly make this company an essential asset to my practice. Their knowledge, integrity, and customer service are second to none."



Dr. Matthew Fien
Periodontist
Deerfield Beach, FL

"I am always impressed with the professional service and knowledge from the team at Salvin Dental. They make it easy for me to recommend them to my colleagues."



Dr. Steve Wallace
Periodontist
Wilmington, NC

"I turn to Salvin when I need answers to my questions. Their knowledge and helpfulness make them a pleasure to do business with."



Dr. Justin Zalewsky
Periodontist
Alexandria, VA

"Salvin has been a business partner for years, because they truly do partner with us. Their products and services are unmatched. Their people care about our practice as much as we do."



Dr. Nancy Newhouse
Periodontist
Independence, MO

"Nobody does it better! Salvin is just what they say they are: 'Everything For Your Implant Practice But The Implants' all wrapped up in the best service in the industry."



Dr. David Lee Hill
Oral Surgeon
Chapel Hill, NC

"I can't think of a single source that offers a more comprehensive quality line of products to support my implant practice than Salvin Dental. They always have what I need."



Dr. James Woodyard
Periodontist
Newburgh, IN

"Knowing that I can expect the right products and great service is why I choose to do business with Salvin. Dealing with them is a first-class experience every time I call."



Dr. James Rutkowski
General Dentist
Clarion, PA

"Honesty and integrity are values that come to mind when working with Salvin Dental. I know that I'm going to be taken care of every time I call."



Dr. Israel Puterman
Periodontist
Chevy Chase, MD

"Getting the Salvin catalog is like being a kid again, flipping through the toy catalog and seeing all the things that I want. The quality of the products is only surpassed by the quality of the people. They truly want your success above all else."



Dr. Paul Fugazzotto
Periodontist
Milton, MA

"Salvin has, in my experience with them over decades, consistently demonstrated a degree of honesty and openness, and a genuine desire to help the clinician, that is unsurpassed in the dental community."



Dr. Jason Souyias
Periodontist
Port Huron, MI

"The team at Salvin is amazing! They know customer service better than anyone in dental. That is why they are so easy to work with."

SALVIN®

Everything For Your Implant Practice But The Implants®

2025 COMPREHENSIVE CATALOG

*We're in business to make your lives easier
and improve your surgical outcomes.*

Call us, we'll take great care of you!



Wishing Bob Salvin A Happy Retirement!

We'd like to thank Bob Salvin for his distinguished leadership, and honor his legacy as he moves to retirement. Words can't adequately express the incredible and lasting impact of Bob's contributions to the dental industry, his accomplishments, and his loyal partnership over more than four decades of exemplary leadership as Founder & CEO of Salvin Dental.

We wish Bob well as he continues to be active in the community, serving on the Board of Trustees of Queens University in Charlotte, spending valuable time with his wonderful wife, Susan, and enjoying his passion for photography and classic cars. He tells us that his new job is learning something new and doing something interesting every day.

William Simmons, Bob Salvin, and Greg Slayton at the 2024 AAP Ice Cream Social Sponsored by Salvin.

800-535-6566 [Toll Free U.S. & Canada] • www.salvin.com



Salvin Dental Specialties, LLC
3450 Latrobe Drive
Charlotte, NC 28211 USA
Ph: 704-442-5400 Fax: 704-442-5424
Email: orders@salvin.com
Hours: M-Th: 8:30-7:00 F: 8:30-6:00
(U.S. Eastern Standard Time)
Quality Management System
ISO 13485:2016 Certified

See It Demonstrated!
Look For This Icon Throughout
Our Catalog



Sustainable Forest Initiative



By supporting the SFI program, Salvin is helping to strengthen forest practices in North America and promote responsible sourcing globally.

SALVIN INSTRUMENT LIFETIME WARRANTY

Salvin Branded Instruments are guaranteed to be free from defects in workmanship and material.

Any Salvin instrument which proves defective in workmanship or material will either be repaired or replaced at our discretion, without charge.

Warranty does not cover routine maintenance, re-tipping of carbide instruments or sharpening of scissors. Individual manufacturer warranties apply to products other than instruments.

Terms Of Sale: Payment terms are Visa, MasterCard, Discover or American Express. Prices in catalog are in U.S. funds and are subject to change due to unforeseen circumstances. Sales and Taxes, if applicable, are customer's responsibility. **Foreign Orders:** Foreign orders may be called in to 704-442-5400, or placed online at www.salvin.com. Duties and import taxes, if any, are customer's responsibility. **Hospital/School Orders:** Institutional orders should be emailed to orders@salvin.com. Terms for hospital/school orders are net 30 with 1.5% per month late fee charged on past due accounts. Customers are responsible for legal and collection fees if required. **Trademarks:** "Salvin", "Salvin Dental Specialties" Logo, "Everything For Your Implant Practice But The Implants", "Salvin Regenerative" Logo, AccessTorq®, AlloSculpt®, Salvin AmnioExcite®, AtrauLux®, AutoTome®, BariOpaque®, CytoSurg®, EasyReach®, OraCOTE®, OraFoam®, OraMEM®, OraMEM-Sustained®, OraMEM+®, OraPLUG®, OraTAPE®, Orastruct®, OsseoFiber®, OsseoFiber-Plus®, PDL-Evator®, Rapid Roll®, Renovix-Plus®, RotoBrush-Titanium®, Salvin AutoTorq®, Salvin Power Point®, SalvinOss®, SalvinOss-B, SalvinOss® Collagen, SalvinTorq®, SurgiPoint® and Xymphony® are registered U.S. and international trademarks of Salvin Dental Specialties, LLC.

Salvin AlloSpark-GF™, Salvin Amplifi™, Enlighten-Pro3™ LED System, Salvin Black Kryptonite™, and Salvin Mesh-Ti™ are trademarks of Salvin Dental Specialties, LLC. Other product brand names may be trademarks or registered trademarks of their respective owners and are mentioned for reference purposes only.

CONTENTS

6

ALLOSPARK-GF
SUPPLEMENTAL ALLOGRAFT
DERIVED PROTEINS

7

ORAGEN
COMPLETE BONE SOLUTION

8

OSSEOFIBER-PLUS
DEMIN CORTICAL FIBERS +
MINERALIZED CANCELLOUS CHIPS

OSSEOFIBER
DEMIN CORTICAL FIBERS

9

PHARMACEUTICALS
AT YOUR FINGERTIPS

10-11

AMNIOEXCITE AMNION / CHORION FULL
THICKNESS PLACENTAL MEMBRANE

12

ALLOSCULPT-3D PLUS
DBM PUTTY WITH
CORTICAL & CANCELLOUS CHIPS

13

SALVINOSS-B
BOVINE XENOGRAFT BONE GRAFT
MATERIAL

14-15

SALVINOSS
XENOGRAFT PARTICULATE

16-17

SALVINOSS COLLAGEN
XENOGRAFT + COLLAGEN PLUG

18-21

LIFENET HEALTH ALLOGRAFT BONE

22-23

ORACELL
DECELLULARIZED DERMIS ALLOGRAFT

24

SALVIN AMPLIFI
LONG TERM MEMBRANE

25

RENOVIX-PLUS
RESORBABLE EXTRA CELLULAR
MATRIX

26

ORAMEM+
RESORBABLE COLLAGEN MEMBRANE

27

XYMPHONY
DRAPEABLE PORCINE PERICARDIUM

28

CYTOSURG
NON-RESORBABLE PTFE MEMBRANE

29

ORAPLUG
RESORBABLE COLLAGEN PLUG

FULFILL
80% CALCIUM PHOSPHATE GRANULES
COMBINED WITH
20% BOVINE COLLAGEN

30

ORATAPE
RESORBABLE COLLAGEN
WOUND DRESSING

ORACOTE
RESORBABLE COLLAGEN
WOUND DRESSING

31

ORASTRUCT
80/20 TCP/HA GRAFT MATERIAL

SALVIN MESH-TI TITANIUM MESH

MESH & WIRE CUTTER

32

ACCELL DBM PUTTY IN A SYRINGE

33

DYNABLAST
DEMIN CORTICAL + MIN CANCELLOUS
CHIPS PASTE OR PUTTY

DYNAGRAFT
DEMIN CORTICAL PASTE OR PUTTY

34-39

SUTURES & WOUND CLOSING

Sutures	34-38
Periacryl Cyanoacrylate	39
Featherlite Suture Removal Scissors	39

40-41

BLEEDING CONTROL

OraFOAM™ Hemostatic Sponges	40
BloodSTOP® Hemostatic Gauze	40
Hemostats	40
Bone Wax	40
HemCon® Dental DressingPro	41
Portable Electrocautery	41

42-45

VITAL SIGN MONITORS

Pulse Ox / BP Monitors	43
Pulse Ox / BP / ECG / CO ₂ Monitors	44-45

46-48

O/R FURNITURE

Implant / Surgical Cart	46
Mayo & IV Stands	47
Table Covers & Mayo Stand Covers	47
Surgical Tables	48
Sedation Chair Pad	48

49-55

DEFIBRILLATOR & SEDATION

Wireless Sedation Stethoscope	49
Defibrillator	49
CO ₂ Nasal Hoods	50
ECG Patches & Wrist Clips	50
CPR Mask & Emergency Airways	50
Nasal Cannula & Oxygen Adapter	51
Anesthesia Instrumentation	51
Emergency Kits	51
IV Catheters & Admin Sets	52-54
Aspirating Anesthetic Syringes	55

56-66

INFECTION CONTROL / PPE

Sterile Drapes, Gowns & Covers	56-60
N95 Masks	61
Towels, Caps & Shoe Covers	61
Patient Skin Prep	61
Cold Packs	62
Scrub Brushes	62
Instrument Care & Cassettes	63-65
Surgical Gloves	66

67-69

CENTRIFUGES

Blood Collection	67
BIO-PRF Centrifuge	68
Variable Speed Centrifuge	69

70-74

PIEZO ULTRASONIC SURGERY & ULTRASONIC SCALERS

Piezotome Ultrasonic Surgery Units	70-71
Piezotome Tips & Procedure Kits	72-73
Newtron Ultrasonic Scaler	74
Air-N-Go Air Polishing	74

75-77

PERI-IMPLANTITIS MANAGEMENT

RotoBrush-Titanium® Brushes	75
Peri-Implantitis Instruments	76-77

78-92

MOTORS, HANDPIECES & DRILLS

Implant Motors	78
Implant Handpieces	79
Irrigation Tubing	80
Bone Saws	81
Drill Stops & Implant Drills	82
Drills, Burs, Trepines & Bur Blocks	83-89
AutoTome® Automatic Osteotome Kit	90
AutoTome® Automatic Extraction Kit	91
Rear-Exhaust High Speed Handpiece	92

93

SALVIN IMPLANT RESCUE KIT

Broken Screw & Failed Implant Removal	93
---	----

94-100

SUCTION & IRRIGATION

Sterile Water & Saline	94
Irrigation	94-95
Suction	96-99
IrriguIDE Adapter	100

101-184

SURGICAL INSTRUMENTS

Mallets & Chisels	101
Atraumatic Extraction	102-119
Bone Curettes	120
Zygomatic Instruments	121
Bone Grafting	122-131
Tenting Screw Kit	132
Fixation Screw Kit	133

Membrane Tack Kits	134-135
Membrane Instruments	136
Sinus Lifts	137-143
Bone Spreading	144-145
Ridge Splitting	146-147
Rongeurs	148-149
Retractors	150-153
Periosteal Elevators	154-156
Scissors	157-159
Needle Holders	160-163
Micro Surgery & Tunneling	164-175
Perio Knives & Surgical Files	176
Perio Curettes & Bone Files	177
Tissue Forceps	178-181
Tissue Punches & Nippers	182-183
Scalpel Handles	184

185

HEADLAMPS

Enlighten-Pro3™ LED Headlamp	185
--	-----

186-190

IMPLANT & COURSE MODELS

Implant Models	186-189
Bone Loss & Course Models	190

191

IMPLANT STABILITY MONITOR

Penguin RFA Implant Stability Monitor	191
---	-----

192-198

DIAGNOSTICS

Spotter Smart Implant Detector	192
Drill Guides	193
Calipers & Gauges	194-195
Photo Mirrors & Contrasters	196-197
Diagnostic Probes	198
Operative Mouth Mirrors	198

199-207

RESTORATIVE & HYGIENE

Implant Buddy Driver Tip Kits	199
EasyReach® Right Angle Wrench	200
AutoTorq® Electric Torque Driver	200
SalvinTorq® Torque Wrench	201
Access-Torq®-Plus RA Torque	202
Abutment Driver Tips	202-203
Abutment Holders	204
Implant Scalers & Curettes	204-205
Crown Grippers	206
Crown Splitters	206
Improv® Implant Cement	207


Complete Index On Pages 208-209



NEW!

Enhance The Regenerative Potential Of Your Choice Of Graft Material!

- With A Natural Blend Of Allograft Derived Growth Factors, AlloSpark-GF Is A Perfect Complement To Any Commercially Available Graft Material
- Contains Nearly 1,000 Unique Proteins, Peptides & Growth Factors Harnessed Through Proprietary Processing Methods
- Provides A Powerhouse Of Osteoinductive, Angiogenic, Proliferative & Chemotactic Properties
- Every Lot Is Tested For Protein Content To Provide Quantified Levels Of Bone Healing Contribution
- Store At Ambient Temperature – 11.0-30.0°C (51.8-86.0°F)



SALVIN® AlloSpark-GF™ Supplemental Allograft Derived Proteins

Size	Product #	Fee	7 or More
0.5 cc	#ALLOSPARK-GF-0.5	240.00	200.00
1.0 cc	#ALLOSPARK-GF-1.0	420.00	350.00
2.5 cc	#ALLOSPARK-GF-2.5	900.00	750.00

OraGEN®

Viable Dental Bone Allograft

The New Gold Standard in Bone Grafting



ORAGEN® IS A COMPLETE BONE GRAFT

Successful Bone Remodeling Requires The Following Three Properties:



OSTEOCONDUCTIVE
Scaffold



OSTEOINDUCTIVE
Signals



OSTEOGENIC
Cells



What Is OraGen®?

OraGen® has three components:

- **Osteogenic:** lineage-committed bone cells that directly participate in bone formation
- **Osteoconductive:** corticocancellous chips that provide a natural scaffold for cells to attach, migrate, and proliferate
- **Osteoinductive:** demineralized bone that exposes natural growth factors within the matrix

How Does OraGen® Work?

- OraGen® contains lineage-committed bone cells and demineralized bone matrix, which can begin producing bone matrix shortly after implantation.¹
- Histological analysis of OraGen® suggests that new bone formation started directly after implantation and continued at a rapid pace.¹

Storage And Shipping

- **Storage and shipping:**
OraGen® is stored and shipped frozen
- **Specialized packaging:**
the product is packaged in a thermal shipper designed to maintain low temperatures during transport
- **Fast delivery:**
to preserve cell viability, OraGen® is typically shipped with next-day delivery options
- **Rapid thaw time:**
specialized packaging allows for rapid thaw time – less than 5 minutes!

When storing OraGen® long-term: store at -70°C or colder. Do not store in liquid phase of liquid nitrogen (LN₂), 12 month shelf life.

Call us for information on special freezers

OraGEN® Viable Dental Bone Allograft

Size	Product #	Fee	4 or More
0.7 cc	#ORAGEN-0.7	650.00	550.00
1.3 cc	#ORAGEN-1.3	795.00	675.00



1. Gianulis E, Wetzell B, Scheunemann D, et al. Characterization of an advanced viable bone allograft with preserved native boneforming cells. Cell and tissue banking. 2022; (E-pub ahead of print). doi: 10.1007/s10561-022-10044-2.

**100% Bone – No Carrier – Handles Like Putty When Hydrated
Facilitates Cellular Proliferation & Cellular Attachment ⁽¹⁾**

**Demineralized Cortical Fiber
PLUS Mineralized Cancellous Chips**

- Demin Cortical Bone Has The Potential To Stimulate New Bone Growth
- Mineralized Cancellous Chips Provide Scaffold For Bone Cell Migration
- For GBR Procedures, Ridge Augmentation, Sinus Lifts & Bony Defects
- Each Lot Is Tested In Vivo For Osteoinductive Potential



Moldable & Easily Shaped
To Pack Into Any Defect



**SALVIN® OsseoFiber-PLUS® Demineralized Cortical Fibers
PLUS Mineralized Cancellous Chips**



Size	Product #	Fee	7 or More
0.5 cc	#OFIBER-PLUS-0.5CC	163.00	144.70
1.0 cc	#OFIBER-PLUS-1.0CC	279.00	250.80
2.5 cc	#OFIBER-PLUS-2.5CC	525.00	475.50


No Refrigeration Required – Store At Room Temperature
Typical Shelf Life: 3 Years

Case Study - Guided Bone Regeneration Using OsseoFiber-Plus®


Case & Photos: Dr. Paul Rosen



Post Surgical Extraction
– Osseous Defect



OsseoFiber-Plus® Has
Putty-Like Handling



Successful GBR
6 Month Post-Op



Implant Placement
Prior To Restoration



“The innovation of making DFDBA in a fiber form has been a tremendous enhancement to this graft. The increase in surface area makes this osseoinductive clot stabilizer a must for one’s regenerative armamentarium. Importantly, this fiber form of DFDBA transforms it into being much more user friendly.”

–Dr. Paul Rosen
Periodontist, Yardley, PA & New York, NY

Demineralized Cortical Fiber

SALVIN® OsseoFiber® Demineralized Cortical Fibers



Size	Product #	Fee	7 or More
0.5 cc	#OFIBER-0.5CC	157.00	139.70
1.0 cc	#OFIBER-1.0CC	267.00	238.80
2.5 cc	#OFIBER-2.5CC	522.00	473.60

No Refrigeration Required – Store At Room Temperature
Typical Shelf Life: 3 Years



1. Rodriguez RU, Kemper N, Breathwaite E, et al. Demineralized bone matrix fibers formable as general and custom 3D printed mold-based implants for promoting bone regeneration. *Biofabrication*. 2016;8(3):035007



For More Info
SCAN HERE



NEW!

PHARMACEUTICALS AT YOUR FINGERTIPS

*We Now Offer A Full Line Of Pharmaceuticals.
Scan The QR Code Or Call Us For More Information.
We'll Take Great Care Of You!*




Orders for pharmaceuticals and the applicable documentation are received, processed, fulfilled, and distributed by PuraGraft and/or Medical Purchasing Solutions, LLC.



Full-Thickness Placental Membrane

- 4X Thicker Than An Amnion Or Chorion Membrane Making It Much Easier To Handle
- Does Not Fold In Upon Itself When Rehydrated

Structure of Salvin AmnioExcite® Full-Thickness Placental Membrane



*AmnioExcite® Includes Both Amnion And Chorion PLUS The Trophoblast Layer
The Trophoblast Layer Contains More Than 50% Of
The Biologically Relevant Proteins In The Placental Membrane (1)*

SALVIN AMNIO EXCITE® Amnio / Chorion Membrane

Size	Product #	Fee	7 or More
10mm x 10mm	#AMNIO-MEM-10X10	99.00	87.00
15mm x 20mm	#AMNIO-MEM-15X20	189.00	174.00
20mm x 30mm	#AMNIO-MEM-20X30	289.00	266.00




1. LifeNet Health®. TR-004-2020 Characterization of the Amnion, Chorion, and Trophoblast Layers of Decellularized and Freeze-Dried Placental Membrane. 2020.




AmnioExcite® Case Studies

Case Study #1 - Socket Preservation Using AmnioExcite® Full-Thickness Placental Membrane


Case & Photos: Dr. Nick Shumaker



Surgical Extraction



*LifeNet Health®
Cortical / Cancellous
Allograft*



*AmnioExcite®
Barrier Membrane*



*3 Week Post-Op –
Excellent Tissue Healing*

Case Study #2- Cemental Tear/Osseous Defect Repair Using AmnioExcite®

Case & Photos: Dr. Nick Shumaker



*Cemental Tear /
Osseous Defect*



*LifeNet Health®
Cortical / Cancellous
Allograft*



*AmnioExcite®
Barrier Membrane*



*6 Month - No Recession
Excellent Aesthetic Result*



"As a user of amnion membranes for over 15 years, I have found AmnioExcite® from Salvin to be a great addition to my practice. AmnioExcite® is much easier to handle and place, but still adapts to adjacent structures such as bone, tooth roots or dental implants. It has the perfect balance of adaptation and handling that has been missing in amnion membranes, and I'm very pleased with the results."

–Dr. Nick Shumaker

Periodontist, Fort Collins, CO

NEW!



Moldable & Easily Shaped To Pack Into Any Defect

DBM Putty With Cortical & Cancellous Chips

- 1-2mm Cortical & Cancellous Chips Provide An Osteoconductive Scaffold For Bone Cell Migration
- Inert Gelatin Carrier Is Easy To Shape & Mold When Hydrated With Blood Or Saline
- The Graft Sets At Body Temperature And Resists Migration Once Implanted
- Terminally Sterilized To Achieve SAL 10⁻⁶ Using Low Temperature, Low Dose Gamma Irradiation
- Each Lot Is Tested In Vivo For Osteoinductive* Potential

SALVIN® AlloSculpt-3D PLUS DBM Putty With Cortical & Cancellous Chips

Size	Product #	Fee	7 or More
1.0 cc	#ALLO-3D-PLUS-1.0	235.00	200.00
2.0 cc	#ALLO-3D-PLUS-2.0	360.00	320.00

No Refrigeration Required - Store At Room Temperature
Product Shelf Life - 2.5 Years



New Manufacturing Efficiency Creates Significantly Reduced Cost!


*These implants were evaluated in a human clinical study and were shown to induce bone formation. Each lot is tested using the athymic nude rat assay to verify osteoinductive potential.

NEW!



Bovine Xenograft Bone Graft Material

- 0.25 - 1.0mm Bovine Granules Facilitate Osteoconduction And New Bone Formation
- Hydrophilic Properties For Rapid Absorption Of Patient's Blood
- Micro & Macro Pores Create Environment For Cell Attachment And New Bone Growth
- Can Be Mixed With Allograft Or Autologous Bone If Desired
- Indicated For Ridge Augmentation, Filling Periodontal Defects, Socket & Ridge Preservation, Sinus Lifts



Excellent Handling Characteristics When Hydrated



Filter On End Of Syringe Keeps Material In Place



Draw Saline Into Syringe For Quick & Easy Delivery

SALVIN OSS-B Bovine Xenograft Bone Graft Material – Vial



Particle Size	Size	Product #	Price	7 or More	14 or More
250-1000 microns	0.54 CC (0.25G)	#SALVINOSS-B-0.50	98.00	87.60	76.00
	1.08 CC (0.5G)	#SALVINOSS-B-1.00	147.00	132.40	118.70
	2.16 CC (1.0G)	#SALVINOSS-B-2.00	254.00	222.00	200.00
	4.32 CC (2.0G)	#SALVINOSS-B-4.00	407.00	386.00	355.00


SALVIN OSS-B Bovine Xenograft Bone Graft Material – Syringe



Particle Size	Size	Product #	Price	7 or More	14 or More
250-1000 microns	0.25 CC	#SALVINOSS-B-SYR-0.25	71.00	66.00	59.50
	0.50 CC	#SALVINOSS-B-SYR-0.50	101.00	90.00	79.50
	1.00 CC	#SALVINOSS-B-SYR-1.00	152.00	132.40	118.70

SalvinOss® Mineralized Cancellous Xenograft

- Indicated For Sinus Lifts, Ridge Augmentation & Guided Bone Regeneration
- Maintains Space & Facilitates Osteoconduction
- 0.25 - 1.0mm High Porosity Particles Become An Integral Part Of Newly Formed Bone
- Radiopaque



Exceptional Handling Characteristics

SALVIN(OSS) Xenograft Bone Graft Material



Particle Size	Size	Product #	Price	7 or More	14 or More
250-1000 microns	0.50 CC (0.25G)	#SALVINOSS-0.50	98.00	87.60	76.00
	1.00 CC (0.5G)	#SALVINOSS-1.00	147.00	132.40	118.70
	2.50 CC (1.25G)	#SALVINOSS-2.50	269.00	244.50	220.60


Download Our White Paper At:
www.salvin.com/Downloads/Guides-Studies.html



SalvinOss® Case Studies

Case Study - Socket Preservation With SalvinOss®

Case & Photos: Dr. Michele Holzinger



Pre-Op



Post-Extraction –
4 Months Healing



Post-Implant –
3 Months Healing




"I trust Salvin Regenerative because they help me produce results that are consistent and reliable for my patients. When I re-enter a site grafted with SalvinOss®, I know exactly what to expect because the results are **predictable**."


–Dr. Michele Holzinger
Periodontist, Middletown, CT

Case Study - Ridge Augmentation Using SalvinOss® Xenograft Bone Material

Case & Photos: Dr. Mitchel Godat



Severely Pneumatized Sinus Cavity



12 Month Post-Op Radiograph



"I have been using SalvinOss® for many years to increase bone volume for space preservation prior to implant placement. SalvinOss® is excellent for sinus lifts, ridge augmentation, and other GBR / GTR procedures with very predictable long-term results."

–Dr. Mitchel Godat
Periodontist, Memphis, TN

Case Study - Guided Tissue Regeneration Using SalvinOss® Xenograft Bone Material

Case & Photos: Dr. Grant King



9-10 mm
Periodontal Defect



8 Month Post-Op –
Significant Pocket Fill



"SalvinOss® handles extremely well, especially when I mix it with blood, making it easy for me use for GTR and GBR. It is a slower resorbing material, holding space and volume more than using allograft alone. Most importantly, SalvinOss® gives me predictable results every time."


–Dr. Grant King
Periodontist, Memphis, TN

SalvinOss® Collagen On Next Page →

SalvinOss® Collagen Xenograft Bone + Bovine Collagen


- 90% SalvinOss® Xenograft Bone Graft Material PLUS 10% Bovine Collagen
- For Socket Preservation, Sinus Lifts, And Periodontal / Furcation Defects
- Plug Shape Helps Make Socket Preservation Fast & Simple
- Retains Shape After Hydrating





Ideal For:


- Anterior Socket Grafting
- Periodontal Defects
- Osseous Defects Caused By Peri-Implantitis



Ideal For:

- Posterior Socket Grafting
- Ridge Augmentation
- Sinus Lifts

Plugs Shown Actual Size



SALVIN(OSS) Collagen Xenograft Bone PLUS Bovine Collagen



Particle Size	Size	Product #	Price	7 or More	14 or More
250-1000 microns	7mm (Diam) x 10mm	#SALVINOSS-COLL-100MG	116.00	108.50	101.50
	10mm (Diam) x 10mm	#SALVINOSS-COLL-250MG	204.00	195.00	184.00

Download Our White Paper At:
www.salvin.com/Downloads/Guides-Studies.html



Scan For



Scan For



Scan For



SalvinOss® Collagen Case Studies

Case Study #1 - Socket Preservation Using SalvinOss® Collagen

Case & Photos: Dr. Tyler Borg



Plug Shape For Easy Socket Preservation And Crestal Sinus Lift

Case Study #2 - Furcation Defect Using SalvinOss® Collagen

Case & Photos: Dr. Tyler Borg



SalvinOss® Collagen Helps When Graft Containment Is Problematic



"SalvinOss® Collagen has become a staple of my practice for regeneration as well as augmentation. The ease of handling coupled with the volumetric stability produces predictability in grafting situations around both teeth and implants alike."

-Dr. Tyler Borg
Periodontist, Northglenn, CO



Over 11 Million Bio-Implants Have Been Processed Using LifeNet Health's® Allowash Technology, With No Disease Transmission



LifeNet Health® Is The World's Largest Full Service Tissue Processor / Organ Procurement Organization

Of Total Donors Screened For Transplantation, Only 3% Are Accepted By LifeNet Health®

All Of LifeNet Health's® Donor Tissue Is Recovered In The United States

Specific Screening Criteria*	LifeNet Health® <small>Saving Lives. Restoring Health. Giving Hope.</small>	AATB <small>(American Association of Tissue Banks)</small>	FDA <small>(U.S. Food & Drug Administration)</small>
Encephalitis (clinically active)	●	●	●
Gonorrhea (within the past year)	●	●	●
Hepatitis B	●	●	●
Hepatitis C	●	●	●
HIV 1/2	●	●	●
Malaria	●	●	●
Meningitis (clinically active)	●	●	●
Sepsis	●	●	●
Syphilis	●	●	●
Transmission spongiform encephalopathy (TSE)	●	●	●
Vaccinia	●	●	●
West Nile Virus (WNV, 4 month exclusion from the time of onset)	●	●	●
Guillain-Barre syndrome (clinically active)	●	●	●
Multiple sclerosis	●		●
Leprosy (Hansen's disease)	●	●	
Clinically significant metabolic bone disease	●	●	
Polyarteritis nodosa	●	●	
Rabies	●	●	
Rheumatoid arthritis (musculoskeletal deferral)	●	●	
Sarcoidosis	●	●	
Systemic lupus erythematosus	●	●	
Systemic mycosis	●	●	
Tuberculosis (clinically active or treated in the past 5 years)	●	●	
Ankylosing spondylitis	●	●	
Autoimmune vasculitis if type is not known	●	●	
Cancer	●	●	
Chagas disease	●	●	
Epstein Barr virus (clinically active)	●	●	
Clostridium difficile infection (defer if occurred during present hospitalization)	●	●	
Endocarditis (clinically active)	●	●	
Herpes (any evidence of active herpes infection is a deferral)	●	●	
Illicit drug use (see Illegal drugs)	●	●	
Peritonitis (clinically active)	●	●	
Rheumatic fever (defer all tissues if clinically active, defer heart valves for any history)	●	●	
Poliomyelitis (clinically active or untreated)	●	●	
Mixed connective tissue disease	●		
Cold agglutinin disease	●		
Pyelonephritis (clinically active)	●		
Reactive arthritis (must be asymptomatic for 6 months for musculoskeletal deferral)	●		
Wegener's granulomatosis	●		

*These Screening Criteria Are Subject To Change

LifeNet Health®, The LifeNet Health® Logo, OrACELL®, OraGRAFT® & Preservon® Are Registered Trademarks Of LifeNet Health®.

LifeNet Health® OraGRAFT®

Mineralized 50/50 Cancellous/Cortical Mix MOST POPULAR



Particle Size	Size	Product #	Fee	7 or More
250-1000 microns	0.25 cc	#CAN-CORT-.25	67.00	59.00
	0.50 cc	#CAN-CORT-.50	98.00	86.60
	1.00 cc	#CAN-CORT-1.00	132.00	117.80
	2.00 cc	#CAN-CORT-2.00	190.00	170.40

Mineralized Cancellous MOST POPULAR



Particle Size	Size	Product #	Fee	7 or More
250-1000 microns	0.25 cc	#MIN-CAN-.25	88.00	77.80
	0.50 cc	#MIN-CAN-.50	101.00	88.70
	1.00 cc	#MIN-CAN-1.00	137.00	121.40
	2.00 cc	#MIN-CAN-2.00	211.00	188.60
1000-2000 microns	0.50 cc	#MIN-CAN-L-.50	101.00	88.70
	1.00 cc	#MIN-CAN-L-1.00	137.00	121.40
	2.00 cc	#MIN-CAN-L-2.00	211.00	188.60

Mineralized Cortical



Particle Size	Size	Product #	Fee	7 or More
250-1000 microns	0.25 cc	#MIN-CORT-.25	48.00	41.30
	0.50 cc	#MIN-CORT-.50	67.00	58.40
	0.70 cc	#MIN-CORT-.75	83.00	72.70
	1.20 cc	#MIN-CORT-1.25	97.00	85.50
	2.50 cc	#MIN-CORT-2.50	143.00	126.80
250-710 microns	0.50 cc	#MIN-CORT-SM-.50	67.00	58.40
	0.70 cc	#MIN-CORT-SM-.75	83.00	72.70
	1.20 cc	#MIN-CORT-SM-1.25	97.00	85.50
	2.50 cc	#MIN-CORT-SM-2.50	143.00	126.80



LifeNet Health® OraGRAFT®

70/30 Mineralized/Demineralized Cortical Mix **MOST POPULAR**

Particle Size	Size	Product #	Fee	7 or More
250-1000 microns	0.50 cc	#MD-7030-.50	92.00	82.60
	1.00 cc	#MD-7030-1.00	134.00	120.50
	2.00 cc	#MD-7030-2.00	189.00	170.70

Demineralized Cortical



Particle Size	Size	Product #	Fee
250-1000 microns	0.25 cc	#DEMIN-CORT-.25	55.00
	0.50 cc	#DEMIN-CORT-.50	73.00
	0.70 cc	#DEMIN-CORT-.75	92.00
	1.20 cc	#DEMIN-CORT-1.25	108.00
	2.50 cc	#DEMIN-CORT-2.50	159.00

Graft Site Suggested Volume Guide

Periodontal Defects		Extraction Sites		Sinus Grafts	
Central Incisor	0.5 cc	Central Incisor	0.7-1.0 cc	Each Sinus	2.5-5.0 cc
Lateral Incisor	0.5 cc	Lateral Incisor	0.5 cc		
Canine	0.5 cc	Canine	0.7-1.0 cc		
1st & 2nd Premolar	0.5 cc	1st & 2nd Premolar	0.7-1.0 cc		
1st & 2nd Molar	0.7-1.0 cc	1st & 2nd Molar	1.2-2.0 cc		
3rd Molar	0.7-1.0 cc	3rd Molar	1.2-2.0 cc		

For optimal safety, all LifeNet Health® allografts distributed by Salvin Dental are gamma irradiated to achieve a Sterility Assurance Level (SAL) 10^{-6} without compromising the inherent biomechanical or biochemical properties of the graft. Typical Shelf Life: 4 Years. Fees are subject to change. Please check current fees at time of placing order.
LifeNet Health® Canadian Registration: CTO #100038

LifeNet Health® OraGRAFT®

Pre-Hydrated Cortical Plate

Ready To Use Out Of The Package

- Preserved Using LifeNet Health’s® Proprietary Preservon® Technology To Maintain A Fully Hydrated State At Ambient Temperature
- Preservon-Treated Grafts Have Been Shown To Have Greater Flexibility Than Freeze-Dried Grafts ⁽¹⁾
- Ideal For “Shell” Technique – Secure Cortical Plate To Host Bone With Fixation Screws To Create A Biologic Container



Plate	Size	Product #	Fee
Cortical Plate	30 mm x 15 mm x 1 mm	#CORT-PLATE-30X15X1	520.00

No Refrigeration Required - Store At Room Temperature
Typical Shelf Life: 5 Years

Pre-Hydrated Ilium Block

Ready To Use Out Of The Package

- Preserved Using LifeNet Health’s® Proprietary Preservon® Technology To Maintain A Fully Hydrated State At Ambient Temperature



Block Size	Size	Product #	Fee
Ilium Block - Small	15 mm x 20 mm	#ILIUM-PRES-SMALL	599.00
Ilium Block - Medium	15 mm x 30 mm	#ILIUM-PRES-MEDIUM	699.00
Ilium Block - Large	20 mm x 40 mm	#ILIUM-PRES-LARGE	899.00

Pre-Hydrated Unicortical Block **MOST POPULAR**

Ready To Use Out Of The Package


- Preserved Using LifeNet Health’s® Proprietary Preservon® Technology To Maintain A Fully Hydrated State At Ambient Temperature



Block Size	Size	Product #	Fee
Unicortical Block - Small	10 mm x 15 mm x 8 mm	#ILIUM-UNICORTICAL-SMALL	339.00
Unicortical Block - Large	10 mm x 25 mm x 8 mm	#ILIUM-UNICORTICAL-LARGE	411.00

Case Study - GBR Using LifeNet Health® Unicortical Block & Renovix-PLUS® Membrane


Case & Photos: Dr. Shaun Young



Large Osseous Defect



LifeNet Health®
Unicortical Block



Renovix-Plus® –
Barrier Membrane



Implant Placed In
Regenerated Bone



“Salvin continues to raise the bar with their regenerative products. The allogeneic blocks have elevated my abilities and given me additional options that I can depend on for ridge augmentation. The Renovix-PLUS® collagen membrane is one of the best on the market due to its versatility as a stiffer membrane, but also very moldable – it’s really the best of both worlds.”

–Dr. Shaun Young
Oral Surgeon, Tampa, FL

1. LifeNet Health. TR 0216 Three-Point Bend Testing of Glycerol-Treated, Freeze-Dried, and Frozen Bone



Decellularized Dermis Allograft

- Excellent For GTR/GBR & Oral Soft Tissue Correction
- Terminally Sterilized For Increased Patient Safety
- Handles Like Autogenous Tissue



**Prehydrated & Ready To Use
Right Out Of The Package**



**No Refrigeration Required –
Store At Ambient Temperature**

OrACELL® Decellularized Dermis Allograft

Size	Thickness	Product #	Fee	7 or More
1.0cm x 1.0cm	0.76 - 1.25mm	#DERMIS-1X1	97.00 EACH	85.70 EACH
1.5cm x 2.0cm	0.76 - 1.25mm	#DERMIS-1.5X2	187.00 EACH	167.50 EACH
1.0cm x 4.0cm	0.76 - 1.25mm	#DERMIS-1X4	214.00 EACH	191.40 EACH
2.0cm x 4.0cm	0.76 - 1.25mm	#DERMIS-2X4	235.00 EACH	210.80 EACH
1.0cm x 1.0cm	1.26 - 1.75mm	#DERMIS-THICK-1X1	97.00 EACH	85.70 EACH
1.5cm x 2.0cm	1.26 - 1.75mm	#DERMIS-THICK-1.5X2	187.00 EACH	167.50 EACH
1.0cm x 4.0cm	1.26 - 1.75mm	#DERMIS-THICK-1X4	214.00 EACH	191.40 EACH
2.0cm x 4.0cm	1.26 - 1.75mm	#DERMIS-THICK-2X4	235.00 EACH	210.80 EACH



OrACELL® Decellularized Dermis Case Studies

Case Study - OrACELL® Decellularized Dermis

Case & Photos: Dr. Matt Ficca



*Pre-Op
Tooth #6 & #7*



6 Year Follow-Up



*Pre-Op
Tooth #10 & #11*



6 Year Follow-Up



“OrACELL® provides an excellent alternative to autogenous tissue for tissue graft procedures. The pre-hydrated samples handle comparative to autogenous or even better. Over the years I have looked for an allograft with consistent quality and thickness, ease of handling, and predictable healing. I want to offer my patients the highest quality allograft on the market. I am confident I've got it covered with OrACELL®.”

–Dr. Matt Ficca
Periodontist, Charlotte, NC

Case Study - OrACELL® Decellularized Dermis

Case & Photos: Dr. Stephen C. Wallace



Pre-Op



OrACELL® Dermis Placed



Suture Closing



3 Month Post-Op



“Based on the results recorded at 3 months, there was significant gain in root coverage, gain in clinical attachment and increase in width of keratinized gingiva for all recession sites treated with OrACELL® Decellularized Dermis Allograft. I have been using OrACELL® for decades with predictable clinical results.”

–Dr. Stephen C. Wallace
Periodontist, Wilmington, NC



Membranes Shown Actual Size



15mm x 20mm



20mm x 30mm



30mm x 40mm

Resorbs In 6-9 Months For Long-Term Barrier Functionability

- Collagen-Fiber Matrix With Intact Intrinsic Crosslinks For Mechanical Strength
- <0.01 Micron Pore Size Allows Nutrients & Growth Factors While Assuring Cell Barrier Function
- Chemically Engineered To Control Stability & Minimize Inflammatory Response
- Significant Tensile Strength For Ease Of Tacking & Suturing



Long Term Membrane


Resorbs In 6 To 9 Months

Size	Product #	Price	7 or More	14 or More
15mm x 20mm	#AMPLIFI-15X20	173.00	153.70	141.60
20mm x 30mm	#AMPLIFI-20X30	216.00	192.40	178.20
30mm x 40mm	#AMPLIFI-30X40	314.00	286.80	266.70



Case Study - Amplifi™ Long-Term Membrane

Case & Photos: Dr. Bradley Ross



Surgical Extraction



Bone Graft Material In Defect



Amplifi™ Sutured In Place




New Bone Formation – Ready For Implant



“The Amplifi™ membrane from Salvin has become a go to biomaterial for me for ridge augmentation procedures. It is extremely easy to manipulate and use, is very resilient, and most importantly it has a long resorption time. The results I have seen with the Amplifi™ membrane have been spectacular.”

—Dr. Bradley Ross
Periodontist, Miami, FL

*Data On File



Our Most Versatile Resorbable Barrier Membrane

- For Bone Regeneration & Soft Tissue Esthetics
- Non-Cross Linked Extracellular Matrix Perfectly Conforms To The Defect When Hydrated
- Contains Type I, III, & IV Collagen
- Excellent Tensile Strength & Tear Resistance
- Trim & Place Either Wet Or Dry
- Resorbs Within 6 Months

Case Study - Renovix-Plus® Guided Healing Collagen

Case & Photos: Dr. Lauralee Nygaard



Aesthetic Zone Surgical Extraction



Renovix-Plus® Over Graft Material



Soft Tissue Aesthetics Preserved



Implant Placed



"Renovix-Plus® has become the membrane I use most for predictable ridge preservation. I use it for ease of handling and preservation of soft tissue architecture. And the bonus of working with team Salvin!"

-Dr. Lauralee Nygaard
Periodontist, Spokane Valley, WA

SALVIN Renovix-Plus® Guided Healing Collagen Membrane

Size	Product #	Price	7 or More	14 or More
15mm x 25mm	#RENOVIX-PLUS-15X25	181.00	162.40	150.60
20mm x 30mm	#RENOVIX-PLUS-20X30	236.00	214.50	199.40
30mm x 40mm	#RENOVIX-PLUS-30X40	307.00	273.70	259.50





NEW!


- *Handles Similarly To Both OraMEM® And OraMEM-Sustained®, With A Resorption Profile That Is In Between.*
- *By Choosing OraMEM®+, You May Be Able To Reduce The Amount Of Inventory You Need On The Shelf.*

Still Want Something That Resorbs A Little Longer? Take A Look At Our Amplifi Membrane On Page 24



Resorbable Collagen Membrane

- *For Use In Guided Tissue Regeneration Procedures*
- *Excellent Handling Properties*
- *Membrane Is Moderately Drapable*
- *Soft Surface Texture*
- *Provides Wound Stabilization*
- *Highly Purified Collagen That Is Biocompatible*
- *Resorbed In 4 To 6 Months*



Scanning Electron Micrograph Of The Dry OraMEM®+ (Cross Section) At 50x Magnification

Membranes Shown Actual Size



15mm x 20mm



20mm x 30mm



30mm x 40mm

OraMEM+

Resorbs In 4 To 6 Months


Size	Product #	Price	7 or More	14 or More
15mm x 20mm	#ORAMEM-PLUS-15X20	169.00	151.60	140.40
20mm x 30mm	#ORAMEM-PLUS-20X30	213.00	191.30	178.00
30mm x 40mm	#ORAMEM-PLUS-30X40	307.00	279.70	260.50





Our Most Drapeable, “No Shape Memory” Barrier Membrane

- Porcine Pericardium Maintains Barrier Function Up To 33 Weeks
- Significant Tensile Strength - Suture Or Tack With Confidence




SALVIN X YMPHONY Contour Adapting Resorbable Membrane

Size	Product #	Price	7 or More	14 or More
15mm x 25mm	#XYMPHONY-15X25	164.00	145.40	133.70
20mm x 30mm	#XYMPHONY-20X30	217.00	193.70	179.40
30mm x 40mm	#XYMPHONY-30X40	296.00	265.70	247.80


1. In vivo animal study data on file.
2. Data on file.





Non-Resorbable PTFE Membrane

- For Socket Preservation When Primary Closure Cannot Be Achieved
- Sufficiently Firm To Protect The Defect – Supple Enough To Drape Over Surrounding Tissue
- Nominal Thickness Of 250 Microns
- Needs To Be Removed In 2-3 Weeks



Space-Making Barrier In The Treatment Of Periodontal Defects

Textured Top Surface Provides Scaffold For Soft Tissue

Smooth Bottom Surface Against Bone For Easy Retrieval

Pore Size Less Than 0.3 Microns



Non-Resorbable PTFE Membrane 1 Membrane Per Box


Size	Product #	Price	10 or More	20 or More
12mm x 24mm	#CYTOSURG-12X24	52.00/Membrane	44.50/Membrane	38.70/Membrane
25mm x 30mm	#CYTOSURG-25X30	89.40/Membrane	77.90/Membrane	68.20/Membrane





Resorbable Bovine Type 1 Collagen Plug

- Use To Stabilize Blood Clots & Preserve Graft Material In Socket Preservation
- Controls Bleeding In Extraction & Biopsy Sites
- Resorbs in 10-14 Days




Resorbs In 10 To 14 Days

Size	Product #	Price	7 or More	14 or More
10mm x 20mm	#ORAPLUG-PKG-10	152.00 BOX OF 10	136.50 Box of 10	125.60 Box of 10



**80% Calcium Phosphate Granules
Combined With 20% Bovine Collagen**

- Simple, Cost Effective Socket Preservation
- Mineral Mimics Natural Bone In Resorption & Remodeling Profile
- Plug Takes The Shape Of Defect When Hydrated To Fill Socket
- Offered In Two Sizes - Five Plugs Per Box



Size	Product #	Price	7 or More	14 or More
10mm x 20mm	#FULFILL-10X20-PKG-5	284.00 BOX OF 5	266.00 Box of 5	251.30 Box of 5
6mm x 25mm	#FULFILL-6X25-PKG-5	284.00 BOX OF 5	266.00 Box of 5	251.30 Box of 5



Resorbable Collagen Tape

- Use To Repair Sinus Membrane Perforations & Stabilize Blood Clots In Minor Oral Wounds
- Resorbs in 10-14 Days

Tape Shown Actual Size - Package Contains 10 Dressings



Resorbs In 10 To 14 Days

Size	Product #	Price	7 or More	14 or More
25mm x 75mm	#ORATAPE-PKG-10	260.00 BOX OF 10	235.70 Box of 10	220.50 Box of 10

Resorbable Collagen Wound Dressings

- For Palatal Donor Sites & To Stabilize Blood Clots In Minor Oral Wounds
- 3X Thicker Than OraTAPE®
- Controls Bleeding & Stabilizes Blood Clots
- Resorbs in 10-14 Days



Cote Shown Actual Size - Package Contains 10 Dressings



Resorbs In 10 To 14 Days

Size	Product #	Price	7 or More	14 or More
20mm x 40mm	#ORACOTE-PKG-10	205.00 BOX OF 10	185.70 Box of 10	167.30 Box of 10

Tricalcium Phosphate + Hydroxyapatite

- Promotes Osteoblast Proliferation For New Host Bone Formation
- Alternative To Allografts Or Xenografts For GBR Procedures
- Slow Resorbing, Excellent For Space Maintenance



SALVIN Orastruct® 80/20 TCP/HA Synthetic Graft Material

Size	Particulate Size	Product #	Price	7 or More	14 or More
0.5cc	0.5 - 1mm	#ORASTRUCT-0.5CC	89.00	76.70	67.50
1.0cc	0.5 - 1mm	#ORASTRUCT-1.0CC	128.00	111.80	100.40
2.5cc	1 - 2mm	#ORASTRUCT-2.5CC	216.00	196.50	181.60
0.5cc In Syringe - Pkg of 5	0.5 - 1mm	#ORASTRUCT-SYR-0.5CC	319.00	290.70	271.50




Salvin Mesh-Ti™ Titanium Mesh



For Large Defects


- 0.2mm Thickness
- 1.48mm (1,480 Micron) Pore Size
- Larger Pore Size For Enhanced Growth Factor Cross Over
- Secure With Salvin Fixation Screws
- Packaged Non-Sterile



Shown Actual Size

50mm x 37mm

#SALVIN-TI-MESH-50X37 116.50 EACH



Shown Actual Size



25mm x 37mm

#SALVIN-TI-MESH-25X37 74.20 EACH

Mesh & Wire Cutter



#MESH-WIRE-CUTTER
..... 199.20 EACH



A Next Generation Demineralized Bone Matrix (DBM)


- Exceptional Handling – Moldable Putty Conforms To The Graft Site
- Each Lot Tested To Ensure Osteoinductive Potential
- Thickens At Body Temperature, For Better Graft Containment
- Packaged Sterile In A Syringe

Evaluation of Osteoinductive Potential

An In Vitro Study Was Conducted To Examine Accessibility To Bioactive Proteins In ABM And DBM¹. The Content Of BMP-2 Was Measured Over Time Using An Enzyme Linked Immunosorbent Assay (ELISA).

Study Results

- Soluble BMP-2 Was Detected In ABM At An Earlier Time Compared To Particulate DBM
- ABM's Open Pore Structure Provides Earlier Accessibility To Bone Proteins



Accell Connexus



Size	Product #	Fee
0.5 cc	#ACCELL-CONN-0.5CC	337.00
1.0 cc	#ACCELL-CONN-1.0CC	514.00
2.5 cc	#ACCELL-CONN-2.5CC	908.00

No Refrigeration Required - Store At Room Temperature
Typical Shelf Life: 6-9 Months

1. Khaliq S, Lollis R, Bell D, Oliver R, Walsh WR, and Ingram R, Evaluation of a Next Generation DBM Putty in a Posterolateral Spinal Fusion Model, (2009) Integra LifeSciences Corporation.



**Demineralized Cortical Bone With Mineralized Cancellous Chips
Paste Or Putty**

- For GBR, Ridge Augmentation, Sinus Lifts & Bony Defects
- Thickens At Body Temp To Keep Graft In Place
- Syringe Option For Easy Placement



DynaBlast Paste

Size	Product #	Fee
0.5 cc	#ORTHO-PASTE-0.5CC	148.00
1.0 cc	#ORTHO-PASTE-1.0CC	206.00
3.0 cc	#ORTHO-PASTE-3.0CC	529.00



DynaBlast Putty

Size	Product #	Fee
1.0 cc	#ORTHO-PUTTY-1CC	252.00
2.5 cc	#ORTHO-PUTTY-2.5CC	505.00
5.0 cc	#ORTHO-PUTTY-5.0CC	894.00

No Refrigeration Required - Store At Room Temperature
Typical Shelf Life: 6-9 Months

**Demineralized Cortical Bone
Paste Or Putty**

- For Osseous Defects
- Thickens At Body Temp To Keep Graft In Place
- Syringe Option For Easy Placement



DynaGraft-D Gel

Size	Product #	Fee
0.5 cc	#DYNA-GEL-0.5CC	148.00
1.0 cc	#DYNA-GEL-1.0CC	206.00
5.0 cc	#DYNA-GEL-5.0CC	840.00



DynaGraft-D Putty

Size	Product #	Fee
1.0 cc	#DYNA-PUTTY-1.0CC	252.00
2.5 cc	#DYNA-PUTTY-2.5CC	505.00
5.0 cc	#DYNA-PUTTY-5.0CC	894.00

No Refrigeration Required - Store At Room Temperature
Typical Shelf Life: 6-9 Months







NEEDLE SELECTION GUIDE

	<p>▽ Reverse Cutting Triangular point with cutting edge on the outer curvature.</p>	
	<p>▽ Premium Reverse Cutting Same as Reverse Cutting, but precision-honed for added sharpness.</p>	
	<p>⊙ Taper Point Round shaft, straight or curved, tapered point with no cutting edge. Special tapered point needle for connective tissue grafts. Ideal for delicate tissue.</p>	
	<p>▽ Taper Cutting Four (4) cutting edges with a taper shaft. Useful in microsurgery for excellent penetration through tough skin.</p>	

The Suture Book By Dr. Lee Silverstein



Beautifully Illustrated Step-By-Step Instructions



The Suture Book: The Definitive Guide To Dental Suturing And Surgical Flap Closure By Dr. Lee Silverstein

- Detailed, Illustrated Step-By-Step Instructions For Dental Suturing
- Characteristics Of Suture Materials & Needle Designs
- Suture Selection Recommendations For Various Procedures

#SUTURE-BOOK 99.00 EACH
6 OR MORE 89.00 EACH



Complete Selection of Sutures
On Pages 35-38



Needle Holders
On Pages 160-163



Non-Resorbable PTFE Sutures

- *Excellent Choice For Bone Grafting & Implant Procedures*
- *Smooth & Flexible – Not Fragile Or Brittle – So Knots Stay Tied*
- *Proprietary Monofilament Is Biologically Inert & Non-Cytotoxic*
- *Longer Length – All PTFE Sutures Now 24" Long*








Black Needle For Higher Visibility & Less Glare
Super Sharp Needle Requires Less Force – Easier Suturing With Less Trauma

“Racetrack” Rather Than Zig-Zag Package



Little-To-No Package Memory For Excellent Handling And Secure Knotting

Ultra-Soft, Ultra-Smooth Suture Designed With Low Surface Porosity

SALVIN CytoSurg Competitor G Competitor O
(2000x Magnification)

	Salvin Item#	Size	Needle	Length	Per Box	10 Or More
	#PTFE-20FS2	2-0	19mm 3/8 Circle - Reverse Cutting	24" (60cm)	115.60 Box 12	102.90 Box 12
MOST POPULAR	#PTFE-30FS2	3-0	19mm 3/8 Circle - Reverse Cutting	24" (60cm)	115.60 Box 12	102.90 Box 12
	#PTFE-30FS3	3-0	16mm 3/8 Circle - Reverse Cutting	24" (60cm)	115.60 Box 12	102.90 Box 12
	#PTFE-40TF	4-0	13mm 1/2 Circle - <u>Taper Point</u>	24" (60cm)	115.60 Box 12	102.90 Box 12
MOST POPULAR	#PTFE-40FS3	4-0	16mm 3/8 Circle - Reverse Cutting	24" (60cm)	115.60 Box 12	102.90 Box 12
	#PTFE-40PS2	4-0	19mm 3/8 Circle - <u>Premium</u> Reverse Cutting	24" (60cm)	136.20 Box 12	121.50 Box 12
	#PTFE-40P3	4-0	13mm 3/8 Circle - <u>Premium</u> Reverse Cutting	24" (60cm)	136.20 Box 12	121.50 Box 12
	#PTFE-50P3	5-0	13mm 3/8 Circle - <u>Premium</u> Reverse Cutting	24" (60cm)	136.20 Box 12	121.50 Box 12
	#PTFE-60P3	6-0	13mm 3/8 Circle - <u>Premium</u> Reverse Cutting	24" (60cm)	136.20 Box 12	121.50 Box 12



SALVIN® SurgiPoint® PGA Polyglycolic Acid Synthetic Absorbable-Braided



- Synthetic Absorbable, Braided Coated Polyglycolic Acid – Similar To Vicryl™ But Greater Tensile Strength
- Unique PolyGlide Coating For Smooth No-Drag Suture Passage Through Tissue

Resorbs In 21-28 Days

	Salvin Item#	Size	Needle	Length	Color	Vicryl™ Comparison	Per Box	10 Or More
 FS-2 (C6) <i>MOST POPULAR</i>	#LV-423B	3-0	3/8 Circle - Reverse Cutting	30" (76cm)	Undyed	J423G	87.50 Box 12	81.70 Box 12
	#LV-422B	4-0	3/8 Circle - Reverse Cutting	30" (76cm)	Undyed	J422G	87.50 Box 12	81.70 Box 12
	#LV-MC397NC	4-0	3/8 Circle - Reverse Cutting	30" (76cm)	Violet		128.60 Box 12	120.10 Box 12
	#LV-421B	5-0	3/8 Circle - Reverse Cutting	30" (76cm)	Undyed	J421G	87.50 Box 12	81.70 Box 12
 FS-3 (C22)	#LV-500B	5-0	3/8 Circle - Reverse Cutting	18" (45cm)	White	J500G	83.90 Box 12	79.10 Box 12
 P-3 (C3) <i>MOST POPULAR</i>	#LV-493B	5-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed	J493G	151.50 Box 12	142.00 Box 12
	#LV-492B	6-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed	J492G	151.50 Box 12	142.00 Box 12
 PS-2 (PC31) <i>MOST POPULAR</i>	#LV-497B	3-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed	J497G	151.50 Box 12	142.00 Box 12
	#LV-496B	4-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed	J496G	151.50 Box 12	142.00 Box 12
	#LV-495B	5-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed	J495G	151.50 Box 12	142.00 Box 12
 SH-1 (T14) <i>MOST POPULAR</i>	#LV-311B	3-0	1/2 Circle - <u>Taper Point</u>	30" (76cm)	Violet	J311G	78.20 Box 12	71.90 Box 12
	#LV-310B	4-0	1/2 Circle - <u>Taper Point</u>	30" (76cm)	Violet	J310G	78.20 Box 12	71.90 Box 12

SALVIN® SurgiPoint® PGA Polyglycolic Acid Synthetic Absorbable-Braided Fast Absorbing



- Good Alternative To Plain Gut Sutures
- Synthetic Absorbable, Braided Coated Polyglycolic Acid – Similar To Vicryl Rapide™ But Greater Tensile Strength
- Unique PolyGlide Coating For Smooth No-Drag Suture Passage Through Tissue

Resorbs In 7-14 Days

	Salvin Item#	Size	Needle	Length	Color	Vicryl™ Comparison	Per Box	10 Or More
 FS-2 (C6) <i>MOST POPULAR</i>	#LVFA-493N	3-0	3/8 Circle - Reverse Cutting	18" (45cm)	Undyed	VR497	102.10 Box 12	95.10 Box 12
	#LVFA-423N	3-0	3/8 Circle - Reverse Cutting	30" (76cm)	Undyed		102.10 Box 12	95.10 Box 12
	#LVFA-492N	4-0	3/8 Circle - Reverse Cutting	18" (45cm)	Undyed	VC422H	102.10 Box 12	95.10 Box 12
	#LVFA-422N	4-0	3/8 Circle - Reverse Cutting	30" (76cm)	Undyed	VR426	102.10 Box 12	95.10 Box 12
 FS-1 (C7)	#LVFA-442N	3-0	3/8 Circle - Reverse Cutting	30" (76cm)	Undyed	VR935	102.10 Box 12	95.10 Box 12
 PS-2 (PC31) <i>MOST POPULAR</i>	#LVFA-497N	3-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed		102.10 Box 12	95.10 Box 12
	#LVFA-496N	4-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed		102.10 Box 12	95.10 Box 12
 P-1 (C1)	#LVFA-490C	5-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed		132.60 Box 12	121.40 Box 12

More Sutures On Next 2 Pages →



SALVIN[®] SurgiPoint[™] Silk *Non-Resorbable*



• Non-Resorbable, Braided, Coated Silk – Excellent Handling & Tying Characteristics

	Salvin Item#	Size	Needle	Length	Ethicon [™] Comparison	Per Box	10 Or More
FS (C8)	#L-770B	2-0	3/8 Circle - Reverse Cutting	18" (45cm)	685-G	57.90 Box 12	53.70 Box 12
FS-2 (C6) MOST POPULAR	#L-784B	3-0	3/8 Circle - Reverse Cutting	18" (45cm)		57.90 Box 12	53.70 Box 12
	#L-781B	4-0	3/8 Circle - Reverse Cutting	18" (45cm)	683G	57.90 Box 12	53.70 Box 12
	#L-754B	5-0	3/8 Circle - Reverse Cutting	18" (45cm)	682G	57.90 Box 12	53.70 Box 12
FS-1 (C7)	#L-786B	3-0	3/8 Circle - Reverse Cutting	18" (45cm)	684G	57.90 Box 12	53.70 Box 12
C-2 (C17)	#L-787B	4-0	3/8 Circle - Reverse Cutting	18" (45cm)	736G	57.90 Box 12	53.70 Box 12
X-1 (C31)	#L-785B	3-0	1/2 Circle - Reverse Cutting	18" (45cm)	632G	57.90 Box 12	53.70 Box 12

SALVIN[®] SurgiPoint[™] Polypropylene *Non-Resorbable*



• Non-Resorbable Blue Monofilament – Inert, Easy To Handle, Passes Easily Through Tissue

	Salvin Item#	Size	Needle	Length	Ethicon [™] Comparison	Per Box	10 Or More
P-3 (C3) MOST POPULAR	#L-1085B	5-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	8698G	168.90 Box 12	156.60 Box 12
	#L-1082B	6-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	8695G	168.90 Box 12	156.60 Box 12
P-1 (C1)	#L-1079B	6-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	8697G	168.90 Box 12	156.60 Box 12
	#L-2468N	7-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	8696G	168.90 Box 12	156.60 Box 12
C-2 (C17)	#L-1011B	6-0	3/8 Circle - Reverse Cutting	10" (25cm)	8889H	168.90 Box 12	156.60 Box 12
FS-2 (C6)	#L-8661B	5-0	3/8 Circle - Reverse Cutting	18" (45cm)	8661G	86.90 Box 12	81.90 Box 12
FS-3 (C22)	#L-8680B	5-0	3/8 Circle - Reverse Cutting	18" (45cm)	8680G	86.90 Box 12	81.90 Box 12

SALVIN[®] SurgiPoint[™] Nylon Monofilament *Non-Resorbable*



• Non-Resorbable Black Monofilament – Passes Easily Through Tissue – Similar To Resolon[™]

	Salvin Item#	Size	Needle	Length	Ethicon [™] Comparison	Per Box	10 Or More
FS-2 (C6) MOST POPULAR	#L-922B	4-0	3/8 Circle - Reverse Cutting	18" (45cm)	662G	67.90 Box 12	62.60 Box 12
	#L-918B	5-0	3/8 Circle - Reverse Cutting	18" (45cm)	661G	67.90 Box 12	62.60 Box 12
P-3 (C3)	#L-698BKC	5-0	3/8 Circle - Black Premium Reverse	18" (45cm)	698G	120.20 Box 12	113.70 Box 12
	#L-1698B	6-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	1698G	120.20 Box 12	113.70 Box 12

SALVIN[®] SurgiPoint[™] Polyviolene-Green Braided Polyester *Non-Resorbable*



• Green Braided Uncoated Polyester – Smooth, Easy To Handle With Good Knot Security

	Salvin Item#	Size	Needle	Length	Ethicon [™] Comparison	Per Box	10 Or More
FS-2 (C6) MOST POPULAR	#L-1180B	4-0	3/8 Circle - Reverse Cutting	18" (45cm)	R683G	57.90 Box 12	53.70 Box 12

SALVIN® SurgiPoint Chromic Gut *Resorbs In 7-14 Days*



- Natural Resorbable Gut – Excellent Suture Except Where There Is Significant Tension Across The Wound

	Salvin Item#	Size	Needle	Length	Ethicon™ Comparison	Per Box	10 Or More
<p>FS-2 (C6) <i>MOST POPULAR</i></p>	#L-560B	3-0	3/8 Circle - Reverse Cutting	27" (68cm)	636G	58.90 Box 12	53.90 Box 12
	#L-559B	4-0	3/8 Circle - Reverse Cutting	27" (68cm)	635G	58.90 Box 12	53.90 Box 12
	#L-525B	3-0	3/8 Circle - Reverse Cutting	18" (45cm)		58.90 Box 12	53.90 Box 12
	#L-558B	4-0	3/8 Circle - Reverse Cutting	18" (45cm)		58.90 Box 12	53.90 Box 12
	#L-557B	5-0	3/8 Circle - Reverse Cutting	27" (68cm)		58.90 Box 12	53.90 Box 12
<p>P-3 (C3) <i>MOST POPULAR</i></p>	#L-595B	4-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	1654G	112.60 Box 12	105.60 Box 12
	#L-1248B	5-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	687G	112.60 Box 12	105.60 Box 12
<p>P-1 (C1)</p>	#1242B-CP	6-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)		185.30 Box 12	173.20 Box 12
<p>SH-1 (T4)</p>	#L-589B	4-0	1/2 Circle - <u>Taper Point</u> For Connective Tissue Grafts	14" (35cm)		72.40 Box 12	66.50 Box 12
<p>RB-1 (T28)</p>	#L-590B	5-0	1/2 Circle - <u>Taper Point</u> For Connective Tissue Grafts	14" (35cm)		58.90 Box 12	53.90 Box 12

SALVIN® SurgiPoint Gut-Plain *Resorbs In 5-10 Days*



- Natural Resorbable Gut – Excellent Suture Except Where There Is Significant Tension Across The Wound

	Salvin Item#	Size	Needle	Length	Ethicon™ Comparison	Per Box	10 Or More
<p>FS-2 (C6) <i>MOST POPULAR</i></p>	#L-554B	3-0	3/8 Circle - Reverse Cutting	27" (68cm)	H822G	58.90 Box 12	53.90 Box 12
	#L-551B	4-0	3/8 Circle - Reverse Cutting	27" (68cm)	H822G	58.90 Box 12	53.90 Box 12
	#L-591B	4-0	3/8 Circle - Reverse Cutting	18" (45cm)		58.90 Box 12	53.90 Box 12
	#L-592B	3-0	3/8 Circle - Reverse Cutting	18" (45cm)		58.90 Box 12	53.90 Box 12
<p>P-3 (C3)</p>	#L-596B	5-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	686G	112.60 Box 12	105.60 Box 12
<p>X-1 (C31)</p>	#L-594B	3-0	1/2 Circle - Reverse Cutting	18" (45cm)	612G	58.90 Box 12	53.90 Box 12
<p>PC-1 (CP1)</p>	#L-1415B	6-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	1916G	189.30 Box 12	177.10 Box 12

SALVIN® SurgiPoint Monoderm™ *Resorbs In 21-28 Days*



- Good Alternative To Chromic Gut Sutures
- Synthetic Resorbable Monofilament Glycolide & E-Caprolactone – Similar to Monocryl™

	Salvin Item#	Size	Needle	Length	Ethicon™ Comparison	Per Box	10 Or More
<p>FS-2 (C6)</p>	#L-388	3-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)		92.60 Box 12	83.50 Box 12
	#L-392C	4-0	3/8 Circle - Black Premium Reverse	18" (45cm)		92.60 Box 12	83.50 Box 12
	#L-Q919	5-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)		92.60 Box 12	83.50 Box 12
<p>FS-3 (C22)</p>	#L-395	4-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)		92.60 Box 12	83.50 Box 12
<p>PS-2 (DGL19) <i>MOST POPULAR</i></p>	L-YD513N	4-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	Y513G	149.70 Box 12	140.70 Box 12
<p>P-3 (C3)</p>	#L-464B	4-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)		156.60 Box 12	147.40 Box 12
	#L-463B	5-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)		156.60 Box 12	147.40 Box 12
	#L-463C	5-0	3/8 Circle - Black Premium Reverse	18" (45cm)	Y463G	156.60 Box 12	147.40 Box 12
	#L-492B	6-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)		156.60 Box 12	147.40 Box 12
	#L-401B	6-0	3/8 Circle - Violet Premium Reverse	18" (45cm)		156.60 Box 12	147.40 Box 12

Suture Names Are Trademarks Of Their Respective Companies.

PeriAcryl® by GluStitch

Purified Cyanoacrylate Dental Adhesive

- Will Adhere To Moist Living Tissue Without Toxic Or Foreign Body Reaction
- Deep Violet Color Provides Visibility For Controlled Placement
- One Year Guaranteed Shelf Life



Single-Use Applicators Kit

Includes: 12 0.2mL Pre-Loaded Applicators

#PERIACRYL-SINGLE	129.20 BOX OF 12
2 OR MORE BOXES	121.50 BOX OF 12



MOST POPULAR

5mL Bottle Multi-Use Kit

Includes: 5mL Bottle, 50 Pipettes & 1 Reusable Administration Tray

#PERIACRYL	125.90 EACH
2 OR MORE	115.10 EACH

High Viscosity PeriAcryl® by GluStitch

- *High Viscosity Helps Control Flow For Precision Application*
- *Nine Times Thicker Than Regular PeriAcryl For Easier Handling*



NEW!

Single-Use Applicators Kit

Includes: 12 0.2mL Pre-Loaded Applicators

#PERIACRYL-HV-SINGLE	133.50 BOX OF 12
2 OR MORE BOXES	127.40 BOX OF 12




5mL Bottle Multi-Use Kit

Includes: 5mL Bottle, 50 Pipettes & 1 Reusable Administration Tray

#PERIACRYL-HV	133.80 EACH
2 OR MORE	124.50 EACH

Micropoint Suture Removal Scissors

Takes The "Wince" Out Of Suture Removal, Even In Sensitive Tissue



Opposing Spring Tips Maintain Force Between Blades For Incredible Cutting Without Pinching Tissue



MOST POPULAR

Featherlite 14cm Curved

Delicate "Castroviejo Style" Scissors

#F-LITE5CURV	335.00 EACH
3 OR MORE	317.30 EACH



Complete Selection Of Suture Scissors On Pages 157-159

OraFOAM™ Absorbable Collagen Hemostatic Sponges



Package Contains 18 Sponges

Absorbable Collagen Hemostatic Sponges

- 0.5" x 1" (1.27cm x 2.54cm)
- Controls Bleeding In 2 To 5 Minutes
- Coherent Sponge Structure Allows Easily Controlled Application

#ORAFOAM-PKG-18	327.00 BOX OF 18
2 OR MORE BOXES	309.30 BOX OF 18

OraFOAM™ Hemostatic Sponge Is Indicated In Surgical Procedures (Other Than Ophthalmological & Urological Surgery) As An Adjunct To Hemostasis When Control Of Bleeding By Ligature Or Conventional Procedures Is Ineffective Or Impractical. Please See Product IFU For Precautions & Adverse Reactions. IFU Available Online At www.salvin.com/Downloads/Guides-Studies.html

BloodSTOP® Absorbent Hemostatic Gauze



BloodSTOP® Absorbent Hemostatic Gauze

- Made Of 100% Natural Cellulose
- Accelerates Coagulation
- Reduces Healing Time & Recurring Bleeding

BloodSTOP® 0.5" x 2"

#BLOODSTOP-IX-0.5X2	195.50 BOX OF 24
2 OR MORE BOXES	185.70 BOX OF 24

BloodSTOP® 2" x 2"


#BLOODSTOP-IX-2X2	345.00 BOX OF 12
2 OR MORE BOXES	329.30 BOX OF 12

Hemostats



Mosquito Hemostat
12.7cm

#MOSQ-CURV-5
..... 71.10 EACH



MOST POPULAR
Kelly Hemostat
Curved - 14cm

#KELLY-CURV-HEMO
..... 71.10 EACH

Bone Wax



- Packaged Sterile In 2.5 Gram Pouch For Easy Dispensing
- Made From Beeswax & Other Natural Materials


#BONE-WAX	
BOX OF 12 TUBES	159.90 BOX
2 OR MORE	152.40 BOX



HemCon® Dental DressingPro

- *Easy To Use – Adhesive Properties Help Secure Dressing To Mucosal Tissue*
- *Safe – Protects The Wound Surface With No Known Contra-Indications*





"I have used HemCon® for many years to rapidly and predictably control bleeding during surgery. I highly recommend it to my students at the Misch Resnik Implant Institute."

–Dr. Randy Resnik
Prosthodontist, Pittsburgh, PA

**Clinical Studies On File*

Package Contains 12 Dressings
Dressings Shown Actual Size



10mm x 12mm
Best For Extraction Sockets



1" x 3" (25mm x 75mm)
Best For Palatal Wound Sites

HemCon® Dental DressingPro 10mm x 12mm

#HEMCON-10X12MM 190.00 BOX OF 12
2 OR MORE BOXES 179.90 BOX OF 12

HemCon® Dental DressingPro 1" x 3" (25mm x75mm)

#HEMCON-1X3INCH 278.00 BOX OF 12
2 OR MORE BOXES 263.70 BOX OF 12

Portable Electrocautery System



Complete Set Includes:

- Handle & Batteries
- 2 Fine Tips & 2 Loop Tips
- 10 Sterile Handle Sheaths

#CAUT-KIT-DEL-1 129.50 SET



Replacement "Fine" Tips

Box Of 10 Single Use "Fine" Tips
With Sterile Handle Sheaths

#CAUT-FINE-H101 89.40 BOX OF 10



Replacement "Loop" Tips

Box Of 10 Single Use "Loop" Tips
With Sterile Handle Sheaths

#CAUT-LOOP-H103 89.40 BOX OF 10

BUYER'S GUIDE - MONITORS

Criticare nGenuity™ 8100
Page 44



- ▶ Pulse Rate
- ▶ O₂ Saturation
- ▶ Auto Blood Pressure
- ▶ Electrocardiogram
- ▶ Expired CO₂

Zoe Medical Nightingale™ PPM3
Page 44



- ▶ Pulse Rate
- ▶ O₂ Saturation
- ▶ Auto Blood Pressure
- ▶ Electrocardiogram
- ▶ Expired CO₂

Mindray ePM 10A
Page 45

MOST POPULAR



- ▶ Pulse Rate
- ▶ O₂ Saturation
- ▶ Auto Blood Pressure
- ▶ Electrocardiogram
- ▶ Expired CO₂

Criticare eQuality™ 506
Page 43



- ▶ Pulse Rate
- ▶ O₂ Saturation
- ▶ Auto Blood Pressure

Mindray VS8
Page 43

MOST POPULAR



- ▶ Pulse Rate
- ▶ O₂ Saturation
- ▶ Auto Blood Pressure

**BUYER'S GUIDE -
FOR MONITORS ON THIS PAGE**

- ▶ Pulse Rate
- ▶ **O₂ Saturation**
- ▶ Auto Blood Pressure

Mindray VS8

MOST POPULAR

- 8.4" Touchscreen Color LCD Display
- Nellcor OxiMax SpO₂
- Built-In Printer
- Three Year Warranty

#MINDRAY-VS8

Call For Lowest Price

Latexfree
PRESSURE CUFF



Built-In Printer



Roll Stand For Mindray VS8

Quick Release Mount & Cable Management

#MINDRAY-ROLL-STAND-VS-EPM

..... **CALL FOR PRICE**



Criticare eQuality™ 506

- **ComfortCuff™** Measures On Inflation For Fast, Comfortable Readings

#CR-506-EQUALITY

Call For Lowest Price

Latexfree
PRESSURE CUFF



Bluetooth Printer Included

Roll Stand For eQuality™ 506 Monitor

Adjusts 37" To 47" High, 20" Diameter Base

#CR-ROLL-506 **CALL FOR PRICE**



BUYER'S GUIDE - FOR MONITORS ON THIS PAGE

▶	Pulse Rate
▶	O ₂ Saturation
▶	Auto Blood Pressure
▶	Electrocardiogram
▶	Expired CO ₂



Zoe Medical Nightingale™ PPM3

- Connectivity With Major Dental EMR Systems
- Masimo ISA CO₂ Analyzer
- Built-In Printer
- Three Year Warranty



WITH CO₂
#ZOE-PPM3-CO2



WITHOUT CO₂
#ZOE-PPM3

Call For Lowest Price

Roll Stand For Zoe Medical Nightingale™ PPM3 Monitor

21" Diameter Base

#ZOE-PPM3-ROLL-STAND . . CALL FOR PRICE



Criticare nGenuity™ 8100

- **ComfortCuff™** Measures On Inflation For Fast, Comfortable Readings
- Three Lead ECG With Freeze & Cascade
- Built-In Printer

WITH CO₂
#CR-NGENUITY-CO2

WITHOUT CO₂
#CR-NGENUITY-8100

Call For Lowest Price

Roll Stand For nGenuity™ 8100 Monitor

Adjusts 37" To 47" High, 20" Diameter Base

#CR-ROLL-8100CALL FOR PRICE

Latexfree
PRESSURE CUFF



	Pulse Rate
	O ₂ Saturation
	Auto Blood Pressure
	Electrocardiogram
	Expired CO ₂

Mindray ePM 10A MOST POPULAR

- 10" Multi-Gesture Touchscreen Display
- Extensive Data Storage Capability For Trend Data, Alarms, And Events
- Integrated Sidestream CO₂
- Built-In Printer
- Three Year Warranty



WITH CO₂
#MINDRAY-EPM-10A-CO2



WITHOUT CO₂
#MINDRAY-EPM-10A

Roll Stand For Mindray ePM 10A Monitor

21" Diameter Base
#MINDRAY-ROLL-STAND-VS-EPM. . . . CALL FOR PRICE



Call For Lowest Price

CO₂ Monitor Accessories



CO₂ Monitor Single Cannula

- Single Use
 - For CO₂ Use Only
- #CANNULA-SINGLE
..... 221.70 BOX OF 50



CO₂ Monitor Divided Cannula

- Single Use
 - Works With Any O₂ System
- #CANNULA-DIVIDED
..... 97.90 BOX OF 25



CO₂ Monitor Sample Line

- Replace Every 20 Uses
 - Connects Cannula & Water Trap
- #SAMPLE-LINE
..... 112.30 BOX OF 50



Zoe Water Trap

330 Working Hours
#ZOE-WATER-TRAP 49.40 EACH



Criticare nGenuity™ Water Trap

- Replace Every 20 Uses
- Water Trap & Filter

#CR-938-FNC 559.00 BOX OF 30



Mindray Water Trap

- For Use Up To One Month

#MINDRAY-EPM10A-CO2-WATER-TRAP
..... 68.00 BOX OF 2

Durable Steel Salvin Implant / Surgical Cart

Accommodates Implant Motors, Piezotomes & Vital Sign Monitors



All Metal Construction

41.5" Working Height For Optimal Visibility Of Equipment Screens


Salvin Implant / Surgical Cart

- Pull Out Side Extension
- Removable Antimicrobial Plastic Top For Easy Cleaning
- No Slam, Soft Closing Drawers
- Removable Ergonomic Side Handles
- 5" Easy Rolling Wheels With Locks
- Power Strip With Six Outlets Attached For Easy & Neat Wire Management
- Foot Pedal Fits In Bottom Space For Mobility & Storage
- 20.5" Wide x 18.5" Deep x 41.5" Tall (Height Includes Wheels)

#SURG-CART-STEEL 1795.00



Spacious Drawers Accommodate Supplies & Wrapped Cassettes



Pull Out Extension Supports Instruments or Equipment

Surfaces & Instruments Shown Undraped For Clarity

Mayo & IV Stands



MOST POPULAR

Heavy Duty Mayo Stand

- Tray Size Is 19 1/8" x 12 5/8" (48.6cm x 32cm)
- Adjusts From 33" To 53" (83cm to 135cm)
- Includes Stainless Steel Tray


#MAYO STAND 333.90 EACH



Center Post Mayo & Implant Motor Stand

- Tray Size Is 19 1/8" x 12 5/8" (48.6cm x 32cm)
- Adjusts From 29 1/2" To 50" (75cm to 127cm)
- Includes Stainless Steel Tray

#MAYO-STAND-CP..... 342.60 EACH




Heavy Duty IV Stand

- Adjusts From 51 1/4" To 94" (130cm to 239cm)
- Chrome Plated Heavy Bar Steel

#IV STAND 179.90 EACH

Sterile Mayo Stand Cover

- Individually Packaged
- Fits Regular Mayo Stands
- 24" x 24" (61 x 61cm)




#MAYO-COVER
CASE OF 24 314.00 CASE
2 OR MORE CASES 302.80 CASE

Latexfree

Sterile Surgical Table Cover

- Individually Packaged
- Fits Surgical Tables
- 53" x 77" (134 x 195cm)



#SURG-TABLE-COVER
CASE OF 20 109.90 CASE
2 OR MORE CASES 102.40 CASE

Latexfree

Over The Patient Tables



MOST POPULAR

"Over The Patient" Surgical Table

- Putty Color Work Top Is 24" x 39" (61cm x 99cm)
- Height Adjusts From 31" To 46" Tall (78.75cm To 116.85cm)
- Dimension Between Legs Is 34" (86.4cm)

#SURG-TABLE-24X39 2050.00 EACH



Stainless Steel "Over The Patient" Surgical Table

- Extra Heavy-Duty Hospital Style Stainless Steel
- Work Top Is 24" x 48" (61cm x 122cm)
- Height Easily Adjusts With Hand Crank From 36" to 56" (91.5cm to 142cm)
- Dimension Between Legs Is 40" (101.6cm)

#SURG-TABLE-STAINLESS 3690.00 EACH

Memory Foam Chair Pad For "Long Appointment" Patient Comfort



MOST POPULAR

- Fits All Dental Chairs - 61" Long x 18" Wide x 3" Thick
- Non-Skid Bottom Won't Slip Or Move On Your Chair
- Sealed Seams & Fluid Proof Top Cover For Infection Control
- Anti-Microbial Foam & Cover

#CHAIR-PAD-SET 749.00 CUSHION & HEADREST
2 OR MORE 696.00 CUSHION & HEADREST

Precordial Amplified Sedation Stethoscope

- Chest Piece With Built-In Microphone Amplifies Lung Sounds While Suppressing Room Noise
- Diagnostic Quality Sound, Less Ambient Noise
- Set Includes:
 - Precordial Amplifier
 - Volume Controlled Earbuds
 - Precordial Pickup Cable With Piezo Chest Piece
 - 2 AAA Batteries
 - Manual & Case



Wireless Bluetooth® Integrated Technology

#PRE-STETH-WIRELESS 758.00 SET

Zoll AED PLUS™ Automated External Defibrillator



- One-Piece, Pre-Connected Electrode With CPR Feedback
- Sensor Tells You If You Are Pressing At The Correct Rate
- Pictures, Text Displays & Voice Prompts Ensure Proper Sequence
- Auto Adjusts Shock Based On Detected Body Mass
- Works With Both Adults & Children
- Includes: Carry Case, Adult Electrodes, Batteries, Prep Kit & Video
- Memory Allows Download Of Rescue Data With Included Software
- Automatically Self-Tests Defibrillator, Battery & Pad Daily
- Uses 5-Year Lithium-123 Duracell Camera Batteries



#ZOLL-AED-PLUS

Call For Lowest Price

Zoll Replacement Electrodes

5 Year Adult Electrodes

#ZOLL-AED-ELECTRODE234.00 SET

2 Year Pedo Electrodes

#ZOLL-AED-PEDO-ELECTR.156.00 SET



CO₂ Capnography Nasal Hoods

- Nasal Hood With Quick Connect To Capnography Monitor Tubing
- For Delivery & Scavenging Of Nitrous Oxide & End-Tidal CO₂ Capnography Monitoring
- Secure Female Luer Lock Sample Line Connection




CO₂ Capnography Nasal Hoods

Unscented Gray - Single Use


#ACC-CO2-HOOD-GRAY 139.70 BOX OF 12

ECG Patches & Wrist Clips



ECG Patches - Wet Gel Foam

- Snap Connectors For ECG Leads
- Wet Gel Foam Electrodes
- Recommended For Use With Zoe ECG Monitors



#CR-LEADLOCK-50

1 BOX 30.90 BOX OF 50
2 OR MORE 28.80 BOX OF 50



ECG Wrist Clips

- Snap Connectors For ECG Leads
- Left Wrist Clip Has 2 Leads & Right Wrist Clip Has 1 Lead Of A 3 Lead ECG
- Use Conductive (#SIGNA-SPRAY) Solution For Best Results

#CR 1166 134.20 PAIR



ECG Patches

- Snap Connectors For ECG Leads
- Packaged In Sets Of 3

#CR-LEADLOCK-20

1 BOX 30.90 BOX OF 20 SETS
2 OR MORE 28.80 BOX OF 20 SETS



Conductive Solution For ECG Wrist Clips

- Use With ECG Snap Leads
- 250mL Bottle

#SIGNA-SPRAY 16.80 EACH

Pocket CPR Mask & Emergency Airways



Single Use Laryngeal Mask Airways (LMA)

- Set Of 3 Allows Rapid Airway Access
- Less Invasive Than An Endotracheal Tube
- Includes #3 (Pedo To Small Adult), #4 (Adult) & #5 (Large Adult)

#LMA-SET-OF-3 75.20 SET OF 3
2 OR MORE 69.50 SET OF 3



Guedel Airways

- Airway, Bite Block & Tongue Depressor
- Coated Plastic Prevents Crimping & Protects Teeth
- 100mm, 90mm, 80mm

#GUEDEL-SET 15.00 SET OF 3
3 OR MORE 13.00 SET OF 3

Nasal Cannula & Oxygen Adapter



MOST POPULAR

Soft-Touch Adult Nasal Cannula

- Curved Tips For Excellent Patient Comfort
- Crush Resistant Tubing
- Individually Packaged - 7' (2.1m) Long

#NASAL-CANNULA

BOX OF 50 83.40 BOX
2 OR MORE BOXES 77.20 BOX

Nitrous System Adapter For Nasal Cannula Tubing

**FOR OXYGEN USE BY ITSELF ONLY
USE YOUR SCAVENGER CIRCUIT WHEN USING NITROUS**

#ACC-CANNULA-ADAPT 77.00 EACH

Anesthesia Instrumentation

Magill Endotracheal Forcep

12cm Shaft Length

#MAGILL 142.90 EACH
3 OR MORE 133.10 EACH



Young Forcep

17.1cm

#YOUNG-FORCEP 225.40 EACH
3 OR MORE 212.10 EACH



Replacement Tips For Young Forcep

Set Of 2

#YOUNG-FORCEP-TIPS 17.50 SET OF 2
3 OR MORE 15.60 SET OF 2

MOST POPULAR

Yankauer Suction Tip

29.8cm

#YANKAUER. 79.30 EACH
3 OR MORE 73.50 EACH



Disposable Sterile Yankauer Tips

- Blue Translucent Plastic Minimizes Visibility Of Fluids
- Flange Tip & Vent

#YANKAUER-DISPOSABLE

BOX OF 50 136.90 BOX
2 OR MORE 127.50 BOX



To Use These Surgical Tips With The High Volume Suction In Your Office, Use Our Tubing Adapter On Page 96

Emergency Kits



MOST POPULAR

Manual Single-Use Resuscitator

- Attached Ventilation Bag
- Connects To Oxygen Outlet For Nasal Cannula

#RESUS-SINGLE 49.40 EACH
2 OR MORE 45.80 EACH



MOST POPULAR



Single-Use Emergency Suction

- Use Metal Aspirator Or Plastic Yankauer
- 1000cc Fluid Capacity

#EMERGENCY-SUCTION 83.90 EACH
2 OR MORE 77.50 EACH

Scan For

**Emergency Suction
Technique Guide**



NEW!

To Our Catalog

UNIVERSAL COLOR CODE SYSTEM – IV CATHETER SIZES

18G – GREEN	20G – PINK	22G – BLUE	24G – YELLOW
-------------	------------	------------	--------------

Terumo SURFLO® IV Catheter

- Clear Flashback Chamber To Confirm Catheter In Vein
- Opaque Catheter



Color	Size	Product #	Price
Green	18g x 1.25"	#IV-CATH-SURFLO-SROX1832CA	94.90 BOX 50
Pink	20g x 1"	#IV-CATH-SURFLO-SROX2025CA	91.90 BOX 50
Pink	20g x 1.25"	#IV-CATH-SURFLO-SROX2032CA	94.90 BOX 50
Blue	22g x 1"	#IV-CATH-SURFLO-SROX2225CA	97.90 BOX 50
Yellow	24g x 0.75"	#IV-CATH-SURFLO-SROX2419CA	91.90 BOX 50

Terumo SurFlash® IV Catheter

- Clear Flashback Chamber To Confirm Catheter In Vein
- Clear Catheter To Show Blood Drawn



Color	Size	Product #	Price
Pink	20g x 1"	#IV-CATH-SURFLASH-SR-FF2025	99.90 BOX 50
Blue	22g x 1"	#IV-CATH-SURFLASH-SR-FF2225	99.90 BOX 50
Yellow	24g x 0.75"	#IV-CATH-SURFLASH-SR-FF2419	99.90 BOX 50

EXEL Safelet IV Catheter

- Clear Flashback Chamber To Confirm Catheter In Vein
- Clear Catheter To Show Blood Drawn




Color	Size	Product #	Price
Green	18g x 1.25"	#IV-CATH-EXEL-26735	49.90 BOX 50
Pink	20g x 1"	#IV-CATH-EXEL-26741	49.90 BOX 50
Pink	20g x 1.25"	#IV-CATH-EXEL-26742	49.90 BOX 50
Blue	22g x 1"	#IV-CATH-EXEL-26746	49.90 BOX 50
Yellow	24g x 0.75"	#IV-CATH-EXEL-26751	49.90 BOX 50

EXEL IV Extension Sets


- 30"
- #IV-SET-EXEL-1730..... 164.50 CASE OF 100



BD Insyte™ Autoguard™ Shielded IV Catheter



Spring Loaded Push Button Mechanism Retracts Needle Into Safety Chamber



Color	Size	Product #	Price
Yellow	24g x 0.75"	#IV-CATH-BD-INSYTE-381412	211.79 BOX 50
Blue	22g x 1"	#IV-CATH-BD-INSYTE-381423	211.79 BOX 50
Pink	20g x 1"	#IV-CATH-BD-INSYTE-381433	211.79 BOX 50
Winged - Yellow	24g x 0.75"	#IV-CATH-BD-INSYTE-WINGED-381512	211.79 BOX 50
Winged - Blue	22g x 1"	#IV-CATH-BD-INSYTE-WINGED-381523	211.79 BOX 50
Winged - Pink	20g x 1"	#IV-CATH-BD-INSYTE-WINGED-381533	211.79 BOX 50

B. Braun Introcath Safety® IV Catheter



Passive Safety Shield – Automatically Covers Needle Tip Upon Needle Withdrawal

Helps Prevent Needlestick Injuries



Color	Size	Product #	Price
Straight – Yellow	24g x 0.75"	#IV-CATH-BB-INTROCAN-4252500-02	160.00 BOX 50
Straight – Blue	22g x 1"	#IV-CATH-BB-INTROCAN-4252519-02	160.00 BOX 50
Straight – Pink	20g x 1"	#IV-CATH-BB-INTROCAN-4252543-02	160.00 BOX 50


B. Braun ULTRASITE® Needle Free Luer Access Device




- Quick, effortless Luer access connection
- High flow promotes rapid fluid delivery
- Microbial barrier design
- Easy to clean and disinfect
- Reduces accidental touch contamination

#IV-VALVE-BB-ULTRASITE-415110. 452.36 CASE OF 100


EXEL IV Administration Sets



Description	Product #	Price
15 Drop/mL Luer Slip – Y Injection Site – 78" Tubing	#IV-SET-EXEL-29081	59.90 CASE 50
60 Drop/mL Luer Slip – Y Injection Site – 78" Tubing	#IV-SET-EXEL-29083	64.90 CASE 50



Description	Product #	Price
15 Drop/mL Luer Lock – Y Injection Site – 78" Tubing	#IV-SET-EXEL-29082	69.90 CASE 50
60 Drop/mL Luer Lock – Y Injection Site – 78" Tubing	#IV-SET-EXEL-29085	69.90 CASE 50



Description	Product #	Price
15 Drop/mL Luer Lock – 78" Tubing	#IV-SET-EXEL-27086	149.90 CASE 50
15 Drop/mL Luer Lock/Luer Slip – 78" Tubing	#IV-SET-EXEL-29086	69.90 CASE 50



Description	Product #	Price
15 Drop/mL Luer Lock/Luer Slip – 2 Y Injection Site – 105" Tubing	#IV-SET-EXEL-29088	73.90 CASE 50

Baxter IV Administration Sets



Description	Product #	Price
10 Drops/mL – 2 Injections Sites Duo Vent – 104" Tubing	#IV-SET-ACT5439	269.90 CASE 48

B. Braun SafeDay™ IV Administration Sets



Description	Product #	Price
15 Drops/mL – 106" Tubing	#IV-SET-SAFEDAY-352636	317.00 CASE 50

Syringes & Needle Capper



MOST POPULAR

Aspirating Syringe CW Type

Works With Standard 1.8cc Cartridges

#SYRINGE CW **48.80 EACH**
3 OR MORE 44.90 EACH



Petite Aspirating Syringe

More Comfort & Control For Smaller Hands

Works With Standard 1.8cc Cartridges
 #SYRINGE-PETITE. **51.10 EACH**
3 OR MORE 46.90 EACH



Aspirating Syringe A Type

Works With Standard 1.8cc Cartridges

#SYRINGE A **45.70 EACH**
3 OR MORE 41.80 EACH



Needle Capper

#NEEDLE CAPPER. **31.20 EACH**
3 OR MORE 28.80 EACH



Intraligamentary Syringe

- Positive Pressure Design Creates Extra Leverage Without Fatigue
- One Complete Trigger Squeeze Delivers 0.2cc
- Uses Standard Carpules

#NTRALIG. **489.00 EACH**



Palatal Anesthesia Device

- Uses Pressure Anesthesia & Gate Control Theory To Reduce Or Eliminate Palatal Injection Pain
- Double-Ended For Right & Left Palatal Injections

#PAD. **79.50 EACH**
3 OR MORE 75.70 EACH

#1090 Salvin Split-Sheet Implant Pack



"Eyes Visible" Patient Drape


MOST POPULAR

#1090 Salvin Split-Sheet Pack

#1090 SOLD IN CASES OF 6

1 CASE 461.70 CASE (76.95 EA)


2 CASES 429.30 CASE (71.55 EA)



BUYER'S GUIDE - SALVIN SPLIT-SHEET PACK


- ▶ "Eyes Visible" Patient Drape
- ▶ One Pack With Two Separate Inner Wraps For Easier Set Up

Doctor/Assistant Inner Wrap



- 2 Large Operating Room Surgical Gowns
- 2 17" x 20" (38 x 51cm) Absorbent Towels
- 2 Bouffant Caps (Doctor & Assistant)
- 2 Face Masks With Ties
- Outer Wrapper Becomes 24" x 24" (61 x 61cm) Table Cover


Patient / Room Setup Inner Wrap



- 1 Implant "Split Sheet" Patient Drape 77" x 108" (195.5 x 274cm)
- 1 Bouffant Cap (For The Patient)
- 1 Cotton O/R Towel
- 1 Mayo Stand Or Back Table Cover
- 1 Oral Surgery Head Drape
- 2 18" x 24" (45.7 x 61cm) Utility Drapes With Tape
- 1 Machine or X-Ray Cover
- 4 4" x 6" (10 x 15.25cm) Light Handle/Chair Button Covers
- 4 Rapid Roll® Covers
- 1 Chairside Bag
- 2 6" (15.25cm) Cotton Tip Applicators
- 1 Surgical Aspirator Tip
- 1 Suction Tip With Out Slots
- 1 Saliva Ejector
- 25 Lint Free 2" x 2" (5 x 5cm) Sponges
- 8 Lint Free 4" x 4" (10 x 10cm) Sponges
- Outer Wrapper Becomes 45" x 45" (114 x 114cm) Table Cover

Latexfree

#1033 Salvin "Eyes Covered" Implant Pack



MOST POPULAR

#1033 Basic Implant Pack

#1033 SOLD IN CASES OF 6

1 CASE 456.00 CASE (76.00 EA)

2 CASES 425.10 CASE (70.85 EA)

- Cooler
- Lighter

*This Implant Pack Contains
2 Aurora™ Gowns*

Scan For

60^{Second} Product Video™

Scan For

60^{Second} Product Video™




Sterile Surgical Gloves
On Page 66

BUYER'S GUIDE - SALVIN "EYES COVERED" PACK


- ▶ "Eyes Covered" Patient Drape
- ▶ **One Pack With Two Separate Inner Wraps For Easier Set Up**

Doctor/Assistant Inner Wrap



- 2 Large Operating Room Surgical Gowns
- 2 Face Masks With Ties
- 2 17" x 20" (38 x 51cm) Absorbent Towels
- Outer Wrapper Becomes 24" x 24" (61 x 61cm) Table Cover
- 2 Bouffant Caps (Doctor & Assistant)


Patient / Room Setup Inner Wrap



- 1 Implant Fenestrated Patient Drape 38 1/2" x 55" (97.5 x 139.5cm)
- 1 Bouffant Cap (For The Patient)
- 1 Cotton O/R Towel
- 1 Mayo Stand Or Back Table Cover
- 3 18" x 24" (45.7 x 61cm) Utility Drapes With Tape
- 1 Machine or X-Ray Cover
- 4 4" x 6" (10 x 15.25cm) Light Handle/Chair Button Covers
- 4 Rapid Roll® Covers
- 1 Chairside Bag
- 2 6" (15.25cm) Cotton Tipped Applicators
- 1 Surgical Aspirator Tip
- 1 Suction Tip With Out Slots
- 1 Saliva Ejector
- 25 Lint Free 2" x 2" (5 x 5cm) Sponges
- 8 Lint Free 4" x 4" (10 x 10cm) Sponges
- Outer Wrapper Becomes 45" x 45" (114 x 114cm) Table Cover

Latexfree

Resnik Institute Implant Surgery Pack



Resnik Institute Implant Surgery Pack

#RESNIK-PACK SOLD IN CASES OF 10
 1 CASE 414.00 CASE (41.40 EA)
 2 CASES 381.00 CASE (38.10 EA)

BUYER'S GUIDE - RESNIK INSTITUTE PACK

▶ "Eyes Covered" Patient Drape

Resnik Institute Implant Surgery Pack

- 1 Implant Fenestrated Patient Drape 38 1/2" x 55" (97.5 x 139.5cm)
- 5 Cotton O/R Towels
- 1 Chair Cover 41" x 85" (104 x 216cm)
- 1 Bouffant Cap
- 2 Rapid Roll® Covers
- 1 20" x 20" (50.8 x 50.8cm) CSR Wrap
- 2 30" x 30" (76.2 x 76.2cm) CSR Wraps
- 4 4" x 6" (10 x 15.25cm) Light Handle/ Chair Button Covers
- 2 Surgical Aspirator Tips
- 1 Suction Tip Without Slots
- 12 Lint Free 4" x 4" (10 x 10cm) Sponges
- 12 Lint Free 2" x 2" (5 x 5cm) Sponges
- Outer Wrapper Becomes 48" x 48" (121.9 x 121.9cm) Table Cover

Latexfree

Salvin E.N.T. Head & Neck Pack



Salvin E.N.T. Head & Neck Pack

#2000 SOLD IN CASES OF 8
 1 CASE 452.40 CASE (56.55 EA)
 2 CASES 417.60 CASE (52.20 EA)

BUYER'S GUIDE - SALVIN E.N.T HEAD & NECK PACK

▶ "Eyes Visible" Patient Drape

Salvin E.N.T. Head & Neck Pack

- 2 Large Operating Room Surgical Gowns
- 2 17" x 20" (38 x 51cm) Absorbent Towels
- 1 Implant "Split Sheet" Patient Drape 77" x 108" (195.5 x 274cm)
- 1 Oral Surgery Head Drape
- 1 Mayo Stand Or Back Table Cover
- 1 18" x 24" (45.7 x 61cm) Utility Drape With Tape
- 1 Chairside Bag
- Outer Wrapper Becomes 54" x 54" (137 x 137cm) Table Cover

Latexfree

- Cooler
- Lighter

This Implant Pack Contains 2 Aurora™ Gowns

Salvin Perio/Minor Surgery Pack

BUYER'S GUIDE - SALVIN PERIO/MINOR SURGERY PACK

Salvin Perio/Minor Surgery Pack


#1034 SOLD IN CASES OF 10
 1 CASE 354.00 CASE (35.40 EA)
 2 CASES 323.00 CASE (32.30 EA)



Sterile Surgical Gloves
 On Page 66

▶ This Pack Includes A Sterile Patient Bib

Salvin Perio/Minor Surgery Pack



- 2 Yellow Isolation Gowns
- 2 17" x 20" (38 x 51cm) Absorbent Towels
- 1 18" x 36" (46 x 91cm) Patient Bib
- 2 Face Masks W/Ties
- 1 Cotton O/R Towel
- 1 Chairside Bag
- 2 18" x 24" (45.7 x 61cm) Utility Drapes With Tape
- 2 4" x 6" (10 x 15.25cm) Light Handle/ Chair Button Covers
- 2 6" (15.25cm) Cotton Tipped Applicators
- 1 Suction Tip Without Slots
- 1 Saliva Ejector
- 25 Lint Free 2" x 2" (5 x 5cm) Sponges
- Outer Wrapper Becomes 36" x 36" (91 x 91cm) Table Cover



Individual Sterile Patient Drapes

MOST POPULAR

Split Sheet "Eyes Visible" Sterile Patient Drape



- Individually Packaged
- 77" x 108" (195.5 x 274cm)

#DRAPE-SPLIT SHEET
 CASE OF 24 490.00 CASE
 2 OR MORE CASES 457.70 CASE

Patient Drape Detail



#DRAPE-SPLIT SHEET


Drape 33 "Eyes Covered" Sterile Patient Drape



- Individually Packaged
- 38 1/2" x 55" (97.5 x 139.5cm)

#DRAPE33
 CASE OF 24 425.00 CASE
 2 OR MORE CASES 403.30 CASE

Patient Drape Detail



#DRAPE 33

Rapid Roll® Covers & Sterile Drapes

SALVIN Rapid Roll®

MOST POPULAR

- Individually Packaged
- 65" (165cm) Long

#RAPID-ROLL
 BOX OF 25 117.80 BOX
 2 OR MORE BOXES 108.50 BOX



Sterile Drapes

- Individually Packaged Poly Lined Paper Drape
- 18" x 26" (45.7 x 66cm)
- Use As Patient Bib, Table Cover, Etc.

#TOWEL-DRAPE-50
 BOX OF 50 58.40 BOX
 2 OR MORE BOXES 55.90 BOX



Salvin Sterile Gown/Towel Packs



Latexfree



Salvin Sterile Gown / Towel Pack - **UNISIZE Fit Gown**

- 1 Sterile O/R Gown (Large-XLarge Fit)
- 1 Absorbent Scrub Towel 15" x 30" (38 x 76cm)

#1001
 CASE OF 24 487.00 CASE
 2 OR MORE CASES 472.20 CASE

Salvin Sterile Gown / Towel Pack - **SMALL Size Gown**

- 1 Sterile O/R Gown (Small-Medium Fit)
- 1 Absorbent Scrub Towel 15" x 30" (38 x 76cm)

#GOWN-SMALL
 CASE OF 30 409.70 CASE
 2 OR MORE CASES 394.20 CASE

Sterile Absorbent Scrub Towels

- 80 Individually Packaged Absorbent Scrub Towels 15" x 30" (38 x 76cm)

#TOWEL-SCRUB-80
 CASE OF 80 161.90 CASE
 2 OR MORE CASES 152.70 CASE

Backhaus Towel Clamp

- Sharp Points
- Use For Cloth Drapes

#BACKHAUS 22.10 EACH
 3 OR MORE 20.10 EACH



Aurora™ Gown

- More Protection
- Cooler
- Lighter

O/R & Isolation Gowns – Non-Sterile



Latexfree

Aurora™ Gown

- More Protection
- Cooler
- Lighter



Latexfree

MOST POPULAR

Disposable Blue O/R Gown (Non-Sterile)

- The Feel Of Cotton Without The Weight
- Breathable Gown With Knitted Cuffs
- Tie Back & Velcro Neck Closure
- Generously Sized Unisize Fit

#OR-GOWN
 CASE OF 50 450.00 CASE
 2 OR MORE CASES 425.40 CASE

Isolation Gown-Yellow (Non-Sterile)

- Elastic Cuff
- Unisize Fit

#ISO-GOWN-50
 CASE OF 50 104.50 CASE
 2 OR MORE CASES 97.10 CASE

Surgical N95 Respirator

- Flat Pouched Design Provides A Cooler Breathing Chamber
- Non-Abrasive Comfort Layer For Moisture Management

Surgical N95 Respirator

- Approved For Use By The FDA And NIOSH
- ASTM F1862 Level 3 Fluid Resistance
- Sold In Boxes Of 50

#N95-MASK-5074.50 BOX OF 50




Towels, Caps & Shoe Covers



Latexfree

Surgeon's Caps

#SURG CAP 100
BOX OF 100 39.50 BOX
2 OR MORE BOXES36.30 BOX



Latexfree

MOST POPULAR

Bouffant Caps


#BOUF CAP 100
BOX OF 10034.70 BOX
2 OR MORE BOXES31.80 BOX



Latexfree

Blue 100% Cotton O/R Towels

16" x 27" (40.6 x 68.6cm)
Single Use
#OR TOWEL 50
CASE OF 50130.90 CASE
2 OR MORE CASES124.70 CASE




Latexfree

Shoe Covers-Unisize

- Fits Over Regular Shoes
- No-Slip Grip On Bottom For Safety In The Office


#SHOE 100
BOX OF 10042.40 BOX
2 OR MORE BOXES39.50 BOX



Latexfree

Surgical Masks W/Ties

Soft Polypropylene Facing For Comfort & Protection
#SURG-MASK-TIE
BOX OF 5025.50 BOX
6 OR MORE BOXES22.80 BOX



Latexfree

MOST POPULAR

Shoe Covers-Sport

- Fits Over Athletic Shoes
- No-Slip Grip On Bottom For Safety In The Office

#SHOE-SPORT-100
BOX OF 10047.90 BOX
2 OR MORE BOXES45.10 BOX

Patient Skin Prep

ChloraPrep® Surgical Skin Prep Individual Swabs

- Chlorhexidine Gluconate/IPA
- Effective Skin Antiseptic
- Will Not Stain Or Irritate Skin
- Rapid Action Against Microorganisms
- Available In US Only

#CHLORA-PREP
BOX OF 2061.70 BOX
2 BOXES58.20 BOX



Latexfree

MOST POPULAR

Patient Prep Individual Swab Sticks

Povidone Iodine (Iodophor)
#SWABSTICK 50
BOX OF 50 22.50 BOX
2 BOXES 20.40 BOX



Latexfree

Cold Therapy



Plastic Outer Covering

Latexfree

MOST POPULAR

Instant Cold Pack

- Squeeze To Activate
- 6" x 8.75" (15.2 x 22.2cm)

#INSTANT-COLD

BOX OF 2451.80 BOX
2 OR MORE BOXES48.20 BOX



"Sweatless" Covering For Comfort

Latexfree

Deluxe Instant Cold Pack

- Squeeze To Activate
- 6.75" x 9" (17.1cm x 22.86cm)

#INSTANT-COLD-DELUXE

BOX OF 2473.90 BOX
2 OR MORE BOXES67.60 BOX



Latexfree

Bilateral Ice Pack

- No Wetness Or Leakage
- Refillable Ice Pack With Velcro Strap

#BILAT-24

CASE OF 24 225.50 CASE
2 OR MORE CASES 215.60 CASE

Sterillium Rub® Surgical Hand Scrub

- Exceeds FDA Requirements For Surgical Hand Antisepsis
- 80% Ethanol For Superior Germ Killing
- Chlorhexidine Gluconate (CHG) Free To Prevent Hand Irritation



Touchless Activation Dispenser

- Mount With Screws (Included) Or Double-Sided Tape
- Uses 4 "D" Batteries (Not Included)

#SCRUB-STERI-DISPENSE176.90 EACH



Sterillium® Rub Surgical Hand Scrub Refills

- Approximately 850 Uses Per Bottle
- Fits Touchless Dispenser On This Page

#SCRUB-STERI-1000. 242.90 PKG OF 2

SPECIAL

Sterillium® Starter Kit

Includes:

- Touchless Activation Dispenser
- 2 1000mL Bottles

**#SCRUB-STERI-KIT
359.00 COMPLETE KIT**

**2 OR MORE STARTER KITS
330.00 COMPLETE KIT**

Scrub Brushes

Iodophor Surgical Scrub Brush/ Sponge

Traditional Hospital Iodine Based Iodophor Solution

#SCRUB-POV-30

BOX OF 30 67.90 BOX
2 BOXES 63.40 BOX



MOST POPULAR

"Clens" Type Surgical Scrub Brush/ Sponge

PCMX "Clens" Solution Is Kinder To Sensitive Hands

#SCRUB-CLEN-30

BOX OF 30 73.60 BOX
2 BOXES 68.40 BOX



Sterile Equipment Covers

**Sterile "Self Stick" MOST POPULAR
Light Handle Covers**

- 20cm x 20cm
- Packaged Sterile In Sets Of 2 (50 Sets Per Case)



#LIGHT-HANDLE-COVER

1 CASE 189.90 CASE
2 OR MORE CASES 180.30 CASE

Sterile Mayo Stand Cover

- Individually Packaged
- Fits Regular Mayo Stands
- 24" x 24" (61 x 61cm)



#MAYO-COVER

CASE OF 24 314.00 CASE
2 OR MORE CASES 304.40 CASE

Latexfree

Sterile Machine Or X-Ray Cover

- Individual Sterile Machine Cover For X-Ray Cassette, Arm Rest, Implant Motor, Etc.



- 24" x 40" (61 x 101.6cm)

#MACHINE-COVER-50

CASE OF 50 214.50 CASE
2 OR MORE CASES 204.40 CASE

Latexfree

Sterile Surgical Table Cover

- Individually Packaged
- Fits Surgical Tables
- 53" x 77" (134 x 195cm)



#SURG-TABLE-COVER

CASE OF 20 109.90 CASE
2 OR MORE CASES 102.40 CASE

Latexfree

Sterilization Wrap For Surgical Kits & Steam Autoclave Tape



18" x 18" Sterilization Wrap

45.7cm x 45.7cm

#INST-WRAP-18 248.30 CASE OF 1000

Latexfree

MOST POPULAR

24" x 24" Sterilization Wrap

61cm x 61cm

#INST-WRAP-24 188.90 CASE OF 500

Latexfree



Steam Autoclave Tape - 1" x 60 YDS

- Changes From Beige To Black To Indicate Proper Sterilization Cycle

#CLAVE-TAPE-1X60 11.60 ROLL
5 OR MORE 9.90 ROLL

Latexfree

Instrument Cleaners



EZ-Zyme Enzyme Instrument Cleaner

- Concentrated - 1 Packet Makes A Gallon
- Dissolves Blood & Tissue
- Instrument Pre-Soak
- Ultrasonic & Evac System Cleaner
- Biodegradable & Neutral PH

#EZ-ZYME

BOX OF 32 PACKETS 96.50 BOX
2 OR MORE 90.50 BOX



Instrument Stain Remover

- Removes Stains & Tarnish
- Will Not Harm Instrument Surfaces

#STAIN 17.90 EACH
2 OR MORE 15.90 EACH

MOST POPULAR

Surgical Instrument Spray Lubricant - 8 oz.



"Instrument Milk" In A Spray

- For Hinged Instruments
- Keeps Scissors Sharp
- Preserves Instrument Life
- Non-Greasy, Anti-Corrosive

#SPRAY LUBE 17.90 EACH
2 OR MORE 15.90 EACH

Stainless Steel Instrument Cassettes



MOST POPULAR

8 Instrument Stainless Steel Cassette


- Silicone Cushions & Locking Crossbar
 - 7.3" x 5.2" x 0.9" (18.5 x 13.3 x 2.3cm)
- #CASS-SS-8 184.60 EACH
 3 OR MORE 175.50 EACH

MOST POPULAR

16 Instrument Stainless Steel Cassette

- Silicone Cushions & Locking Crossbar
 - 10.7" x 7.3" x 0.9" (27.3 x 18.5 x 2.3cm)
- #CASS-SS-16 219.30 EACH
 3 OR MORE 207.50 EACH

Plastic Instrument Cassettes



10" x 6" Cassette

- "Finger Tip" Mat
 - Autoclave, Dry Heat, Chemical Sterilization
 - 25cm x 15.3cm - 0.75" (1.9cm) Height
- #CASSETTE-10X6 135.90 EACH
 3 OR MORE 129.10 EACH



Ideal For Large Instruments

MOST POPULAR

10" x 6" Deep Cassette

- "Finger Tip" Mats
 - Autoclave, Dry Heat, Chemical Sterilization
 - 25cm x 15.3cm - 1.5" (3.7cm) Height
- #CASSETTE-DEEP-BASE 162.90 EACH
 3 OR MORE 152.80 EACH

Sterilization Pouches For Instrument Cassettes



12" x 15" Self-Sealing Instrument Cassette Sterilization Pouches

- Fits All Salvin Instrument Cassettes
 - Fits Other Brand Cassettes Up To 8" x 12" (20.3cm x 30.5cm)
 - External & Internal Indicators With See-Through Packaging For Visibility
- #POUCH-12X15-STER 117.80 BOX OF 200
 2 OR MORE 111.30 BOX OF 200

Latexfree

Salvin “Where Is It?” Small Parts Cassettes

Perfect For Busy Practices – No More Looking For Small Parts!



For Larger Parts



For Smaller Parts



“Where Is It?” Stainless Steel Cassette - Large

- 7.1cm Wide x 5.3cm Deep x 2.3cm Tall


#WHERE-IS-IT-LARGE 47.40 EACH
3 OR MORE 43.70 EACH

MOST POPULAR

“Where Is It?” Stainless Steel Cassette - Small

- 7.1cm Wide x 3.7cm Deep x 1.4cm Tall

#WHERE-IS-IT-SMALL 39.60 EACH
3 OR MORE 36.90 EACH



Includes Removable Silicone Mat



MOST POPULAR

“Where Is It?” - Plastic Cassette

- Removable “Finger Tip” Mat
- 7.3cm Wide x 4cm Deep x 1.9cm Tall

#WHERE-IS-IT-PLASTIC 30.50 EACH
3 OR MORE 28.10 EACH

SPECIAL

Set Of 3 Salvin “Where Is It?” Cassettes

#WHERE-IS-IT-SET-3 **.87.60 SET**

Small Parts Organizing System

- One 10”x 6” Deep Cassette & Nine “Where-Is-It” Plastic Cassettes
- Autoclave, Dry Heat, Chemical Sterilization

#CASS-ORG-KIT **385.00 SET**



Bur Blocks On Page 89

Signature™ Latex Essential Surgical Gloves



Synthetic Polymer Coating For Exceptional Damp Hand Donning

- Beaded Cuff For Easy Donning Without Roll-Back
- Proprietary Coating For Easy Damp/Dry Donning

Glove Guide	
▶	Powder Free
▶	Anatomically Shaped
▶	Sterile
▶	Latex

Size	Product Number	
6	#GL-MSG-5960	
6 1/2	#GL-MSG-5965	
7	#GL-MSG-5970	
7 1/2	#GL-MSG-5975	
8	#GL-MSG-5980	BOX OF 50 PAIRS 129.90 BOX
8 1/2	#GL-MSG-5985	2 OR MORE BOXES 121.60 BOX

SensiCare® - Synthetic Surgical Gloves

MOST POPULAR



Fit & Feel Of Latex But Without The Latex Proteins & Allergens

Pure Aloe Vera Coating Is Applied To Interior Glove Surface. Body Heat Activates The Aloe, Coating & Protecting Your Hands.

- Beaded Cuff For Easy Donning Without Roll-Back
- Unique Styrene-Isoprene-Styrene (SIS) Synthetic
- Proprietary Coating For Easy Damp/Dry Donning
- Bisque Finish For A Better Grip

Glove Guide	
▶	Powder Free
▶	Anatomically Shaped
▶	Sterile
▶	Latex Free Synthetic
▶	Aloe Vera

Size	Product Number	
6	#GL-MDS-377060	
6 1/2	#GL-MDS-377065	
7	#GL-MDS-377070	
7 1/2	#GL-MDS-377075	
8	#GL-MDS-377080	BOX OF 25 PAIRS 132.40 BOX
8 1/2	#GL-MDS-377085	2 OR MORE BOXES 124.10 BOX

SensiCare® PI Micro - Synthetic Surgical Gloves



15% Thinner Than Standard Gloves To Enhance Tactile Sensitivity

- Beaded Cuff For Easy Donning Without Roll-Back
- Proprietary Coating For Easy Damp/Dry Donning

Glove Guide	
▶	Powder Free
▶	Anatomically Shaped
▶	Sterile
▶	Latex Free Synthetic

Size	Product Number	
6	#GL-MSG-9660	
6 1/2	#GL-MSG-9665	
7	#GL-MSG-9670	
7 1/2	#GL-MSG-9675	
8	#GL-MSG-9680	BOX OF 50 PAIRS 229.60 BOX
8 1/2	#GL-MSG-9685	2 OR MORE BOXES 218.10 BOX

See For Yourself! Please Call For Samples Of Any Of Our Gloves

Blood Collection



Basic Blood Collection Setup

Vacutainer® Needle 21 Ga x 1" (2.54cm)
#VACUTAIN-NEEDLE 73.60 BOX OF 100

Vacutainer® Needle Holder
#VACUTAIN-HOLDER. 51.00 BOX OF 50

IV Adapter
#BIO-PRF-IV-CONNECT 204.00 BOX OF 100

Tourniquet-Latex Free
#VACUTAIN-TOURNIQUET 37.10 BOX OF 25

Vacutainer® Tube-Red - 10mL - No Additives
#VACUTAIN-RED-TUBE 147.00 PKG OF 100

Vacutainer® Tube-Yellow - 10mL (8.5mL Draw) Citrate Additive
#VACU-TUBE-YELLOW 169.90 PKG OF 100

10cc Sterile Syringe
#SYRINGE-10CC. 70.90 PKG OF 100

3cc Sterile Syringe
#SYRINGE-STERILE-3CC. 66.10 PKG OF 100

Blunt Catheter Sterile 16 Ga x 3.5" (9cm)
#BLUNT-CATH-16G-9CM 395.00 BOX OF 100

Blunt Catheter Sterile 18 Ga x 1" (2.54cm)
#BLUNT-CATH-18G-2.5CM. 124.90 BOX OF 100

Blood Collection Tube Rack
#VACU-TUBE-RACK 38.60 EACH

3cc Sterile Syringe With 21 Ga x 1.5" (3.81cm) Catheter
#SYRINGE-3CC-21G 62.50 BOX OF 100

Butterfly To Vacutainer® Adapter 21 Ga x .75" (1.9cm)
#BUTTERFLY-21G 105.80 BOX OF 50

Rutkowski "Bio-Compress" Instrument

Closes to 1mm Thickness To Create Buffy Coat Membrane
#BIO-COMPRESS 201.20 EACH

Vacutainer® is a registered trademark of Becton, Dickinson and Company

BIO-PRF Centrifuge




BIO-PRF Horizontal Centrifuge

- Up To 4X More Platelets And Leukocytes
- Higher Growth Factor Release When Compared To Fixed-Angle Centrifugation
- Better Horizontal Layer Separation Between PRF-Clot And Red Blood Cells
- More Even Distribution Of Regenerative Cells Throughout The PRF-Clot
- Less Cell Damage / Death On The Tube Walls Produced On Fixed-Angle Centrifugation
- Less RBC Accumulation

Includes:

- BIO-PRF Centrifuge
- BIO-PRF Box
- Set Of 6 BIO-PRF Instruments
- BIO-PRF Bowl - For Sticky Bone
- BIO-PRF Tray - For Membranes
- BIO-PRF Mixing Bowl
- BIO-PRF 24 Tube Holders
- BIO-PRF Latex Free Tourniquet
- BIO-PRF Cassette
- Box Of 100 BIO-PRF Red Tubes
- Box Of 100 BIO-PRF Blue Tubes
- Box Of 24 BIO-PRF Butterfly Needles

#BIO-PRF-COMP-KIT 3450.00



BIO-PRF Red Tubes

#BIO-PRF-RED-TUBES
.....135.00 BOX OF 100

BIO-PRF Blue Tubes

#BIO-PRF-BLUE-TUBES
.....115.00 BOX OF 100

BIO-HEAT System



BIO-HEAT System To Extend PRF Resorption

- Extend PRF Resorption From 2-3 Weeks To 4-6 Months

Indicated For:

- Barrier Membrane For GBR Procedures
- Covering The Lateral Window During Sinus Lifts
- Root Coverage
- Mixing With Bone Grafting Materials
- Insertion Into Perio Pockets Non-Surgically For Stability

#BIO-HEAT-KIT 2050.00

Centrifuge



Variable Speed Centrifuge

- Adjustable Speed Range Of 500rpm - 4,000rpm
- Compact Size 10.2" x 11.8" x 8" (26cm x 30cm x 20cm)
- Very Quiet Sealed Motor
- 1 Year Parts & Labor Warranty
- 110 Volts

#CENTRIFUGE-VARISPEED..... 1475.00



"Rutkowski" Blood Collection Starter Package

Disposables Package For Approximately 50 Patients Includes:

- Box Of 50 Vacutainer® Needle Holders
- Box Of 25 Latex Free Tourniquets
- 1 Blood Collection Tube Rack
- Pkg Of 100 Red Top Vacutainer® Tubes
- Pkg Of 100 Yellow Top Vacutainer® Tubes
- Pkg Of 100 3cc Sterile Syringes
- Box Of 50 Butterfly To Vacutainer® Adapters
- PRF Box

#BLOOD-COLL-RUTKOWSKI
..... 860.00/PKG FOR APPROX 50 PATIENTS

Disposable Cost Approximately 17.00 Per Patient

PRF Box




PRF Box

- Form Fibrin Membrane Or Plug Shape From Clot
- Weighted Top Creates 1mm Thick Membrane - No Force Needed

#PRF-BOX.....343.00



Form Fibrin Membrane Shape From Clot



Form Fibrin Plug Shape From Clot

Piezotome Cube

Safely Cuts Bone While Preserving Soft Tissue

Scan For
60th Second Product Video™



*Power Adjusts Dynamically For Effortless Cuts
Will Not Damage Soft Tissue*

Automatically
Reduces Power
In Softer Bone

Automatically
Increases Power
In Denser Bone



BUYER'S GUIDE	
▶	Console With Irrigation
▶	LED Piezo Handpiece
▶	Touch Screen Display

PROCEDURE GUIDE	
▶	Atraumatic Extractions
▶	Block Grafting
▶	Osteoplasty
▶	Lateral Sinus Lift
▶	Crestal Sinus Lift
▶	Ridge Splitting
▶	Osteotomies
▶	Crown Lengthening



Call For Package Specials



"There is nothing like the Piezotome for lifting the sinus membrane, harvesting autogenous bone & removing lingual tori while protecting the soft tissue in the floor of the mouth."

-Dr. John Russo
Periodontist, Sarasota, FL



"The use of the Piezotome has been the most revolutionary development in extractions in my career. Most extractions are completed in a matter of minutes with NO loss of alveolar bone."

-Dr. Michael Sonick
Periodontist, Fairfield, CT

Implant Center 2



**Piezotome + Implant Motor
Combination**

Scan For
**60⁺ Second
Product
Video™**



*Motor Fits Your Existing Implant Handpieces
Cuts Bone While Preserving Soft Tissue
LED Handpiece For Excellent Visibility
60 Watts Of Power For Faster Cuts*

BUYER'S GUIDE	
▶ Console With Irrigation	▶ LED Implant Handpiece
▶ LED Piezo Handpiece	▶ 100 to 40,000 RPM
▶ Touch Screen Display	

PROCEDURE GUIDE	
▶ Implant Placement	
▶ Atraumatic Extractions	
▶ Block Grafting	
▶ Osteoplasty	
▶ Lateral Sinus Lift	
▶ Crestal Sinus Lift	
▶ Ridge Splitting	
▶ Osteotomies	
▶ Crown Lengthening	



Call For Package Specials

Piezotome Accessories



Autoclavable Tubing With Pump

Includes 10 Disposable Spikes

#PIEZO-IRR-AC 205.00 EACH

Disposable Irrigation Spikes

Pkg Of 10 Spikes

#PIEZO-IRR-SPIKE 60.00 PKG OF 10



*Salvin Surgical
Cart On Page 46*


Tips Shown By Procedure On Next Two Pages ➔

Atraumatic Extraction Kit *MOST POPULAR*

For Piezotome 2 or Cube

#PIEZO-2-EXTRACTION-KIT

WITH CASSETTE. 1833.00 KIT



LC1 Tip



LC2 Tip



LC2L Tip



LC2R Tip



Ninja Tip


Piezotome 2 or Cube	#PIEZO2-TIP-LC1	#PIEZO2-TIP-LC2	#PIEZO2-TIP-LC2L	#PIEZO2-TIP-LC2R	#PIEZO2-TIP-NINJA
	305.00 EACH	305.00 EACH	305.00 EACH	305.00 EACH	305.00 EACH

Sinus Lift Kit *MOST POPULAR*

For Piezotome 2 or Cube

#PIEZO-2-SL-TIP-KIT

WITH CASSETTE. 1500.00 KIT



SL1 Diamond Tip



SL2 Diamond Tip



SL3 Tip



SL4 Tip



SL5 Tip


Piezotome 2 or Cube	#PIEZO2-TIP-SL1	#PIEZO2-TIP-SL2	#PIEZO2-TIP-SL3	#PIEZO2-TIP-SL4	#PIEZO2-TIP-SL5
	305.00 EACH	305.00 EACH	305.00 EACH	305.00 EACH	305.00 EACH

Bone Surgery Kit

For Piezotome 2 or Cube

#PIEZO-2-BS-TIP-KIT

WITH CASSETTE. 1833.00 KIT



BS1-Slim Tip



BS2L Tip



BS2R Tip



BS4 Tip



BS5 Tip



BS6 Tip


Piezotome 2 or Cube	#PIEZO2-TIP-BS1-SLIM	#PIEZO2-TIP-BS2L	#PIEZO2-TIP-BS2R	#PIEZO2-TIP-BS4	#PIEZO2-TIP-BS5	#PIEZO2-TIP-BS6
	305.00 EACH	305.00 EACH	305.00 EACH	305.00 EACH	305.00 EACH	305.00 EACH

Essential Kit

For Piezotome 2 or Cube

#PIEZO2-ESSENTIAL-KIT

WITH CASSETTE. 1833.00 KIT



BS1-Slim Tip



BS4 Tip



SL1 Diamond Tip



SL2 Diamond Tip



SL3 Tip



LC2 Tip

Piezotome 2 or Cube	#PIEZO2-TIP-BS1-SLIM	#PIEZO2-TIP-BS4	#PIEZO2-TIP-SL1	#PIEZO2-TIP-SL2	#PIEZO2-TIP-SL3	#PIEZO2-TIP-LC2
	305.00 EACH	305.00 EACH	305.00 EACH	305.00 EACH	305.00 EACH	305.00 EACH

Crest Splitting Kit **MOST POPULAR**

For Piezotome 2 or Cube

#PIEZO-2-CS-TIP-KIT

WITH CASSETTE. 1833.00 KIT



CS1 Tip
0.55mm



CS2 Tip
0.85mm



CS3 Tip
0.5mm - For Side Cuts



CS4 Tip
1.8mm



CS5 Tip
2.75mm



CS6 Tip
3.75mm

Piezotome 2 or Cube	#PIEZO2-TIP-CS1	#PIEZO2-TIP-CS2	#PIEZO2-TIP-CS3	#PIEZO2-TIP-CS4	#PIEZO2-TIP-CS5	#PIEZO2-TIP-CS6
	305.00 EACH	305.00 EACH	305.00 EACH	305.00 EACH	305.00 EACH	305.00 EACH

Osteotomy Kit

For Piezotome 2 or Cube

#PIEZO-2-INTRALIFT-KIT

WITH CASSETTE. 1500.00 KIT



TKW1 Tip



TKW2 Tip



TKW3 Tip



TKW4 Tip



TKW5 Tip

Piezotome 2 or Cube	#PIEZO2-TIP-TKW1	#PIEZO2-TIP-TKW2	#PIEZO2-TIP-TKW3	#PIEZO2-TIP-TKW4	#PIEZO2-TIP-TKW5
	305.00 EACH	305.00 EACH	305.00 EACH	305.00 EACH	305.00 EACH

Crown Lengthening Kit

For Piezotome 2 or Cube

#PIEZO-2-CE-TIP-KIT

WITH CASSETTE. 1225.00 KIT



CE1 Tip



CE2 Tip



CE3 Tip



BS6 Tip

Piezotome 2 or Cube	#PIEZO2-TIP-CE1	#PIEZO2-TIP-CE2	#PIEZO2-TIP-CE3	#PIEZO2-TIP-BS6
	305.00 EACH	305.00 EACH	305.00 EACH	305.00 EACH

**For Accelerated
Orthodontic
Bone Cuts**



Rounded Piezo Tip

#PIEZO2-TIP-BS1-ROUND
305.00 EACH

Newtron® P5 XS BLED Ultrasonic Scaler

- Perform Faster, More Comfortable Treatments With 80+ Unique Tips
- Auto-Responsive Handpiece & Precision Tips Ensure Optimal Protection For Teeth & Soft Tissue
- Smooth, Linear Vibrations Reduce Gingival Impact & Preserve Cementum

NEW!
To Our Catalog

SPECIAL

Newtron® P5 XS BLED Package With Intro Scaling Set

Includes:

- Newtron® P5 XS BLED Unit
- One LED Handpiece
- Newtron® Intro Scaling Kit With 3 Tips
- Autoclavable Cassette

#NEWTRON-F62117

Call For Lowest Price



Newtron® Intro Scaling Kit

- Includes: Tips – H3, H4L, H4R With 3 Blue Wrenches
- Autoclavable Cassette

#NEWTRON-INTRO-KIT

..... **442.00**



Implant Protect Tip Kit

- Includes: Tips – IP1, IP2L, IP2R, IP3L, IP3R With Universal Wrench
- Autoclavable Cassette

#PIEZO-F02120 **736.00**

Air-N-Go® Air Polishing



Air-N-Go® Easy

- Versatile Nozzles & Tips Provide Effective Plaque, Stain, And Implant Cleaning In Both Supragingival & Subgingival Applications
- Gentle, Comfortable Treatments: Supra And Perio Powders Offer Pain-Free Cleaning
- Easy Setup & Use: Quick Start-Up – Connects Directly To High-Speed Handpieces

Call For Lowest Price



Air-N-Go® Classic Air Polishing Powder Variety Pack

- Box Of Four 250gm Bottles
- Variety Pack Contains One Bottle of Each Flavor: Cola, Raspberry, Peppermint, Lemon
- Also Available Per Flavor In Packs Of Four 250gm Bottles

#AIR-N-GO-F10255 **142.00 BOX OF 4**

**SALVIN®
RotoBrush-Titanium®
For Exposed Implant Threads**

- Cleans Infected Implant Surface Without Grinding The Threads
- Maintains Structure & Shape Of The Implant
- Gentler On Implant Surface Than Diamond Burs & Ultrasonics
- Start At 300rpm Working Up To 1,200rpm As Needed



MOST POPULAR

RotoBrush-Titanium® For Vertical Access

Single Use
#ROTOBRUSH-POCKET 41.60 EACH
3 OR MORE 36.70 EACH



RotoBrush-Titanium® For Posterior Access

Single Use
#ROTOBRUSH-SHORT-POCKET 41.60 EACH
3 OR MORE 36.70 EACH



RotoBrush-Titanium® For Wide Area

Single Use
#ROTOBRUSH-OMEGA 41.60 EACH
3 OR MORE 36.70 EACH



RotoBrush-Titanium® Short For Wide Area

Single Use
#ROTOBRUSH-DELTA-S 41.60 EACH
3 OR MORE 36.70 EACH



MOST POPULAR

RotoBrush-Titanium® For Non-Flap Case

Single Use
#ROTOBRUSH-NANO 41.60 EACH
3 OR MORE 36.70 EACH



MOST POPULAR


RotoBrush-Titanium® For Flap Open Case

Single Use
#ROTOBRUSH-2D-LONG 41.60 EACH
3 OR MORE 36.70 EACH

Scan For
60 Second
Product
Video™



*Fits Latch Type
Implant Handpiece*



*Gently Cleans
Implant Surface*

SPECIAL

Salvin RotoBrush-Titanium® Kit



MOST POPULAR

Includes:

- 2 RotoBrush-Titanium® Pocket
- 2 RotoBrush-Titanium® Nano
- 2 RotoBrush-Titanium® 2D-Long
- **FREE** Autoclavable Bur Block

#ROTOBRUSH-KIT
209.00 KIT WITH BUR BLOCK

Peri-Implantitis Tungsten Carbide Bur Kit

SPECIAL

Peri-Implantitis Bur Kit

- Tungsten Carbide Burs Clean & Smooth Titanium Surfaces
- 32mm Long FG Shanks To Reach Difficult To Access Implants
- Different Shapes & Sizes For Variety Of Individual Cases
- Includes: 4 Egg Shaped Burs & 4 Flame Shaped Burs
- Fits Friction Grip Handpieces

#PERIIMPLANTITIS-BURS **143.15 KIT W/BLOCK**



PERI-IMPLANTITIS NON-SURGICAL TREATMENT – IMPLA-SAFE PROTOCOL

1 Probe The Sulcus & Assess For Calculus Using Impla-Safe Probe To Avoid Tissue Seal Rupture



If Calculus Present



Remove Calculus Using Impla-Safe Soft Curette (Only On Abutments – Never On Implant Threads)

2 Use Impla-Safe Cement Explorer For Detection Of Residual Cement



If Residual Cement Is Present



Use Impla-Safe Cement Remover To Chip Off Cement

3 Subgingival Irrigation With Saline

4 Shoe-Shine Implant Surface – Use Impla-Safe Ti Micropliers To Safely Push Gauze Soaked In Saline In The Sulcus



“Impla-Safe is the first fully implant-specific instrument set for peri-implant maintenance, and for the prevention and non-surgical management of peri-implant inflammation. From the flexible probe to the soft titanium alloy instruments for implant-specific factors, such as excess cement, this set addresses the iatrogenic damage caused to implants by dental instruments and scaling using an evidence-based workflow. Should be integral in the maintenance of patients with dental implants.”

Dr. George Kotsakis
Periodontist, Hoboken, NJ



“These exceptionally well crafted implant site specific instruments for detecting and removing excess cement are a brilliant innovation and must have for all implant surgically based clinicians.”

Dr. Michael A. Pikos
Oral Surgeon, Palm Harbor, FL



“The Impla-Safe set is a game-changing, revolutionary implant-specific solution!”

Dr. Jennifer Hirsch Doobrow
Periodontist, Cullman, AL



“I’m excited to see an implant specific solution with instruments that adapt to the shape of implant-supported restorations and can help detect excess cement during and after cementation.”

Dr. Konstantinos Chochlidakis
Prosthodontist, Rochester, NY



“The Impla-Safe Protocol and instrument set addresses what has been missing in peri-implantitis treatment: an evidence-based, implant-specific solution for our patients.”


Dr. Eduardo R. Lorenzana
Periodontist, San Antonio, TX

SPECIAL

Impla-Safe™ Peri-Implant Management Set

INCLUDES:

- Impla-Safe™ Probe
 - Impla-Safe™ Mirror Handle
 - Box of 6 Impla-Safe™ Mirrors
 - Impla-Safe™ Titanium Curette
 - Impla-Safe™ Cement Explorer™
 - Impla-Safe™ Cement Remover
 - Impla-Safe™ Titanium Micropliers
 - Instrument Cassette
- #IMPLA-SAFE-INSTR-SET**
.....**694.00 SET**





*For Prevention And Non-Surgical Management
Of Peri-Implantitis*

Peri-Implant Management Set

Designed By Dr. George Kotsakis



Depth Marked
At 1-15mm

**Impla-Safe™
Perio Probe**

Titanium
Flexible
Ultra Thin
#IMPLA-SAFE-PROBE
.....99.00 EACH



**Impla-Safe™
Mirror #4**

Double Sided
#IMPLA-SAFE-MIRROR
..... 99.00 BOX OF 6

**Impla-Safe™
Mirror Handle**

#MIRROR-HNDL-BLK
.....65.00 EACH



**Impla-Safe™
Soft Tissue Curette 3/4**

Titanium
For Supragingival Calculus
#IMPLA-SAFE-SOFT-TI-CU
.....99.00 EACH



**Impla-Safe™
Cement Explorer™**

Titanium
High Modulus Of Elasticity For
Optimal Tactile Sensitivity
#IMPLA-SAFE-EXPLORER
.....99.00 EACH



**Impla-Safe™
Cement Remover**

Titanium
Micro-Serrated For Efficient
Cement Removal
#IMPLA-SAFE-CEM-REMOVER
.....99.00 EACH



**Impla-Safe™
Titanium Micropliers -
Curved**

Titanium
#IMPLA-SAFE-MICROPLIERS
.....129.00 EACH

Salvin Motor Systems For Implantology & Oral Surgery



Single Motor System



BUYER'S GUIDE - IMPLANT MOTORS

▶	One Motor Unit
▶	300 to 40,000 RPM Brushless Motor
▶	Pre-Selectable Torque Levels 10 -- 70Ncm
▶	Memory Backup For Program Changes
▶	Autoclavable Motor & Cord
▶	Internal & External Irrigation Pump

SPECIAL

Single Motor – 1 Handpiece Package

Includes:
Single Motor Console & 1 Motor, 1 Handpiece

MOST POPULAR

With 20:1 Handpiece (25 To 2,000 RPM)
#MOTOR-SINGLE-20-1..... **5190.00**

With 32:1 Handpiece (15 To 1,250 RPM)
#MOTOR-SINGLE-32-1..... **5190.00**

With 16:1 Handpiece (30 to 2,500 RPM)
#MOTOR-SINGLE-16-1..... **5190.00**

Dual Motor System



BUYER'S GUIDE - IMPLANT MOTORS

▶	Two Motor Unit
▶	300 to 50,000 RPM Brushless Motor
▶	Pre-Selectable Torque Levels 10 -- 70Ncm
▶	Memory Backup For Program Changes
▶	Autoclavable Motor & Cord
▶	Internal & External Irrigation Pump

SPECIAL

Dual Motor – 1 Handpiece Package

Includes:
Dual Motor Console & 2 Motor, 1 Handpiece

With 20:1 Handpiece (25 To 2,500 RPM)
#MOTOR-DUAL-20-1..... **6490.00**

With 32:1 Handpiece (15 To 1,550 RPM)
#MOTOR-DUAL-32-1..... **6490.00**

With 16:1 Handpiece (30 to 3,125 RPM)
#MOTOR-DUAL-16-1..... **6490.00**

Irrigation Tubing On Page 80

Implant & Oral Surgery Handpieces



32:1 Reduction Push Button Latch Contra Angle

**FOR DRILLING, TAPPING & PLACEMENT
RECOMMENDED SPEED FOR MOST IMPLANT SYSTEMS**

- 10 To 1,000 RPM For Drilling, Tapping & Placement
- Internal & External Irrigation
- Autoclave Or Chemiclave Up To 275° F (135° C)

#NOU-32-1 **1099.00 EACH**



16:1 Reduction Push Button Latch Contra Angle

**ORDER THIS HANDPIECE FOR IMPLANT SYSTEMS
REQUIRING DRILLING SPEEDS EXCEEDING 1,250 RPM**

- 20 To 2,400 RPM For Drilling & Placement
- Internal & External Irrigation
- Autoclave Or Chemiclave Up To 275° F (135° C)

#NOU-16-1 **1099.00 EACH**



MOST POPULAR

20:1 Reduction Push Button Latch Contra Angle

- 15 To 1,700 RPM
- Internal & External Irrigation
- Autoclave Or Chemiclave Up To 275° F (135° C)
- Works With X-Guide

#NOU-20-1 **1195.00 EACH**



MOST POPULAR

1:1 40,000 RPM Straight Surgical

For Bone & Jaw Surgery

- 300 To 40,000 RPM
- External Irrigation
- Autoclave Or Chemiclave Up To 275° F (135° C)

#NOU-1-1-STRAIGHT **1085.00 EACH**



1:1 40,000 RPM Angled Surgical

For Bone & Jaw Surgery

- 300 To 40,000 RPM
- External Irrigation
- Autoclave Or Chemiclave Up To 275° F (135° C)

#NOU-1-1-ANGLED **1195.00 EACH**



1:1 Push Button Latch Contra Angle

- 300 To 40,000 RPM
- Internal & External Irrigation
- Autoclave Or Chemiclave Up To 275° F (135° C)

#NOU-1-1-LATCH **799.00 EACH**




1:2 Straight Surgical (Multiplies RPM X 2)

For Bone & Jaw Surgery

- 1,000 To 80,000 RPM
- External Irrigation
- Autoclave Or Chemiclave Up To 275° F (135° C)

#NSK-1-2-HANDPIECE **1445.00 EACH**

**Handpieces On This Page Fit All
Standard "E" Type Implant Motors**



Disposable Irrigation Tubing Fits Nouvag & Salvin Implant Motor Systems

- Includes Tubing From Solution Bag To Handpiece
- Works With Internal, External Irrigation, Or Both
- Packaged Sterile For One-Time Use

#NOU-IRR-TUBE-DISP198.60 PKG OF 10
 2 OR MORE183.10 PKG OF 10



Nouvag & Salvin Motor Tubing Irrigation Needle

#NOU-IRR-NEEDLE.....54.00 EACH



Fits NOUVAG & SALVIN Single & Dual Motors

Disposable Irrigation Tubing Fits W&H Implant Motor Systems

- Includes Tubing From Solution Bag To Handpiece
- Works With Internal, External Irrigation, Or Both
- Packaged Sterile For One-Time Use


#WH-IRR-TUBE-DISP198.60 PKG OF 10
 2 OR MORE183.10 PKG OF 10



Fits W&H Motors Except Elcomed 100 Or 200 System

Motor & Handpiece Manufacturer's Names Are Trademarks Of Their Respective Companies.

Bone Saws Fit Your Existing Implant Motor



Saws Fit "E" Type Implant Motors

- Easily Cut Channels & Bone Blocks
- Very Useful For Ridge Osteoplasty
- Universal "E" Type Coupling Fits All Implant Motors
- Includes External Irrigation & Tubing
- Maximum 15,000 Oscillations Per Minute



Sagittal Bone Saw

Easily Harvests Blocks From The Ramus Or Mandible

#SAG-SAW 1696.00 EACH
(Blades Sold Separately)



MOST POPULAR



4mm Blade-Medium Shank

#SAG-BLAD-4-M

96.50 PKG OF 6



6mm Round Blade-Medium Shank

#SAG-BLAD-6R-M

96.50 PKG OF 6



6mm Blade-Medium Shank

#SAG-BLAD-6-M

96.50 PKG OF 6



6mm Blade-Long Shank

#SAG-BLAD-6-L

96.50 PKG OF 6



10mm Blade-Long Shank

#SAG-BLAD-10-L

96.50 PKG OF 6

(Saw Blades Shown Actual Size)

Oscillating Bone Saw

#OSCIL-SAW 1696.00 EACH
(Blades Sold Separately)



MOST POPULAR



6mm Blade

#OSCIL-BLAD-6

82.40 EACH



10mm Blade

#OSCIL-BLAD-10

82.40 EACH



5mm Blade


#OSCIL-BLAD-5

82.40 EACH

(Saw Blades Shown Actual Size)

Reciprocating Bone Saw

MOST POPULAR
#RECIP-SAW 1696.00 EACH
(Blades Sold Separately)



MOST POPULAR



11mm Cutting Length Blade

#RECIP-BLAD-10

96.50 PKG OF 6



14mm Cutting Length Blade

#RECIP-BLAD-15

96.50 PKG OF 6



17mm Cutting Length Blade

#RECIP-BLAD-17

96.50 PKG OF 6



22mm Cutting Length Blade

#RECIP-BLAD-21


96.50 PKG OF 6

(Saw Blades Shown Actual Size)

Salvin Universal Drill Stop Kit



*Universal Drill Stop Kit Fits
The Implant Drills You Already Have!*



*With The Mount On The Drill,
Measure From The Tip Of The
Drill To The Mount*



*Select & Mount
The Appropriate Length
Drill Stop*



*Measure With Drill Stop
Attached To Confirm
Desired Drill Depth*



*Includes Autoclavable Block,
13 Stops & 2 Mounts*

BUYER'S GUIDE




For Standard Drills Without Collars
#DRILL-STOP-KIT 679.00 KIT
(Fits Drills Less Than 4.0mm)



For Drills With Wide Collars
#DRILL-STOP-WIDE-KIT 679.00 KIT
(Fits Drills Less Than 4.0mm)

Salvin Tri-Spade Intermediate Implant Drills Save You Money



2.00mm Implant Drill #IMPL-DRILL-2.00 99.60 EACH



2.75mm Implant Drill #IMPL-DRILL-2.75 99.60 EACH



3.30mm Implant Drill #IMPL-DRILL-3.30 99.60 EACH



3.50mm Implant Drill #IMPL-DRILL-3.50 99.60 EACH

SPECIAL 4 Or More Intermediate Drills (YOUR ASSORTMENT) **.89.90 EACH**

**Salvin
POWER POINT® Pilot Drills**

Ultra Sharp For Precise Re-Angulation In Any Bone Density

- Easy To Read Depth Markings
- Chatter-Free




MOST POPULAR

**Salvin POWER POINT®
Ultra Sharp Implant Pilot Drills**

1.3mm Diameter - Overall Length 33mm

#PILOT-DRILL 155.20 SET OF 2
2 OR MORE 144.30 SET OF 2




MOST POPULAR

**Salvin POWER POINT®
Ultra Sharp Implant Pilot Drills-Long**

1.3mm Diameter - Overall Length 43mm

#PILOT-DRILL-LONG 155.20 SET OF 2
2 OR MORE 144.30 SET OF 2



Drill Extender & Pilot Drills



**Implant Drill Extender -
Extends 16mm**

- Fits Latch Type Implant Drills
- #DRILL-EXTENDER . . 77.70 EACH



2mm Trepine Pilot Drill

- Harvest Bone As You Drill Primary Osteotomy
- #TREP-2 115.50 EACH



MOST POPULAR

Lance Pilot Drill

- 27mm Length
- #IMPL-DRILL-LANCE . . 60.90 EACH
2 OR MORE 55.40 EACH



Implant Pilot Drill - Round

- 22mm Length #6 Round
- #IMPL-DRILL-ROUND . . 18.90 EACH
5 OR MORE 16.50 EACH

SPECIAL

"Hole In One" Pilot Drill Kit

MOST POPULAR

Includes:

- 2 Salvin **POWER POINT®** Pilot Drills - 43mm
- 2 Salvin **POWER POINT®** Pilot Drills - 33mm
- 1 Lindemann Side Cutting Bur - 27mm
- 1 Round Implant Pilot Drill
- 1 Stainless Steel Bur Block

#HOLE-IN-ONE-KIT
. **399.00 KIT W/BLOCK**



This Kit Fits Latch-Type Handpieces

Meisinger Diamond Saw Discs



7mm Diameter Diamond Saw Disc

- 0.3mm Disc Thickness
 - Fits Latch Type Handpieces
- #SAW-DISC-231DC 79.75 EACH
3 OR MORE 76.80 EACH



10mm Diameter Diamond Saw Disc

- 0.3mm Disc Thickness
 - Fits Latch Type Handpieces
- #SAW-DISC-231DC/100 83.15 EACH
3 OR MORE 79.80 EACH



13mm Diameter Diamond Saw Disc

- 0.3mm Disc Thickness
 - Fits Latch Type Handpieces
- #SAW-DISC-231DC/130 96.60 EACH
3 OR MORE 90.80 EACH



Saw Disc Kit For Latch Type Handpieces

- 7mm Saw Disc
 - 10mm Saw Discs (2)
 - 13mm Saw Disc
- #SAW-DISC-KIT 306.35 KIT

MOST POPULAR

Oral Surgery Burs



Lindemann Side Cutting Bur-SU

- Fits Friction Grip Handpieces
- 11mm Head Length

#LINDEMANN-SU
 INDIVIDUALLY82.80 EACH
 2 OR MORE77.10 EACH



Lindemann Side Cutting Bur Non-Irrigated - 27mm Long

MOST POPULAR

#LINDEMANN
 INDIVIDUALLY82.80 EACH
 2 OR MORE77.10 EACH



Lindemann Side Cutting Bur Non-Irrigated - 32mm Long

MOST POPULAR

#LINDEMANN-LONG
 INDIVIDUALLY82.80 EACH
 2 OR MORE77.10 EACH



Lindemann Side Cutting Bur-HP

- Fits Straight Handpieces
- 44mm Long

#LINDEMANN-HP
 INDIVIDUALLY82.80 EACH
 2 OR MORE77.10 EACH



700 XXXL Implant Carbides

- Fits Friction Grip Handpiece
- 30mm Long

#700 XXXL
 SET OF 5 88.90 SET OF 5
 2 OR MORE SETS . . . 82.20 SET OF 5



700 XXL Implant Carbides

MOST POPULAR

- Fits Friction Grip Handpiece
- 25mm Long

#700 XXL
 INDIVIDUALLY16.60 EACH
 12 OR MORE14.50 EACH



Pikos Scoring Bur

MOST POPULAR

- #8 Round Scoring Bur
- Depth Marked At 15mm
- Fits Straight Handpieces
- 44mm Long

#PIKOS-SCORINGBUR
 INDIVIDUALLY7.50 EACH
 10 OR MORE5.90 EACH



Pikos Osteotomy Bur

MOST POPULAR

- 702-L Fits Straight Handpieces
- 44mm Long

#PIKOS-OSTEO-BUR-S
 INDIVIDUALLY7.50 EACH
 10 OR MORE5.90 EACH

- 703 Fits Straight Handpieces
- 44mm Long

#PIKOS-OSTEO-BUR-L
 INDIVIDUALLY7.50 EACH
 10 OR MORE5.90 EACH



Zygomatic Trough Bur

- 2.3mm Pear Shape
- Fits Straight Handpieces
- 44mm Long

#MEIS-HM77MF-023-HP
 INDIVIDUALLY 29.90 EACH
 10 OR MORE 28.35 EACH

SPECIAL


Pikos Block Graft Kit

MOST POPULAR

Includes:

- 2 Pikos 702-L Osteotomy Burs
- 2 Pikos #8 Round Scoring Burs
- 1 Pikos Large Round Block Contouring Bur
- 1 0.8mm Diameter Perforation Drill
- 1 Stainless Steel Bur Block

#PIKOS-BUR-KIT 153.20 KIT W/BLOCK



Special Purpose Drills & Burs



Pikos Block Contouring Bur

MOST POPULAR

- Fits Latch Type Handpieces
 - 25mm Long
- #PIKOS-BLOCK-CONTOUR-RA**
INDIVIDUALLY **49.90 EACH**
2 OR MORE **44.70 EACH**
- Fits Straight Handpieces
 - 44mm Long
- #PIKOS-BLOCK-CONTOUR**
INDIVIDUALLY **49.90 EACH**
2 OR MORE **44.70 EACH**



0.8mm Bone Perforation Decorticating Drill

- Fits Latch Type Handpieces
 - 27mm Long
- #PERF-DRILL-0.8-RA**
INDIVIDUALLY **88.30 EACH**
3 OR MORE **78.40 EACH**
- Fits Straight Handpieces
 - 43mm Long
- #PERF-DRILL-0.8-HP**
INDIVIDUALLY **88.30 EACH**
3 OR MORE **78.40 EACH**



Sinus Scoring #8 Diamond

MOST POPULAR

- Fits Latch Type Handpieces
 - 22mm Long
- #SINUS-SCORE-DIAM-RA**
INDIVIDUALLY **18.75 EACH**
10 OR MORE **16.40 EACH**
- Fits Straight Handpieces
 - 44mm Long
- #SINUS-SCORE-DIAM-HP**
INDIVIDUALLY **18.75 EACH**
10 OR MORE **16.40 EACH**




Sectioning Bur

- 703L Surgical Length Friction Grip Bur
 - 25mm Long
- #SECTION-BUR-703FG**
INDIVIDUALLY **7.35 EACH**
2 OR MORE **6.30 EACH**



Ridge Contouring Bone Bur

- Contour & Reshape Bone Anatomy
 - Fits Latch Type Handpiece
 - 25mm Long
- #RIDGE-CONTOUR**
INDIVIDUALLY **139.90 EACH**
2 OR MORE **132.50 EACH**



1.2mm Step Drill Urban Decorticating Bur

- Fits Latch Type Handpieces
 - 21mm Long
- #MEIS-2035-012-RA**
PACK OF 2 **103.90 PKG/2**
2 OR MORE PACKS . . . **94.60 PKG/2**



Neumeier Root Planing Bur-FG

- Fits Friction Grip Handpieces
 - 20mm Long
- #NEUMEYER-FG**
INDIVIDUALLY **17.65 EACH**
10 OR MORE **15.50 EACH**



Wire Passing Bur

- Fits Straight Handpieces
 - 1.5mm Diameter
 - 5cm (2") Long
- #WIRE-PASS-1.5**
INDIVIDUALLY **105.50 EACH**
2 OR MORE **97.90 EACH**



MetalKut™ Carbides

- Our MetalKut™ Carbides Cut Through Metal Like Butter!
 - Fits Friction Grip Handpiece
 - 20mm Long
- #MKUT-PLUS**
INDIVIDUALLY **13.25 EACH**
12 OR MORE **11.30 EACH**

SPECIAL

Rasmussen Biologic Shaping Bur Kit

Includes:

- 2 1.6mm Coarse Flat End Taper Friction Grip Burs
- 2 1.6mm Fine Flat End Taper Friction Grip Burs
- 1 1.0mm Round End Needle Tip Diamond
- 1 2.3mm Round Friction Grip Diamond
- 1 2.3mm Fine Football Shaped Friction Grip Bur
- 1 Stainless Steel Bur Block

#RASMUSSEN-BUR-KIT 108.70 KIT W/BLOCK



SPECIAL

Peri-Implantitis Bur Kit

- Tungsten Carbide Burs Clean & Smooth Titanium Surfaces
- 32mm Long FG Shanks To Reach Difficult To Access Implants
- Different Shapes & Sizes For Variety Of Individual Cases
- Fits Friction Grip Handpieces

Includes:

- 4 Egg Shaped Burs
- 4 Flame Shaped Burs

#PERIIMPLANTITIS-BURS **143.15 KIT W/BLOCK**



SPECIAL

P.E.T. Partial Extraction Therapy Bur Kit - FG

FG Shanks

Fits Friction Grip Handpieces

Includes:

- Short 1.2mm Diameter Bone Cutter
- Long 1.2mm Diameter Bone Cutter
- Coarse 1.0mm Round Diamond Bur
- 1.8mm Round Prep Bur
- 2.3mm Round Prep Bur
- 2.9mm Round Prep Bur
- Stainless Steel Bur Block

#PET-BUR-KIT **143.15 KIT W/BLOCK**



SPECIAL

P.E.T. Partial Extraction Therapy Oral Surgery Bur Kit - HP

HP Shanks

Fits Straight Handpieces

Includes:

- 701L Bur HP
- 1703 Bur HP
- 1.8mm Round Diamond HP
- 3.1mm Round Bur HP
- Lindemann Bur HP
- 2.3mm Pointed Football Diamond HP
- Stainless Steel Bur Block

#PET-OS-BUR-KIT **143.15 KIT W/BLOCK**



SPECIAL

Degranulation Bur Kit

- Multiuse Burs (Typically 100 Procedures)
- Designed For Quick & Effective Bone Debridement
- Removes Remnants Of Tissue Attached To Bone
- Fits Latch Type Implant Handpieces

Includes:

- 1.0mm Degranulation Bur
- 1.8mm Degranulation Bur
- 2.5mm Degranulation Bur
- 3.0mm Degranulation Bur
- 3.5mm Degranulation Bur



#DEGRANULATION-KIT **235.00 KIT W/BLOCK**

SPECIAL

Alveoplasty Bur Kit

- Smooth & Re-Contour Bone Following Tooth Extraction To Provide Ideal Platform For Denture Preparation Or Full Arch Restoration
- Prepare Bone For Implant Placement

Includes:

- 2 Lindemann Bone Cutters (Fits Friction Grip Handpieces)
- 2 Initial Starter Drills (Fits Friction Grip Handpieces)
- 1 Round Diamond (Fits Straight Handpieces)
- 5 Carbide Alveoplasty Burs With Round Ends To Protect Tissue (Fits Straight Handpieces)



#MEIS-ALV18 **289.25 KIT W/BLOCK**

SPECIAL

Crown Removal Bur Kit

- For Cutting Crowns & Bridges
- Includes Burs For Every Type Of Restorative Material
- Recommended Cutting Angle Of 45°
- Optimum Speed Of 160,000rpm
- Fits Friction Grip Handpieces

#MEIS-BC04 **145.00 KIT W/BLOCK**



**SALVIN®
RotoBrush-Titanium®**
For Exposed Implant Threads
On Page 75









Salvin Long Trepine Drills

For Failed Implant Retrieval & Bone Harvesting



Internally Irrigated

-  2.0mm Inside Diameter x 22mm Barrel
#TREP-H-2 115.50 EACH
-  3.5mm Inside Diameter x 22mm Barrel
#TREP-H-3.5 115.50 EACH
-  4.25mm Inside Diameter x 22mm Barrel
#TREP-H-4 115.50 EACH
-  5.0mm Inside Diameter x 22mm Barrel
#TREP-H-5 115.50 EACH
-  6.0mm Inside Diameter x 22mm Barrel
#TREP-H-6 115.50 EACH
-  8.0mm Inside Diameter x 22mm Barrel
#TREP-H-8 115.50 EACH

SPECIAL

Long Trepine Kit **MOST POPULAR**



Internally Irrigated

Set Of 6 Long Trephines With **FREE** Bur Block

- Fits Latch Type Handpieces
- Marked At 8, 10, 13 & 15mm

#TREP-H-KIT
..... 579.00 KIT W/BLOCK

Salvin Short Trepine Drills

For Bone Harvesting



-  2.0mm Inside Diameter x 10mm Barrel
#TREP-SHORT-2 115.50 EACH
-  4.0mm Inside Diameter x 10mm Barrel
#TREP-SHORT-4 115.50 EACH
-  6.0mm Inside Diameter x 10mm Barrel
#TREP-SHORT-6 115.50 EACH
-  8.0mm Inside Diameter x 10mm Barrel
#TREP-SHORT-8 115.50 EACH
-  10.0mm Inside Diameter x 10mm Barrel
#TREP-SHORT-10 115.50 EACH
-  12.0mm Inside Diameter x 10mm Barrel
#TREP-SHORT-12 115.50 EACH

SPECIAL

Short Trepine Kit



Set Of 6 Short Trephines With **FREE** Bur Block

- Fits Latch Type Handpieces
- Marked At 2, 4, 6, 8 & 10mm

#TREP-SHORT-KIT
..... 579.00 KIT W/BLOCK

Autoclavable Bur Blocks & Baskets



MOST POPULAR

6 Bur Block - Latch & Friction Grip

#BUR-BLOCK-6-RA.....	58.00 EACH
3 OR MORE	53.90 EACH



MOST POPULAR

6 Bur Block - Straight Handpiece

BUR-BLOCK-6-HP.....	58.00 EACH
3 OR MORE	53.90 EACH



8 Bur Block - Latch & Friction Grip

#BUR-BLOCK-8-RA.....	58.00 EACH
3 OR MORE	53.90 EACH



12 Bur Block - Latch & Friction Grip

#BUR-BLOCK-12-RA.....	58.00 EACH
3 OR MORE	53.90 EACH



Autoclavable Bur/Parts Basket

46mm Diameter


#BUR-BASKET	72.60 EACH
3 OR MORE	69.10 EACH



12 Compartment Autoclavable Parts Cassette

21cm x 6cm x 2cm

#CASS-SS-PARTS-12.....	99.90 EACH
3 OR MORE	93.80 EACH




Instrument Cassettes On Page 64-65


SALVIN
AutoTome®-Osteotome
Automatic Osteotome Kit



- Automatic Osteotome For Crestal Sinus Lifts & Maxillary Bone Spreading
- Ensures Precision And Saves Time
- Enhanced Patient Safety And Comfort
- One-Hand Operation For Visibility & Ease Of Use
- Creates An Osteotomy For Most Major Implant Systems



Atraumatic, Controlled Bone Preparation



Gentle Auto Tapping Easily Advances Osteotomy While Spreading Bone




Fits Your Existing Implant Motor



Set Your Motor To 12,000-15,000 RPM Translates To 120-150 Gentle Taps Per Minute (Max 15,000 RPM)

SPECIAL

AutoTome® Osteotome Starter Kit With Offset Tips




Includes:

- 1 AutoTome® Handpiece
- 1 Manual Handpiece For Use With Mallet When Desired
- 8 Osteotome Tips (4 Offset Concave, 4 Offset Convex) (2.0-2.8mm, 2.5-3.3mm, 3.0-3.9mm, 3.5-4.5mm Diameters)
- 1 Sterilization Cassette

#AUTOTOME-OST-STARTER-KIT **2290.00 KIT**

SPECIAL

AutoTome® Combination Osteotome & Extraction Kit



Includes:

- 1 AutoTome® Handpiece
- 1 Manual Handpiece
- 16 Osteotome Tips
- Sterilization Cassette
- 6 Periotome Tips
- 1 Autoclavable Tip Block


#AUTOTOME-COMP-KIT **3690.00 KIT**




Conventional Osteotomes On Pages 138-139 & 145

SALVIN
AutoTome®-Extract
Automatic Extraction Kit

- Automatic Periotome For Root Elevation & Extraction
- Ensures Precision And Saves Time
- Enhanced Patient Safety And Comfort
- One-Hand Operation For Visibility & Ease Of Use
- Gentle Tapping Facilitates Atraumatic Controlled Extractions



Bone Preserving
Atraumatic Extraction



Gentle Auto Tapping
Easily Facilitates
Controlled Extractions



Fits Your
Existing
Implant Motor

SPECIAL



AutoTome® Extraction Kit With Handpiece

Includes:

- 1 AutoTome® Handpiece
- 6 Periotome Tips
- 1 Autoclavable Tip Block

#AUTOTOME-EXTRACT. **1790.00 KIT**

SPECIAL



AutoTome® Extraction Tip Kit

(If You Already Have An AutoTome® Handpiece)


Includes:

- 6 Periotome Tips
- 1 Autoclavable Tip Block

#AUTOTOME-EXTRAC-TIPS. **850.00 KIT**

Replacement AutoTome® Extraction Tips

Straight Narrow - 2.1mm	#AUTOTOME-EXT-TIP-1	144.00 EACH
Straight Wide - 2.6mm	#AUTOTOME-EXT-TIP-2	144.00 EACH
Angled Vertical - 2.1mm	#AUTOTOME-EXT-TIP-3	144.00 EACH
Angled Horizontal - 2.1mm	#AUTOTOME-EXT-TIP-4	144.00 EACH
Straight Narrow Concave - 2mm	#AUTOTOME-EXT-TIP-5	144.00 EACH
Straight Medium Concave - 3mm	#AUTOTOME-EXT-TIP-6	144.00 EACH



Conventional Periotomes On Pages 102-105

“Rear-Exhaust” 45° High Speed Surgical Handpiece

Creates Maximum Access & Visibility On 3rd Molars

- No Air, Just Water In Surgical Field - Air Exits At Rear Of Handpiece For Safety
- High Speed: 380,000 - 450,000rpm
- High Torque: 21W Of Cutting Power



Push-Button Bur Release Uses Standard Friction Grip Burs



45° High Speed Handpiece With Swivel 4 Hole Connector

New Version With Higher Torque

#SURG-45-PKG. 1099.00 EACH

Handpiece Without Connector

New Version With Higher Torque

#SURG-45. 899.00 EACH



💡 Fiber Optic 45° High Speed Handpiece With Swivel 5 Hole Fiber Optic Connector

New Version With Higher Torque

#SURG-45-FO-PKG. 1375.00 EACH

Handpiece Without Connector

New Version With Higher Torque

#SURG-45-FO. 1120.00 EACH



💡 Fiber Optic Fits Kavo® MULTiflex® Connector

#SURG-45-FO-K 1285.00 EACH

Lubricant For 45° High Speed Handpiece



Pana Spray Lubricant - 10.5 oz.

Available In US Only

#PANA-SPRAY 71.80 EACH

**Remove Failed Implants Or
Broken Abutment Screws**



SPECIAL

Implant Rescue Kit

#IMPLANT-RESCUE-KIT **1785.00 KIT**

SPECIAL

Narrow Implant Broken Screw Kit

For Broken Screw Removal On Small Diameter & Bone Level Implants

#IMPLANT-RESCUE-NARROW **835.00 KIT**

Broken Abutment Screw Removal

- Helps To Remove Broken Screws Without Damaging The Implant
- Includes Centering Guides To Fit Most Implant Systems



Place Guide On Implant & Hold With Stabilizing Handle.

Set Motor To **REVERSE** At 1,000-1,250RPM & 50-70Ncm Torque.



Drill In **REVERSE** Using "Up & Down" Motion To Prepare 1-2mm Deep Dimple Into Top Of Broken Screw.

Use Irrigation To Minimize Heat.



Suction To Remove Metal Shavings.

Then Insert Tap Into Implant Handpiece

Set Motor To **REVERSE** At 70-80RPM & 50-70Ncm Torque.



Insert Tap Into The 1-2mm Dimple In The Top Of The Broken Screw.

Use The Tap In **REVERSE** To Remove Broken Screw.



Broken Screw Removed.

Failed Implant Removal

- Removes Failed Implants Without Compromising Surrounding Bone Or Damaging Neighboring Structures
- Works With Any Threaded Implant System
- Kit Includes Small, Medium And Large Fixture Removers



Insert Fixture Remover Screw Into Hand Adapter.



Turn Counter-Clockwise Into The Implant.



Use Wrench Counter-Clockwise Until The Implant Becomes Mobile.



Remove Failed Implant.



For Difficult Cases: Trephine Top 3-4mm Of Bone.

Sterile Water



Sterile Water-Bottles

- Bottles For Pouring
- 500 mL

#H2O-500-BOT-18
 CASE OF 18 110.90 CASE

Please Order Carefully. No Returns On Opened Cases.

Sterile Saline



Sterile Saline-250 mL

- 0.9% Sodium Chloride Bags
- 250 mL

#NACL-250-BAG-24
 CASE OF 24 169.90 CASE



MOST POPULAR

Sterile Saline-500 mL

- 0.9% Sodium Chloride Bags
- 500 mL

#NACL-500-BAG-24
 CASE OF 24 175.90 CASE



Sterile Saline-1000 mL

- 0.9% Sodium Chloride Bags
- 1000 mL

#NACL-1000-BAG-12
 CASE OF 12 85.90 CASE



Sterile Saline-10 mL Syringes

- Convenient For Mixing Graft Materials
- 0.9% Sodium Chloride
- 10 mL

#NACL-10-SYRINGE
 CASE OF 30 56.60 CASE

Please Order Carefully. No Returns On Opened Cases.

Sterile D5W (5% Dextrose Solution)



Sterile D5W

5% Dextrose Bags 500 mL

#D5W-500-BAG-24
 CASE OF 24 157.90 CASE

Please Order Carefully. No Returns On Opened Cases.

Irrigating Syringes & Catheters



Sterile Irrigating Syringes

- 60cc Luer Lock Syringes
- Sterile - Individually Packaged

#SYRINGE-STERILE-60CC
 BOX OF 30 73.90 BOX
 2 OR MORE BOXES 69.30 BOX



Sterile Irrigating Catheters

- Sterile Blunt Catheters
- 18 Ga x 1" (2.54cm)
- Individually Packaged

#BLUNT-CATH-18G-2.5CM
 BOX OF 100 124.90 BOX
 2 OR MORE BOXES 114.80 BOX



Monoject #412 Irrigating Syringes

- 12cc Disposable Syringes
- Curved Tip
- Packaged Non-Sterile

#MONOJECT-412
 BOX OF 50 59.50 BOX
 2 OR MORE BOXES 55.60 BOX

External Irrigation Wand Systems

- Available With Or Without Suction
- Uses Pump Up Pressure Cuff
- Includes Autoclavable Silicone Irrigation Tubing & IV Bag Striker Insert

Scan For
**Suction Irrigator
Technique Guide**



MOST POPULAR

Wand Irrigator - Without Suction

Complete Kit
#IRR-SET 209.00 SET



Wand Irrigator - With Suction

Complete Kit Including Built In Suction Aspirator
Order Suction Tubing Separately
#SUCT-IRR-SET 311.00 SET



Replacement Irrigation Tubing

#SILICONE-IRR-TUBING
. 22.10 FOR 7FT



Wand Irrigator With Suction Connects
To 1/4" Disposable Suction Tubing Or
Salvin Autoclavable Suction Tubing
On Pages 96-97

Stainless Steel Bowls



Solution Bowl

- Heavy Duty Stainless Steel (20 oz)
- 5" x 2 3/8" (12.7 x 6.1cm)

#SOLN BOWL 5 34.10 EACH
2 OR MORE 31.60 EACH



MOST POPULAR

Prep Cup

- Heavy Duty Stainless Steel (6 oz)
- 3 1/2" x 2" (8.9 x 5.1cm)

#PREP CUP 23.90 EACH
2 OR MORE 21.70 EACH

Autoclavable Tubing Adapter

Use Autoclavable Surgical Tubing With Your Regular High Volume Suction



Latexfree

Connects To Your Regular Dental Suction



MOST POPULAR

Autoclavable Tubing Adapter

- Autoclavable Hundreds Of Times
- Non-Crimping "U" Adapter Connects To Regular Dental Suction
- 5'(152.4cm) x 1/4"(0.6cm) ID Silicone Surgical Tubing
- Fits All Aspirator Tips

#TUBING ADAPT 49.40 EACH
2 OR MORE 47.80 EACH

● 1/4"(0.6cm) Inside Diameter

Autoclavable Silicone Surgical Suction Tubing

Remains Flexible...Autoclavable Hundreds Of Times



Latexfree

MOST POPULAR


Silicone Tubing By The Foot

1/4"(0.6cm) Internal Diameter
#TUBING. 4.65 FOOT
50 FT OR MORE 4.10 FOOT



End To End Connector

#TUBE CONNECT. . . . 16.90 EACH
3 OR MORE 15.50 EACH



MOST POPULAR

Autoclavable Cut Off Clamp

#CUT. 10.90 EACH
3 OR MORE 9.50 EACH



Connects High Volume Surgical Aspirators To Surgical Tubing

High Volume Adapter

Connects High Volume Surgical Aspirators To Surgical Tubing
#HVE NECK ADAP . . . 48.60 EACH

Autoclavable Suction Cut Off Valves

- *Smooth On/Off Lever Action*
- *Easy To Remove & Clean*



MOST POPULAR

Flip Saliva Ejector Valve

- Fits 7/32" Saliva Ejector Tubing
 - Autoclavable Silicone Head
 - Includes Cleaning Brush
- #FLIP-SE 69.00 EACH



Flip-D Cut Off Valve

- Fits Standard 1/2" Tubing
 - Includes Cleaning Brush
- #FLIP-D 69.00 EACH



Flip-A Cut Off Valve

- Fits Standard 1/2" Tubing
 - Includes Cleaning Brush
- #FLIP-A 69.00 EACH

**SALVIN
Rapid
Roll®**



MOST POPULAR

- Easy To Use Cover Ring Slides Over Long Hoses In Seconds
- Individually Packaged - 65" (165cm) Long

#RAPID-ROLL
 BOX OF 25..... 117.80 BOX
 2 OR MORE BOXES 108.50 BOX



Latexfree

Disposable Sterile Suction Tubing & Yankauer Suction Tips



Latexfree

MOST POPULAR

6' Sterile Suction Tubing

- 6' (182.9cm) x 1/4" (0.6cm) Inside Diameter
- Individually Packaged

#SUCTION-TUBE-6
 BOX OF 50 185.90 BOX (2.99 EACH)
 2 OR MORE 175.30 BOX (2.81 EACH)



Disposable Sterile Yankauer Tips

- Blue Translucent Plastic Minimizes Patient Visibility Of Fluids
- Flange Tip & Vent

#YANKAUER-DISPOSABLE
 BOX OF 50 136.90 BOX
 2 OR MORE 126.60 BOX

10' Sterile Suction Tubing

- 10' (304.8cm) x 1/4" (0.6cm) Inside Diameter
- Individually Packaged

#SUCTION-TUBE-10
 BOX OF 50 196.90 BOX (3.19 EACH)
 2 OR MORE 185.30 BOX (2.99 EACH)










MOST POPULAR

Surgical Tubing To Dental Suction Adapter


Adapts Surgical Tubing To Standard Perio & General Dentistry High Volume Suction

#TUBETIP 29.90 EACH

Surgical Aspirators







 <p>1.5mm</p>	<p>1.5mm Surgical Aspirator #15P1A 41.30 EACH</p>
 <p>2.5mm</p>	<p>2.5mm Surgical Aspirator #15P2A 41.30 EACH</p>
 <p>3.0mm</p>	<p>MOST POPULAR 3mm Surgical Aspirator #15P3A 41.30 EACH</p>
 <p>4.0mm</p>	<p>MOST POPULAR 4mm Surgical Aspirator #16RSA 41.30 EACH <i>Recommended Cleaning Brushes: Medium & Small</i></p>
 <p>2.0mm</p>	<p>2mm Modified Frazier #FR2A 53.70 EACH</p>
 <p>3.0mm</p>	<p>3mm Modified Frazier #FR3A 53.70 EACH</p>
 <p>4.0mm</p>	<p>MOST POPULAR 4mm Modified Frazier #FR4A 53.70 EACH <i>Recommended Cleaning Brushes: Small & Medium</i></p>

Control Suction By
Covering Relief Hole



**Aspirators On This Page Fit
Surgical Tubing
Or
Standard High Volume Suction**

Autoclavable Cleaning Brushes

	 <p>#BRUSH-MEDIUM 39.90 DOZEN</p>
	 <p>#BRUSH-SMALL 39.90 DOZEN MOST POPULAR</p>
	 <p>#BRUSH-X-SMALL 39.90 DOZEN</p>

Silicone Tip Aspirator Kit

- Autoclavable Silicone Tips Prevent Metal Contact With Implant
- Tips Are Easily Modified With Scissors Or Knife



MOST POPULAR

Silicone Tip Aspirator Kit

Stainless Steel Aspirator &
1 Trimmable Silicone Tip

#46YST 41.30 EACH


Silicone Tips

Set Of 3 Trimmable Silicone Tips

#YST 33.00 SET OF 3

*Recommended Cleaning Brushes:
Medium*

No-Clog Surgical Aspirators



● 2.5mm
or
● 3.0mm

Scan For
**No-Clog Aspirator
Technique Guide**



No-Clog Aspirator

Built-In Clog Ejector
Sliding Button Forward Ejects
The Clog

(2.5MM TIP)

#B2 NO CLOG 135.00 EACH

MOST POPULAR

(3.0MM TIP)

#B3 NO CLOG 135.00 EACH

*Recommended Cleaning Brushes:
Small & Medium*

Micro Aspirators For Endo/Root Tip Suctioning



● 1.0mm
● 2.0mm

1mm Micro/Endo Aspirator

Bendable Tubing
Includes Cleaning Stylet

#02BA 96.00 EACH

MOST POPULAR

2mm Micro/Endo Aspirator

Bendable Tubing
Includes Cleaning Stylet

#03BA 96.00 EACH

*Recommended Cleaning Brushes:
X-Small*

Micro/Endo Suction Handle



MOST POPULAR

Micro/Endo Suction Handle

Connects To Saliva Ejector Suction
Includes 10 Non-Sterile
Blunt Catheters

#SUCT-LL 53.90 EACH

Sterile Luer Lock Catheters

#BLUNT-CATH-18G-2.5CM
. 124.90 BOX OF 100

Non-Sterile Luer Lock Catheters

#SUCT-LL-TIP-14GA 41.30 BOX OF 100
#SUCT-LL-TIP-18GA 41.30 BOX OF 100
#SUCT-LL-TIP-20GA 41.30 BOX OF 100





IrriguIDE

ADAPTER

Developed by Dr. Thomas Jackson

NEW!

- *Compatible With Any Guided Dental Implant System*
- *Ensures Continuous Direct Water Flow To The Osteotomy Site*
- *Minimizes The Risk Of Thermal Damage To Surrounding Tissues*
- *Improves Visibility During Implant Procedures By Eliminating Obstructions Within The Surgical Field*
- *Eliminates The Need To Stop And Start The Drilling Process For Irrigation*



SPECIAL

IrriguIDE Adapter Starter Set

Includes:

- IrriguIDE Adapters - 5 pack
- IrriguIDE Seating Bur
- IrriguIDE Trimming Bur
- Cyanoacrylate glue
- Assembly Instructions

#IRRIGUIDE-STARTER-KIT
250.00



Mallets & Chisels



Cottle Curved Chisels

Depths Marked At
6, 9, 12, 14 & 17mm

MOST POPULAR

#COTTLE-4 (4MM WIDE)

#COTTLE-6 (6MM WIDE)

INDIVIDUALLY
..... 133.90 EACH
3 OR MORE
..... 126.40 EACH



Sheehan Osteotomes

Depths Marked At
6, 9, 12, 14 & 17mm


#SHEE 2 (2MM WIDE)

MOST POPULAR

#SHEE 4 (4MM WIDE)

#SHEE 6 (6MM WIDE)

INDIVIDUALLY
..... 133.90 EACH
3 OR MORE
..... 126.40 EACH



Mead Mallet

#MEAD
..... 167.80 EACH



TG Perio Chisel

#PERIO-CHISEL-TG
..... 106.50 EACH
3 OR MORE
..... 100.90 EACH



MOST POPULAR

TG-O Perio Chisel

#PERIO-CHISEL-TGO
..... 106.50 EACH
3 OR MORE
..... 100.90 EACH

SALVIN
PDL-Evator®
Precise Tip Luxating Elevator

Combines The Benefit Of A Periosteal Elevator & Elevator In One Instrument



MOST POPULAR **PDL-Evator® Straight - 2.5mm Tip**
2.5mm Narrow Tip Is Thin Enough To Insert In Tight Spaces
Between The Root And The Socket Wall
#PDL-STRAIGHT..... 150.70 EACH



MOST POPULAR **PDL-Evator® Curved**
#PDL-CURVED 150.70 EACH



MOST POPULAR **PDL-Evator® Bayonet**
#PDL-BAYONET 150.70 EACH



PDL-Evator® Angled "In" Small
#PDL-IN-SMALL 150.70 EACH



PDL-Evator® Angled "Out" Small
#PDL-OUT-SMALL 150.70 EACH



PDL-Evator® Angled "In" Large
#PDL-IN-LARGE 150.70 EACH



PDL-Evator® Angled "Out" Large
#PDL-OUT-LARGE..... 150.70 EACH

SPECIAL



SALVIN
PDL-Evator®
Precise Tip Luxating Elevator

MOST POPULAR

Includes:

- Set Of 7 Precise Tip Luxating Elevators
- **FREE** Stainless Steel Cassette

#PDL-PKG
1035.00 SET OF 7 W/CASSETTE

“Interchangeable Tip” Periotomes



MOST POPULAR

**Periotome Handle With
5 Interchangeable Tips**

#PERIOTOME-5-TIPS286.20 SET



*Use BUR-BLOCK-6-HP On Page 89
To Easily Store Periotome Tips*

Thin Straight Replacement Tip

#PERIOTOME-TIP-1 128.00 PKG OF 3



Wide Straight Replacement Tip

#PERIOTOME-TIP-2 128.00 PKG OF 3



Thin Angled Replacement Tip

#PERIOTOME-TIP-3 128.00 PKG OF 3



Flat Angled Replacement Tip

#PERIOTOME-TIP-4 128.00 PKG OF 3



Flat Wide Angled Replacement Tip

#PERIOTOME-TIP-5 128.00 PKG OF 3



Third Molar Elevators



Third Molar Elevator - Left

#ELEV-THIRD-MOLAR-L. 108.40 EACH



Third Molar Elevator - Right

#ELEV-THIRD-MOLAR-R. 108.40 EACH

Periotomes



MOST POPULAR

Straight Periotome With Fixed Straight Tip


#PERIOTOME-STR
 104.90 EACH
 3 OR MORE
 99.20 EACH



Straight Periotome With Fixed Rounded Tip

#PERIOTOME-STR-ROUND
 104.90 EACH
 3 OR MORE
 99.20 EACH

Serrated Tip



Straight Periotome With Serrated Tip


#PERIOTOME-SERR-STR
 125.10 EACH
 3 OR MORE
 117.80 EACH



Lane Elevator

Designed By Dr. Bill Lane

#LANE-ELEVATOR
 103.70 EACH
 3 OR MORE
 96.90 EACH



Periosteal-Tome™

Periosteal Elevator/
 Periotome Combo

#PERIOSTEAL-TOME
 108.20 EACH
 3 OR MORE
 101.90 EACH

Double-Ended Periostomes



MOST POPULAR


Double Ended Periostome #1

#PERIOTOME-DE-1
 102.70 EACH
 3 OR MORE
 95.10 EACH



Double Ended Periostome #2

#PERIOTOME-DE-2
 102.70 EACH
 3 OR MORE
 95.10 EACH



Serrated Tips

Serrated Tip Periostome #1

#PERIOTOME-SERR-1
 125.10 EACH
 3 OR MORE
 117.80 EACH



Serrated Tips

Serrated Tip Periostome #2

#PERIOTOME-SERR-2
 125.10 EACH
 3 OR MORE
 117.80 EACH

Serrated Tips

Serrated Tips

SALVIN®
AtrauLux®
Ergonomic Elevators



Handle Fits Comfortably In Palm
For Increased Control



MOST POPULAR
AtrauLux® Elevator
2mm Straight

#ATRAU-ELEV-2-STR
.....94.60 EACH

SALVIN
**Ergonomic
Handle**



MOST POPULAR
AtrauLux® Elevator
3mm Straight

#ATRAU-ELEV-3-STR
.....94.60 EACH

SALVIN
**Ergonomic
Handle**



AtrauLux® Elevator
4mm Straight

#ATRAU-ELEV-4-STR
.....94.60 EACH

SALVIN
**Ergonomic
Handle**



AtrauLux® Elevator
2mm Curved

#ATRAU-ELEV-2-CUR
.....94.60 EACH

SALVIN
**Ergonomic
Handle**



AtrauLux® Elevator
3mm Curved

#ATRAU-ELEV-3-CUR
.....94.60 EACH

SALVIN
**Ergonomic
Handle**



AtrauLux® Elevator
4mm Curved

#ATRAU-ELEV-4-CUR
.....94.60 EACH

SALVIN
**Ergonomic
Handle**

SPECIAL

Set Of 6
Salvin AtrauLux®
Elevators

MOST POPULAR

#ATRAU-ELEV-SET-6
499.00 SET OF 6

Serrated Elevators



Serrated Elevator 2.2mm

#DEL-SERR-2
..... 102.80 EACH



MOST POPULAR
Serrated Elevator 3mm

#DEL-SERR-3
..... 102.80 EACH



Serrated Elevator 4mm

#DEL-SERR-4
..... 102.80 EACH



Serrated Elevator 4mm - Left

#DEL-SERR-4-L
..... 102.80 EACH



Serrated Elevator 4mm - Right

#DEL-SERR-4-R
..... 102.80 EACH

SPECIAL	MOST POPULAR	Set Of 5 Bayonet Elevators #DEL-SERR-SETOF5 429.00 SET OF 5
----------------	---------------------	--

Atraumatic Bayonet Elevators



MOST POPULAR
Bayonet Elevator Straight

#BAYONET-ELEV-1
.... 102.80 EACH



Bayonet Elevator Offset "In" - Small

#BAYONET-ELEV-2
.... 102.80 EACH



Bayonet Elevator Offset "Out" - Small

#BAYONET-ELEV-3
.... 102.80 EACH



Bernard Bayonet Elevator Offset - Left

#BERNARD-L
..... 102.80 EACH



Bernard Bayonet Elevator Offset - Right

#BERNARD-R
..... 102.80 EACH

SPECIAL	MOST POPULAR	Set Of 5 Bayonet Elevators #BAYONET-SETOF5 429.00 SET OF 5
----------------	---------------------	---

Bone Preserving Atraumatic Luxating Elevators



3mm Luxating Elevator Curved

#LUX-3-CURVED
..... 96.30 EACH



MOST POPULAR
3mm Luxating Elevator Straight

#LUX-3-STRAIGHT
..... 96.30 EACH



5mm Luxating Elevator Curved

#LUX-5-CURVED
..... 96.30 EACH



5mm Luxating Elevator Straight

#LUX-5-STRAIGHT
..... 96.30 EACH

SPECIAL

Set Of 4 Luxating Elevators



MOST POPULAR

#LUX-SETOF4
336.00 SET OF 4

"Anatomically Grooved" Luxating Elevators



MOST POPULAR
2mm Grooved Luxating Elevator

#LUX-GROOVE-2
..... 96.30 EACH



3mm Grooved Luxating Elevator

#LUX-GROOVE-3
..... 96.30 EACH



4mm Grooved Luxating Elevator

#LUX-GROOVE-4
..... 96.30 EACH



5mm Grooved Luxating Elevator

#LUX-GROOVE-5
..... 96.30 EACH

SPECIAL

Set Of 4 Grooved Elevators



#LUX-GROOVE-SETOF4
336.00 SET OF 4

Luxating Elevators With Periotope Tip



MOST POPULAR

3mm Luxating Elevator/Periotope

#LUX-TOME-3
..... 96.30 EACH



4mm Luxating Elevator/Periotope

#LUX-TOME-4
..... 96.30 EACH



5mm Luxating Elevator/Periotope

#LUX-TOME-5
..... 96.30 EACH

SPECIAL

**Set Of 3
Luxating Elevators**



#LUX-TOME-SETOF3
252.00 SET OF 3

Thin Handle Luxating Elevators



MOST POPULAR

2.5mm Thin Handle Luxating Elevator

#LUX-THIN-2.5
..... 96.30 EACH



3.5mm Thin Handle Luxating Elevator

#LUX-THIN-3.5
..... 96.30 EACH

SPECIAL

**Set Of 2
Luxating Elevators**



#LUX-THIN-SETOF2
162.00 SET OF 2

Salvin Twist Elevators

- Twist Design Provides Improved Vertical Movement Of The Tooth Over Conventional Elevators
- Sharp Tips Engage Tooth & Provide Deeper Subgingival Access



MOST POPULAR

**Twist Elevator
Straight Left**

#TWIST-ELEV-STR-L
..... 113.40 EACH



MOST POPULAR

**Twist Elevator
Straight Right**

#TWIST-ELEV-STR-R
..... 113.40 EACH



**Twist Elevator
Angled Left**

#TWIST-ELEV-ANG-L
..... 113.40 EACH



**Twist Elevator
Angled Right**

#TWIST-ELEV-ANG-R
..... 113.40 EACH

SPECIAL

**Set Of 4 Salvin
Twist Elevators**



**#TWIST-ELEV-SET-4
384.00 SET OF 4**

Elevators



**Potts Elevator-
Left**

#POTTS-L
..... 92.70 EACH




**Potts Elevator-
Right**

#POTTS-R
..... 92.70 EACH



**Potts "T-Handle"
Elevator-Left**

#POTTS-T-L
..... 98.80 EACH



**Potts "T-Handle"
Elevator-Right**

#POTTS-T-R
..... 98.80 EACH



Coupland Gouge Elevator

Useful For Removal Of Teeth
With Amputated Crowns

#COUPLAND-1C... 98.40 EACH

Elevators



Your Choice Of Elevators

INDIVIDUALLY . . . 93.90 EACH
3 OR MORE 88.50 EACH

SALVIN®
AtrauLux®
Extraction Forceps


- Tapered Tips Adapt Exactly To The Tooth With Excellent Subgingival Contact
- Lightweight, Ergonomic Design Fits Both Right & Left Handed Practitioners



SPECIAL **MOST POPULAR**


Set Of 9 Salvin AtrauLux® Forceps

#ATRAU-FORCEP-SET-9 . . . 2150.00 SET OF 9




AtrauLux® Forcep Cow Horn

Serrated Beak
Lower Atraumatic Extraction
#ATRAU-FORCEP-COWHORN . . 284.70 EACH




AtrauLux® Forcep #1

Serrated Beak
Upper Anterior Atraumatic Extraction
#ATRAU-FORCEP-1 284.70 EACH




AtrauLux® Forcep #2

Serrated Beak
Lower Anterior Atraumatic Extraction
#ATRAU-FORCEP-2 284.70 EACH




AtrauLux® Forcep #3

Serrated Beak
Upper Premolar Atraumatic Extraction
#ATRAU-FORCEP-3 284.70 EACH




AtrauLux® Forcep #4

Serrated Beak
Lower Premolar Atraumatic Extraction
#ATRAU-FORCEP-4 284.70 EACH




AtrauLux® Forcep #5

Lower Root Atraumatic Extraction
#ATRAU-FORCEP-5 284.70 EACH




AtrauLux® Forcep #6

Upper Root Atraumatic Extraction
#ATRAU-FORCEP-6 284.70 EACH



AtrauLux® Forcep #7

Serrated Beak
Lower Molar Atraumatic Extraction
#ATRAU-FORCEP-7 284.70 EACH



AtrauLux® Forcep #8

Serrated Beak
Lower Atraumatic Extraction
#ATRAU-FORCEP-8 284.70 EACH

DEEP GRIP Extended Tip Length Serrated Forceps


BUYER'S GUIDE	
▶	Superior Grip & Tactile Feel Even With Limited Contact
▶	Less Time & Force Required, Even For Difficult Extractions

*Tips Marked At
1mm - 8mm Depths*



SPECIAL

MOST POPULAR



Set Of 4 DEEP GRIP Forceps

#DGA-SET-4 895.00 SET OF 4






DEEP GRIP
Extended Tip
Serrated Forcep




DGA-34


Serrated Beak
Upper Atraumatic Extraction

#DGA-34..... 284.70 EACH





DEEP GRIP
Extended Tip
Serrated Forcep





MOST POPULAR

DGA-35


Serrated Beak
Upper Atraumatic Extraction

#DGA-35..... 284.70 EACH






DEEP GRIP
Extended Tip
Serrated Forcep




DGA-36


Serrated Beak
Lower Atraumatic Extraction

#DGA-36..... 284.70 EACH





DEEP GRIP
Extended Tip
Serrated Forcep




MOST POPULAR

DGA-37

Serrated Beak
Lower Atraumatic Extraction

#DGA-37..... 284.70 EACH

Tungsten Carbide Tipped Root Tip Forceps




MOST POPULAR

Root Tip-Lower

Carbide Serrated Beak
Lower Root Tips & Molars

#ROOT CU293.40 EACH



TUNGSTEN CARBIDE


MOST POPULAR

Root Tip-Upper

Carbide Serrated Beak
Upper Root Tips & Molars

#ROOT TC.293.40 EACH


German Stainless Steel Extraction Forceps



DEF 150AS

Serrated Beak
Upper Anteriors, Split Beaks


#DEF150AS.212.70 EACH



DEF 151AS

Serrated Beak
Lower Anteriors, Split Beaks


#DEF151AS.212.70 EACH



DEF 222

Lower Molars
Excellent Access To 3rd Molars

#DEF-222212.70 EACH




MOST POPULAR

DEF 23-Lower

Lower Cow Horn
1st & 2nd Lower Molars


#DEF 23212.70 EACH



DEF 23-Upper

Upper Cow Horn
1st & 2nd Upper Molars

#DEF-23-UPPER212.70 EACH



Wrap Around Extraction Forcep


Right Angle "Wrap Around" Design
Provides Superior Access

Serrated Beak
#DEF-68-WRAP 212.70 EACH

Tungsten Carbide Tip Extraction Forceps – Better Grip & Tactile Feel – Less Trauma

BUYER'S GUIDE

- ▶ Superior Grip & Tactile Feel Even With Limited Contact
- ▶ Less Time & Force Required, Even For Difficult Extractions



MOST POPULAR

DEF 150 TC

Carbide Serrated Beak
Upper Incisors, Cuspids,
Bicuspid & Roots, Universal

#DEF 150 TC 293.40 EACH



MOST POPULAR

DEF 151 TC

Carbide Serrated Beak
Lower Incisors, Cuspids,
Bicuspid & Roots, Universal

#DEF 151 TC 293.40 EACH



DEF 15TD TC

Carbide Serrated Beak
Upper Incisors, Cuspids


#DEF 15TD TC . . . 293.40 EACH



DEF 17 TC

Carbide Serrated Beak
1st & 2nd Lower Molars,
Universal


#DEF 17 TC 293.40 EACH



DEF 53L TC

Carbide Serrated Beak
1st & 2nd Upper Molars, Left


#DEF 53L TC 293.40 EACH



DEF 53R TC

Carbide Serrated Beak
1st & 2nd Upper Molars, Right

#DEF 53R TC 293.40 EACH



DEF 65 TC

Carbide Serrated Beak
Upper Incisors,
Roots & Fragments

#DEF 65 TC 293.40 EACH



MOST POPULAR

DEF 69 TC

Carbide Serrated Beak
Upper & Lower
Fragments & Roots

#DEF 69 TC 293.40 EACH



DEF 150S TC Pedo Forcep

Carbide Serrated Beak
Pedo, Upper Primary
Teeth & Roots, Universal

#DEF 150S TC . . . 293.40 EACH



DEF 151S TC Pedo Forcep

Carbide Serrated Beak
Pedo, Lower Primary
Teeth & Roots, Universal

#DEF 151S TC . . . 293.40 EACH

Bayonet Atraumatic Extraction Forceps



MOST POPULAR

Bernard Forcep 153

Upper Atraumatic Extraction

#BERNARD-153 222.30 EACH



MOST POPULAR

Bernard Forcep 146

Lower Atraumatic Extraction

#BERNARD-146 222.30 EACH

"Fish-Hook" Atraumatic Extraction Root Forceps



"Fish Hook"
Serrations
For Better
Traction



**45° Angled Root Forcep
"Fish-Hook" Serrations**

Upper Atraumatic Extraction

#FHA-45 279.00 EACH



**90° Angled Root Forcep
"Fish-Hook" Serrations**

Lower Atraumatic Extraction

#FHA-90 279.00 EACH

SPECIAL

#FISH-HOOK-SET-2 499.00 SET OF 2

Physics Forceps® Standard Series

- Predictable & Efficient Bone Preservation Extractions In Less Than Four Minutes
- Elimination Of Root Tip Fractures
- No Need For Elevating Or Laying A Flap
- Very Little Operator Movement Or Force Necessary



Photo: Dr. Randy Resnik



Upper Left Forcep

#PHYSICS-FORC-UP-L
..... 630.00 EACH



Upper Anterior Forcep

#PHYSICS-FORC-UP-ANT
..... 630.00 EACH



Upper Right Forcep

#PHYSICS-FORC-UP-R
..... 630.00 EACH



Lower Universal Forcep

#PHYSICS-FORC-LOWER-U
..... 630.00 EACH



Bumper Guards

#PHYSICS-GUARDS-24
..... 26.60 EACH

SPECIAL

Set Of 4 Physics Forceps®

Includes:

- Upper Left Forcep
- Upper Right Forcep
- Upper Anterior Forcep
- Lower Universal Forcep
- 24 Bumper Guards

#PHYSICS-FORCEPS-SET4 **2009.00 SET**



Partial Extraction Therapy

SPECIAL

P.E.T. Partial Extraction Therapy Bur Kit

Fits Friction Grip Handpieces

Includes:

- Short 1.2mm Diameter Bone Cutter
- Long 1.2mm Diameter Bone Cutter
- Coarse 1.0mm Round Diamond Bur
- 1.8mm Round Prep Bur
- 2.3mm Round Prep Bur
- 2.9mm Round Prep Bur
- Stainless Steel Bur Block

#PET-BUR-KIT **143.15 KIT W/BLOCK**



Forceps & Root Tip Picks



Tungsten Carbide Tip
20mm Tip Length
For Deep Access

TUNGSTEN CARBIDE



MOST POPULAR

TC Coated Tip
Root Fragment Forcep

Extended Tips For Excellent
Access & Gripping

#FRAGMENT-FORCEP-TC 243.50 EACH
2 OR MORE 230.10 EACH



Warwick Root Tip
Elevator - Left

#WARWICK-L
. 80.50 EACH



Warwick Root Tip
Elevator - Right

#WARWICK-R
. 80.50 EACH

SPECIAL

**Set Of
Warwick Root Tip
Elevators**

**#WARWICK-SET
134.60 SET OF 2**



Heidbrink H1

Single-Ended
#HEID 1
. 73.80 EACH



MOST POPULAR

Heidbrink H2

Single-Ended
#HEID 2
. 73.80 EACH



MOST POPULAR

Heidbrink H3

Single-Ended
#HEID 3
. 73.80 EACH




Heidbrink H2-H3

Double-Ended
#HEID 2/3
. 73.80 EACH



Meisinger Benex Pro Atraumatic Extraction Kit


- *Screw Action Lifting Makes Thin Buccal Plate Extractions Easy!*
- *Excellent For Use In The Aesthetic Zone*




"I use the Benex Extractor on a daily basis in my office for difficult extractions to save time and preserve bone for implant placement."
-Dr. John Russo
 Periodontist, Sarasota, FL



Drill Hole In Root



Insert Extraction Screw



Turn Knob To Gently & Safely Extract Root

SPECIAL

Salvin Benex Atraumatic Kit

MOST POPULAR

Includes:

- Benex Extraction Handle
- Drills, Extraction Screws, Pull Wires
- Organizer Holds Both Parts & Handle
- Support Tray To Stabilize Handle If Required

#BENEX-PRO2445.80 COMPLETE KIT



Replacement Parts For Benex Pro Atraumatic Extraction Kit



Narrow Extraction Screw

- 1.3mm x 10mm
#MEIS-BE001-060 119.35 EACH
- 1.3mm x 16mm
#MEIS-BE001-230-070 119.35 EACH



Wide Extraction Screw

- 2.1mm x 10mm
#MEIS-BE001-230-080 119.35 EACH
- 2.1mm x 16mm
#MEIS-BE001-230-085 119.35 EACH



Traction String

- #MEIS-BE001-230020 143.35 EACH




**Salvin
POWER POINT®
Pilot Drills
On Page 83**

Serrated Lucas Bone Currettes

Serrated Tips Help Debride Extraction Sockets & Create Bleeding Points For Faster Healing


Designed By Dr. Jin Kim



Serrated Lucas Bone Currettes

INDIVIDUALLY	81.30 EACH
3 OR MORE	75.50 EACH

Non-Serrated Lucas Currettes



Non-Serrated Lucas Bone Currettes

INDIVIDUALLY	76.80 EACH
3 OR MORE	71.20 EACH

ZAGA Zygomatic Surgical Instruments

Designed By Dr. Carlos Aparicio – Creator Of The ZAGA Method



“It is my pleasure to share this great set of instruments designed to sequentially facilitate the placement of zygomatic implants following the ZAGA philosophy.”


Dr. Carlos Aparicio
Periodontist, Barcelona, Spain



**Salvin ZAGA 1 -
Zygoma Dissector**

Detaches The Periosteum Beyond The Angle Formed By Temporal & Frontal Anatomy Of The Zygoma To Prepare For The Anchorage Of The Retractor

#ZAGA-1
..... 133.00 EACH



**Salvin ZAGA 2 -
Zygoma Retractor**

- Increases Field Of Vision
- Directs The Drilling
- Prevents Perforation Of The Orbit
- Protects Soft Tissue While Drilling Superior Cortex

#ZAGA-2
..... 209.00 EACH



**Salvin ZAGA 3 -
Maxilla Dissector**

Detaches The Periosteum From Posterior Part Of The Maxillary Wall & Detaches The Floor Of The Nose


#ZAGA-3
..... 133.00 EACH



**Salvin ZAGA 4 -
Posterior Dissector**

Finishes Detachment Of The Periosteum In More Posterior Superior Sectors

#ZAGA-4
..... 133.00 EACH



**Salvin ZAGA 5 -
Nasal Protector**

Protects The Mucosa Of The Nose In Case Of Perforation Of The Nasal Floor

#ZAGA-5
..... 133.00 EACH



**Zygomatic
Trough Bur**

- 2.3mm Pear Shape
- Fits Straight Handpieces
- 44mm Long

#MEIS-HM77MF-023-HP
INDIVIDUALLY . . . 29.90 EACH
10 OR MORE . . . 28.35 EACH

SPECIAL

Salvin ZAGA Instrument Set

Includes:

- Zygoma Dissector
- Zygoma Retractor
- Maxilla Dissector
- Posterior Dissector
- Nasal Protector
- Plus **FREE** Salvin 10" x 6" Deep Instrument Cassette

#ZAGA-INSTR-KIT
..... **723.00 SET**



Graft Material Packers - Titanium

MOST POPULAR

Titanium Graft Packer / Placement Site Depth Gauge

- Titanium For Biocompatibility
- Excellent In Difficult Access Sites
- 3.3mm & 4.0mm Diameter Tips
- Marked At 8, 11, 13, 15 & 18mm

#GRAFT-PACKER-TI 132.40 EACH




MOST POPULAR

Serrated Titanium Tip Graft Packer-Small

- Excellent In Difficult Access Sites
- 3mm & 4mm Diameter Tips


#GRAFT-PACK-SERR-S 132.40 EACH



Serrated Titanium Tip Graft Packer-Large

- Excellent In Difficult Access Sites
- 6mm & 8mm Diameter Tips

#GRAFT-PACK-SERR-L 132.40 EACH



Graft Material Packers - Stainless Steel

Graft Packer

- Excellent In Difficult Access Sites
- 6mm & 8mm Diameter Smooth Tips


#GRAFT-PACKER 111.70 EACH



Pikos Graft Packer-Angled

- Concave Tips
- Excellent For Use In Sinus
- 4mm & 6mm Diameter Tips

#PIKOS-PACKER-ANGLED 126.90 EACH



Curved Graft Packers For Pockets Around Implants - Titanium

MOST POPULAR

**Curved Titanium Pocket Packer
1.75mm Radius**

3mm x 1mm Tip

#POCKET-PACKER-1.75

..... 102.90 EACH



**Curved Titanium Pocket Packer
2.5mm Radius**

4mm x 1mm Tip

#POCKET-PACKER-2.5

..... 102.90 EACH



Graft Material Transfer Instruments

**Graft Material Spoon With
4mm Serrated Graft Packer**

Concave Spoon Tip For Larger Grafts

#BONE-SPOON-PACK-4 122.10 EACH



MOST POPULAR

**Graft Material Spoon With
3mm Serrated Graft Packer**

Concave Spoon Tip For Larger Grafts

#BONE-SPOON-PACKER 122.10 EACH



**Sinus Floor Elevator With
3mm Graft Packer**

- Marked At 7, 10, 12, 14, 16 & 19mm
- Designed By Dr. Richard Martin

#MARTIN-PLUGGER. 122.10 EACH



Stainless Steel Graft Material Syringes

MOST POPULAR

5mm Diameter Stainless Steel Syringe-Curved


#GRAFTSYRINGE-CURVED 109.20 EACH

3.6mm Diameter Stainless Steel Syringe-Curved

#GRAFT-SYR-S-CURV 109.20 EACH

2.8mm Diameter Stainless Steel Syringe-Curved

#GRAFT-SYR-2.8-CU 109.20 EACH



Autoclavable Plastic Syringes For Smaller Grafts & Perio Defects

Perio Defect Graft Syringes

- Disassembles For Easy Cleaning
- Autoclavable, Chemiclavable

#HA SYRINGE 64.90 SET OF 2
2 OR MORE SETS 59.20 SET OF 2



Holds Approximately 0.5 cc
Sold As A Set Of 2 Syringes



Single Use Plastic Bone Graft Delivery Syringes

Bone Graft Delivery Syringes

- Packaged Sterile
- Single Use
- Blunt 5mm Wide Tip For Ease Of Delivery

#BONE-SYRINGE-1CC 49.90 PACK OF 10
2 OR MORE PACKS 47.40 PACK OF 10



Holds 1.0 cc
Sold As A Pack Of 10 Syringes



Stainless Steel Graft Material Syringes




3.5mm Diameter Ring Handle Syringe-Curved
 #GRAFT-SYRINGE-3.5. 206.30 EACH



MOST POPULAR

4.5mm Diameter Ring Handle Syringe-Curved
 #GRAFT-SYRINGE-4.5. 206.30 EACH



6.0mm Diameter Ring Handle Syringe-Curved
 #GRAFT-SYRINGE-6.0. 206.30 EACH



Graft Material Mixing Dishes



MOST POPULAR

Graft Material Mixing Cups

- Deep Dish With Rounded Bottom
- Holds 4 To 6 Grams Of Hydrated Material
- Autoclavable, Chemiclavable
- Will Not Crack During Sterilization

#GRAFT MIX. 53.40 SET OF 2
 2 OR MORE SETS 49.20 SET OF 2



Stainless Steel Weighted Graft Material Dish

- Deep Dish With Weighted Bottom To Prevent Tipping
- Holds 20 Grams Of Hydrated Material
- Autoclavable, Chemiclavable

#GRAFT-DISH 58.40 EACH
 2 OR MORE 54.40 EACH

OsseoScape™ Angled Single Use Bone Harvester



MOST POPULAR



- Angled Lightweight Ergonomic Handle
- Collection Chamber Holds Approximately 1cc
- Autoclave Before Use

#OSSEO-SCRAPE 179.00 BOX OF 5
(Hex Driver Tip & Thumb Knob Sold Separately)



*For A New Sharp Edge,
Use A 1.25mm Hex Driver Tip
To Turn The Blade 90°*



.05" (1.25mm) Hex Long - 12mm

#HTLT 58.90 EACH



Thumb Knob

#TK 79.90 EACH

Sterile Single Use Bone Harvester



- Angled Lightweight Ergonomic Handle
- Collection Chamber Holds Approximately 1cc
- Packaged Sterile

#BONE-SCRAPER-STERILE 50.40 EACH
(Hex Driver Tip & Thumb Knob Sold Separately)



Angled Tip For Easy Access




*Collection Chamber Holds
Approximately 1cc*



*For A New Sharp Edge,
Use A 1.25mm Hex Driver Tip
To Turn The Blade 90°*

Salvin Stainless Steel Bone Scraper



Stainless Steel Bone Scraper

- Clean & Autoclave For Repeated Use
- No Blade To Change

#BONE-SCRAPER 184.90 EACH

Stainless Steel Bone Scraper

- Collects And Stores Approximately 0.5cc Of Harvested Bone
- Blade Can Reattach For Unlimited Volume Per Case
- Low Cost – Replace Blade, Not Entire Grafters

#BONE-SCRAPE-SS. 180.00 EACH

NEW!



Bone Scrape SS Blade

- Easily Interchangeable Blade

#BONE-SCRAPE-SS-BLADE 44.00 EACH



Thin Profile Head Provides Easy Access



Replaceable Blade For Economy Of Use



Collects Approximately 0.5cc Of Bone


ACM – Autobone Chip Maker

- Collect 1cc Of Autogenous Bone In 10-20 Seconds
- Stopper Prevents Over Drilling & Helps Contain Bone Chips

NEW!



Easily Operated By Mounting On A Contra-Angle Handpiece



- Drilling Speed: 300rpm
- Available In Four Diameters: 4.0mm, 4.5mm, 5.0mm, 6.0mm
- Includes: One Drill, Seven Stoppers

#ACM-4.0	225.00
#ACM-4.5	225.00
#ACM-5.0	225.00
#ACM-6.0	225.00



Osseous Collection Trap

- Stop Suctioning Usable Bone!
- Simplify Bone Harvesting With Our Osseous Trap

Scan For **60th Second Product Video™**



Scan For **Osseous Collection Trap Technique Guide**



MOST POPULAR



- Uses Any Aspirator That Fits Into A High Volume Suction
- Fits Regular High Volume Suction As Well As Silicone Or Latex Surgical Suction Tubing



3mm Stainless Steel Aspirator Tip

#15P3A	41.30 EACH
--------	-------	------------

Osseous Coagulum Trap With 3 Single Use Filters

- Use "In Line" With Your Suction
- Osseous Material Is Collected Directly Into A Filter
- Can Be Chemiclaved, Autoclaved, Or Dry Heat Sterilized

#OCT	132.00 EACH
------	-------	-------------



Set Of 3 Replacement Single Use Filters

• Supplied Non-Sterile – Autoclave Before Using	
#OCT FILTERS26.90 PKG OF 3
3 OR MORE SETS25.30 PKG OF 3

Urban Master Core Kit

Simple And Safe Extraction Of Bone Cylinders For Bone Augmentation

Developed By Dr. Istvan Urban

- Depth Limited Trepines
- Depth Markings For Additional Depth Cont
- Black Coating To Prevent Glare

Master Core Kit

- 6 7mm Trepines Of 3mm-8mm Lengths
- 5 5mm Trepines Of 4mm-8mm Lengths
- 3 3mm Trepines Of 6mm-8mm Lengths
- 1 3mm Trepine Of 2mm Length
- 1 5mm Trepine Of 2mm Length
- 1 7mm Trepine Of 2mm Length
- 1 Drill Extender

#MEIS-BMC00 1220.15 KIT



Master Core Kit - Basic

- 5 7mm Trepines Of 3mm-7mm Lengths
- 5 5mm Trepines Of 3mm-7mm Lengths
- 1 5mm Trepine Of 2mm Length
- 1 7mm Trepine Of 2mm Length
- 1 Drill Extender

#MEIS-BMCBA 1087.15 KIT



Autogenous Bone Grafting Bur Kits



Bone Chip Maker Bur Kit With Stops

4 Bone Harvesting Trepines (3.4mm, 3.9mm, 4.2mm, 5.2mm), 4 Trepine Stops (5mm, 7mm, 10mm, 11.5mm), 1 Bone Remover & 1 Latch To Square Adapter

#BONE-CHIP-MAKER329.00 KIT



Bone Collector Bur Kit

5 Bone Collecting Drills (2.8mm, 3.3mm, 3.8mm, 4.3mm, 4.8mm), 5 Drill Stops (6mm, 8mm, 10mm, 11.5mm, 13mm)

#BONE-COLLECTOR-KIT329.00 KIT

Tissue Rakes & Flap Holding Forceps

For Greater Flap Release & Bone Graft Site Access



Separate Fibers For Increased Flap Release



Separate Fibers For Increased Flap Release



Retract Flap For Easy Access

Tissue Rake #1

Designed By Dr. Hom-Lay Wang
For Increased Flap Extension

#TISSUE-RAKE-199.90 EACH

Tissue Rake #2

Designed By Dr. Hom-Lay Wang
For Increased Flap Extension


#TISSUE-RAKE-299.90 EACH

Flap Holding Forcep


For Easy Bone Graft Site Access

#FLAP-FORCEP.....105.90 EACH

Bone Block Clamps



Screw Fits Through Notched Tip
Photos: Dr. Yvan Poitras



Anterior Cortical Block Clamp

- Ratchet Lock Holds Block In Place
- 19.7cm Long

#BONE-BLOCK-CLAMP359.50 EACH



Tightens To Hold Block In Place



Adjustable Straight Cortical Block Clamp

- Adjustable Width To Hold Block In Place
- 10cm Long

#BONE-CLAMP-ADJ.....299.00 EACH


Hand-Held Bone Morselizer

MOST POPULAR

Bone Morselizer

- Collected Bone Chips Ready To Use For Augmentation
- No Loss Of Bone Material
- Easy To Clean
- Adjustable Screw Allows You To Achieve Preferred Particle Size

#BONE-MORSELIZER 609.00 EACH




Bone Mill For Fast & Easy Crushing Of Cortical Bone

Stainless Steel Bone Mill

- Easy Assembly & Cleaning
- No Loss Of Bone Material
- Bone Collected In Easy To Use Dish
- Collected Bone Chips Ready To Use For Bone Augmentation

#BONE-MILL-PLUS 1925.00 EACH



Titanium Bone Cruncher™

Bone Cruncher™

1. Place Harvested Bone In Well Of Instrument
2. Cover Bone With Top, Serrated Side Down
3. Use Mallet To Reduce To Desired Particle Size

#BONE-CRUNCHER 339.00 EACH

Mead Mallet

#MEAD 167.80 EACH



Salvin Titanium **Tenting** Screws





Photo: Dr. Nick Shumaker

5mm Diameter Head Provides Broad Surface Area During GBR Tenting Procedures



Screws Have 4mm Of Thread For Stabilization In Bone

Available With 3, 4 or 5mm Unthreaded Surface For Each Unique Case



Tenting

Dental Latch Friction Fit Screwdriver #FIX-DRIVER-LATCH. 107.50 EACH


Drill Bit-Dental Latch (1.1mm) #FIX-DRILL-RA 107.50 EACH

1.5mm Diameter Tenting Screws - 5mm Diameter Head

1.5mm x 7mm Tenting Screw	#TENT-1.5X7	53.70 EACH
1.5mm x 8mm Tenting Screw	#TENT-1.5X8	53.70 EACH
1.5mm x 9mm Tenting Screw	#TENT-1.5X9	53.70 EACH

SPECIAL

Salvin Tenting Screw Kit



Includes:

- 12 Tenting Screws (Your Assortment)
- Bone Fixation Screwdriver
- Latch Screwdriver
- Latch Drill Bit
- Autoclavable Organizer Block

REGULAR PRICE OF COMPONENTS IS 1648.40

SAVE 573.40

#TENT-KIT-W/DRIVER. 1087.00 KIT

*Screwdriver Tips Are Proprietary To Manufacturer's Screw System And Are Not Interchangeable
The Salvin Fixation Screwdriver Works With Both Our Fixation Screws And Our Tenting Screws*

Salvin Titanium Fixation Screws






Photo: Dr. Yvan Poitras


Use To Secure Autogenous or Allograft Bone Blocks In The Mandible Or Maxilla

Dental Latch Friction Fit Screwdriver 
#FIX-DRIVER-LATCH..... 107.50 EACH

Drill Bit-Dental Latch (1.1mm) 
#FIX-DRILL-RA 107.50 EACH

Lag Drill Bit-Dental Latch (1.5mm) 
Oversize Drill Bit For Drilling Enlarged Hole In Graft
#FIX-DRILL-LAG 107.50 EACH


Drill Bit-Straight Handpiece (1.1mm)
#FIX-DRILL-HP 107.50 EACH



Lag Drill Bit-Straight Handpiece (1.5mm)
Oversize Drill Bit For Drilling Enlarged Hole In Graft
#FIX-LAG-DRILL-HP 107.50 EACH



Scan For




Locking Tip Driver Keeps Screw Secure Until You're Ready To Release It



Block Fixation

1.5mm Diameter Screws - 3mm Diameter Head		
1.5mm x 4mm Bone Screw	#FIX-1.5X4-PT	53.70 EACH
1.5mm x 6mm Bone Screw	#FIX-1.5X6-PT	53.70 EACH
1.5mm x 8mm Bone Screw	#FIX-1.5X8-PT	53.70 EACH
1.5mm x 10mm Bone Screw	#FIX-1.5X10-PT	53.70 EACH
1.5mm x 13mm Bone Screw	#FIX-1.5X13-PT	53.70 EACH
1.5mm x 15mm Bone Screw	#FIX-1.5X15-PT	53.70 EACH



Block Fixation

2.0mm Diameter Emergency Screws - 3mm Diameter Head		
2.0mm x 4mm Bone Screw	#FIX-2.0X4-PT	53.70 EACH
2.0mm x 6mm Bone Screw	#FIX-2.0X6-PT	53.70 EACH
2.0mm x 8mm Bone Screw	#FIX-2.0X8-PT	53.70 EACH
2.0mm x 10mm Bone Screw	#FIX-2.0X10-PT	53.70 EACH
2.0mm x 13mm Bone Screw	#FIX-2.0X13-PT	53.70 EACH
2.0mm x 15mm Bone Screw	#FIX-2.0X15-PT	53.70 EACH

SPECIAL

Salvin Block Fixation Screw Kit



or



Includes:

- 18 Fixation Screws (Your Assortment)
- Latch Screwdriver
- Bone Fixation Screwdriver
- Drill Bit (Latch or Straight HP)
- Lag Bit (Latch or Straight HP)
- Autoclavable Organizer Block

REGULAR PRICE OF COMPONENTS IS 2078.10
SAVE 648.00

#FIX-KIT-SPECIAL **1430.00 KIT**

Salvin Titanium Bone Tacks



Scan For

Tack System
Technique Guide



- Available In 3mm & 5mm Lengths
- Packaged 5 Tacks To A Vial, Non-Sterile
- Head Of Tack Is 2.5mm Diameter

3MM LENGTH
#TACK-3-PT 71.90 PKG OF 5

5MM LENGTH
#TACK-5-PT 71.90 PKG OF 5

Use Tacks To Secure Membrane Over Surgical Site

Salvin Tack Kits Available With Either:

Mallet & Placement Instrument



Load Tack On Instrument



Lightly Tap To Place Tack

Or Spring Loaded PowerTack™ Delivery Handle



Pull End To Activate PowerTack™ Delivery Handle



Load Tack On Tip Of Handle



Push Button To Place Tack

Meisinger "Pin Control" Membrane Fixation Kit

- **Fixation Of Resorbable Membranes, Non-Resorbable Membranes & Titanium Mesh**
- **Hybrid Screw & Pin Easily Removed With Hand Screwdriver**

Designed By Dr. Istvan Urban

Includes:

- 2 1.8mm Latch-Type Tri-Spade Pilot Drills
- 2 0.8mm Latch-Type Twist Perforation Drills
- 2 0.6mm Latch-Type Twist Perforation Drills
- 1 Placement Instrument
- 1 Fixation Handle
- 1 Hand Screwdriver



Basic Kit With 10 Pins

#MEIS-BMPBA 1175.00 KIT

Master Kit With 34 Pins

#MEIS-BMP00 1685.00 KIT

Replacement Pins - 3.5mm Length

10 Pins, Non-Sterile

#MEIS-MP10..... 228.50 PKG OF 10

SPECIAL



Scan For
60th Second Product Video™



Salvin Membrane Tack Kit With Mallet & Placement Instrument

MOST POPULAR

Includes:

- 25 Tacks
- Autoclavable Tack Block
- Tack Placement Instrument
- Tack Mallet
- Tack Drill For Dense Cortical Bone

REGULAR PRICE IS 925.20

SAVE 226.20

#TACK-KIT-TI 699.00 KIT



For The Salvin Membrane Tack Kit With Mallet



Titanium Tack Forcep

#TACK-FORCEP 137.50 EACH



Optional Offset Tack Instrument

#TACK-INSTR-OFFSET 207.40 EACH



Replacement Tack Pilot Drill

#TACK-DRILL 68.90 EACH

SPECIAL

Scan For
60th Second Product Video™



Salvin Membrane Tack Kit With PowerTack™ Handle

Includes:

- 20 Tacks
- Autoclavable Tack Block
- PowerTack™ Delivery Handle
- Tack Forcep

REGULAR PRICE IS 1361.10

SAVE 275.10

#POWERTACK-KIT 1086.00 KIT



If You Already Have Our Regular Tack Kit

Includes:

PowerTack™ Delivery Handle, Autoclavable Tack Block

#POWERTACK 945.00 HANDLE & BLOCK

Membrane Placement Instruments



**Membrane Forcep
With Ring Tip**

Non-Serrated For Membranes

#MEMB-RING-FORC

.....121.90 EACH

MOST POPULAR

**Membrane Forcep
"Paddle Tip"**

Non-Serrated For Membranes

#MEMB-PADDLE-FORCEP

.....90.20 EACH



**Membrane Forcep - Long
"Paddle Tip"**

Non-Serrated For Membranes

#MEMB-PADDLE-LONG

.....90.20 EACH



MOST POPULAR

SuperCut Iris Curved 11cm

Ideal For Delicate Membranes

#SC IRIS CU

.....142.90 EACH



MOST POPULAR

Membrane Placement Instrument

Micro-Point Tip For Precise Placement & Manipulation Of Membranes

#MEMB-PLACER.....83.30 EACH



Membrane Positioner

Stabilize & Condense

Designed By Dr. Larry Towning

#MEMB-CONDENSER128.30 EACH

Comprehensive Lateral & Crestal Sinus Kit

Minimizes Risk Of Membrane Perforation



Complete Crestal & Lateral Sinus Kit

#SINUS-COMP-KIT 1595.00 EACH



*Diamond Tip
Minimizes Risk
Of Perforation*



*Trephinated Bur
Cuts & Detaches
Bony Wall*



*Dome Shape
Minimizes Risk
Of Perforation*



*Irrigation Lifts
Membrane*



*Syringe Holder Helps Control
Irrigation Flow*



6 Crestal Drills



8 Drill Stops



4 Crestal Burs



2 Lateral
Trephinated Drills



2 Lateral Dome
Shaped Burs



*Irrigation Syringe
With Tubing*



8 Irrigation Tips



Syringe Holder



1 Sinus Depth Gauge/Packer

**Replacement Irrigation Syringe
With Tubing**

#SINUS-MAS2-01-09 9.40 EACH



Lateral Reamer 6.5mm x 2mm

#SINUS-REAM-6.5X2 203.00 EACH



Lateral Reamer 6.5mm x 3mm

#SINUS-REAM-6.5X3 203.00 EACH




Concave Tip Sinus Lift Osteotomes

A conventional osteotomy is made to the floor of the sinus. The Sinus Osteotome is inserted into the osteotomy & gently tapped apically to separate & lift the cortical bone at the sinus floor 2mm-4mm.



MOST POPULAR
Offset Sinus Lift Set

- 5 Offset Osteotomes With Adjustable Stops
 - FREE Stainless Steel Cassette
- #OST-SINUS-OFF-PKG 756.00 SET



Straight Sinus Lift Set

- 5 Straight Osteotomes With Adjustable Stops
 - FREE Stainless Steel Cassette
- #OST-SINUS-STR-PKG 756.00 SET





Straight + Offset Sinus Lift Set

- 5 Straight Osteotomes With Adjustable Stops
 - 5 Offset Osteotomes With Adjustable Stops
 - FREE Stainless Steel Cassette
- #OST-SINUS-BOTH-PKG 1359.00 SET

Extended Length Osteotomes For Surgical Guides

Designed By Dr. Jack Krauser





*NobelGuide™ is a
Trademark of NobelBiocare*



Extended Length Osteotome Set

- 5 Offset Osteotomes With Adjustable Stops
- FREE Stainless Steel Cassette

#GUIDE-TOME-PKG 761.00 SET



SALVIN
AutoTome®
Automatic Osteotome Kit
On Page 90


Meisinger Crestal Sinus Lift Kit

- *Depth Stops Prevent Puncturing Of The Sinus Membrane*
- *Concave Tipped Drills*
- *Crestal Drills Collect Bone Chips*
- *Can Be Used With Any Implant System*



Complete Crestal Sinus Lift Kit

6 Crestal Drills, 11 Drill Stops, 4 Pilot Drills, 1 Cylinder Bur, 2 Condensers & 1 Depth Gauge
#CRESTAL-LIFT-CONTROL 1950.00 KIT



Pilot Drill

Lift Membrane

Collect Bone

Depth Gauge

Place Graft



Basic Crestal Sinus Lift Kit 2 MOST POPULAR

4 Crestal Drills (3.3, 3.6, 3.8, 4.1mm), 7 Drill Stops (3-9mm), 1 Pilot Drill & 1 Bone Graft Condenser
#MEIS-BCLBA2997.25 KIT



Depth Gauge For Crestal Sinus Lift Kit

12mm Depth Gauge
#MEIS-CL021154.70 EACH

Internal Sinus Manipulation Set

Minimally Invasive Sinus Floor Augmentation

Designed By Dr. Jason M. Yamada



Small Rounded Head Bone Graft Spreader & Packer
#ISM-1 131.20 EACH



Large Rounded Head Bone Graft Spreader & Packer
#ISM-2 131.20 EACH



Sinus Lift Instrument -
 Keep Sharp Working End Against Bony Floor
#ISM-3 131.20 EACH



Sinus Lift Instrument -
 Keep Sharp Working End Against Bony Floor
#ISM-4 131.20 EACH



Membrane Perforation Repair Instrument
#ISM-5 131.20 EACH

SPECIAL

Internal Sinus Manipulation Set




Set Of 5 Instruments With
FREE Stainless Steel Cassette
#ISM-KIT
580.00 SET

Kerrison Rongeurs For Widening The Sinus Window



Low Profile Head
For Visibility
& Safe Access



Kerrison Rongeur - 4mm Bite

10cm Long Beak

#KERRISON 630.00 EACH



MOST POPULAR

Kerrison Rongeur - 3mm Bite

8cm Long Beak

#KERRISON-SMALL 630.00 EACH

Serrated Tip Sinus Window Bone Tampers

To Precisely Separate Sinus Window From Surrounding Bone

MOST POPULAR

4mm Sinus Window / Bone Tamper

#BONE-TAMPER-4MM. 84.50 EACH



7mm Sinus Window / Bone Tamper

#BONE-TAMPER-7MM. 84.50 EACH



SPECIAL **Set Of 2 Bone Tampers**
 #BONE-TAMPER-SET **149.00 SET OF 2**

Sinus Scoring Diamonds



Sinus Scoring Diamond-Latch Type

Fits Latch Type Implant Handpieces

#SINUS-SCORE-DIAM-RA

INDIVIDUALLY 18.75 EACH

10 OR MORE 16.40 EACH



Sinus Scoring Diamond-HP Type


Fits Straight Implant & Oral Surgery Handpieces

#SINUS-SCORE-DIAM-HP


INDIVIDUALLY 18.75 EACH

10 OR MORE 16.40 EACH

Sinus Lift Instruments



Sinus Lift Instruments



SPECIAL

Salvin Basic Sinus Set

MOST POPULAR



Basic Set Includes:

- 4mm Bone Tamper
 - Small Serrated Graft Packer
 - 4.5mm Graft Syringe
 - Salvin Sinus Curette #1
 - Salvin Sinus Curette #5
 - Salvin Sinus Curette #6
 - Salvin Sinus Freer
 - Salvin Sinus Freer Curved
- (Or Your Choice Of Any 5 Curettes)
- Plus **FREE** Cassette

**#SINUS-KIT-WINDOW
844.00 SET**



Bone Expander Kit



- Alternative To Osteotomes For The Expansion & Condensing Of The Atrophic Mandible & Maxilla
- Use With Ratchet Wrench Or Low Speed Handpiece
- Improve Stability, Maintain Bone Density & Increase Fixation
- Can Be Used With Implants Up To 16mm Long
- Contains 5 Bone Expanders, 1 Ratchet Wrench, 1 Pilot Drill, 1 Ratchet Wrench Connector, 1 Latch-Type Connector, 1 Thumb Knob, 1 Saw Disc

#BONE-EXPANDER-KIT-2737.50 KIT



PHARMACEUTICALS
AT YOUR FINGERTIPS

NEW!



*We Now Offer A Full Line Of Pharmaceuticals.
Scan The QR Code Or Call Us For More Information.
We'll Take Great Care Of You!*

For More Info
SCAN HERE



Orders for pharmaceuticals and the applicable documentation are received, processed, fulfilled, and distributed by PuraGraft and/or Medical Purchasing Solutions, LLC.


Bone Spreading Osteotomes For Implant Placement In Soft Bone



MOST POPULAR

Offset Bone Spreading Set

- 6 Offset Spreaders With Adjustable Stops
 - FREE Stainless Steel Cassette
- #OST-SPREAD-OFF-PKG. . . . 899.00 SET**



Straight + Offset Bone Spreading Set


- 6 Offset Spreaders With Adjustable Stops
 - 6 Straight Spreaders With Adjustable Stops
 - FREE Stainless Steel Cassette
- #OST-SPREAD-BOTH-PKG 1683.00 SET**



Straight Bone Spreading Set

- 6 Straight Spreaders With Adjustable Stops
 - FREE Stainless Steel Cassette
- #OST-SPREAD-STR-PKG. . . . 899.00 SET**



Double Sided Ridge Splitting / Expanding Chisels - Offset



Ridge Splitting Offset Chisel Set

- Set Of 4 Ridge Splitting Offset Chisels
 - Stainless Steel Cassette
- #RIDGE-SPLIT-OFF-SET 576.00 SET

Double Sided Ridge Splitting / Expanding Chisels - Straight




MOST POPULAR

Ridge Splitting Chisel Set


- Set Of 4 Ridge Splitting Chisels
 - Stainless Steel Cassette
- #RIDGE-SPLIT-SET 576.00 SET

Single Bevel Ridge Splitting / Expanding Chisels



Meisinger Mandibular Bone Spreading Kit

Horizontal Widening Up To 5mm




MOST POPULAR

Comprehensive Ridge Splitting Kit

3 Spreaders, 5 Burs, 3 Saw Discs, 2 Thumb Drivers & 1 Ratchet
#CREST-CONTROL 1888.65 KIT

Meisinger Maxillary Bone Spreading Kit

Minimally Invasive Alternative To Osteotomes



MOST POPULAR

Comprehensive Bone Spreading Kit

- Controlled Dilatation Of Horizontally Resorbed Bone Up To 4mm
- 6 Spreaders, 7 Burs, 2 Saw Discs, 1 Adapter & 1 Ratchet

#SPLIT-CONTROL 1340.00 KIT

Basic Bone Spreading System

- Controlled Dilatation Of Horizontally Resorbed Bone Up To 4mm
- 3 Spreaders, 2 Burs, 1 Saw Disc, 2 Adapters & 1 Ratchet

#SPLIT-CONTROL-BASIC 840.00 KIT

Splitting / Spreading Combination Kit


Mandibular Ridge Splitting + Maxillary Bone Spreading Kit

- Comprehensive Ridge Splitting System Components
- Comprehensive Bone Spreading System Components
- Large Organizing Cassette
- Plus **FREE** Motor Driver Attachments For Spreaders

TOTAL PRICE OF COMPONENTS IS 2705.00 **SAVE 200.00**
#MEIS-CREST-SPLIT 2554.40 KIT



Rongeurs



MOST POPULAR

Friedman Rongeur

#FRIED
..... 299.00 EACH



**Micro Friedman
Right Angle Rongeur**


#MICRO FRIED R
..... 299.00 EACH



MOST POPULAR

**Micro Friedman
Rongeur**

#MICRO FRIED
..... 299.00 EACH



Tip Laser
Marked At
1-10mm

MOST POPULAR

Blumenthal Rongeur

#BLUM
..... 299.00 EACH



**Blumenthal Rongeur
Right Angle**


#BLUM-RA
..... 299.00 EACH



MOST POPULAR

**Mini Blumenthal
Rongeur**

#MINI BLUM
..... 299.00 EACH



**Mead 1-A Oral Surgery
Rongeur**

#MEAD-1-A RONG
..... 299.00 EACH

Kerrison Rongeurs For Widening The Sinus Window



Kerrison Rongeur - 4mm Bite

10cm Long Beak

#KERRISON 630.00 EACH



MOST POPULAR

Kerrison Rongeur - 3mm Bite

8cm Long Beak

#KERRISON-SMALL 630.00 EACH



Low Profile Head
For Visibility
& Safe Access



Double Action Rongeurs - Increased Leverage & Effortless Cutting



3.5mm Bite


DOUBLE ACTION

MOST POPULAR

**Beyer Rongeur
17.8cm Long**

Double Action For Increased Leverage

#BEYER. 456.00 EACH



5.0mm Bite

DOUBLE ACTION

**Ruskin Rongeur
17.8cm Long**

Double Action For Increased Leverage

#RUSKIN. 456.00 EACH

Retractors



MOST POPULAR

Implant Retractor

Use As Upper Or Lower Vestibular Retractor

#IMPLANT-RETRACT
..... 191.90 EACH

Implant Retractor Ergonomic Handle

Use As Upper Or Lower Vestibular Retractor

#IMPLANT-RETRACT-ERGO
..... 247.50 EACH



Molt Mouth Gag

#MOLT GAG 246.50 EACH

Molt Mouth Gag - Small

#MOLT-GAG-SMALL
..... 246.50 EACH



Replacement Guards For Molt Mouth Gag

Autoclavable Silicone Rubber

#MOLT-GAG-GUARDS 14.90 SET OF 2



MOST POPULAR

Plastic Cheek Retractor Single Handle Autoclavable

#CHEEK-SINGLE
SET OF 2 37.90 SET
2 OR MORE SETS 34.90 SET



Plastic Cheek Retractor Double Ended Autoclavable

#CHEEK
SET OF 2 56.90 SET
2 OR MORE SETS 53.10 SET



Orringer Retractor

Spring Steel Flexes For Comfortable Fit

#ORRINGER 88.70 EACH
2 OR MORE 82.90 EACH



Orringer Retractor-Small

Spring Steel Flexes For Comfortable Fit

#ORRINGER-SMALL . . 88.70 EACH
2 OR MORE 82.90 EACH



Retractors



**Cawood-Minnesota
Retractor**

#CAWOOD-RETRACT
INDIVIDUALLY 42.60 EACH
2 OR MORE 38.80 EACH



MOST POPULAR
Minnesota Retractor

#MINN
INDIVIDUALLY 42.60 EACH
2 OR MORE 38.80 EACH

**SALVIN
Comfort
Handle**



**Minnesota Retractor
Ergonomic Handle**

#MINN-ERGO
INDIVIDUALLY 140.20 EACH
2 OR MORE 129.90 EACH



**Dog Bone Sinus Retractor
- 18cm**

Sinus Flap Retractor
#DOG-BONE 65.90 EACH



MOST POPULAR
Columbia Retractor

Excellent For Surgery & Photography
#COL RETR 53.20 EACH
2 OR MORE 49.30 EACH



Austin Retractor

#AUSTIN-RETRACTOR 69.70 EACH



*Easily Bent By Hand
To Desired Shape*



Bendable Retractor Stainless Steel

20cm Long x 18mm Wide
#RETRACTOR-BENDABLE 57.90 EACH

Retractors



MOST POPULAR

Weider Tongue Retractor

44mm Wide At Tip

#WEIDER 81.30 EACH
2 OR MORE 75.40 EACH

Weider Tongue Retractor - Small

31mm Wide At Tip

#WEIDER-SMALL 81.30 EACH
2 OR MORE 75.40 EACH

Salvin Zygoma Retractor

- Increases Field Of Vision
 - Directs The Drilling
 - Prevents Perforation Of The Orbit
 - Protects Soft Tissue While Drilling Superior Cortex
- #ZAGA-2 209.00 EACH




Photo: Dr. Thomas Balshi

SALVIN
Comfort
Handle





Photo: Dr. Michael Pikos

SALVIN
Comfort
Handle

MOST POPULAR



Laster 3rd Molar Retractor

#LASTER-RETRACT
INDIVIDUALLY 79.80 EACH
2 OR MORE 75.60 EACH

Balshi Zygomatic Retractor

Designed By Dr. Thomas Balshi
#BALSHI-ZYGO-RETRACT 193.90 EACH

Pikos Ramus Retractor

Designed By Dr. Michael Pikos
16mm Wide At Tip
#PIKOS-RETRACTOR 193.90 EACH

Retractors



MOST POPULAR

Henahan Retractor Large

19.5cm

- #HENAHAN-LARGE 99.90 EACH
- #3 OR MORE 94.70 EACH

Seldin 23

- #SELDIN 99.90 EACH
- #3 OR MORE 94.70 EACH

Mini Seldin

Round Handle For Comfort

- #MINI-SELDIN 99.90 EACH
- #3 OR MORE 94.70 EACH

Misch Sinus Retractor

Designed By Dr. Carl Misch

- #MISCH-SINUS-RETRACT 97.90 EACH



MOST POPULAR

Dual View Photo Mirror

View Buccal, Occlusal & Lingual Surfaces At Once

- #PHOTO-DUALVIEW 105.90 EACH



Latexfree

MOST POPULAR

Autoclavable Silicone Mouth Props

Autoclavable To 275°F Or 134°C
Dry Heat To 350°F Or 160°C

- #MOUTH-PROP-CHILD
SET OF 2 40.90 SET
2 OR MORE SETS 37.70 SET
- #MOUTH-PROP-ADULT
SET OF 2 40.90 SET
2 OR MORE SETS 37.70 SET

Periosteal Elevators



MOST POPULAR

Molt 2/4

#MOLT 2/4
 84.80 EACH
 3 OR MORE
 80.10 EACH




Molt 2 Angled Double-Ended

#MOLT-2-DE
 84.80 EACH
 3 OR MORE
 80.10 EACH



Molt 4 Angled Double-Ended

#MOLT-4-DE
 84.80 EACH
 3 OR MORE
 80.10 EACH



Periosteal-Tome™

Periosteal Elevator/
 Periotome Combo
 #PERIOSTEAL-TOME
 108.20 EACH
 3 OR MORE
 102.60 EACH



MOST POPULAR

Prichard Elevator

#PRICHARD
84.80 EACH
 3 OR MORE
80.10 EACH



MOST POPULAR

Seldin 23

#SELDIN
 103.80 EACH
 3 OR MORE
 97.50 EACH



Mini Seldin

Round Handle
 For Comfort
 #MINI-SELDIN
 103.80 EACH
 3 OR MORE
 97.50 EACH



Freer

#FREER
 84.80 EACH
 3 OR MORE
 80.10 EACH



Tip Marked
 at 1 - 10mm
 For Spacing

Periosteal Elevators



MOST POPULAR

Buser

#BUSER
.....84.80 EACH
3 OR MORE
.....80.10 EACH




Goldman-Fox Elevator

17cm
#GOLDMAN-FOX-ELEV
.....84.80 EACH
3 OR MORE
.....80.10 EACH



Hirschfeld #20

#HIRSCHFELD-20
.....84.80 EACH
3 OR MORE
.....80.10 EACH



MOST POPULAR

Woodson

#WOODSON
.....84.80 EACH
3 OR MORE
.....80.10 EACH



Quinn-Titanium

Gold Titanium Tips For
Use Around Implants
#Q-ELEV-TI
.....105.90 EACH
3 OR MORE
.....98.40 EACH



MOST POPULAR

Quinn

#Q-ELEV
.....84.80 EACH
3 OR MORE
.....80.10 EACH



24G Elevator

#EP24G
.....84.80 EACH
3 OR MORE
.....80.10 EACH



Benque

#BENQUE
.....84.80 EACH
3 OR MORE
.....80.10 EACH

Periosteal Elevators



Molt 4

#MOLT-4
 94.90 EACH
 3 OR MORE
 89.90 EACH



Molt 2

#MOLT-2
 94.90 EACH
 3 OR MORE
 89.90 EACH



West

#WEST
 85.60 EACH
 3 OR MORE
 80.80 EACH



Henahan Retractor Large

19.5cm
 #HENAHAN-LARGE
 99.90 EACH
 3 OR MORE
 94.70 EACH



MOST POPULAR

Molt 9

#MOLT 9
 85.80 EACH
 3 OR MORE
 80.90 EACH



Molt 9 With Suture Hole

#MOLT 9-SUTURE
 85.80 EACH
 3 OR MORE
 80.90 EACH



Molt 9 With 15 Aspirating Holes


#MOLT 9-ASP
 85.80 EACH
 3 OR MORE
 80.90 EACH



SuperCut Serrated Scissors Cut Effortlessly Through Tissue Or Suture

BUYER'S GUIDE - SUPERCUT SCISSORS

- ▶ **Razor Sharp Specially Machined Blades Mesh Perfectly To Provide Superior Shearing Force At The Tips**
- ▶ **Micro Serrations On Lower Blade Grip Tissue For Slip-Free Cutting**



SUPERCUT

**Goldman Fox Scissors
13.2cm Curved**

#SC GF SCISS142.90 EACH
3 OR MORE135.10 EACH



Double Curved

MOST POPULAR

SUPERCUT

**LaGrange Scissors
11.5cm Double Curved**

#SC-LAGRANGE . . .142.90 EACH
3 OR MORE135.10 EACH




MOST POPULAR

SUPERCUT

**Iris Scissors
11cm Curved**

#SC IRIS CU142.90 EACH
3 OR MORE135.10 EACH



SUPERCUT

**Iris Scissors
11cm Straight**

#SC IRIS ST142.90 EACH
3 OR MORE135.10 EACH

SALVIN Black Kryptonite™ Ceramic Coated Serrated Surgical Scissors

BUYER'S GUIDE - SALVIN BLACK KRYPTONITE™ SCISSORS

- ▶ **Long Lasting, Razor Sharp Precise Cutting**
- ▶ **Lasts 4-6 Times Longer Than Conventional Scissors**
- ▶ **Non-Glare Finish For Better Tissue Visibility**



**SALVIN
Black
Kryptonite™**

**Goldman Fox Scissors
12.5cm Curved**

#Kryp-GF-SCISS . . . 259.90 EACH
3 OR MORE 247.30 EACH



Double Curved

**SALVIN
Black
Kryptonite™**

**LaGrange Scissors
11.5cm Double Curved**

#Kryp-LAGRANGE . . .259.90 EACH
3 OR MORE 247.30 EACH



MOST POPULAR

**SALVIN
Black
Kryptonite™**

**Iris Scissors
11.5cm Curved**

#Kryp-IRIS-CU259.90 EACH
3 OR MORE247.30 EACH



MOST POPULAR

**SALVIN
Black
Kryptonite™**

**Castroviejo Micro Scissors
12cm Curved**

#Kryp-CASTRO420.00 EACH
3 OR MORE391.80 EACH

Surgical & Suture Scissors



Goldman Fox 12.7cm

Serrated For Slip-Free Cutting
 #GF SCISS. 109.20 EACH
 3 OR MORE 104.60 EACH



Double Curved

MOST POPULAR

LaGrange 11.5cm

Serrated For Slip-Free Cutting
 #LAGRANGE. 109.20 EACH
 3 OR MORE 104.60 EACH



Iris Scissors Straight 11.5cm

#IRIS-SCISS. 76.30 EACH
 3 OR MORE 72.20 EACH



MOST POPULAR

Iris Scissors Curved 11.5cm

#IRIS-SCISS-CURVED 76.30 EACH
 3 OR MORE 72.20EACH



Long Suture Scissors 15cm

Hooked End To Lift Suture
 #LONG 143.90 EACH
 3 OR MORE 135.60 EACH



Littauer Angled Suture Scissors 11.5cm

Hooked End To Lift Suture
 #LITTAUER-ANGLED . 115.20 EACH
 3 OR MORE 108.90 EACH



MOST POPULAR

Littauer Suture Scissors 11.5cm

Hooked End To Lift Suture
 #LITTAUER 91.40 EACH
 3 OR MORE 87.70 EACH



Metzenbaum-Mayo Dissecting Scissors 14.5cm

Blunt Tip To Dissect Tissue
 #METZ-5. 104.80 EACH
 3 OR MORE 97.60 EACH

Surgical & Suture Scissors



**Locklin Scissors Angled
15.5cm**

Serrated For Slip-Free Cutting
#LOCKLIN 142.90 EACH
3 OR MORE 135.80 EACH



**Wagner Scissors Angled
11.5cm**

#WAGNER-SCISS-ANGLED
. 119.00 EACH
3 OR MORE 113.40 EACH



Kelly Curved 15.9cm

Serrated For Slip-Free Cutting
#KELLY CURV 108.40 EACH
3 OR MORE 100.90 EACH



MOST POPULAR
Dean 17.2cm

Serrated For Slip-Free Cutting
#DEAN 142.90 EACH
3 OR MORE 135.80 EACH

Featherlite Micropoint Scissors



*Takes The "Wince" Out Of Suture Removal,
Even In Sensitive Tissue*



MOST POPULAR

Featherlite 14cm Curved

#F-LITE5CURV 335.00 EACH
3 OR MORE 317.30 EACH



*Opposing Spring Tips Maintain Force
Between Blades For Incredible
Cutting Without Pinching Tissue*

Suture Scissors / Suture Forcep Combination Instrument

Cut & Remove Suture With One Instrument



#SUT-REM-FORCEP 335.00 EACH
3 OR MORE 317.30 EACH



*Engage Suture On Right
Side Of Knot...*

*...Cut, Squeeze & Pull
To Remove Suture...*

**Tungsten Carbide Tipped Needle Holders
For Larger Needle 3-0 & 4-0 Sutures**



Designed For 3-0 To 4-0 Sutures

Accommodates Larger Needles Without Damaging The Locking Clasp & Instrument Tip



Serrated Tips

TUNGSTEN CARBIDE

FOR 3-0 TO 4-0 SUTURES

MOST POPULAR

**Castroviejo Needle Holder
Accommodates 3-0 To 4-0 Sutures**

Spring Handle For Precise Control
14cm

#CASTRO-3.0-4.0 467.00 EACH



Serrated Tips

TUNGSTEN CARBIDE

FOR 3-0 TO 4-0 SUTURES

**Jacobson Needle Holder
Accommodates 3-0 To 4-0 Sutures**

Spring Handle For Precise Control
17.8cm

#JACOBSON-3.0-4.0 551.00 EACH

Tungsten Carbide Tipped Needle Holder With Side Cutting Scissors

Scissor On Side For Faster Suturing Without Changing Instruments




TUNGSTEN CARBIDE


FOR 2-0 TO 4-0 SUTURES



Hold Suture Taut In Slot...



...Slide Down To Knot...



...Squeeze To Cut,
Leaving A 2mm Tag

Side Cutting Mayo Hegar 6"

Carbide Tipped With Scissors
15.3cm

#SIDE CUT MAYO 339.00 EACH

3 OR MORE 322.60 EACH

Tungsten Carbide Tipped Castroviejo & Micro Needle Holders




BUYER'S GUIDE - CARBIDE TIPPED NEEDLE HOLDERS	
▶	Superior Grip Holds Needle Without Slipping
▶	Less Hand Fatigue & Cramping During Long Surgeries
▶	Last 4-6 Times Longer Than Standard Instruments
▶	Can Be Economically Re-Tipped



Castroviejo Needle Holder

Spring Handle For Precise Control - 14cm
 #CASTRO 467.00 EACH
 #CASTRO-CURVED. 467.00 EACH



Jacobson Needle Holder

Spring Handle For Precise Control - 17.8cm
 #JACOBSON 551.00 EACH
 #JACOBSON-CURVED 551.00 EACH



ThumbLok Snagless Needle Holder

17.8cm - Straight Tips
 #NH-ROUND-7-TCL-R. 551.00 EACH



Round Jacobson Needle Holder

Spring Handle For Precise Control - 17.8cm
 #JACOBSON-ROUND. 551.00 EACH



Micro Jacobson Needle Holder

17.8cm - Straight Tips
 #MICRO-JACOBSON-STR. 551.00 EACH
 17.8cm - Curved Tips
 #MICRO-JACOBSON-CUR 551.00 EACH



Micro Castro Needle Holder

13.3cm - Straight Tips
 #MICRO-CASTRO-STR 551.00 EACH
 13.3cm - Curved Tips
 #MICRO-CASTRO-CUR. 551.00 EACH

Tungsten Carbide Tipped Needle Holders



Serrated Tips
TUNGSTEN CARBIDE
 FOR 3-0 TO 5-0 SUTURES

Crile Wood Needle Holder - Curved

15.3cm
 #CRILE-CURVED 202.70 EACH
 3 OR MORE 193.50 EACH



Serrated Tips
TUNGSTEN CARBIDE
 FOR 3-0 TO 5-0 SUTURES

MOST POPULAR
Crile Wood Needle Holder

15.3cm
 #CRILE 202.70 EACH
 3 OR MORE 193.50 EACH



Serrated Tips
TUNGSTEN CARBIDE
 FOR 2-0 TO 4-0 SUTURES

MOST POPULAR
Mayo Hegar Needle Holder

15.3cm
 #MAYO HEGAR 202.70 EACH
 3 OR MORE 193.50 EACH



Serrated Tips
TUNGSTEN CARBIDE
 FOR 3-0 TO 5-0 SUTURES

Crile Wood Needle Holder - Long

17.8cm
 For Suturing Extracted 3rd Molar Sockets
 #CRILE-LONG 202.70 EACH
 3 OR MORE 193.50 EACH



Serrated Tips
TUNGSTEN CARBIDE
 FOR 3-0 TO 5-0 SUTURES

Mathieu Needle Holder

"Squeeze Through" Quick Release
 12.7cm
 #MATHIEU 276.90 EACH
 3 OR MORE 264.70 EACH



Serrated Tips
TUNGSTEN CARBIDE
 FOR 5-0 TO 7-0 SUTURES

SALVIN
Comfort
Handle

**Castroviejo Needle Holder
 Straight Handle**

12.7cm
 #CASTRO-STR-HANDLE 514.00 EACH



Sutures On Pages 34-38

Tungsten Carbide Tipped Needle Holders




Serrated Tips
TUNGSTEN CARBIDE
FOR 4-0 TO 6-0 SUTURES

MOST POPULAR

Par Perio Needle Holder

12.7cm
 Delicate Serrated
 #PAR 199.80 EACH
 3 OR MORE 187.30 EACH



Serrated Tips
TUNGSTEN CARBIDE
FOR 4-0 TO 6-0 SUTURES

Derf Needle Holder

12.1cm
 #DERF 199.80 EACH
 3 OR MORE 187.30 EACH



Serrated Tips
TUNGSTEN CARBIDE
FOR 3-0 TO 5-0 SUTURES

Olsen Hegar Scissors/Needle Holder

14cm
 #OH 5 235.00 EACH
 3 OR MORE 225.50 EACH



Serrated Tips
TUNGSTEN CARBIDE
FOR 3-0 TO 5-0 SUTURES

Olsen Hegar Scissors/Needle Holder

16.5cm
 #OH 6 235.00 EACH
 3 OR MORE 225.50 EACH

Salvin Non-Glare Featherweight Microsurgery Instruments
Custom Designed For Periodontal Microsurgery



- ▶ Light & Long
- ▶ Smooth Tips
- ▶ Rounded Handles
- ▶ Non-Locking
- ▶ Black Non-Glare Coating



"The Salvin microsurgical instruments exhibit a semicircular grip to enhance handling and haptic feedback; these instruments are light in weight to prevent muscle fatigue and have been designed without locking mechanisms to prevent undesirable traumatic instrument manipulation. The dark surface coating minimizes glaring when working under coaxial illumination. Once you use these finely crafted instruments it will be hard to work with any other set!"

Dr. Diego Velasquez
 Periodontist, Fenton, MI



"The Salvin microsurgery set combines the features needed in high quality surgical instrumentation. The non-reflective black coating effectively eliminates glare during surgeries, allowing for better visibility and precision. The lightweight design enhances our tactile control, and the smooth tips are conducive for manipulation of micro sutures. Salvin listened to what we needed and delivered!"

Dr. Irene Marron-Tarrazzi
 Periodontist, Miami, FL



"When performing perio plastic surgery under magnification and illumination, using specialized instrumentation is critical to my success. The Salvin microsurgery instruments check all of the requirements I need: Light, long, non-reflective, rounded handles and manufactured with the highest quality. This is exactly what I use in my practice."

Dr. John Russo
 Periodontist, Sarasota, FL



Non-Glare Micro Tunneling - 1.2mm

#MICRO-BLK-TUNNEL-1.2
 107.80 EACH

Non-Glare Micro Tunneling - 2.0mm

#MICRO-BLK-TUNNEL-2.0
 107.80 EACH

Non-Glare Micro Tunneling - 2.8mm

#MICRO-BLK-TUNNEL-2.8
 107.80 EACH

Double-Ended Micro Tunnel Elevators Help Reflect & Elevate Tissue From Recipient Site In Precise Micro Flaps



Non-Glare Lucas 85 Bone Curette

Sharp Tips Efficiently Debride Sockets
 #LUCAS-85-BLK
 92.20 EACH



Non-Glare Prichard Elevator

Double-Ended Micro Tunnel Elevators Help Reflect & Elevate Tissue From Recipient Site In Precise Micro Flaps
 #PRICHARD-BLK
 87.70 EACH



Non-Glare Perio Probe

#PERIO-PROBE-BLK
 102.80 EACH



Non-Glare Papilla Elevator

- Double-Ended Elevator
- Rounded Tip For Blunt Tissue Dissection

#PAPILLA-ELEV-BLK
 102.80 EACH

SPECIAL



Salvin Non-Glare Featherweight Microsurgery Set

Includes:

- Non-Glare Micro Tunneling - 2.0mm
 - Non-Glare Micro Tunneling - 2.8mm
 - Non-Glare Tissue Forcep - Curved OR
 - Non-Glare Tissue Forcep - Straight
 - Non-Glare Castroviejo
 - Non-Glare Micro Needle Holder - Curved
 - Non-Glare Micro Scissor - Curved
 - Non-Glare Lucas 85 Bone Curette
 - Non-Glare Prichard Elevator
 - Non-Glare Perio Probe
 - Non-Glare Papilla Elevator
 - Non-Glare Micro Blade Handle
 - Non-Glare Micro Mouth Mirror - 5mm
 - Non-Glare Mirror Handle
 - Plus **FREE** Instrument Cassette
- #NON-GLARE-MICRO-KIT**
.....**.2390.00 SET**



**Non-Glare Castroviejo
Micro Needle Holder -
Curved**

- Light (<20 grams)
And Long (18cm)
For Easy Handling Of
Micro Sutures
 - Smooth Tips Grasp
Delicate Sutures, While
The Rounded Handles
Allow For Precise Control
- #MICRO-BLK-CASTRO-CUR**
..... **528.00 EACH**



**Non-Glare
Micro Scissor - Curved**

- Sharp, Slightly Curved,
Long Scissors For Fine
Resections Of Soft Tissues &
Trimming Of Mucosal Grafts
- #MICRO-BLK-SCISS-CUR**
..... **475.00 EACH**



**Non-Glare
Mirror Handle**

- #MIRROR-HNDL-BLK**
..... **67.40 EACH**



**Non-Glare
Micro Mouth Mirror - 5mm**

- #MICRO-MIRR-5-BLK**
..... **151.00 EACH**



**Non-Glare
Micro Mouth Mirror
- 3x6mm**

- #MICRO-MIRR-3X6-BLK**
..... **151.00 EACH**




**Non-Glare
Tissue Forcep -
Curved**

- #MICRO-BLK-FORCEP-CUR**
..... **475.00 EACH**



**Non-Glare
Tissue Forcep -
Straight**

- #MICRO-BLK-FORCEP-STR**
..... **475.00 EACH**



**Non-Glare
Suture Pliers**

- #MICRO-BLK-SUTURE-PLIER**
..... **424.00 EACH**



**Non-Glare
Micro Blade Handle**

- Fits Micro "Beaver-Type"
Scalpel Blades
 - Collet Locks Blade In Place
- #MICRO-BLK-BLADE-HNDL**
..... **140.00 EACH**

Kern Tunneling Instruments

Designed By Dr. Ricardo Kern

- *Minimally Invasive Periodontal Microsurgery*
- *Avoids Need For Large Flap*
- *Anatomical Silicone Handles With Stainless Steel Tips*



"Precise and delicate, these instruments have been designed to ensure that you have no limitations in performing to the best of your ability. In this kit you'll find the sophistication you need to explore new surgical possibilities and maximize the performance of current mucogingival reconstruction techniques!"

Dr. Ricardo Kern
Periodontist, Brazil

SPECIAL



Kern Tunneling Set

NEW!

Includes:

- Kern Tunneler #1
- Kern Tunneler #2
- Kern Tunneler #3
- Kern Tunneler #4
- Kern Tunneler #5
- 1 Cassette

#KERN-TUNNELING-KIT
795.00



Tunnelling Blades And Handles



Aluminum Tunnel Blade Scalpel Handle

#TUNNEL-BLADE-HNDL-A
..... 79.00 EACH



Titanium Tunnel Blade Scalpel Handle

#TUNNEL-BLADE-HNDL-TI
..... 109.00 EACH



Stingray Tunnel Blade 1

- "Back Action" Greatly Facilitates The Division Of Muscle Fibers In The Mucosal Or Periosteal Plane.
 - Pack Of 5 – Individually Wrapped Sterile
- #TUNNEL-BLADE-STINGRAY 69.00 PACK OF 5



Stingray Tunnel Blade 2

- Micro-Blade For Access To Sulcular Tunneling. Small & Delicate, Perfect For Thin & Slender Phenotypes Even In The Narrow Contours Of The Lower Incisors.
 - Pack Of 5 – Individually Wrapped Sterile
- #TUNNEL-BLADE-TB2 69.00 PACK OF 5

Bendable Miniature & Micro Scalpel Blades

Uncompromised Sharpness & Cutting Ability


Blades Are Sterile & Individually Wrapped



Bendable Blades

Use As Straight Blade Or Bend With Needle Holder To Desired Angle (Up To 45°)

Bendable Miniature Surgical Blades

-  Mini Blade #64
#BLAD-MINI-64 68.90 BOX OF 12
-  Mini Blade #67
#BLAD-MINI-67 68.90 BOX OF 12
-  Mini Blade #63 - *Full Radius Edge*
#BLAD-MINI-63 68.90 BOX OF 6
-  Mini Blade #69 - *Full Radius Edge*
#BLAD-MINI-69 68.90 BOX OF 6

MOST POPULAR

Bendable Micro-Miniature Surgical Blades

-  Micro Blade #64
#BLAD-MICRO-64 53.30 BOX OF 6
-  Nordland Micro Blade #6900 - *Full Radius Edge*
#BLAD-NORDLAND-69 101.90 BOX OF 6

MOST POPULAR

MOST POPULAR

Blades Shown Actual Size



Blade Handle For Miniature & Micro Surgical Blades

#BLAD-HANDLE 76.30 EACH



Miniature & Micro Surgical Blade Handle

#TVA-BLAD-HANDLE 141.90 EACH

Standard Scalpel Blades



Carbon Steel Blades For Regular Scalpel Handles

Extremely Sharp, Sterile & Individually Wrapped

	1 Box	2 Or More Boxes
#BLAD 12	BOX OF 100 67.30	BOX 63.10
#BLAD 12B	BOX OF 100 67.30	BOX 63.10
#BLAD 15	BOX OF 100 67.30	BOX 63.10
#BLAD 15C	BOX OF 100 96.50	BOX 91.30

Blades May Be Assorted For "2 Or More Box" Pricing

Tunneling & Vestibular Approach Instrument Set



Perio Probe

#TVA-PERIO-PROBE
..... 53.40 EACH



Perio Elevator #1

#TVA-PERIO-ELEV-1
..... 110.60 EACH



Perio Elevator #2

#TVA-PERIO-ELEV-2
..... 110.60 EACH




Perio Elevator #3

#TVA-PERIO-ELEV-3
..... 110.60 EACH



Tunnel Knife #1

#TVA-TUNNEL-KNIFE-1
..... 110.60 EACH



Tunnel Knife #2

#TVA-TUNNEL-KNIFE-2
..... 110.60 EACH



Gingivectomy Knife

#TVA-GINGI-KNIFE
..... 110.60 EACH



Flap Prep Instrument #1

#TVA-FLAP-PREP-1
..... 110.60 EACH



Flap Prep Instrument #2

#TVA-FLAP-PREP-2
..... 110.60 EACH



Perio Chisel

#TVA-PERIO-CHISEL
..... 110.60 EACH

Tunneling & Vestibular Approach Instrument Set



Tissue Plugger

#TVA-TISSUE-PLUGGER
..... 110.60 EACH



Tunnel Spoon #1

#TVA-TUNNEL-SPOON-1
..... 110.60 EACH



Tunnel Spoon #2

#TVA-TUNNEL-SPOON-2
..... 110.60 EACH



Tunnel Spoon #3

#TVA-TUNNEL-SPOON-3
..... 110.60 EACH




Micro Blade Handle

#TVA-BLAD-HANDLE
..... 141.90 EACH

SPECIAL

Tunneling & Vestibular Approach Instrument Set

MOST POPULAR



- Perio Probe
 - Perio Elevator #1
 - Perio Elevator #2
 - Perio Elevator #3
 - Tunnel Knife #1
 - Tunnel Knife #2
 - Gingivectomy Knife
 - Flap Prep Instrument #1
 - Flap Prep Instrument #2
 - Perio Chisel
 - Tissue Plugger
 - Tunnel Spoon #1
 - Tunnel Spoon #2
 - Tunnel Spoon #3
 - Micro Blade Handle
 - **FREE** Stainless Steel Cassette
- #TVA-TUNNEL-KIT **1710.00 SET**

Tunneling & Soft Tissue Instruments



Tunneling Elevator #1

#TUNNEL-ELEV-1
.....94.90 EACH



MOST POPULAR

Tunneling Elevator #2

#TUNNEL-ELEV-2
.....94.90 EACH



Tunneling Elevator #3

#TUNNEL-ELEV-3
.....94.90 EACH



PD Miller
Tunneling Elevator

#PDMILLER-ELEV-1
.....94.90 EACH




Micro Forceps

#FORCEP-MICRO
.....118.90 EACH



Tissue Augmenter

#TISSUE-AUGMENT
.....118.90 EACH



MOST POPULAR

Corn Pliers - Angled

Suture Forcep With
"Pass-Through" Hole
15.4cm
#CORN-PLIER
.....94.80 EACH



MOST POPULAR

Suture Pliers

"Delicate Procedure"
15.8cm Suture Forcep
1.6mm Ring Tip
#SUTURE-PLIER
.....154.20 EACH

Orban & Hirschfeld For Tunneling



MOST POPULAR

Orban Knife

#ORBAN
.....85.40 EACH
3 OR MORE
.....80.90 EACH



MOST POPULAR

Orban 1/2 Short

#ORBAN-1/2-SHORT
.....85.40 EACH
3 OR MORE
.....80.90 EACH



Orban 1/2 Long

#ORBAN-1/2-LONG
.....85.40 EACH
3 OR MORE
.....80.90 EACH



Hirschfeld 16/17 Tunneling File

Rounded Tip Prevents Puncturing Tissue
#HIRSCHFELD-16/17-FIL . . . 88.90 EACH
3 OR MORE 83.30 EACH

SPECIAL

Salvin Tunneling Basic Instrument Set #2

MOST POPULAR



Includes:

- Tunneling Forcep
 - Micro Jacobson Curved N/H
 - Tunneling Elevator #1
 - Tunneling Elevator #2
 - Tunneling Elevator #3
 - Orban 1/2 Short
 - Orban 1/2 Long
 - Hirschfeld 16/17 Tunneling File
 - PD Miller Tunneling Elevator
 - Younger Good 7/8
 - Suture Pliers
 - **FREE** Stainless Steel Cassette
- #TUNNEL-KIT-2. **1725.00 SET**

Micro Surgery Instruments



MOST POPULAR

Rhodes Chisel

17.2cm
Back Action
#RHODES
..... 104.80 EACH



MOST POPULAR

Micro Elevator #1

17.2cm
#MICRO-ELEV-1
..... 104.80 EACH



Micro Elevator #2

17.2cm
#MICRO-ELEV-2
..... 104.80 EACH



Micro Elevator #3

17.2cm
#MICRO-ELEV-3
..... 104.80 EACH



Micro Elevator #4

17.2cm
#MICRO-ELEV-4
..... 104.80 EACH



MOST POPULAR

Micro Chisel #1

17.2cm
#MICRO-CHISEL-1
..... 104.80 EACH



Micro Chisel #2

17.2cm
#MICRO-CHISEL-2
..... 104.80 EACH

SPECIAL

Salvin Micro Surgery Set

MOST POPULAR



- Micro Castro Needle Holder Curved
- Micro Adson Tissue Forcep
- Ross Scissors
- Rhodes Chisel
- Micro Chisel #1
- Micro Elevator #1
- Micro Elevator #3
- **FREE** Stainless Steel Cassette



#MICRO-SURG-KIT-1
1182.00 SET

Pikos / Ordonez Soft Tissue Manipulation Instrument Set

Ideal For Tunneling & Pocket Flap Procedures

Collaboratively Designed By Dr. Mike Pikos & Dr. Alvaro Ordonez



Pikos/Ordonez #1
Tunneling Instrument
#SOFTIMAK-1
..... 97.90 EACH



Pikos/Ordonez #2
Used For Pocket Flaps
#SOFTIMAK-2
..... 97.90 EACH



Pikos/Ordonez #3
Used For Pocket Flaps
#SOFTIMAK-3
..... 97.90 EACH



Pikos/Ordonez #5
Used For Pocket Flaps
#SOFTIMAK-5
..... 97.90 EACH



Pikos/Ordonez #6
Used For Pocket Flaps
#SOFTIMAK-6
..... 97.90 EACH



Pikos/Ordonez #8
Used For Reverse Flap
#SOFTIMAK-8
..... 97.90 EACH



Pikos/Ordonez #7
Thin Modified Orban
#SOFTIMAK-7
..... 97.90 EACH



Pikos/Ordonez #4
Thin Microperiosteal
#SOFTIMAK-4
..... 97.90 EACH



Pikos/Ordonez Tissue Forcep
TC Coated
#SOFTIMAK-TF-1
..... 124.20 EACH

SPECIAL **Pikos/Ordonez Soft Tissue Manipulation Set**

MOST POPULAR




Set Of 9 Instruments With **FREE** Stainless Steel Cassette
#SOFTIMAK-KIT
842.00 SET


Tunneling Forceps

For Pulling Tissue Graft Through Tunnel

- Micro Thin For Uninterrupted Passage Through Tunnel
- 75° Angle For All Areas & Applications
- Pressure Limiting Spring Prevents Damage To Graft
- Diamond Coated For Slip Resistance



Guide Through Tunnel To Grasp Tissue Graft



Close Forceps & Pull Graft Through Tunnel



MOST POPULAR

Tunnel Forcep Straight

#TUNNEL-FORCEP 400.00 EACH



Tunnel Forcep Curved Left

#TUNNEL-FORCEP-CU-L 400.00 EACH



Tunnel Forcep Curved Right

#TUNNEL-FORCEP-CU-R 400.00 EACH

SPECIAL Set Of 3 Tunnel Forceps
 #TUNNEL-FORCEP-SET-3 927.00 SET OF 3

Micro Needle Holders



Smooth Tips
 TUNGSTEN CARBIDE
 FOR 6-0 TO 10-0 SUTURES

SALVIN Comfort Handle

Micro Castro Needle Holder

- 13.3cm - Straight Tips
 #MICRO-CASTRO-STR 529.00 EACH
- 13.3cm - Curved Tips
 #MICRO-CASTRO-CUR 529.00 EACH



Serrated Tips
 TUNGSTEN CARBIDE
 FOR 5-0 TO 7-0 SUTURES

SALVIN Comfort Handle

Micro Needle Holder

- 13.3cm - Straight Tips
 #MICRO-N-HOLDER-STR 529.00 EACH
- 13.3cm - Curved Tip
 #MICRO-N-HOLDER 529.00 EACH



MOST POPULAR
Smooth Tips
 TUNGSTEN CARBIDE
 FOR 6-0 TO 10-0 SUTURES

SALVIN Comfort Handle

Micro Jacobson Needle Holder

- 17.8cm - Straight Tips
 #MICRO-JACOBSON-STR 529.00 EACH
- 17.8cm - Curved Tips
 #MICRO-JACOBSON-CUR 529.00 EACH

Micro Surgery Instruments



**Ross Scissors
Curved**

15cm
#ROSS.....379.40 EACH
3 OR MORE360.50 EACH



MOST POPULAR

**Castroviejo Micro Scissors
Curved**

15cm
#MICRO-CASTRO-SCISS
.....509.00 EACH



**Micro Goldman Fox Scissors
Curved**

11.5cm
#MICRO-GF-SCISS...146.70 EACH
3 OR MORE140.30 EACH



*Extra Thin
Micro Tip*

**Extra Thin Micro
Iris Scissors - Curved**

11cm
#IRIS-THIN-CURVED...123.50 EACH
3 OR MORE115.30 EACH



**Micro Adson 1x2
Tissue Forcep**

1x2 Teeth - 0.5mm Tips
12.1cm
#MICRO-ADSON-1X2
.....93.90 EACH
3 OR MORE89.20 EACH



MOST POPULAR

**Micro Adson Serrated
Tissue Forcep**

0.5mm Serrated Jaws
12.1cm
#MICRO-ADSON
.....93.90 EACH
3 OR MORE89.20 EACH



**Bishop Harmon
Micro Tissue Forcep**

1x2 Teeth - 0.3mm Tips
8.6cm
#BISHOP-HARMON
.....146.70 EACH
3 OR MORE140.30 EACH

Periodontal Knives & Surgical Files



MOST POPULAR

Orban Knife

#ORBAN
85.40 EACH
 3 OR MORE
80.90 EACH



MOST POPULAR

Orban 1/2 Short

17.2cm
 #ORBAN-1/2-SHORT
85.40 EACH
 3 OR MORE
80.90 EACH



Orban 1/2 Long

Extended Length For
 Tunneling Procedures
 17.2cm
 #ORBAN-1/2-LONG
85.40 EACH
 3 OR MORE
80.90 EACH



Kirkland Knife

#KIRKLAND
85.40 EACH
 3 OR MORE
80.90 EACH



**1S/2S Sugarman
 Periodontal File**

#SUGARMAN-1S/2S
88.90 EACH



MOST POPULAR

**3S/4S Sugarman
 Periodontal File**

#SUGARMAN-3S/4S
88.90 EACH



**Hirschfeld 16/17
 Tunneling File**

Rounded Tip
 Prevents Puncturing Tissue
 During Tunneling Procedures
 #HIRSCHFELD-16/17-FIL
88.90 EACH

Perio Curettes



MOST POPULAR

Prichard 1/2 Universal Perio Curette


#PRICHARD-1/2 47.40 EACH
3 OR MORE 43.90 EACH



Kramer #3 Universal Perio Curette

#KRAMER-3 47.40 EACH
3 OR MORE 43.90 EACH

Bone Files



Bone File #10

#FIL10 86.70 EACH



Bone File #12

#FIL12 86.70 EACH



MOST POPULAR

Bone File #12CA

#FIL-12CA 86.70 EACH

Tissue Forceps



MOST POPULAR

Adson Long 1x2 Tissue Forceps

1x2 Teeth
15cm

#ADSON-LONG-1X2
..... 76.40 EACH
3 OR MORE
..... 71.20 EACH




MOST POPULAR

Adson Long Serrated Tissue Forceps

Serrated Teeth
15cm

#ADSON-LONG-SERR
..... 76.40 EACH
3 OR MORE
..... 71.20 EACH



Adson Curved Serrated Tissue Forceps

Serrated Teeth
12cm

#ADSON-SERR-CURV
..... 76.40 EACH
3 OR MORE
..... 71.20 EACH



Debakey Atraumatic Tissue Forcep - Straight

Serrated Teeth
16cm

#DEBAKEY-STRAIGHT
..... 134.00 EACH
3 OR MORE
..... 128.10 EACH



Debakey Atraumatic Tissue Forcep - Curved

Serrated Teeth
16cm

#DEBAKEY-CURVED
..... 134.00 EACH
3 OR MORE
..... 128.10 EACH



Chalazion Mucocele Forcep

12mm x 16mm Ring
8.5cm Length

#CHALAZION
..... 183.10 EACH
3 OR MORE
..... 174.50 EACH

Tissue Forceps



MOST POPULAR

**Adson 1x2
Tissue Forceps**

1x2 Teeth
12cm

#ADSON 1X2
.....71.90 EACH
3 OR MORE
.....67.70 EACH



MOST POPULAR

**Adson 1x2 Curved
Tissue Forceps**

1x2 Teeth
12cm

#ADSON 1X2 CUR
.....76.50 EACH
3 OR MORE
.....71.30 EACH



MOST POPULAR

**Micro Adson
Serrated
Tissue Forcep**

0.5mm Serrated
12cm

#MICRO-ADSON
.....93.90 EACH
3 OR MORE
.....89.30 EACH



**Adson 7x7
Tissue Forceps**

7x7 Teeth
12cm

#ADSON-7X7
.....94.50 EACH
3 OR MORE
.....89.90 EACH



**Adson 7x7 Curved
Tissue Forceps**

7x7 Teeth
12cm

#ADSON-7X7-CURVED
.....94.50 EACH
3 OR MORE
.....89.90 EACH



**Dressing/Tissue
Forcep**

Serrated Tip
15.2cm

#DRESS-FORCEP-6
.....61.30 EACH
3 OR MORE
.....56.60 EACH



Cotton Plier

Serrated Tip
15cm

#COTPLI
.....40.40 EACH
3 OR MORE
.....37.90 EACH

Tissue Forceps



Gerald 1x2

Delicate Tissue Forceps
1x2 Teeth
18cm

#GERALD 1X2
.....92.40 EACH
3 OR MORE
.....87.40 EACH



Gerald Serrated

Delicate Tissue Forceps
Serrated Teeth
18cm

#GERALD
.....92.40 EACH
3 OR MORE
.....87.40 EACH



MOST POPULAR

Gerald Serrated Straight

Delicate Tissue Forceps
Serrated Teeth
18cm

#GERALD-STRAIGHT
.....92.40 EACH
3 OR MORE
.....87.40 EACH




MOST POPULAR

Gerald 1x2 Teeth Straight

Delicate Tissue Forceps
1x2 Teeth
18cm

#GERALD-1X2-STRAIGHT
.....92.40 EACH
3 OR MORE
.....87.40 EACH



MOST POPULAR

Russian Long

Tissue Forceps
Serrated Tip
19.7cm

#RUSS-LONG
.....90.50 EACH
3 OR MORE
.....85.40 EACH



Russian

Tissue Forceps
Serrated Tip
15.2cm

#RUSS
.....81.50 EACH
3 OR MORE
.....77.10 EACH

Tissue Forceps



Adson Titanium Curved Tissue Forceps

Serrated Teeth
12cm
#ADSON-TI-CURVED
..... 254.90 EACH
3 OR MORE
..... 242.50 EACH



MOST POPULAR

Adson Carbide Tipped Tissue Forcep

Serrated Carbide Tips For Atraumatic Grasp Of Tissue
12cm
#ADSON TC
..... 178.00 EACH
3 OR MORE
..... 171.30 EACH



TC Gripper Carbide Tipped Tissue Forcep

Serrated Carbide Tips For Atraumatic Grasp Of Tissue
14.5cm
#TC GRIPPER
..... 189.00 EACH
3 OR MORE
..... 179.80 EACH



TC Gripper Long Carbide Tipped Tissue Forcep

Serrated Carbide Tips For Atraumatic Grasp Of Tissue
18cm
#TC-GRIPPER-LONG
..... 189.00 EACH
3 OR MORE
..... 179.80 EACH


Specialty Suturing Forceps



MOST POPULAR

Corn Pliers - Angled

Suture Forcep With "Pass-Through" Hole
15.4cm
#CORN-PLIER
..... 94.80 EACH
3 OR MORE
..... 89.90 EACH



NEW!

V Forceps

Serrated Teeth Designed By Dr. Richard Martin
#V-FORCEP-BLK
..... 92.50 EACH
3 OR MORE
..... 87.70 EACH



MOST POPULAR

Suture Pliers

Delicate Suture Forcep With 1.6mm Ring Tip
15.8cm
#SUTURE-PLIER
..... 154.20 EACH
3 OR MORE
..... 146.80 EACH

Super Sharp Disposable Tissue Punches

Sterile Disposable Tissue Punches

- Razor Sharp Cutting Edge
- Individually Packaged Sterile - Box Of 50
- Can Be Used In Single Tooth Sites

MOST POPULAR


MOST POPULAR

● 3.5mm

● 4.0mm

● 5.0mm

● 6.0mm



3.5MM PUNCHES-BOX OF 50
 4.0MM PUNCHES-BOX OF 50
 5.0MM PUNCHES-BOX OF 50
 6.0MM PUNCHES-BOX OF 50

	1 Box	2 Boxes
#PUNCH-3.5	177.80 BOX	165.90 BOX
#PUNCH-4.0	177.80 BOX	165.90 BOX
#PUNCH-5.0	177.80 BOX	165.90 BOX
#PUNCH-6.0	177.80 BOX	165.90 BOX

Short Sterile Disposable Tissue Punches

- Razor Sharp Cutting Edge
- Individually Packaged Sterile - Box Of 50
- Designed For Edentulous Areas - Not For Single Tooth Sites



3.5MM SHORT PUNCHES-BOX OF 50
 4.0MM SHORT PUNCHES-BOX OF 50
 5.0MM SHORT PUNCHES-BOX OF 50
 6.0MM SHORT PUNCHES-BOX OF 50

	1 Box	2 Boxes
#PUNCH-SHORT-3.5	177.80 BOX	165.90 BOX
#PUNCH-SHORT-4.0	177.80 BOX	165.90 BOX
#PUNCH-SHORT-5.0	177.80 BOX	165.90 BOX
#PUNCH-SHORT-6.0	177.80 BOX	165.90 BOX

Disposable Punches Are One Size Per Box
 Boxes May Be Combined For "2 or More Box" Pricing

Stainless Steel Multi-Use Rotary Tissue Punches



3.0MM ROTARY PUNCH
 4.0MM ROTARY PUNCH
 5.0MM ROTARY PUNCH
 6.0MM ROTARY PUNCH
 7.0MM ROTARY PUNCH

#PUNCH-3-ROTARY	129.00 EACH
#PUNCH-4-ROTARY	129.00 EACH
#PUNCH-5-ROTARY	129.00 EACH
#PUNCH-6-ROTARY	129.00 EACH
#PUNCH-7-ROTARY	129.00 EACH



Use With Latch Type Implant Handpiece


SPECIAL

Rotary Tissue Punch Kit MOST POPULAR

Set Of 5 Punches With Autoclavable Block
 #PUNCH-ROTARY-KIT
 569.00 SET OF 5 W/BLOCK



Stainless Steel Tissue Punches




Stainless Steel Tissue Punches - Offset

	Individually	3 or More
4.0MM OFFSET PUNCH	#PUNCH-SS-4-RA. 107.10 EACH	101.90 EACH
5.0MM OFFSET PUNCH	#PUNCH-SS-5-RA. 107.10 EACH	101.90 EACH
6.0MM OFFSET PUNCH	#PUNCH-SS-6-RA. 107.10 EACH	101.90 EACH
8.0MM OFFSET PUNCH	#PUNCH-SS-8-RA. 107.10 EACH	101.90 EACH

Stainless Steel Punches May Be Assorted For "3 or More" Pricing

Tissue Flap Punch



**Garvar Tissue Flap Punch
Implant Exposing Forcep**

Designed By Dr. Lanny Garvar

- Use To Contour Tissue Around Implant
- 5mm Diameter Ring

#FLAP-PUNCH 299.00 EACH

Tissue Punch Forcep



Tissue Punch Forcep

To Remove Tissue After Using Punch


#TISSUE-PUNCH-FORC. 80.90 EACH
3 OR MORE 76.50 EACH



Start With Your Choice Of Tissue Punch → Remove Tissue With Forcep

Scalpel Handles

Blade Locks At 90°, 60°, 30°, Straight, 330°, 300°, & 270° Angles



Tighten / Loosen To Lock / Unlock



MOST POPULAR

6 Way Adjustable Angled Scalpel Handle

#SCALP-6-WAY-ADJ
 255.00 EACH
 2 OR MORE
 241.90 EACH

MOST POPULAR

Siegel Round Scalpel Handle

For Precise Control When Rotating Angle Of Incision

#SIEGEL 93.30 EACH
 3 OR MORE 85.90 EACH



Siegel Round Angled Scalpel Handle

For Precise Control When Rotating Angle Of Incision

#SIEGEL-ANGLED..... 93.30 EACH
 3 OR MORE 85.90 EACH




MOST POPULAR

#3 Scalpel Handle

#SCALP..... 28.90 EACH
 3 OR MORE .. 26.30 EACH



Push Lever Forward To Remove Blade



Safety Ejector Scalpel Handle

Safe & Easy Blade Removal

#SCALP-SAFETY..... 222.90 EACH



Paquette Blade Holder

Palatal, Gingival & Connective Tissue Grafts
 Uses Standard Double-Edged Razor Blade
 15cm Length

#SCALP-PAQUETTE
 229.00 EACH



Double Bladed Scalpel Handle

1.5mm Blade Separation
 Uses #12 or #15 Blades

#SCALP-DUAL-1.5
 185.00 EACH



**SALVIN
ENLIGHTEN-PRO3™
LED SYSTEM**

- Continuously Variable Adjustable Intensity **MOST POPULAR**
- 100,000 Hour LED Life
- Dual Lithium Polymer Battery Pack
 - 6.5 Hour Capacity On High
 - 14 Hour Capacity On Medium
 - 50 Hour Capacity On Low
- Spot Size Adjusts To 3.5" At 18" Working Distance
- Total Weight: 11oz (311g)
- Optimized Weight Balance And Comfort
- 2 Year Warranty



Includes Batteries & Charger



Includes Anti-Curing Filter



Use With Your Loupes



Includes Carrying Case

SPECIAL

**Enlighten-Pro3™
LED Headlamp System**

Includes:

- Enlighten-PRO3™ LED Headlamp
- Secure Carrying Case
- Anti-Curing Filter
- 2 Lithium Batteries, Charger & Cable

#HL-PRO3-PKG..... 1999.00

- **Patented Technology**
- **Reduces LED Temperature**
- **No Cables Or Belt Packs**
- **100,000 LUX At 10" (25cm)**
- **15,000 LUX At 18" (46cm)**
- **Preferred Working Distance For Most Doctors**
*(Important To Measure At Your Working Length
As Luminance Decreases With Distance)*

Enlighten-Pro3™ Accessories



Extra Enlighten-Pro3™ Lithium Batteries - Set Of 2

#HL-PRO3-BAT (SET OF 2)..... 305.00



Optional Enlighten-Pro3™ Four-Bay Charging Station

- Charge Up To Four Batteries At Once
- Charge Status Indicators For Each Battery
- Batteries Not Included

#HL-PRO-DOC 555.00

Implant Model Package Specials

SPECIAL

Bar Overdenture Model + Crown & Bridge Model

MOST POPULAR



MOD D4



MOD-C7

Implant Model Special #1

Regular Price For Package
Is 876.00 - **Save 125.00**

#MOD-SPECIAL-D4&C7
751.00 FOR BOTH MODELS

SPECIAL

Locator™ Model + Crown & Bridge Model

MOST POPULAR



MOD LOC-2



MOD-C7

Implant Model Special #2

Regular Price For Package
Is 795.00 - **Save 120.00**

#MOD-SPECIAL-LOC2&C7
675.00 FOR BOTH MODELS

SPECIAL

Bar Overdenture + Locator™ + Crown & Bridge + "On 4" Models

MOST POPULAR



MOD D4



MOD LOC-2



MOD-C7



MOD-COMplete-4

Implant Model Special #3

Regular Price For Package
Is 1655.00 - **Save 244.00**

#MOD-SPECIAL-SETOF4
1411.00 FOR ALL 4 MODELS

SPECIAL

Set Of Both Mandible & Maxilla Bone Loss Models

MOST POPULAR



MOD MM



MOD MMM

Set Of Both Mandible & Maxilla Bone Loss Models

#MOD-SET
253.00 FOR BOTH MODELS

Crown & Bridge Implant Models



**3 Unit Bridge On 2 Implants
Single Crown On 1 Implant**
#MOD C2 324.00



MOST POPULAR
**3 Unit Bridge On 3 Implants
Single Crown On 1 Implant
3 Unit Natural Tooth Bridge**
#MOD-C7 451.00



**Single Posterior Implant Crown
Single Anterior Implant Crown
3 Unit Natural Tooth Bridge**
#MOD-C8 451.00



6 Implant High Water Bridge
#MOD HWB-6 351.00



6 Implant Full Arch Cemented
#MOD C4 526.00



MOST POPULAR
Complete "On 4" Restoration
#MOD-COMplete-4 435.00



Complete "On 6" Restoration
#MOD-COMplete-6 451.00



**Most of Our Models Can Be Customized
With Your Practice Or Study Club Name
Please Call For Quantity Pricing**

More Models On Next Pages →

Hader Bar Implant Models



MOST POPULAR

4 Implant Hader Clip Bar
Lower Overdenture
#MOD D4..... 425.00



6 Implant Hader Clip Bar
Upper Overdenture
#MOD-D6-UPPER..... 475.00

"O" Ring Abutment Implant Models



2 Implant "O" Ring
Lower Overdenture
#MOD D2-O..... 279.00



MOST POPULAR

4 Implant "O" Ring
Lower Overdenture
#MOD D4-O..... 309.00



*Most of Our Models Can Be Customized
With Your Practice Or Study Club Name
Please Call For Quantity Pricing*

Locator™ Abutment Implant Models



MOST POPULAR

**4 Implant Locator™ Abutment
Upper Overdenture**
#MOD-LOC-4-UPPER 434.00



**4 Implant Locator™ Abutment
Lower Overdenture**
#MOD-LOC-4 434.00



**3 Implant Locator™ Abutment
Lower Overdenture**
#MOD-LOC-3 405.00



MOST POPULAR

**2 Implant Locator™ Abutment
Lower Overdenture**
#MOD-LOC-2 344.00

Locator™ Is A Trademark of Zest Anchors, Inc.



**Most of Our Models Can Be Customized
With Your Practice Or Study Club Name
Please Call For Quantity Pricing**

Bone Loss Models



MOST POPULAR

Mandible Bone Loss Model Set

- 4 Mandibles Show Progressive Bone Loss
- Shows Benefits Of Augmentation & Consequences Of No Treatment

#MOD MM 137.00 EACH

Maxilla Bone Loss Model Set

- 4 Maxillas Show Progressive Bone Loss
- Shows Benefits Of Augmentation & Consequences Of No Treatment

#MOD MMM 165.00 EACH

Atrophied Ridge Mandible

Atrophied Ridge Mandible Shows Vertical & Horizontal Defects

#MAND-RIDGE-ATROPH 58.00 EACH

SPECIAL Set Of Both Mandible & Maxilla Bone Loss Models **MOST POPULAR**

#MOD-SET 253.00 FOR BOTH

Teaching Models For Hands-On Courses



MOST POPULAR

Suturable Implant Drilling Mandibles

Similar In Density To Bone With Suturable Rubber Surface

#MAND SUTURABLE 69.40 EACH
 12 OR MORE 66.70 EACH
 24 OR MORE 63.90 EACH



MOST POPULAR

Implant Drilling Mandibles

Similar In Density To Bone

#MAND TEACH 24.00 EACH
 12 OR MORE 22.20 EACH
 50 OR MORE 20.30 EACH



Drillable Sinus Lift Teaching Model

Similar In Density To Bone With Rubber Sinus Membrane

#SINUS TEACH 122.00 EACH
 12 OR MORE 113.30 EACH

Call For Study Club & Implant Course Volume Pricing

Penguin^{RFA} Implant Stability Monitor (Resonance Frequency Analysis)

- *Easy-To-Use, Reliable & Affordable Way To Measure Implant Stability*
- *Cordless & Lightweight*
- *Provides Accurate ISQ Measurement*



Install The MultiTipeg™ Into The Implant And The Measurement Is Made In Seconds.



Penguin^{RFA} II

- **Includes:** Stand-Up Wireless Charging Base, MultiTipeg™ Driver, And ISQ Test Peg
- Rubberized Grip For Easy Handling
- 2 Year Warranty

#PENGUIN-22490.00



MultiTipeg™ Multi-Use Sensing Devices

- Available For All Major Implant Systems
- Made Of Durable, Tissue Friendly Titanium
- Laser Marked With Type Number
- ISQ Standard Calibrated™
- Autoclavable Up To 20 Times

CALL TO CONFIRM WHICH PEGS YOU NEED



NEW!

Penguin Sterile Covers (20 Pack)

- Durable, Disposable Covers Designed For Use With The Penguin II RFA ISQ Monitors

#PENGUIN-SLEEVE-20 44.00 PACK OF 20

Spotter Smart Implant Detector



Accurately Locates Covered Dental Implants



- Easily Finds Buried Titanium Implants
- Saves Time And Minimizes Recovery Period



Spotter Smart Implant Detector

Includes One Sterilizable Detachable Sensor (Recommended For 10 Uses)

#SPOTTER 755.00



Spotter Sensor Replacements

Package of 5

#SPOTTER-SENSORS-5 140.00 PACKAGE OF 5

Surgical Rulers & Markers



MOST POPULAR

Surgical Ruler

Stainless Steel 6" (15.25cm) Surgical Ruler

#SURGRULER
 INDIVIDUALLY 15.90 EACH
 2 OR MORE 13.90 EACH



Sterile Marking Pens

Box of 10 Individually Packaged Markers

#TISSUE-MARK
 BOX OF 10 53.90 BOX
 2 OR MORE BOXES 50.30 BOX

X-Ray Distortion Markers




X-Ray Distortion Markers

- Precise 5mm Stainless Balls Calibrate For Distortion On Periapical & Panoramic Films
- Incorporate Markers Into A Plastic Stent When Taking Radiograph, Then Use 5mm Diameter As A Reference When Taking Measurements On Your Film

#RADIO MARK 33.50 BAG OF 50

Salvin Implant Drill Guide Set


- *Helps Provide Proper Spacing, Angulation & Parallelism*
- *Autoclavable*



Regular Drill Guide Set Includes:

- Set Of 3 Color Coded Drill Guides & Autoclavable Cassette

#IMPL-GUIDE-SETOF3 207.00 SET OF 3 WITH CASSETTE
US Patent #6869283B2



“Regular Implant” Guide Protocol

1. Select Color Coded Guide Based On Implant Size. Position Guide Against Adjacent Tooth At A Right Angle To The Jaw

2. Use Perio Probe To Locate Crest Of Ridge, Then Drill Pilot Hole With 2.0mm Drill Parallel To The Guide Pin

3. Place Guide Pin Into Osteotomy & Align Guide Over Ridge

4. Drill Second Pilot Hole Through Metal Ring, Use Paralleling Pin To Confirm Parallel Alignment With Guide Pin

Salvin Paralleling Guide



MOST POPULAR

Easy Alignment Of Implants In Edentulous Mandibles

- 1.9mm Diameter Guide Pin
- 7mm Length Guide Pin

#PARALLEL-GUIDE 134.00 EACH

Salvin Paralleling Guide Pins

Paralleling Guide Pins
1.9mm Diameter For 2mm Pilot Drill
#GUIDE-PIN-PKG-4. . . . 103.90 SET OF 4 PINS
#GUIDE-PIN-PKG-8. . . . 199.00 SET OF 8 PINS
Includes **FREE** Organizing Block

SPECIAL

Salvin Guide Pin
Set Of 8 Guide Pins With **FREE** Sterilization / Organizing Block
#GUIDE-PIN-PKG-8
199.00 KIT W/BLOCK

Calipers & Gauges




Flexes For Increased Tactile Feel While Probing

SALVIN Ergonomic Handle

Flexible Tip Implant Depth Gauge

Marked At 6, 8, 10, 12, 14 & 16mm

#DEPTH-FLEX 111.90 EACH



6mm Wide Flat Handle For Quick Reference Measurements


Ball Tip End Marked At 7-20mm

MOST POPULAR

"Ball Tip" Implant Depth Gauge - Titanium

- Ball Tip Probe Marked At 7mm To 20mm
- Flat End Of Handle Is 6mm Wide For Reference Measurements


#DEPTH-BALL-TIP 99.90 EACH



Electronic Digital Caliper

- 0-150mm Range
- For Extra-Oral Use Only

#CALIPER-DIGITAL 99.90 EACH



Short Tips For Access To Posterior

Thumb Control Locks Tips In Place




MOST POPULAR

Sliding Caliper

- Facilitates Positioning of Implants
- Measures Both Inside & Outside Dimensions
- Sliding Depth Measurement Probe
- 0-90mm Range

#CALIPER-SLIDING 299.00 EACH



Boley Gauge Caliper

0-100mm Range

#BOLEY 127.80 EACH

Calipers



MOST POPULAR

Ridge Mapping Caliper

- For Mapping Bone Width
- Measures 0-25mm Range

#RIDGE-CALIPER379.00 EACH



MOST POPULAR

Castroviejo Caliper Short

- Measures 0-40mm Range
- 9cm Long

#CASTRO-CAL341.50 EACH



Castroviejo Caliper Short - 30° Angled

- Measures 0-40mm Range
- 9cm Long

#CASTRO-CAL-ANGLED-S 341.50 EACH



Socket Wall Caliper

- Measures Flap & Socket Wall Thickness
- Measures 0-10mm Range

#CALIPER-SOCKET56.90 EACH



Castroviejo Caliper Long

- Measures 0-40mm Range
- 17cm Long

#CASTRO-CAL-LONG439.00 EACH



Castroviejo Caliper Long - 30° Angled

- Measures 0-40mm Range
- 17cm Long

#CASTRO-CAL-ANGLED439.00 EACH

Salvin Bendable "Magic" Mirror Handle For Glass Photo Mirrors



MOST POPULAR



Handle Securely Holds Glass Photo Mirrors

Easily Bends & Stays In Place For Your Desired Angle

Bendable Silicone Handle For Glass Photo Mirrors

Securely Holds Photography Mirror

#PHOTO-MIRR-MAGIC-HND	77.90 EACH
2 OR MORE	73.50 EACH

Dual View Photography Mirror



Dual View Photo Mirror

View Buccal, Occlusal & Lingual Surfaces With One Mirror

#PHOTO-DUALVIEW	105.90 EACH
-----------------------	-------------

Photo Contrasters

Set Of 2 Contrasters Provide Better Visualization Of Shade & Detail

MOST POPULAR

For Effective Communication With Patient & Laboratory
Black Background For Better Visualization Of Color, Detail,
Shade Variations & Translucencies



Without Contraster



With Contraster

Photos: Dr. Ross W. Nash




Large & Small Contraster Set

#CONTRASTER-STR-2	66.50 SET
2 OR MORE SETS	62.80 SET

Operative Mirrors On Page 198

Front Surface Float Glass For Reflected Image With No Distortion

- *Mirror Surface On Both Sides*
- *Sterilization Up To 428°F (220°C)*



MOST POPULAR

MOST POPULAR

Pedo Occlusal

60mm Wide
85mm Long

#PHOTO OCC PED
INDIVIDUALLY
..... 79.90 EACH
2 OR MORE
..... 75.50 EACH

Adult Occlusal

68mm Wide
100mm Long

#PHOTO OCC AD
INDIVIDUALLY
..... 79.90 EACH
2 OR MORE
..... 75.50 EACH

Adult Occlusal XL

68mm Wide
135mm Long

#PHOTO OCC XL
INDIVIDUALLY
..... 101.80 EACH
2 OR MORE
..... 94.90 EACH

Buccal

57mm Wide
152mm Long

#PHOTO BUCCAL
INDIVIDUALLY
..... 79.90 EACH
2 OR MORE
..... 75.50 EACH


Buccal Narrow

41mm Wide
152mm Long

#PHOTO-BUCCAL-2
INDIVIDUALLY
..... 79.90 EACH
2 OR MORE
..... 75.50 EACH

Autoclavable Photo Mirror Cassette

- Silicone Cushions & Locking Crossbar
 - Holds Set Of 4 Mirrors:
Pedo Occlusal, Adult Occlusal, Buccal Narrow, Buccal
 - 29 x 18.3 x 3.7cm (11.5" x 7.2" x 1.5")
- #CASS-SS-MIR-4..... 169.70 EACH



SPECIAL

Mirror, Cassette & "Magic Handle" Set

MOST POPULAR



Includes:

- Pedo Occlusal Mirror
- Adult Occlusal Mirror
- Buccal Mirror
- Buccal Narrow Mirror
- Bendable Silicone "Magic" Mirror Handle
- Stainless Steel Cassette



#PHOTO-KIT-4-MAGIC
479.00 COMPLETE SET

Diagnostic Probes



**Titanium
UNC Probe**

#PROBE-UNC-TITANIUM
..... 53.90 EACH
3 OR MORE
..... 52.30 EACH




**Titanium
CP12 Probe**

#PROBE-CP12-TITANIUM
..... 53.90 EACH
3 OR MORE
..... 52.30 EACH



**Stainless Steel
UNC Probe**

#PROBE-UNC
..... 37.70 EACH
3 OR MORE
..... 34.30 EACH



**PerioWise Plastic
Perio Probes**

- Set Of 6 Perio Probes
- Autoclavable Plastic
Will Not Scratch Implants

#PERIO-PROBE-6
..... 112.90 SET OF 6
2 OR MORE SETS
..... 106.80 SET OF 6

Operative Mouth Mirrors & Handles



MOST POPULAR
**Mirror Handle
Cone Socket
Lightweight**

#MIRR-HNDL
..... 21.10 EACH
3 OR MORE
..... 18.80 EACH



**SALVIN
Ergonomic
Handle**

**Mirror Handle
Cone Socket
Lightweight
Anatomical**

#MIRR-HNDL-A
..... 43.90 EACH
3 OR MORE
..... 40.60 EACH



**#4 Cone Socket
Front Surface Mirror**

Rhodium-Coated To
Reduce Ghost Imaging

#MIRR-4-CS
..... 54.20 BOX OF 12
5 OR MORE BOXES
..... 50.60 BOX OF 12



**#5 Cone Socket
Front Surface Mirror**

Rhodium-Coated To
Reduce Ghost Imaging

#MIRR-5-CS
..... 54.20 BOX OF 12
5 OR MORE BOXES
..... 50.60 BOX OF 12

Implant Buddy

NEW!



**Implant Buddy
Square-Type Dental Implant Driver Set**

- Service 99% Of Implant Cases
- Prosthetic Connections Included:
 - Straight Screwdriver
 - SCS Screwdriver
 - UniGrip Driver
 - ASC Angled Access Driver
 - Hex 1.3mm/.050 in.
 - Hex 1.2mm/.048 in.
 - Hex 0.9mm/.035 in.
- Short & Long Drivers For Each Prosthetic Screw Connection
Short Drivers: Total Length 23mm; Shaft Length 16mm
Long Drivers: Total Length 39mm; Shaft Length 32mm
- Tapered And Color Coded Drivers
- Smart Thumb Knob
- Autoclavable
- Magnetic Base

#IMPLANT-BUDDY-SET-SQUARE 575.00 EACH



**Implant Buddy LT
Latch-Type Dental Implant Driver Set**

- Prosthetic Connections Included:
 - Straight Screwdriver
 - Multi-Unit Driver
 - SCS Screwdriver
 - ASC Angled Access Driver
 - UniGrip Driver
 - Hex 1.3mm/.050 in.
 - Hex 1.2mm/.048 in.
 - Hex 0.9mm/.035 in.
- Seven Latch Drivers:
Total Length 24mm; Shaft Length 10.5mm
Multi-Unit Driver: 20mm
- Tapered And Color Coded Drivers
- Autoclavable
- Magnetic Base

#IMPLANT-BUDDY-SET-LATCH. 389.00 EACH

DRIVER / PROSTHETIC CONNECTION		SYSTEM
—	Straight Diver	Straight Systems
★	SCS Diver	Straumann®
★	UniGrip Driver	Nobel BioCare®
⬡	Hex 1.3mm/.050 in.	Lifecore®, BioHorizons®, Astra®, Zimmer®
⬡	Hex 1.2mm/.048 in.	Biomet 3i®, MegaGen®, Hiossen®
⬡	Hex 0.9mm/.035 in.	Friadent®, Ankylos®, Biomet 3i® Cover Screw
▲	Angled Driver	Dynamic Systems
■	Square 1.3mm/.050 in.	Gold Screw/Nobel BioCare®

SALVIN
EASYReach®
RIGHT ANGLE PROSTHETIC WRENCH

Scan For
60⁺ Second Product Video™



- Excellent Ergonomics With Great Access In Difficult To Reach Areas
- Secures Driver Tip For Patient Safety

- ▶ Right Angle
- ▶ Latch Type

SALVIN
Ergonomic Handle



Finger Tip Rotation Easily Tightens Or Loosens Screws

Fits "Latch" Driver Tips



MOST POPULAR

EasyReach® Right Angle Prosthetic Wrench

Autoclavable

#EASY-REACH 573.00 EACH

SALVIN®
AUTO TORQ®
Automatic Electric Torque Driver

- Adjustable Torque From 10Ncm To 35Ncm
- Adjustable Speed From 30rpm To 50rpm
- Includes Reverse Function
- Excellent Ergonomics With Great Access In Difficult To Reach Areas
- Secures Driver Tip For Patient Safety

- ▶ Torque Wrench
- ▶ Adjustable Torque
- ▶ Right Angle
- ▶ Latch Type

Now Torques Up To 35Ncm



One Button Operation To Achieve Desired Torque

Fits "Latch" Driver Tips



AutoTorq® Automatic Electric Torque Driver

#AUTOTORQ-35 1500.00 EACH

SalvinTorq®
Adjustable Torque Wrench

- ▶ Torque Wrench
- ▶ Adjustable Torque
- ▶ Square Connection



Nine Torque Settings:
15, 20, 25, 30, 35, 40,
45, 50 & 60 Ncm

MOST POPULAR

Adjustable Torque Wrench

#SALVIN-TORQ 609.00 EACH

Use With Driver Tips On Page 203



Lubrication Hole
At 15Ncm Setting

SPECIAL SalvinTorq® Package

SalvinTorq® Wrench, Thumb Knob & Cassette
With Your Choice Of Four Driver Tips

#SALVINTORQ-KIT-4 **879.00 KIT**

SalvinTorq® Wrench, Thumb Knob & Cassette
With Your Choice Of Six Driver Tips

#SALVINTORQ-KIT-6 **922.00 KIT**

Please Specify Implant System & We'll Make Sure That You Have The Correct Driver Tips



Salvin Torque Wrenches

- ▶ Torque Wrench
- ▶ Square Connection



-  20 Newton Centimeters
#TORQ 20.229.90 EACH
-  30 Newton Centimeters
#TORQ 30.229.90 EACH
-  35 Newton Centimeters
#TORQ 35.229.90 EACH

Use These Wrenches With Driver Tips On Page 203

SALVIN® ACCESS TORQ®-PLUS

Right Angle Variable Torque Driver

Perfect For All-on-4® Procedures

	Torque Wrench
	Adjustable Torque
	Right Angle
	Latch Type



MOST POPULAR

AccessTorq®-PLUS Right Angle Variable Torque Driver

Autoclavable

Lubricate With Pana-Spray (Page 92) After Each Use & Before Sterilization

#ACCESSTORQ-PLUS 999.00 EACH

Rotation Easily Torques Or Loosens Screws

Seven Torque Settings:
10, 15, 20, 25, 30, 32 & 35 Ncm

Fits "Latch" Driver Tips

All-on-4® is a trademark of NobelBiocare®


"Latch Type" Prosthetic Driver Tips

	.048" (1.22mm) Tip Hex Driver #HN-LATCH-DRIVER 58.40 EACH
	.05" (1.25mm) Tip Hex Driver #HT-LATCH-DRIVER 58.40 EACH
	.05" (1.25mm) Tip Hex Driver - Short #HT-LATCH-DRIVER-S 58.40 EACH
	.05" (1.25mm) Tip Square Driver #CER-LATCH-DRIVER 58.40 EACH
	PROC Tip "Star" Driver #PROC-LATCH-DRIVER 80.90 EACH
	PROC Tip "Star" Driver - Short #PROC-LATCH-DRIVER-S 80.90 EACH
	LOC Triangle Driver #TRIRA 58.40 EACH
	ITI Driver #ITI-LATCH-DRIVER 73.60 EACH
	ITI Driver - Short #ITI-LATCH-DRIVER-S 73.60 EACH

"Latch" Tip To Square Driver Adapter












Use Your Latch Type Driver Tip With These Salvin Torque Wrenches

#LATCH-ADAPT 46.90 EACH



Square Connection Abutment Driver Tips - Fit Salvin Square Connection Torque Wrenches

Specify Implant Systems When Ordering – We'll Make Sure You Get The Right Driver Tips!

 <p>.048" (1.22mm) Hex - 6mm #HNT 58.90 EACH</p> <p>.048" (1.22mm) Hex Long - 12mm #HNTL 58.90 EACH</p> <p>For NobelBiocare, 3i, Keystone, MegaGen, Hiossen, Dentium</p>	 <p>.05" (1.25mm) Hex - 6mm #HTT 58.90 EACH</p> <p>.05" (1.25mm) Hex Long - 12mm #HTLT 58.90 EACH</p> <p>For Zimmer, Astra, BioHorizons, Intra-Lock, MIS, SteriOss, Sybron, Implant Direct, PerioSeal</p>	 <p>PROC Driver - 6mm #PROC-S 76.90 EACH</p> <p>PROC Driver Long - 12mm #PROC-L 76.90 EACH</p> <p>"Star" Driver For NobelBiocare Active & Unigrip Screw, Neoss</p>
 <p>Screw Driver - 6mm #SN 58.90 EACH</p> <p>Screw Driver Long - 12mm #SNL 58.90 EACH</p> <p>Flat Blade Screwdriver</p>	 <p>.035" (0.9mm) Hex - 6mm #HC 58.90 EACH</p> <p>.035" (0.9mm) Hex Long - 12mm #HCL 58.90 EACH</p> <p>For NobelBiocare, 3i, Keystone, BioHorizons, Friadent, PerioSeal</p>	 <p>.05" (1.25mm) Square - 6mm #CER-S 58.90 EACH</p> <p>.05" (1.25mm) Square Long - 12mm #CER-L 57.40 EACH</p> <p>For 3i</p>
 <p>ITI Driver - 6mm #ITIS 76.90 EACH</p> <p>ITI Driver Long - 12mm #ITIL 76.90 EACH</p> <p>"ITI" Driver For Straumann SCS</p>	 <p>LOC Triangle Driver - 6mm #TRIS 57.40 EACH</p> <p>LOC Triangle Driver Long - 12mm #TRIL 57.40 EACH</p> <p>For Locator™ Abutment</p>	 <p>.07" (1.77mm) Hex - 6mm #HTSI 57.40 EACH</p> <p>For Paragon TSI Abutment</p>
 <p>.08" (2.03mm) Female Hex - 6mm #HFD 57.40 EACH</p> <p>For NobelBiocare, 3i</p>	 <p>TG Driver #TG 90.30 EACH</p> <p>For 3i TG Abutment, Straumann SCS</p>	

Abutment Driver Adapters & Handles

Use Torque Wrench For Final Seating

Torque Wrench Adapter For NobelBiocare




#TORQ-ADAP-NOBEL89.50 EACH
Allows Use Of All Salvin Driver Tips
On This Page With Your Wrench

Torque Wrench Adapter For Straumann



#TORQ-ADAP-ITI89.50 EACH
Allows Use Of All Salvin Driver Tips
On This Page With Your Wrench

Square Connection



Straight Driver Tip Handle

#SIH 169.30 EACH

Reversible Ratchet

#RATCHET 162.90 EACH



Thumb Knob

Convenient For Starting Screws
#TK 79.90 EACH

Implant & Abutment Names Are Trademarks Of Their Respective Companies

Abutment & Prosthetic Parts Holders



MOST POPULAR

Hebel-Melton Abutment Holder

Self Tensioning Handle Keeps Parts Secure During Pick Up & Transfer

#HEBEL-MELTON-TI 161.20 EACH
3 OR MORE 154.40 EACH



Titanium Parts Forcep

For Slip-Free Grasp Of Prosthetic Parts
13cm

#TI-PARTS-GRIPPER 192.20 EACH
3 OR MORE 183.00 EACH



TC Gripper Forcep - Carbide Tipped

For Slip-Free Grasp Of Prosthetic Parts
14.5cm

#TC GRIPPER 189.00 EACH
3 OR MORE 179.70 EACH

Premier® Autoclavable Graphite Implant Scalers

MOST POPULAR

Universal Implant Scaler - Columbia 4L/4R

- Safe For Use Around Implants
- Resharpen With Clean Sharpening Stone
- Autoclavable Up To 140°C

#IMPL-SCALER-UNIV 205.70 SET OF 5
2 OR MORE 195.80 SET OF 5



Implant-Prophy+™ Plastic Hygiene Instruments

Plastic Implant Scalers With Sharpening Stone

- Special Plastic Will Not Scratch Abutment Surface
- Easily Resharpened With Sharpening Stone
- Autoclavable Up To 190°C (375°F)

Set Includes:
1 Gracey 5/6
1 Gracey 11/12
1 Gracey 13/14
1 Columbia 13/14U
Sharpening Stone

#PROPHY-INSTRSET 140.80 SET OF 4



Titanium Tip Implant Curettes

Langer 1/2 Titanium Curette

Angled Like A Gracey But
With Two Cutting Edges

#LANGER-1/2 75.50 EACH

Titanium Tipped



Langer 3/4 Titanium Curette

Angled Like A Gracey But
With Two Cutting Edges

#LANGER-3/4 75.50 EACH

Titanium Tipped



Langer 5/6 Titanium Curette

Angled Like A Gracey But
With Two Cutting Edges

#LANGER-5/6 75.50 EACH

Titanium Tipped



SPECIAL Set Of 3 Langers #TI-LANG-SET-3 **179.00 SET OF 3**

H6/H7 Titanium Curette

#TI-H6/H7 75.50 EACH

Titanium Tipped



204 SD Titanium Curette

#TI-204SD 75.50 EACH

Titanium Tipped



SPECIAL H6/H7 & 204 SD #TI-H6-204-SET-2 **127.40 SET OF 2**

SPECIAL **5 Instrument Titanium Curette Set**

MOST POPULAR

Includes:

- Titanium Langer 1/2
- Titanium Langer 3/4
- Titanium Langer 5/6
- Titanium H6/H7
- Titanium 204 SD
- Stainless Steel Cassette

#TI-CURETTE-SET-5
429.00 SET OF 5 W/CASSETTE

Crown Grippers



Carbide Tipped Crown Gripper

Slip-Free Gripping Of Crowns, Temporaries & Prosthetic Parts

#CROWN-GRIP245.00 EACH
2 OR MORE235.60 EACH



MOST POPULAR

Rubber Tipped Crown Plier



Includes 10g Of Silicon Carbide Gripping Powder

#CROWN-PLIER-II 325.00 EACH
2 OR MORE 316.50 EACH

Replacement Silicon Carbide Gripping Powder

10g Container

#CROWN-PLIER-POWDER 45.90 EACH

Replacement Crown Plier Tips



#CROWN-PLIER-TIPS-II 33.70 SET OF 2
2 OR MORE SETS 31.30 SET OF 2



Carbide Tipped Locking Crown Gripper

Slip-Free Gripping Of Crowns, Temporaries & Prosthetic Parts

#CROWN-GRIP-LOCKING245.00 EACH
2 OR MORE235.60 EACH

Crown Splitters

- *Slit The Crown With A Diamond*
- *Easily Pry Crown Apart*
- *Reinforced Tip*



Crown Splitter – Straight

#CROWN-SPLIT-STR 124.40 EACH



MOST POPULAR

Crown Splitter – Offset



#CROWN-SPLIT 154.00 EACH

Improv® Temporary Implant Cement For Implants & Natural Teeth

- *90 Seconds Working Time*
- *Final Set In 3.5 Minutes*



Temporary Implant Cement

- 90 Seconds Working Time
- Compressive Strength: 57 MPa
- Water Sorption: 1.314mg/cm³
- Final Set In 3.5 Minutes
- Film Thickness: <20 Microns
- Shelf-Life: 1.5 Years
- Contains:
 - 2 - 6mL Dual-Barrel Syringes
 - 12 - Mixing Tips
 - 1 - Mixing Pad

#IMPROV-CEMENT. 119.00 BOX
3 OR MORE 109.90 BOX



Extra Mixing Tips

Package Of 30
#IMPL-CEM-TIPS 51.20 BOX

45° Surgical Handpiece	92	Cortical Bone Fibers, OsseoFiber®	8	Horizontal Centrifuge	69
Abutment Holders	204	Covers, Light Handle	41	Hygiene, Implant	204-205
AccessTork®-Plus Torque Wrench	202	Covers, Machine & Mayo Stand	63	Ice Packs	62
ACM Autobone Chip Maker	128	Covers, Shoe	61	Impla-Safe Peri-Implant Management Set	76
AED	49	Crest-Control, Meisinger	147	Implant Abutment Holders	204
Air-N-Go® Air Polishing	74	Crestal Sinus Lift Kit	137, 140	Implant Brush, Titanium	75
AlloSculpt®-3D Plus Bone Putty	12	Criticare Monitors	43-44	Implant Buddy Implant Driver Set	199
AlloSpark-GF™ Growth Factor	6	Crown Grippers	206	Implant Center Piezo + Motor	71
Alveoplasty Burs	87	Crown Removal Burs	87	Implant Detector, Spotter	192
AmnioExcite® Amnion/Chorion Membrane	10-11	Crown Removers	206	Implant Equipment Cart	46
Amplifi® Long Term Membrane	24	Curettes, Perio	177	Implant Cement	207
Anesthetic Syringes	55	Curettes, Titanium	205	Implant Drill Stops	82
Aspirators, Surgical	98-99	Cyanoacrylate	39	Implant Drills	82-83
Atraulux® Elevators	106	CytoSurg® PTFE Sutures	35	Implant Hygiene	204-205
Atraulux® Extraction Forceps	112	CytoSurg® PTFE Membrane	28	Implant Positioning Guides	193
Atraumatic Extraction	102-119	Defibrillators	49	Implant Models	186-189
Autoclave Tape	63	Degranulation Burs	87	Implant Motor Stands	47
Automatic Membrane Tack System	135	Depth Gauges	194	Implant Motors & Handpieces	78-79
AutoTome® Extraction Kit	91	Dermis Allograft	22-23	Implant Packs	56-59
AutoTome® Osteotome Kit	90	Diagnostic Instruments	191-198	Implant Removal Kit	93
AutoTorq® Electric Torque Driver	200	Distortion Markers, X-Ray	192	Implant Rescue Kit	93
Balshi Zygomatic Retractor	152	Drapes, Patient	56-59	Implant Scalers	204-205
Bendable Scalpel Blades	167	Drill Extender	83	Implant Stability, Penguin	191
Benex Extraction System	119	Drill Guides	193	Improv Implant Cement	207
BIO-PRF Horizontal Centrifuge	68	Drill Stops	82	Infection Control	56-66
Bleeding Control	39-41	Drills & Burs	82-88	Instrument Cassettes	64
Blood Collection	67-69	Driver Tips	199, 202-203	Instrument Cleaner	63
Blood Pressure Monitors	42-45	DynaGraft-D and DynaBlast Bone Putty	33	Instrument Wrap & Pouches	63-64
BloodStop® Hemostatic Gauze	40	EasyReach® Prosthetic Wrench	200	Irrigation Needles & Syringes	94
Bone Chip Maker	129	ECG Monitors	42-45	Irrigation Tubing	80, 96
Bone Chisels	101	ECG Patches & Wrist Clips	50	Irrigation, External	94-95
Bone Collector Burs	128-129	Electric Torque Driver	200	Irrigation Syringes	94
Bone Curettes, Lucas	120	Electrocautery	41	Isolation Gowns	60
Bone Expander Kits	144	Elevators	102-111	IV Stand	47
Bone Fibers, OsseoFiber®	8	Elevators, PDL	102	IV Catheters & Admin Sets	52-54
Bone Files	177	Elevators, Periosteal	154-156	Kern Tunneling Instruments	166
Bone Graft Packers	122-123	Emergency Equipment	49-51	Knives, Perio	176
Bone Graft Syringes	124-125	Enlighten® LED Headlamp	185	Laryngeal Tubes	50
Bone Mills & Bone Morselizer	131	External Irrigation	94-95	Laster 3rd Molar Retractor	152
Bone Putty, AlloSculpt®-3D Plus	12	Extraction Forceps	112-118	LED Headlight	185
Bone Putty, OsseoFiber®	8	Extraction, Atraumatic	102-119	LifeNet Health® Allograft Bone	18-21
Bone Saws	81	Featherlite Scissors	39, 159	LifeNet Health® Dermis Allograft	22-23
Bone Scrapers	126-127	Files, Bone	177	Light Handle Covers	41
Bone Fixation Screws	133	Files, Periodontal	176	Lindemann Burs	84
Bone Spreaders	144-147	Fixation Screws	133	Live Cells Bone Allograft	7
Bone Trap	128	Flap Holding Forcep	183	Locator Implant Models	189
Bone Wax	40	Forceps, Membrane	136	Long Term Membrane, Amplifi®	24
Bouffant Caps	61	Forceps, Extraction	112-118	Lucas Bone Curettes	120
Bovine Xenograft, SalvinOss®-B	13	Forceps, Root Tip	114, 118	Luxating Elevators	106-109
Broken Screw Remover Kit	93	Forceps, Suture	181	Machine Covers	41
Bur Blocks	89	Forceps, Tissue	178-181	Mallet, Mead	101
Burs & Drills	82-88	Fragment Forceps	118	Mandibles, Drilling	190
Calipers	194-195	Furniture, O/R	46-48	Marking Pens, Sterile	192
Cannulas, Nasal	51	Gloves, Sterile	66	Masks, N95	61
Capnography Monitors & Accessories	42-45	Gowns, Sterile & Non-Sterile	60	Masks, Surgical	61
Cassettes, Instrument	64	Graft Mixing Dishes	125	Master Core Kit	129
Castroviejo Needle Holders	160-162	Graft Packers	122-123	Mayo Stand Covers	41
Cautery Systems	41	Graft Syringes	124-125	Mayo Stands	47
Cement, Implant	207	Growth Factor, AlloSpark-GF™	6	Mead Mallet	101
Centrifuges	68-69	Guide Pins	193	Meisinger Bone Master Core Kit	129
Chair Pads	48	Handpiece Lubricant	92	Meisinger Pin Control	134
Cheek Retractors	150-153	Handpiece, Third Molar	92	Membranes	10-11, 24-28
Chisels, Bone	101	Handpieces, Implant	79	Membrane Forceps	136
Clamps, Suction	96	Hand Scrubs	62	Membrane Tacks	134-135
CO ₂ Monitors	42-45	Headlamp	185	Mesh, Titanium	31
CO ₂ Nasal Hoods	50	HemCon® Dental DressingPro	41	Micro Needle Holders	161, 165, 174
Cold Packs	62	Hemostatic Sponges & Gauze	40-41	Micro Scalpel Blades	167
Collagen Plugs	29	Hemostats	40	Micro Scissors	159, 165, 175
Contrasters, Photography	196	High Speed Handpieces	92	Micro Surgery	164-175

Mindray Vital Sign Monitors	43, 45	Probes, Plastic & Stainless	198	Suction Cut-Off Valves	96
Mirror Handles	198	Prosthetic Wrench, Right Angle	200-202	Suction Irrigator	95
Mirrors	196-198	Prosthetic Parts Drivers	199-203	Suction Tips	98-99
Models, Implant	186-189	PTFE Membrane, CytoSurg®	28	Suction Tubing	96-97
Molt Mouth Gag	150	PTFE Sutures	35	Suction, Emergency	51
Monitors, Vital Sign	42-45	Pulse Oximeters	42-45	Supercut Scissors	157
Motors, Implant	78	Punches, Tissue	182-183	Surgeon's Caps	61
Motor Stands	47	Radiographic Marking Balls	192	Surgical Aspirators	98-99
Mouth Props	153	Ramus Retractors	152	Surgical Gloves	66
N95 Respirator Masks	61	Rapid Roll® Covers	97	Surgical Guide Osteotomes	139
NaCl Sterile Saline	94	Renovix-Plus® Collagen Membrane	25	Surgical Handpieces	79, 92
Nasal Cannula	93	Rescue Kit, Implants & Broken Screws	93	Surgical Motors	78
Nasal Hoods, CO ₂	50	Resnik Implant Pack	58	Surgical Rulers	192
Needle Holders	160-163	Resonance Frequency Analysis (RFA)	191	Surgical Tables	48
Needles, Irrigation	94	Respirator Masks, N95	61	Sussman Implant Positioning Guides	193
Neumeyer Root Planing Bur	85	Retractors, Surgical	150-153	Suture Book	34
Newtron® Ultrasonic Scaler	74	Ridge Contouring Bur	85	Suture Forceps & Pliers	181
No-Clog Aspirators	99	Ridge Mapping Caliper	195	Suture Scissors	39, 157-159
Non-Glare Microsurgery Kit	164-165	Right Angle Prosthetic Wrench	200-202	Sutures	34-38
Nordland Micro Blade	167	Ridge Splitting Instruments	144-147	Synthetic Graft Material, Orastruct®	31
O/R Gowns	60	Rongeurs	148-149	Syringes, Anesthetic	55
O/R Towels	59-60	Root Planing Bur, Neumeyer	85	Syringes, Bone Grafting	124-125
OrACELL® Dermis Allograft	22-23	Root Tip Forceps	114, 118	Syringes, Irrigation	94
OraFoam™ Hemostatic Sponges	40	Rotary Tissue Punches	182	Table Covers	63
OraMEM®+ Collagen Membrane	26	RotoBrush-Titanium™	75	Tables, Surgical	48
OraPLUG® Collagen Plugs	29	Rulers, Surgical	192	Tacks, Membrane	134-135
Orastruct® Synthetic Graft Material	31	Saline, Sterile	94	Tenting Screws	132
OraTAPE & OraCOTE Wound Dressings	30	SalvinMesh-Ti™ Titanium Mesh	31	Third Molar Handpiece	92
Orban Knife	176	SalvinOss® Porcine Xenograft	14-15	Third Molar Retractor	152
OsseoFiber® Allograft Bone Putty	10	SalvinOss®-B Bovine Xenograft	13	Tissue Forceps	178-181
OsseoFiber-PLUS® Allograft Bone Putty	11	SalvinOss® Collagen Xeno + Bovine	16-17	Tissue Forceps, Micro	175
OsseoScrape™ Bone Scraper	126	Saw Discs	83	Tissue Marking Pens	192
Osseous Trap	128	Scalers, Implant	204-205	Tissue Punches	182-183
Osteotomes, Bone Spreading	145	Scalpel Blades & Handles	184	Tissue Rakes	130
Osteotomes, Sinus	138-139	Scissors, Micro	159, 165, 175	Tissue Scissors	157-159
Packs, Sterile	56-59	Scissors, Supercut	157	Titanium Implant Brush	75
Paralleling Guides	193	Scissors, Suture	39, 157-159	Titanium Mesh	31
Parts Cassettes	65	Scissors, Tissue	157-159	Tongue Retractor	152
Patient Drapes	56-59	Scoring Bur, Sinus	85	Torque Wrenches	199-202
PDL-Evators®	68	Screws, Fixation	133	Towel Clamps	60
P.E.T Partial Extraction Therapy Burs	102	Scrub Brushes	62	Towels, O/R	59-60
Penguin RFA Implant Stability Analysis	191	Scrub Towels	61	Trephines	88
Peri-Implant Management Set	76-77	Sectioning Bur	85	Tubing Covers	97
Peri-Implantitis Brushes & Burs	75	Sedation Monitors	42-45	Tubing, Handpiece Irrigation	80, 96
Periacryl	39	Serrated Elevators	107	Tubing, Suction	96-97
Perio Curettes	177	Serrated Periostomes	105	Tunneling Forceps	174
Perio Knives	176	Shoe Covers	61	Tunneling Instruments	164-175
Perio Probes, Plastic & Stainless	198	Side Cutting Needle Holder	160	Twist Elevators	110
Periodontal Surgical Files	176	Silverstein Suture Book	34	Ultrasonic Scaler, Newtron®	74
Periosteal Elevators	154-156	Sinus Crestal Kit	137, 140	Ultrasonic Surgery, Piezotome	70-73
Periostomes	102-105	Sinus Lateral Burs	137, 141	Urban Master Core Kit	129
Pharmaceuticals	6, 144	Sinus Lift Instruments	137-143	Vacutainer®	67
Photography Contrasters	165	Sinus Osteotomes	138-139	Vital Sign Monitors	42-45
Photography Mirror Handle	196	Sinus Rongeur	141, 149	Water, Sterile	94
Photography Mirrors	196-197	Sinus Scoring Diamond	85	"Where Is It?" Parts Cassettes	65
Physics Forceps	117	Socket Caliper	195	Wire Cutter	31
Piezotome Ultrasonic Surgery	70-73	Soft Tissue Manipulation Kit	173	Wrenches, Torque	199-202
Pikos Burs	84-85	Solution Bowl	95	Wrist Clips, ECG	50
Pikos Ramus Retractor	152	Split-Control, Meisinger	147	X-Ray Distortion Markers	192
Pilot Drills	83	Sponges, Collagen	40	Xenograft, Bovine, SalvinOss®-B	13
Pin-Control, Meisinger	134	Spotter Smart Implant Detector	192	Xenograft, Porcine, SalvinOss®	14-15
Plugs	29	Stands, Implant Motor	47	Xeno + Bovine, SalvinOss® Collagen	16-17
Pneumatic Crown Remover	206	Stands, Mayo	47	Xymphony No Memory Membrane	27
Pocket Packers	123	Sterile Packs	56-59	Yankauer Suction	51, 97
Porcine Xenograft, SalvinOss®	14-15	Sterile Saline & Water	94	ZAGA Zygomatic Instrument Kit	121
PowerTack™ Automatic Tack System	135	Sterilization Cassettes	64	Zoe Vital Sign Monitor	44
Precordial Stethoscope	49	Sterilization Wrap & Pouches	63-64	Zygomatic Instruments	121
Prep Cup	95	Sterillium Rub Surgical Hand Scrub	62		
PRF Box	69	Stethoscope, Precordial	49		



Toll Free 800.535.6566
www.salvin.com

Salvin Dental Specialties, LLC
3450 Latrobe Drive
Charlotte, NC 28211 U.S.A.

CHANGE SERVICE REQUESTED

PRSRRT STD
U.S. POSTAGE
PAID
PPCO

**Whether You Are New To Implant Dentistry Or Performing
Advanced Procedures, We've Got The Solutions You Need To
Make Your Life Easier And Your Practice More Productive.**

SALVIN®

*Everything For Your Implant Practice
But The Implants®*