



SALVIN[®] 2026 CATALOG

**Everything For Your Implant Practice
But The Implants[®]**



Dr. Michael A. Pikos
Oral Surgeon
Palm Harbor, FL

"Salvin Dental is an important resource for our practice and our course participants. They are great people to work with and won't let you make a mistake."



Dr. Scott Froum
Periodontist
New York, NY

"Not only does Salvin have a great product line, but their support of academics and advancement of dental technology is unparalleled."



Dr. Hom-Lay Wang
Periodontist
Ann Arbor, MI

"Salvin Dental has an excellent product line. Their commitment to quality and service has been a valuable resource for both my practice and my students."



Dr. Shaun Young
Oral Surgeon
Tampa, FL

"Salvin continues to be a significant portion of our regenerative portfolio due to its consistency. They have always proven to be at the top when it comes to grafting materials and our practice continues to see great results."



Dr. Michael Sonick
Periodontist
Fairfield, CT

"There's a rare, consistent quality to Salvin products and service that's hard to find anywhere else. I use the Salvin catalog confidently for the implant courses and seminars I teach."



Dr. Alvaro Ordonez
General Dentist
Coral Gables, FL

"The team at Salvin Dental is always looking out for my best interest. Relationships like that are hard to find and impossible to replace."



Dr. Dwayne Karateew
Periodontist/Prosthodontist
Chicago, IL

"Every time I place an order with Salvin Dental my expectations for quality and service are exceeded. I consider them to be an invaluable partner in my practice and an asset to the profession."



Dr. Randy Resnik
Prosthodontist
Pittsburgh, PA

"The entire Salvin team is always willing to go above and beyond. I know I can count on them to get me exactly what I need, when I need it."



Dr. James Wilson
Periodontist
Tampa Bay, FL

"For over 20 years Salvin has been my most valuable corporate partner. They are customer satisfaction oriented and have a comprehensive portfolio. Their culture frowns on high-pressure salesmanship making them seem more like friends than sales reps."



Dr. Scott Ganz
Prosthodontist
Fort Lee, NJ

"As implant dentistry has evolved, so has Salvin Dental as a valuable resource to my practice. They are always well informed, and provide a personal touch that is very much appreciated by me and my staff."



Dr. Dmitriy Ivanov
Periodontist,
Redondo Beach, CA

"We have been using Salvin Dental products in our periodontal practices for over 10 years with the highest success. The quality of their products is exceptional, and they consistently deliver outstanding results for our patients."



Dr. Mia Geisinger
Periodontist
Birmingham, AL

"As an educator and a researcher, I am consistently impressed by Salvin Dental's support for residents and new practitioners. They have been an amazing resource for our practice and our students!"



Dr. Brian Jackson
General Dentist
Utica, NY

"I take pride in being able to say that I have been with Salvin Dental Specialties since its very beginning. Today, Salvin is a critical player on my team. I could not imagine performing implant surgery without Salvin."



Dr. Jaime Lozada
Prosthodontist
Loma Linda, CA

"The Salvin catalog is a snapshot of the very best products for the up-to-date implant practice. They always have the newest technology, and going the extra mile is the standard procedure at Salvin."



Dr. Pat Allen
Periodontist
Dallas, TX

"It's the quality people behind the quality products that sets Salvin Dental apart from the rest. Their knowledgeable, helpful, and friendly staff makes working with them a very pleasant experience."



Dr. Christopher Richardson
Periodontist
Richmond, VA

"Salvin Dental Specialties is so much more than a dental supply company. I consider Salvin and their team members to be an essential and collaborative partner in the care I provide my patients."



Dr. Russ Baer
General Dentist
Chicago, IL

"There are three things you need in Implant Dentistry: a well trained dentist, a good implant system, and your Salvin catalog. Salvin's products and customer support are an irreplaceable part of my implant practice."



Dr. Fadi Hasan
Periodontist,
McLean, VA

"What truly sets Salvin apart is their support team and management, who provide a personalized, one-of-a-kind support experience at any time of day or night. Their commitment to excellence makes them an invaluable partner in our daily practice."



Dr. Tim Kosinski
General Dentist
Bingham Farms, MI

"Working with the Salvin team for the past 30 years, I know that I can count on them for top of the line products and world class service. They are a valuable resource for my practice."



Dr. Gordon Christensen
Prosthodontist
Provo, UT

"I consider Salvin Dental to be a premium company in providing implant accessory supplies. The Salvin catalog is not only a resource for my practice, but also for the doctors who attend our courses."



Dr. Jennifer Hirsch Doobrow
Periodontist
Cullman, AL

"Loyalty, integrity, and excellence are what Salvin Dental embodies. Their commitment to impeccable service and outstanding quality is why they are a leader in the dental implant community."



Dr. George Kotsakis
Periodontist
Hoboken, NJ

"It's always a pleasure to work with Salvin Dental because it's always a partnership rather than a business transaction. They genuinely care and always deliver. I would highly recommend Salvin."



Dr. Natalie Wong
Prosthodontist
Toronto, Ontario, Canada

"Delivering excellent quality instrumentation with great service is what you can expect from Salvin. Professional, knowledgeable, and courteous, I consider them an integral part of my dental implant team."



Dr. Bradley McAllister
Periodontist
Tualatin OR

"Salvin has proven to be an innovator time and again. They are my go to company for everything but the implants."



Dr. Burt Melton
Prosthodontist
Albuquerque, NM

"Salvin Dental has all the things you need for a day-to-day implant practice as well as all the special things you need from time to time. Salvin's products are always available for you when you need them in a hurry."



Dr. Bach Le
Oral Surgeon
Whittier, CA

"Salvin Dental has one of the most professional sales teams in the dental industry. They know that my time is valuable and always deliver more than they promise."



Dr. Lee Silverstein
Periodontist
Atlanta, GA

"When Salvin says everything for the implant practice but the implants, they mean it! Salvin is our prime resource for expertise and incredible service. They are a doctor's true friend in our industry."



Dr. Joseph Cristoforo
Periodontist
Madison, WI

"Knowledge. Service. Commitment. Salvin Consistently provides vast product knowledge, unparalleled service, and commitment to enhancing the clinical experience for both the provider and patient."



Dr. Todd Engel
General Dentist
Charlotte, NC

"Salvin Dental consistently delivers great products and excellent service. They have become an important resource for my practice."



Dr. Peter Nordland
Periodontist
La Jolla, CA

"The best aspect of Salvin Dental is the personal attention and knowledge that all of the staff members have whenever I call. I am always greeted with a kind, courteous, knowledgeable representative."



Dr. George Romanos
Periodontist
Stony Brook, NY

"The Salvin Group is the best partner of my practice due to the immediate customer service and professionalism. Their products are of high quality and a joy to work with them."



Dr. Carlos Rubio
General Dentist
Algodones, Mexico

"Salvin has been a major part in my professional growth as a dentist and a teacher. I have always been able to count on them to have the best quality material, as well as the service. Thanks Salvin for being there along the way."



Dr. Richard Rasmussen
Periodontist
Tampa, FL

"Salvin has a great team with one thing in mind...To make my practice and me their first priority, every time I call. I know that when I call Salvin, I've called the right place."



Dr. Sayward E. Duggan
Periodontist
Yorktown, VA

"Salvin consistently continues to set the bar high in the dental industry with excellent products, competitive prices, and above all superior client services. Thank you Salvin for your loyal professional support."



Dr. Tyler Borg
Periodontist
Northglenn, CO

"The team at Salvin Dental offers consistent support, continued product innovation, and top notch customer service. All in the pursuit of the ideal doctor and patient experience alike."



Dr. Robert Heller
Prosthodontist
Lewis Center, OH

"I confidently recommend Salvin Dental to our Midwest Institute students. Their knowledge and helpfulness make them a pleasure to do business with."



Dr. Ross Chavkin
Oral Surgeon
Melville, NY

"We rely on Salvin regenerative products because they consistently deliver predictable, high-quality outcomes. Their materials support reliable bone and soft-tissue healing, allowing us to provide patients with the strongest foundation for long-term oral health."



Dr. Craig M. Misch
Oral Surgeon/Prosthodontist
Sarasota, FL

"Salvin has a thorough product selection with great personalized customer service. I keep going back to Salvin for whatever I need for my implant practice and highly recommend them to my colleagues."



Dr. Rick Ferguson
General Dentist
Pembroke Pines, FL

"Dealing with Salvin is a first-class experience. They have a reputation for always putting their customers first. They provide great service, value, and expertise to my practice and the students at our course."



Dr. Jack Krauser
Periodontist
Boca Raton, FL

"The Salvin catalog is perhaps the most valuable tool in implant dentistry. Their unique products and excellent service are a tremendous resource for my private practice and for our students at the school."



Dr. Jason Stoner
Periodontist
Dublin, OH

"Every time I call Salvin, I can count on speaking with a knowledgeable person with the answers I need. They know that my time is valuable and always deliver what they promise."



Dr. Adam S. Foleck
General Dentist
Norfolk, VA

"Salvin has provided expertise and experience necessary to take my implant practice to the next level. They do everything they can to ensure my success. Salvin helps me equip my offices to offer the highest level of treatment to my patients."



Dr. John Russo
Periodontist
Sarasota, FL

"Salvin Dental is a great resource for my practice and especially for the doctors who take my courses. Salvin has what you need whether you're getting started or doing advanced procedures."



Dr. Michael Pruett
General Dentist
Augusta, GA

"I love the way that Salvin does business. Their willingness to go the extra mile and their commitment to serving makes it a pleasure to deal with them."



Dr. Nick Shumaker
Periodontist
Fort Collins, CO

"When it comes to quality products and great customer service, nobody does it better than Salvin!"



Dr. Fred Sakamoto
Periodontist
Blacklick, OH

"Salvin Dental runs a first class operation with excellent customer service and products. We consider Salvin as part of our team with periodontal, implant and related surgical cases."



Dr. Kanyon Keeney
Oral Surgeon
Mechanicsville, VA

"Salvin Dental has one of the most professional sales teams in the dental industry. They know that my time is valuable and always deliver more than they promise."



Dr. Lauralee Nygaard
Periodontist
Spokane Valley, WA

"When Salvin says everything for the implant practice but the implants, they mean it! Salvin is our prime resource for expertise and incredible service. They are a doctor's true friend in our industry."



Dr. Etienne Lacrampe
Periodontist
Walnut Creek, CA

"I use Salvin for all my biologics. There is no one I trust more to provide the best results for both my patients and my referring dentists."



Dr. Richard Martin
Oral Surgeon
Lewisville, TX

"When you work with Salvin Dental, you do not feel like you are working with a company, you feel like you are working with your family. Their attention to detail and reliability is unmatched in the dental field."



Dr. Ryan Voelkert
Periodontist
Greenville, SC

"I have been working with Salvin for over 20 years. The quality of their products, prices and level of service are unmatched. There simply is not a company I trust more! It's hard to imagine our surgical/implant practice without them."



Dr. Avi Schetritt
Periodontist
Miami, FL

"Great products, great service, and great people! Their professionalism is what sets Salvin apart. They deliver every time!"



Dr. Bradley Ross
Periodontist
Miami, FL

"Salvin provides top quality surgical armamentarium and biomaterials, but their people truly make this company an essential asset to my practice. Their knowledge, integrity, and customer service are second to none."



Dr. Matthew Fien
Periodontist
Deerfield Beach, FL

"I am always impressed with the professional service and knowledge from the team at Salvin Dental. They make it easy for me to recommend them to my colleagues."



Dr. Steve Wallace
Periodontist
Wilmington, NC

"I turn to Salvin when I need answers to my questions. Their knowledge and helpfulness make them a pleasure to do business with."



Dr. Justin Zalewsky
Periodontist
Alexandria, VA

"Salvin has been a business partner for years, because they truly do partner with us. Their products and services are unmatched. Their people care about our practice as much as we do."



Dr. Nancy Newhouse
Periodontist
Independence, MO

"Nobody does it better! Salvin is just what they say they are: 'Everything For Your Implant Practice But The Implants' all wrapped up in the best service in the industry."



Dr. David Lee Hill
Oral Surgeon
Chapel Hill, NC

"I can't think of a single source that offers a more comprehensive quality line of products to support my implant practice than Salvin Dental. They always have what I need."



Dr. James Woodyard
Periodontist
Newburgh, IN

"Knowing that I can expect the right products and great service is why I choose to do business with Salvin. Dealing with them is a first-class experience every time I call."



Dr. James Rutkowski
General Dentist
Clarion, PA

"Honesty and integrity are values that come to mind when working with Salvin Dental. I know that I'm going to be taken care of every time I call."



Dr. Israel Puterman
Periodontist
Chevy Chase, MD

"Getting the Salvin catalog is like being a kid again, flipping through the toy catalog and seeing all the things that I want. The quality of the products is only surpassed by the quality of the people. They truly want your success above all else."



Dr. Paul Fugazzotto
Periodontist
Milton, MA

"Salvin has, in my experience with them over decades, consistently demonstrated a degree of honesty and openness, and a genuine desire to help the clinician, that is unsurpassed in the dental community."



Dr. Jason Souyias
Periodontist
Port Huron, MI

"The team at Salvin is amazing! They know customer service better than anyone in dental. That is why they are so easy to work with."

SALVIN®

2026 COMPREHENSIVE CATALOG

We're in business to make your lives easier
and improve your surgical outcomes.
Call us, we'll take great care of you!

OUR VALUES



Deliver Surgical Excellence

Provide the highest-quality supplies and biomaterials that empower dental surgeons to achieve predictable outcomes.



Share Knowledge & Expertise

Be more than a supplier, be a trusted source, guiding our customers with deep industry knowledge and expert-driven solutions.



Provide Exceptional Service & Reliability

Support every customer with fast, personal and dependable service – because when it comes to surgery, every detail matters.



Build Lasting Relationships & Trust

Earn customer loyalty through integrity, reliability and a commitment to being the most trusted name in dental surgery.

800-535-6566 [Toll Free U.S. & Canada] • www.salvin.com



Salvin Dental Specialists, LLC
3450 Latrobe Drive
Charlotte, NC 28211 USA
Ph: 704-442-5400 Fax: 704-442-5424
Email: orders@salvin.com
Hours: M-Th: 8:30-7:00 F: 8:30-6:00
(U.S. Eastern Standard Time)

Quality Management System
ISO 13485:2016 Certified

See It Demonstrated!
Look For This Icon Throughout
Our Catalog



Sustainable Forest Initiative



By supporting the SFI program, Salvin is helping to strengthen forest practices in North America and promote responsible sourcing globally.

SALVIN INSTRUMENT LIFETIME WARRANTY

Salvin Branded Instruments are guaranteed to be free from defects in workmanship and material.

Any Salvin instrument which proves defective in workmanship or material will either be repaired or replaced at our discretion, without charge.

Warranty does not cover routine maintenance, re-tipping of carbide instruments or sharpening of scissors. Individual manufacturer warranties apply to products other than instruments.

Terms Of Sale: Payment terms are Visa, MasterCard, Discover or American Express. Prices in catalog are in U.S. funds and are subject to change due to unforeseen circumstances. Sales and Taxes, if applicable, are customer's responsibility. **International Orders:** International orders may be called in to 704-442-5400, or placed online at www.salvin.com. Duties and import taxes, if any, are customer's responsibility. **Hospital/School Orders:** Institutional orders should be emailed to orders@salvin.com. Terms for hospital/school orders are net 30 with 1.5% per month late fee charged on past due accounts. Customers are responsible for legal and collection fees if required. **Trademarks:** "Salvin", "Salvin Dental Specialists" Logo, "Everything For Your Implant Practice But The Implants", "Salvin Regenerative" Logo, AccessTorq®, AlloSculpt®, Salvin AmnioExcite®, AtrauLux®, AutoTome®, CytoSurg®, EasyReach®, OraCOTE®, OraFoam®, OraMEM®, OraMEM⁺, OraPLUG®, OraTAPE®, Orastruct®, OsseoFiber®, OsseoFiber-Plus®, PDL-Evator®, Rapid Roll®, Renovix-Plus®, RotoBrush-Titanium®, Salvin AutoTorq®, Salvin Power Point®, SalvinOss®, SalvinOss^{-B}, SalvinOss® Collagen, SalvinTorq®, SurgiPoint® and Xymphony® are registered U.S. and international trademarks of Salvin Dental Specialists, LLC.

Salvin AlloSpark-GF™, Salvin Amplifi™, Enlighten-Pro3™ LED System, Salvin Black Kryptonite™, and Salvin Mesh-Ti™ are trademarks of Salvin Dental Specialists, LLC. Other product brand names may be trademarks or registered trademarks of their respective owners and are mentioned for reference purposes only.

CONTENTS

6-7

ORAGEN
COMPLETE BONE SOLUTION

8

ALLOSPARK-GF
SUPPLEMENTAL ALLOGRAFT
DERIVED PROTEINS

9

FIBERPUTTY ALLOGRAFT PUTTY
IN A SYRINGE

10

ALLOSCULPT-3D PLUS
DBM PUTTY WITH
CORTICAL & CANCELLOUS CHIPS

11

OSSEOFIBER-PLUS
DEMIN CORTICAL FIBERS +
MINERALIZED CANCELLOUS CHIPS

OSSEOFIBER
DEMIN CORTICAL FIBERS

12-13

AMNIOEXCITE AMNION / CHORION
FULL THICKNESS PLACENTAL
MEMBRANE

14-17

LIFENET HEALTH ALLOGRAFT BONE

18-19

ORACELL
DECELLULARIZED DERMIS
ALLOGRAFT

20-21

SALVINOSS
XENOGRAFT PARTICULATE

22-23

SALVINOSS COLLAGEN
XENOGRAFT + COLLAGEN PLUG

24

SALVINOSS-B
BOVINE XENOGRAFT BONE GRAFT
MATERIAL

25

ORASTRUCT
80/20 TCP/HA GRAFT MATERIAL

SALVIN MESH-TI TITANIUM MESH

MESH & WIRE CUTTER

26

ORATAPE
RESORBABLE COLLAGEN
WOUND DRESSING

ORACOTE
RESORBABLE COLLAGEN
WOUND DRESSING

27

ORAPLUG
RESORBABLE COLLAGEN PLUG

FULFILL - 80% CALCIUM PHOSPHATE
GRANULES COMBINED WITH
20% BOVINE COLLAGEN

28

ACCELL DBM PUTTY IN A SYRINGE

29

DYNABLAST
DEMIN CORTICAL + MIN CANCEL-
LOUS CHIPS PASTE OR PUTTY

DYNAGRAFT
DEMIN CORTICAL PASTE OR PUTTY

30

SALVIN AMPLIFI
LONG TERM MEMBRANE

31

RENOVIX-PLUS
RESORBABLE EXTRA CELLULAR
MATRIX

32

ORAMEM+
RESORBABLE COLLAGEN MEMBRANE

33

XYMPHONY
DRAPEABLE PORCINE PERICARDIUM

34

CYTOSURG
NON-RESORBABLE PTFE MEMBRANE

35-40

SUTURES & WOUND CLOSING

Sutures 35-39
Periacryl Cyanoacrylate 40
Featherlite Suture Removal Scissors . . 40

41-42

BLEEDING CONTROL

HemCon® Dental DressingPro 41
Portable Electrocautery 41
OraFOAM™ Hemostatic Sponges . . . 42
BloodSTOP® Hemostatic Gauze 42
Hemostats 42
Bone Wax 42

43

ORASOO[®]THE ORAL HYDROGEL WOUND DRESSING

PERIODT GEL BIOFILM

44-46

VITAL SIGN MONITORS

Pulse Ox / BP Monitors	44
Wireless Sedation Stethoscope	45
Pulse Ox / BP / ECG / CO ₂ Monitors	45-46

47-53

DEFIBRILLATOR & SEDATION

Defibrillator	47
Nasal Cannula & Oxygen Adapter	48
Anesthesia Instrumentation	48
Emergency Kits	48
CO ₂ Nasal Hoods	49
ECG Patches & Wrist Clips	49
Emergency Airways	49
IV Catheters & Admin Sets	50-52
Aspirating Anesthetic Syringes	53

54-56

O/R FURNITURE

Implant / Surgical Cart	54
Mayo & IV Stands	55
Table Covers & Mayo Stand Covers	55
Surgical Tables	56
Sedation Chair Pad	56

57

HEADLAMPS

Enlighten-Pro3 [™] LED Headlamp	57
--	----

58-68

INFECTION CONTROL / PPE

Sterile Drapes, Gowns & Covers	58-62
N95 Masks	63
Towels, Caps & Shoe Covers	63
Patient Skin Prep	63
Cold Packs	64
Scrub Brushes	64
Instrument Care & Cassettes	65-67
Surgical Gloves	68

69-71

CENTRIFUGES

Blood Collection	69
BIO-PRF Centrifuge	70
Variable Speed Centrifuge	71

72-76

PIEZO ULTRASONIC SURGERY & ULTRASONIC SCALERS

Piezotome Ultrasonic Surgery Units	72-73
Piezotome Tips & Procedure Kits	74-75
Newtron Ultrasonic Scaler	76
Air-N-Go Air Polishing	76

77-79

PERI-IMPLANTITIS MANAGEMENT

RotoBrush-Titanium [®] Brushes	77
Peri-Implantitis Instruments	78-79

80-92

MOTORS, HANDPIECES & DRILLS

Implant Motors	80
Implant Handpieces	81
Irrigation Tubing	82
Implant Drills	82
Bone Saws	83
Drills, Burs, Trephines & Bur Blocks	83-89
AutoTome [®] Automatic Osteotome Kit	90
AutoTome [®] Automatic Extraction Kit	91
Rear-Exhaust High Speed Handpiece	92

93

SALVIN IMPLANT RESCUE KIT

Broken Screw & Failed Implant Removal	93
---	----

94-99

SUCTION & IRRIGATION

IrriGUIDE Adapter	94
Irrigation	94-95
Suction	96-99

100-183

SURGICAL INSTRUMENTS

Atraumatic Extraction	100-116
Mallets & Chisels	117
Bone Curettes	118
Zygomatic Instruments	119
Bone Grafting	120-129
Tenting Screw Kit	130
Fixation Screw Kit	131

Membrane Tack Kits	132-133
Membrane Instruments	134
Sinus Lifts	135-141
Ridge Splitting	142-143
Bone Spreading	144-145
Rongeurs	146-147
Retractors	148-151
Periosteal Elevators	152-154
Scissors	155-157
Needle Holders	158-161
Micro Surgery & Tunneling	162-173
Perio Knives & Surgical Files	174
Perio Curettes & Bone Files	175
Tissue Forceps	176-180
Tissue Punches & Nippers	181
Scalpel Blades & Handles	182-183

184-190

IMPLANT & COURSE MODELS

Implant Models	184-187
Bone Loss & Course Models	188

189

IMPLANT STABILITY MONITOR

Penguin RFA Implant Stability Monitor	189
---	-----

190-196

DIAGNOSTICS

Spotter Smart Implant Detector	190
Drill Guides	191
Calipers & Gauges	192-193
Photo Mirrors & Contrasters	194-195
Diagnostic Probes	196
Operative Mouth Mirrors	196

197-207

RESTORATIVE & HYGIENE

Implant Buddy Driver Tip Kits	197
EasyReach [®] Right Angle Wrench	198
AutoTorq [®] Electric Torque Driver	198
SalvinTorq [®] Torque Wrench	199
Access-Torq [®] -Plus RA Torque	200
Abutment Driver Tips	200-201
Implant Curettes	202
American Eagle Scalers & Curettes	203
Abutment Holders	204
Implant Scalers	204
Crown Grippers	205
Crown Splitters	205
Improv [®] Implant Cement	206
Preat Partnership	207

Complete Index On Pages 208-209

OraGEN®

Cellular Bone Allograft

A Predictable Graft For The Most Challenging Cases^{1,2}

ORAGEN® IS A COMPLETE BONE GRAFT^{1,2}

Successful Bone Remodeling Requires The Following Three Properties:



OSTEOCONDUCTIVITY
Scaffold



OSTEOINDUCTIVITY
Signals



OSTEOGENICITY
Cells



What Is OraGen®?

OraGen® has three components^{1,2}:

- **Osteogenic:** lineage-committed bone cells that directly participate in bone formation
- **Osteoconductive:** corticocancellous chips that provide a natural scaffold for cells to attach, migrate, and proliferate
- **Osteoinductive:** demineralized bone that exposes natural growth factors within the matrix

How Does OraGen® Work?

- OraGen® contains lineage-committed bone cells and demineralized bone matrix, which can begin producing bone matrix shortly after implantation.³⁻⁵
- Histological analysis of OraGen® suggests that new bone formation started directly after implantation and continued at a rapid pace.³⁻⁵

Storage And Shipping

- **Storage and shipping:**
OraGen® is stored and shipped frozen
- **Specialized packaging:**
the product is packaged in a thermal shipper designed to maintain low temperatures during transport
- **Fast delivery:**
to preserve cell viability, OraGen® is typically shipped with next-day delivery options
- **Rapid thaw time:**
specialized packaging allows for rapid thaw time – less than 5 minutes!

When storing OraGen® long-term:
store at -70°C or colder. Do not store in liquid phase of liquid nitrogen (LN₂), 12 month shelf life.

Call us for information on special freezers

OraGEN® Cellular Bone Allograft

Size	Product #	Fee	4 or More
0.7 cc	#ORAGEN-0.7	650.00	550.00
1.3 cc	#ORAGEN-1.3	795.00	675.00

1. LifeNet Health, PQ13-033 and ES 22-019

2. LifeNet Health, ES-17-011

3. Are Cellular Bone Matrix Allografts a Viable Option for Mandibular Tissue Engineering and Reconstruction. Marshall, J. S., Davis, S. S., Jones, L., & Kushner, G. M. (2024). Journal of Oral and Maxillofacial Surgery. <https://doi.org/10.1016/j.joms.2024.05.040>.

4. Tentindo, G. (a), Brown K, Gianulis E (b), Softic D (b), Qui X (b). Histology of ViviGen® Mandibular Clinical Explant Reveals Extensive Bone Formation in a Patient with Chronic Osteomyelitis. (a) Oral Facial Surgery Institute, 621 S New Ballas Rd, Saint Louis, MO, USA; (b) LifeNet Health, 1864 Concert Dr., Virginia Beach, VA, USA

5. Gianulis et al. Cell and Tissue Banking (2022) E-pub ahead of print. Doi: 10.1007/s10561-022-10044-2 Khan SN, Cammisa FP Jr., Sandhu HS, Diwan AD, Girardi FP, Lane JM. The biology of bone grafting. J Am Acad Orthop Surg. 2005; 13(1):77-86. <https://doi.org/10.1097/00005131-189809000-0002>.

OraGen® Case Study

Case Study - Lateral Sinus Elevation Using 100% OraGen®

Case & Photos: Dr. Edgard El Chaar

Patient & Diagnosis:

A 65-year-old male, missing second premolar and with one first molar that is supporting a three-unit bridge. Recurrent cavities were detected that required removal of restoration. **The site in question has less than 1mm in crestal bone available.**

Previous Treatment:

None other than the extraction and making the three unit bridge.

Intervention:


A lateral sinus elevation using OraGen 100% placed.

Outcome:

An adequate amount of bone regeneration occurred for the placement of an implant in the position of the second premolar. **The patient was able to be successfully restored with three individual crowns despite having <1mm of crestal bone prior to treatment.**



Lack of the crestal bone as seen in radiograph



In these three clinical views, we can appreciate the lateral window and the placement of 100% OraGen



Immediately post grafting



One year post implant placement



"OraGen® has expanded the potential of Guided Bone Regeneration. It has provided us clinicians more possibilities for patients who were faced with limited options that would have caused them a lot of morbidities."

—Dr. Edgard El Chaar
Periodontist, New York, NY






Enhance The Regenerative Potential Of Your Choice Of Graft Material!

- With A Natural Blend Of Allograft Derived Growth Factors, AlloSpark-GF Is A Perfect Complement To Any Commercially Available Graft Material
- Contains Nearly 1,000 Unique Proteins, Peptides & Growth Factors Harnessed Through Proprietary Processing Methods
- Provides A Powerhouse Of Osteoinductive, Angiogenic, Proliferative & Chemotactic Properties
- Every Lot Is Tested For Protein Content To Provide Quantified Levels Of Bone Healing Contribution
- Store At Ambient Temperature – 11.0-30.0°C (51.8-86.0°F)

**A High Concentration
Of Growth Factors Include,
But Are Not Limited To:
BMP-2, PDGF, TGF-1 & VEGF**

Concentration amount of various growth factors (ng/g)



SALVIN® AlloSpark-GF™ Supplemental Allograft Derived Proteins

Size	Product #	Fee	7 or More
0.5 cc	#ALLOSPARK-GF-0.5	240.00	200.00
1.0 cc	#ALLOSPARK-GF-1.0	420.00	350.00
2.5 cc	#ALLOSPARK-GF-2.5	900.00	750.00





Moldable and cohesive making it easy to place; maintains volume in irregular defects



Patented fiber process – provides a natural scaffold to support optimal bone growth^{3,4}



Ready-to-use, pre-filled syringe – no graft preparation or hydration prior to placement

100% Allograft Putty – No Carriers, Fillers Or Binders

- 100% human allograft composition significantly reduces slumping to maintain bone volume.
- Patented fiber process creates a moldable, cohesive putty that maintains its integrity when exposed to the oral environment.
- Cortical Bone Fibers with a proven 4.1x higher OI (Osteoinductive) score than Cortical Particulate provide a scaffold to support optimal bone growth.^{1,2}
- Ideal handling characteristics for irregular defects or cases where bone particulate alone might not hold.

Indicated for:

Sinus Lifts, Ridge Augmentations, Exposed Implant Threads

SALVIN FiberPutty™ 100% Allograft Putty In A Syringe

Size	Product #	Fee	7 or More
0.5 cc	FIBERPUTTY-0.5CC	139.00	118.00
1.0 cc	FIBERPUTTY-1.0CC	199.00	169.00
2.5 cc	FIBERPUTTY-2.5CC	399.00	339.00

Store at ambient temperature – 3-year shelf life

1. DCI Data on file.
2. DCI Data on file.
3. US Patent 11,927,393 Fiber Slurry Tray and Process
4. US Patent D977,672 Tray and Baffles Design





*Moldable & Easily Shaped
To Pack Into Any Defect*

DBM Putty With Cortical & Cancellous Chips

- 1-2mm Cortical & Cancellous Chips Provide An Osteoconductive Scaffold For Bone Cell Migration
- Inert Gelatin Carrier Is Easy To Shape & Mold When Hydrated With Blood Or Saline
- The Graft Sets At Body Temperature And Resists Migration Once Implanted
- Terminally Sterilized To Achieve SAL 10^{-6} Using Low Temperature, Low Dose Gamma Irradiation
- Each Lot Is Tested In Vivo For Osteoinductive* Potential

SALVIN® AlloSculpt-3D PLUS DBM Putty With Cortical & Cancellous Chips

Size	Product #	Fee	7 or More
1.0 cc	#ALLO-3D-PLUS-1.0	242.00	206.00
2.0 cc	#ALLO-3D-PLUS-2.0	371.00	330.00

No Refrigeration Required - Store At Room Temperature
Product Shelf Life - 2.5 Years



***New Manufacturing Efficiency Creates
Significantly Reduced Cost!***

*These implants were evaluated in a human clinical study and were shown to induce bone formation. Each lot is tested using the athymic nude rat assay to verify osteoinductive potential.

**100% Bone – No Carrier – Handles Like Putty When Hydrated
Facilitates Cellular Proliferation & Cellular Attachment¹**

Demineralized Cortical Fiber PLUS Mineralized Cancellous Chips

- Demin Cortical Bone Has The Potential To Stimulate New Bone Growth
- Mineralized Cancellous Chips Provide Scaffold For Bone Cell Migration
- For GBR Procedures, Ridge Augmentation, Sinus Lifts & Bony Defects
- Each Lot Is Tested In Vivo For Osteoinductive Potential



*Moldable & Easily Shaped
To Pack Into Any Defect*



SALVIN[®] OsseoFiber-PLUS[®] Demineralized Cortical Fibers **PLUS** Mineralized Cancellous Chips

Size	Product #	Fee	7 or More
0.5 cc	#OFIBER-PLUS-0.5CC	170.00	150.70
1.0 cc	#OFIBER-PLUS-1.0CC	290.00	260.80
2.5 cc	#OFIBER-PLUS-2.5CC	545.00	494.50

No Refrigeration Required – Store At Room Temperature
Typical Shelf Life: 3 Years



Case Study - Guided Bone Regeneration Using OsseoFiber-Plus[®]


Case & Photos: Dr. Paul Rosen



*Post Surgical Extraction
– Osseous Defect*



*OsseoFiber-Plus[®] Has
Putty-Like Handling*



*Successful GBR
6 Month Post-Op*



*Implant Placement
Prior To Restoration*



"The innovation of making DFDBA in a fiber form has been a tremendous enhancement to this graft. The increase in surface area makes this osseoinductive clot stabilizer a must for one's regenerative armamentarium. Importantly, this fiber form of DFDBA transforms it into being much more user friendly."

–Dr. Paul Rosen

Periodontist, Yardley, PA & New York, NY

Demineralized Cortical Fiber

SALVIN[®] OsseoFiber[®] Demineralized Cortical Fibers



Size	Product #	Fee	7 or More
0.5 cc	#OFIBER-0.5CC	163.00	144.70
1.0 cc	#OFIBER-1.0CC	277.00	247.80
2.5 cc	#OFIBER-2.5CC	542.00	492.60

No Refrigeration Required – Store At Room Temperature
Typical Shelf Life: 3 Years



1. Rodriguez RU, Kemper N, Breathwaite E, et al. Demineralized bone matrix fibers formable as general and custom 3D printed mold-based implants for promoting bone regeneration. Biofabrication. 2016;8(3):035007






Full-Thickness Placental Membrane

- 4x Thicker Than An Amnion Or Chorion Membrane Making It Much Easier To Handle
- Does Not Fold In Upon Itself When Rehydrated

Structure of Salvin AmnioExcite® Full-Thickness Placental Membrane



AmnioExcite® Includes Both Amnion And Chorion PLUS The Trophoblast Layer
The Trophoblast Layer Contains More Than 50% Of
The Biologically Relevant Proteins In The Placental Membrane¹

SALVIN® AMNIOEXCITE® Amnio / Chorion Membrane

Size	Product #	Fee	7 or More
10mm x 10mm	#AMNIO-MEM-10X10	104.00	91.00
15mm x 20mm	#AMNIO-MEM-15X20	198.00	182.00
20mm x 30mm	#AMNIO-MEM-20X30	299.00	275.00


1. LifeNet Health®. TR-004-2020 Characterization of the Amnion, Chorion, and Trophoblast Layers of Decellularized and Freeze-Dried Placental Membrane. 2020.



AmnioExcite® Case Studies

Case Study #1 - Socket Preservation Using AmnioExcite® Full-Thickness Placental Membrane


Case & Photos: Dr. Nick Shumaker



Surgical Extraction



*LifeNet Health®
Cortical / Cancellous
Allograft*



*AmnioExcite®
Barrier Membrane*



*3 Week Post-Op –
Excellent Tissue Healing*

Case Study #2- Cemental Tear/Osseous Defect Repair Using AmnioExcite®

Case & Photos: Dr. Nick Shumaker



*Cemental Tear /
Osseous Defect*



*LifeNet Health®
Cortical / Cancellous
Allograft*



*AmnioExcite®
Barrier Membrane*



*6 Month - No Recession
Excellent Aesthetic Result*




"As a user of amnion membranes for over 15 years, I have found AmnioExcite® from Salvin to be a great addition to my practice. AmnioExcite® is much easier to handle and place, but still adapts to adjacent structures such as bone, tooth roots or dental implants. It has the perfect balance of adaptation and handling that has been missing in amnion membranes, and I'm very pleased with the results."

—Dr. Nick Shumaker
Periodontist, Fort Collins, CO



**Over 11 Million Bio-Implants Have Been Processed
Using LifeNet Health's® Allowash Technology,
With No Disease Transmission**



*LifeNet Health® Is The World's Largest Full Service
Tissue Processor / Organ Procurement Organization*

*Of Total Donors Screened For Transplantation,
Only 3% Are Accepted By LifeNet Health®*

*All Of LifeNet Health's® Donor Tissue
Is Recovered In The United States*

Specific Screening Criteria*	LifeNet Health® Saving Lives. Restoring Health. Giving Hope.	AATB (American Association of Tissue Banks)	FDA (U.S. Food & Drug Administration)
Encephalitis (clinically active)	•	•	•
Gonorrhea (within the past year)	•	•	•
Hepatitis B	•	•	•
Hepatitis C	•	•	•
HIV 1/2	•	•	•
Malaria	•	•	•
Meningitis (clinically active)	•	•	•
Sepsis	•	•	•
Syphilis	•	•	•
Transmission spongiform encephalopathy (TSE)	•	•	•
Vaccinia	•	•	•
West Nile Virus (WNV, 4 month exclusion from the time of onset)	•	•	•
Guillain-Barre syndrome (clinically active)	•	•	•
Multiple sclerosis	•	•	•
Leprosy (Hansen's disease)	•	•	•
Clinically significant metabolic bone disease	•	•	•
Polyarteritis nodosa	•	•	•
Rabies	•	•	•
Rheumatoid arthritis (musculoskeletal deferral)	•	•	•
Sarcoidosis	•	•	•
Systemic lupus erythematosus	•	•	•
Systemic mycosis	•	•	•
Tuberculosis (clinically active or treated in the past 5 years)	•	•	•
Ankylosing spondylitis	•	•	•
Autoimmune vasculitis if type is not known	•	•	•
Cancer	•	•	•
Chagas disease	•	•	•
Epstein Barr virus (clinically active)	•	•	•
Clostridium difficile infection (defer if occurred during present hospitalization)	•	•	•
Endocarditis (clinically active)	•	•	•
Herpes (any evidence of active herpes infection is a deferral)	•	•	•
Illicit drug use (see Illegal drugs)	•	•	•
Peritonitis (clinically active)	•	•	•
Rheumatic fever (defer all tissues if clinically active, defer heart valves for any history)	•	•	•
Poliomyelitis (clinically active or untreated)	•	•	•
Mixed connective tissue disease	•	•	•
Cold agglutinin disease	•	•	•
Pyelonephritis (clinically active)	•	•	•
Reactive arthritis (must be asymptomatic for 6 months for musculoskeletal deferral)	•	•	•
Wegener's granulomatosis	•	•	•

**These Screening Criteria Are Subject To Change*

LifeNet Health®, The LifeNet Health® Logo, OrACELL®, OraGRAFT® & Preservon® Are Registered Trademarks Of LifeNet Health®.

LifeNet Health® OraGRAFT®

Mineralized 50/50 Cancellous/Cortical Mix **MOST POPULAR**



Particle Size	Size	Product #	Fee	7 or More
250-1000 microns	0.25 cc	#CAN-CORT-.25	70.00	61.00
	0.50 cc	#CAN-CORT-.50	102.00	89.60
	1.00 cc	#CAN-CORT-1.00	138.00	122.80
	2.00 cc	#CAN-CORT-2.00	199.00	178.40

Mineralized Cancellous **MOST POPULAR**



Particle Size	Size	Product #	Fee	7 or More
250-1000 microns	0.25 cc	#MIN-CAN-.25	92.00	80.80
	0.50 cc	#MIN-CAN-.50	106.00	92.70
	1.00 cc	#MIN-CAN-1.00	143.00	126.40
	2.00 cc	#MIN-CAN-2.00	220.00	196.60
1000-2000 microns	0.50 cc	#MIN-CAN-L-.50	106.00	92.70
	1.00 cc	#MIN-CAN-L-1.00	143.00	126.40
	2.00 cc	#MIN-CAN-L-2.00	220.00	196.60

Mineralized Cortical



Particle Size	Size	Product #	Fee	7 or More
250-1000 microns	0.25 cc	#MIN-CORT-.25	50.00	42.30
	0.50 cc	#MIN-CORT-.50	70.00	60.40
	0.70 cc	#MIN-CORT-.75	87.00	75.70
	1.20 cc	#MIN-CORT-1.25	101.00	88.50
	2.50 cc	#MIN-CORT-2.50	150.00	132.80
250-710 microns	0.50 cc	#MIN-CORT-SM-.50	70.00	60.40
	0.70 cc	#MIN-CORT-SM-.75	87.00	75.70
	1.20 cc	#MIN-CORT-SM-1.25	101.00	88.50
	2.50 cc	#MIN-CORT-SM-2.50	150.00	132.80



LifeNet Health® OraGRAFT®

70/30 Mineralized/Demineralized Cortical Mix **MOST POPULAR**

Particle Size	Size	Product #	Fee	7 or More
250-1000 microns	0.50 cc	#MD-7030-.50	97.00	86.60
	1.00 cc	#MD-7030-1.00	140.00	125.50
	2.00 cc	#MD-7030-2.00	198.00	178.70

Demineralized Cortical



Particle Size	Size	Product #	Fee
250-1000 microns	0.25 cc	#DEMIN-CORT-.25	58.00
	0.50 cc	#DEMIN-CORT-.50	77.00
	0.70 cc	#DEMIN-CORT-.75	97.00
	1.20 cc	#DEMIN-CORT-1.25	113.00
	2.50 cc	#DEMIN-CORT-2.50	167.00

Graft Site Suggested Volume Guide

Periodontal Defects

Central Incisor	0.5 cc
Lateral Incisor	0.5 cc
Canine	0.5 cc
1st & 2nd Premolar	0.5 cc
1st & 2nd Molar	0.7-1.0 cc
3rd Molar	0.7-1.0 cc

Extraction Sites

Central Incisor	0.7-1.0 cc
Lateral Incisor	0.5 cc
Canine	0.7-1.0 cc
1st & 2nd Premolar	0.7-1.0 cc
1st & 2nd Molar	1.2-2.0 cc
3rd Molar	1.2-2.0 cc

Sinus Grafts

Each Sinus	2.5-5.0 cc
------------	------------

For optimal safety, all LifeNet Health® allografts distributed by Salvin Dental are gamma irradiated to achieve a Sterility Assurance Level (SAL) 10^{-6} without compromising the inherent biomechanical or biochemical properties of the graft. Typical Shelf Life: 4 Years. Fees are subject to change. Please check current fees at time of placing order.

LifeNet Health® Canadian Registration: CTO #100038

LifeNet Health® OraGRAFT®

Pre-Hydrated Cortical Plate

Ready To Use Out Of The Package

- Preserved Using LifeNet Health's® Proprietary Preservon® Technology To Maintain A Fully Hydrated State At Ambient Temperature
- Preservon-Treated Grafts Have Been Shown To Have Greater Flexibility Than Freeze-Dried Grafts¹
- Ideal For "Shell" Technique – Secure Cortical Plate To Host Bone With Fixation Screws To Create A Biologic Container




Plate	Size	Product #	Fee
Cortical Plate	30 mm x 15 mm x 1 mm	#CORT-PLATE-30X15X1	545.00

No Refrigeration Required - Store At Room Temperature
Typical Shelf Life: 5 Years

Pre-Hydrated Ilium Block

Ready To Use Out Of The Package

- Preserved Using LifeNet Health's® Proprietary Preservon® Technology To Maintain A Fully Hydrated State At Ambient Temperature



Block Size	Size	Product #	Fee
Ilium Block - Small	15 mm x 20 mm	#ILIUM-PRES-SMALL	627.00
Ilium Block - Medium	15 mm x 30 mm	#ILIUM-PRES-MEDIUM	730.00
Ilium Block - Large	20 mm x 40 mm	#ILIUM-PRES-LARGE	935.00

Pre-Hydrated Unicortical Block

MOST POPULAR

Ready To Use Out Of The Package


- Preserved Using LifeNet Health's® Proprietary Preservon® Technology To Maintain A Fully Hydrated State At Ambient Temperature



Block Size	Size	Product #	Fee
Unicortical Block - Small	10 mm x 15 mm x 8 mm	#ILIUM-UNICORTICAL-SMALL	355.00
Unicortical Block - Large	10 mm x 25 mm x 8 mm	#ILIUM-UNICORTICAL-LARGE	411.00

Case Study - GBR Using LifeNet Health® Unicortical Block & Renovix-PLUS® Membrane

Case & Photos: Dr. Shaun Young



Large Osseous Defect



*LifeNet Health®
Unicortical Block*



*Renovix-Plus® –
Barrier Membrane*




*Implant Placed In
Regenerated Bone*



"Salvin continues to raise the bar with their regenerative products. The allogeneic blocks have elevated my abilities and given me additional options that I can depend on for ridge augmentation. The Renovix-PLUS® collagen membrane is one of the best on the market due to its versatility as a stiffer membrane, but also very moldable – it's really the best of both worlds."

–Dr. Shaun Young
Oral Surgeon, Tampa, FL

1. LifeNet Health. TR 0216 Three-Point Bend Testing of Glycerol-Treated, Freeze-Dried, and Frozen Bone



Decellularized Dermis Allograft

- Excellent For GTR/GBR & Oral Soft Tissue Correction
- Terminally Sterilized For Increased Patient Safety
- Handles Like Autogenous Tissue



**Prehydrated & Ready To Use
Right Out Of The Package**



**No Refrigeration Required –
Store At Ambient Temperature**

OrACELL® Decellularized Dermis Allograft

Size	Thickness	Product #	Fee	7 or More
1.0cm x 1.0cm	0.76 - 1.25mm	#DERMIS-1X1	101.00 EACH	88.70 EACH
1.5cm x 2.0cm	0.76 - 1.25mm	#DERMIS-1.5X2	196.00 EACH	175.50 EACH
1.0cm x 4.0cm	0.76 - 1.25mm	#DERMIS-1X4	224.00 EACH	200.40 EACH
2.0cm x 4.0cm	0.76 - 1.25mm	#DERMIS-2X4	246.00 EACH	220.80 EACH
1.0cm x 1.0cm	1.26 - 1.75mm	#DERMIS-THICK-1X1	101.00 EACH	88.70 EACH
1.5cm x 2.0cm	1.26 - 1.75mm	#DERMIS-THICK-1.5X2	196.00 EACH	175.50 EACH
1.0cm x 4.0cm	1.26 - 1.75mm	#DERMIS-THICK-1X4	224.00 EACH	200.40 EACH
2.0cm x 4.0cm	1.26 - 1.75mm	#DERMIS-THICK-2X4	246.00 EACH	220.80 EACH



OrACELL® Decellularized Dermis Case Studies

Case Study - OrACELL® Decellularized Dermis

Case & Photos: Dr. Matt Ficca



*Pre-Op
Tooth #6 & #7*



6 Year Follow-Up



*Pre-Op
Tooth #10 & #11*



6 Year Follow-Up



"OrACELL® provides an excellent alternative to autogenous tissue for tissue graft procedures. The pre-hydrated samples handle comparative to autogenous or even better. Over the years I have looked for an allograft with consistent quality and thickness, ease of handling, and predictable healing. I want to offer my patients the highest quality allograft on the market. I am confident I've got it covered with OrACELL®."

—Dr. Matt Ficca
Periodontist, Charlotte, NC

Case Study - OrACELL® Decellularized Dermis

Case & Photos: Dr. Stephen C. Wallace



Pre-Op



OrACELL® Dermis Placed



Suture Closing



3 Month Post-Op



"Based on the results recorded at 3 months, there was significant gain in root coverage, gain in clinical attachment and increase in width of keratinized gingiva for all recession sites treated with OrACELL® Decellularized Dermis Allograft. I have been using OrACELL® for decades with predictable clinical results."

—Dr. Stephen C. Wallace
Periodontist, Wilmington, NC

SalvinOss® Mineralized Cancellous Xenograft

- Indicated For Sinus Lifts, Ridge Augmentation & Guided Bone Regeneration
- Maintains Space & Facilitates Osteoconduction
- 0.25 - 1.0mm High Porosity Particles Become An Integral Part Of Newly Formed Bone
- Radiopaque



Exceptional Handling Characteristics



SALVIN(OSS)® Xenograft Bone Graft Material



Particle Size	Size	Product #	Price	7 or More	14 or More
250-1000 microns	0.50 CC (0.25G)	#SALVINOSS-0.50	102.00	90.60	78.00
	1.00 CC (0.5G)	#SALVINOSS-1.00	153.00	137.40	122.70
	2.50 CC (1.25G)	#SALVINOSS-2.50	280.00	254.50	229.60


Download Our White Paper At:
www.salvin.com/Downloads/Guides-Studies.html



SalvinOss® Case Studies

Case Study - Socket Preservation With SalvinOss®

Case & Photos: Dr. Michele Holzinger



Pre-Op



*Post-Extraction –
4 Months Healing*



*Post-Implant –
3 Months Healing*



"I trust Salvin Regenerative because they help me produce results that are consistent and reliable for my patients. When I re-enter a site grafted with SalvinOss®, I know exactly what to expect because the results are predictable."

–Dr. Michele Holzinger
Periodontist, Middletown, CT

Case Study - Ridge Augmentation Using SalvinOss® Xenograft Bone Material

Case & Photos: Dr. Mitchel Godat



Severely Pneumatized Sinus Cavity



12 Month Post-Op Radiograph



"I have been using SalvinOss® for many years to increase bone volume for space preservation prior to implant placement. SalvinOss® is excellent for sinus lifts, ridge augmentation, and other GBR / GTR procedures with very predictable long-term results."

–Dr. Mitchel Godat
Periodontist, Memphis, TN

Case Study - Guided Tissue Regeneration Using SalvinOss® Xenograft Bone Material

Case & Photos: Dr. Grant King



*9-10 mm
Periodontal Defect*



*8 Month Post-Op –
Significant Pocket Fill*




"SalvinOss® handles extremely well, especially when I mix it with blood, making it easy for me use for GTR and GBR. It is a slower resorbing material, holding space and volume more than using allograft alone. Most importantly, SalvinOss® gives me predictable results every time."

–Dr. Grant King
Periodontist, Memphis, TN


SalvinOss® Collagen Xenograft Bone + Bovine Collagen

- 90% SalvinOss® Xenograft Bone Graft Material PLUS 10% Bovine Collagen
- For Socket Preservation, Sinus Lifts, And Periodontal / Furcation Defects
- Plug Shape Helps Make Socket Preservation Fast & Simple
- Retains Shape After Hydrating



Ideal For:

- Anterior Socket Grafting
- Periodontal Defects
- Osseous Defects Caused By Peri-Implantitis




Ideal For:


- Posterior Socket Grafting
- Ridge Augmentation
- Sinus Lifts

Plugs Shown Actual Size

Ridge Height



Ridge Width



■ SalvinOss® Collagen ■ Xenograft B

SalvinOss® Collagen Demonstrated Greater Ridge Height & Width Compared To Xenograft B Throughout Our 24-Week Study (Data On File)

SALVIN(OSS) Collagen Xenograft Bone PLUS Bovine Collagen



Particle Size	Size	Product #	Price	7 or More	14 or More
250-1000 microns	7mm (Diam) x 10mm	#SALVINOSS-COLL-100MG	121.00	112.50	104.50
	10mm (Diam) x 10mm	#SALVINOSS-COLL-250MG	212.00	202.00	190.00

Download Our White Paper At:
www.salvin.com/Downloads/Guides-Studies.html



SalvinOss® Collagen Case Studies

Case Study #1 - Socket Preservation Using SalvinOss® Collagen

Case & Photos: Dr. Tyler Borg



Plug Shape For Easy Socket Preservation And Crestal Sinus Lift

Case Study #2 - Furcation Defect Using SalvinOss® Collagen

Case & Photos: Dr. Tyler Borg



SalvinOss® Collagen Helps When Graft Containment Is Problematic



"SalvinOss® Collagen has become a staple of my practice for regeneration as well as augmentation. The ease of handling coupled with the volumetric stability produces predictability in grafting situations around both teeth and implants alike."


—Dr. Tyler Borg

Periodontist, Northglenn, CO



Bovine Xenograft Bone Graft Material

- 0.25 - 1.0mm Bovine Granules Facilitate Osteoconduction And New Bone Formation
- Hydrophilic Properties For Rapid Absorption Of Patient's Blood
- Micro & Macro Pores Create Environment For Cell Attachment And New Bone Growth
- Can Be Mixed With Allograft Or Autologous Bone If Desired
- Indicated For Ridge Augmentation, Filling Periodontal Defects, Socket & Ridge Preservation, Sinus Lifts



Excellent Handling Characteristics When Hydrated



Filter On End Of Syringe Keeps Material In Place



Draw Saline Into Syringe For Quick & Easy Delivery

SALVIN OSS-B Bovine Xenograft Bone Graft Material – Vial



Particle Size	Size	Product #	Price	7 or More	14 or More
250-1000 microns	0.54 CC (0.25G)	#SALVINOSS-B-0.50	102.00	90.60	78.00
	1.08 CC (0.5G)	#SALVINOSS-B-1.00	153.00	137.40	122.70
	2.16 CC (1.0G)	#SALVINOSS-B-2.00	264.00	231.00	208.00
	4.32 CC (2.0G)	#SALVINOSS-B-4.00	422.00	400.00	368.00

SALVIN OSS-B Bovine Xenograft Bone Graft Material – Syringe



Particle Size	Size	Product #	Price	7 or More	14 or More
250-1000 microns	0.25 CC	#SALVINOSS-B-SYR-0.25	74.00	68.00	60.50
	0.50 CC	#SALVINOSS-B-SYR-0.50	105.00	93.00	81.50
	1.00 CC	#SALVINOSS-B-SYR-1.00	158.00	137.40	122.70

Tricalcium Phosphate + Hydroxyapatite

- Promotes Osteoblast Proliferation For New Host Bone Formation
- Alternative To Allografts Or Xenografts For GBR Procedures
- Slow Resorbing, Excellent For Space Maintenance



SALVIN Orastruct 80/20 TCP/HA Synthetic Graft Material

Size	Particulate Size	Product #	Price	7 or More	14 or More
0.5cc	0.5 - 1mm	#ORASTRUCT-0.5CC	93.00	79.70	69.50
1.0cc	0.5 - 1mm	#ORASTRUCT-1.0CC	133.00	115.80	103.40
2.5cc	1 - 2mm	#ORASTRUCT-2.5CC	224.00	203.50	187.60
0.5cc In Syringe - Pkg of 5	0.5 - 1mm	#ORASTRUCT-SYR-0.5CC	329.00	299.70	279.50




Salvin Mesh-Ti™ Titanium Mesh



For Large Defects


- 0.2mm Thickness
- 1.48mm (1,480 Micron) Pore Size
- Larger Pore Size For Enhanced Growth Factor Cross Over
- Secure With Salvin Fixation Screws
- Packaged Non-Sterile



Shown Actual Size

50mm x 37mm

#SALVIN-TI-MESH-50X37 121.00 EACH



Shown Actual Size

25mm x 37mm

#SALVIN-TI-MESH-25X37 77.20 EACH

Mesh & Wire Cutter

#MESH-WIRE-CUTTER
..... 207.20 EACH





Resorbable Collagen Tape

- Use To Repair Sinus Membrane Perforations & Stabilize Blood Clots In Minor Oral Wounds
- Resorbs in 10-14 Days

Tape Shown Actual Size - Package Contains 10 Dressings




Resorbs In 10 To 14 Days

Size	Product #	Price	7 or More	14 or More
		270.00	244.70	228.50
25mm x 75mm	#ORATAPE-PKG-10	BOX OF 10	Box of 10	Box of 10



Resorbable Collagen Wound Dressings

- For Palatal Donor Sites & To Stabilize Blood Clots In Minor Oral Wounds
- 3X Thicker Than OraTAPE®
- Controls Bleeding & Stabilizes Blood Clots
- Resorbs in 10-14 Days



Cote Shown Actual Size - Package Contains 10 Dressings



Resorbs In 10 To 14 Days

Size	Product #	Price	7 or More	14 or More
		213.00	192.70	173.30
20mm x 40mm	#ORACOTE-PKG-10	BOX OF 10	Box of 10	Box of 10



Resorbable Bovine Type 1 Collagen Plug

- Use To Stabilize Blood Clots & Preserve Graft Material In Socket Preservation
- Controls Bleeding In Extraction & Biopsy Sites
- Resorbs in 10-14 Days

Sponges Shown Actual Size - Package Contains 10 Sponges



Resorbs In 10 To 14 Days




Size	Product #	Price	7 or More	14 or More
10mm x 20mm	#ORAPLUG-PKG-10	152.00 BOX OF 10	136.50 Box of 10	125.60 Box of 10




80% Calcium Phosphate Granules Combined With 20% Bovine Collagen

- Simple, Cost Effective Socket Preservation
- Mineral Mimics Natural Bone In Resorption & Remodeling Profile
- Plug Takes The Shape Of Defect When Hydrated To Fill Socket
- Offered In Two Sizes - Five Plugs Per Box

Composite Plugs Shown Actual Size - Package Contains 5 Composite Plugs
10mm x 20mm



Size	Product #	Price	7 or More	14 or More
10mm x 20mm	#FULFILL-10X20-PKG-5	295.00 BOX OF 5	276.00 Box of 5	260.30 Box of 5
6mm x 25mm	#FULFILL-6X25-PKG-5	295.00 BOX OF 5	276.00 Box of 5	260.30 Box of 5



A Next Generation Demineralized Bone Matrix (DBM)


- Exceptional Handling – Moldable Putty Conforms To The Graft Site
- Each Lot Tested To Ensure Osteoinductive Potential
- Thickens At Body Temperature, For Better Graft Containment
- Packaged Sterile In A Syringe

Evaluation of Osteoinductive Potential

An In Vitro Study Was Conducted To Examine Accessibility To Bioactive Proteins In ABM And DBM¹. The Content Of BMP-2 Was Measured Over Time Using An Enzyme Linked Immunosorbent Assay (ELISA).

Study Results

- **Soluble BMP-2 Was Detected In ABM At An Earlier Time Compared To Particulate DBM**
- **ABM's Open Pore Structure Provides Earlier Accessibility To Bone Proteins**



Accell Connexus



Size	Product #	Fee
0.5 cc	#ACCELL-CONN-0.5CC	350.00
1.0 cc	#ACCELL-CONN-1.0CC	533.00
2.5 cc	#ACCELL-CONN-2.5CC	941.00

No Refrigeration Required - Store At Room Temperature
Typical Shelf Life: 6-9 Months

1. Khaliq S, Lollis R, Bell D, Oliver R, Walsh WR, and Ingram R, Evaluation of a Next Generation DBM Putty in a Posterolateral Spinal Fusion Model, (2009) Integra LifeSciences Corporation.



Demineralized Cortical Bone With Mineralized Cancellous Chips Paste Or Putty

- For GBR, Ridge Augmentation, Sinus Lifts & Bony Defects
- Thickens At Body Temp To Keep Graft In Place
- Syringe Option For Easy Placement



DynaBlast Paste



Size	Product #	Fee
0.5 cc	#ORTHO-PASTE-0.5CC	154.00
1.0 cc	#ORTHO-PASTE-1.0CC	214.00
3.0 cc	#ORTHO-PASTE-3.0CC	548.00

DynaBlast Putty



Size	Product #	Fee
1.0 cc	#ORTHO-PUTTY-1CC	262.00
2.5 cc	#ORTHO-PUTTY-2.5CC	524.00
5.0 cc	#ORTHO-PUTTY-5.0CC	926.00

No Refrigeration Required - Store At Room Temperature
Typical Shelf Life: 6-9 Months

Demineralized Cortical Bone Paste Or Putty

- For Osseous Defects
- Thickens At Body Temp To Keep Graft In Place
- Syringe Option For Easy Placement



DynaGraft-D Gel



Size	Product #	Fee
0.5 cc	#DYNA-GEL-0.5CC	154.00
1.0 cc	#DYNA-GEL-1.0CC	214.00
5.0 cc	#DYNA-GEL-5.0CC	871.00

DynaGraft-D Putty



Size	Product #	Fee
1.0 cc	#DYNA-PUTTY-1.0CC	262.00
2.5 cc	#DYNA-PUTTY-2.5CC	524.00
5.0 cc	#DYNA-PUTTY-5.0CC	926.00

No Refrigeration Required - Store At Room Temperature
Typical Shelf Life: 6-9 Months



Resorbs In 6-9 Months For Long-Term Barrier Functionability

- Collagen-Fiber Matrix With Intact Intrinsic Crosslinks For Mechanical Strength
- <0.01 Micron Pore Size Allows Nutrients & Growth Factors While Assuring Cell Barrier Function
- Chemically Engineered To Control Stability & Minimize Inflammatory Response
- Significant Tensile Strength For Ease Of Tacking & Suturing



Long Term Membrane


Resorbs In 6 To 9 Months

Size	Product #	Price	7 or More	14 or More
15mm x 20mm	#AMPLIFI-15X20	180.00	159.70	146.60
20mm x 30mm	#AMPLIFI-20X30	224.00	199.40	184.20
30mm x 40mm	#AMPLIFI-30X40	326.00	297.80	276.70

Membranes Shown
Actual Size



15mm x 20mm



20mm x 30mm



30mm x 40mm



Case Study - Amplifi™ Long-Term Membrane

Case & Photos: Dr. Bradley Ross



Surgical Extraction



*Bone Graft Material
In Defect*



*Amplifi™ Sutured
In Place*



*New Bone Formation –
Ready For Implant*



"The Amplifi™ membrane from Salvin has become a go to biomaterial for me for ridge augmentation procedures. It is extremely easy to manipulate and use, is very resilient, and most importantly it has a long resorption time. The results I have seen with the Amplifi™ membrane have been spectacular."

—Dr. Bradley Ross
Periodontist, Miami, FL

*Data On File

Toll Free: 800-535-6566 [U.S. & Canada]

© 2026 Salvin Dental Specialties, LLC. All Rights Reserved.



Our Most Versatile Resorbable Barrier Membrane

- For Bone Regeneration & Soft Tissue Esthetics
- Non-Cross Linked Extracellular Matrix Perfectly Conforms To The Defect When Hydrated
- Contains Type I, III, & IV Collagen
- Excellent Tensile Strength & Tear Resistance
- Trim & Place Either Wet Or Dry
- Resorbs Within 6 Months

SALVIN Renovix-Plus® Guided Healing Collagen Membrane

Size	Product #	Price	7 or More	14 or More
15mm x 25mm	#RENOVIX-PLUS-15X25	188.00	168.40	155.60
20mm x 30mm	#RENOVIX-PLUS-20X30	245.00	222.50	206.40
30mm x 40mm	#RENOVIX-PLUS-30X40	319.00	284.70	269.50



Case Study - Renovix-Plus® Guided Healing Collagen

Case & Photos: Dr. Lauralee Nygaard



*Aesthetic Zone
Surgical Extraction*



*Renovix-Plus®
Over Graft Material*



*Soft Tissue Aesthetics
Preserved*



Implant Placed



"Renovix-Plus® has become the membrane I use most for predictable ridge preservation. I use it for ease of handling and preservation of soft tissue architecture. And the bonus of working with team Salvin!"

—Dr. Lauralee Nygaard

Periodontist, Spokane Valley, WA




- *Handles Similarly To Both OraMEM® And OraMEM-Sustained®, With A Resorption Profile That Is In Between.*
- *By Choosing OraMEM®+, You May Be Able To Reduce The Amount Of Inventory You Need On The Shelf.*

*Still Want Something That Resorbs A Little Longer?
Take A Look At Our Amplifi Membrane On Page 30*




Resorbable Collagen Membrane

- For Use In Guided Tissue Regeneration Procedures
- Excellent Handling Properties
- Membrane Is Moderately Drapable
- Soft Surface Texture
- Provides Wound Stabilization
- Highly Purified Collagen That Is Biocompatible
- Resorbed In 4 To 6 Months




Scanning Electron Micrograph Of
The Dry OraMEM®+ (Cross Section)
At 50x Magnification


Membranes Shown
Actual Size



15mm x 20mm



20mm x 30mm



30mm x 40mm

OraMEM+

Resorbs In 4 To 6 Months

Size	Product #	Price	7 or More	14 or More
15mm x 20mm	#ORAMEM-PLUS-15X20	175.00	156.60	144.40
20mm x 30mm	#ORAMEM-PLUS-20X30	221.00	198.30	184.00
30mm x 40mm	#ORAMEM-PLUS-30X40	319.00	290.70	270.50






Membranes Shown
Actual Size



15mm x 25mm
25mm Dimension
For Effective
Buccal - Lingual Coverage




20mm x 30mm



30mm x 40mm

Our Most Drapeable, “No Shape Memory” Barrier Membrane

- Porcine Pericardium Maintains Barrier Function Up To 33 Weeks
- Significant Tensile Strength - Suture Or Tack With Confidence



SALVIN **Xymphony** Contour Adapting Resorbable Membrane

Size	Product #	Price	7 or More	14 or More
15mm x 25mm	#XYMPHONY-15X25	169.00	149.40	136.70
20mm x 30mm	#XYMPHONY-20X30	224.00	199.70	184.40
30mm x 40mm	#XYMPHONY-30X40	305.00	273.70	254.80

1. In vivo animal study data on file.
2. Data on file.





Membranes Shown
Actual Size



12mm x 24mm



25mm x 30mm

Non-Resorbable PTFE Membrane

- For Socket Preservation When Primary Closure Cannot Be Achieved
- Sufficiently Firm To Protect The Defect – Supple Enough To Drape Over Surrounding Tissue
- Nominal Thickness Of 250 Microns
- Needs To Be Removed In 2-3 Weeks



Space-Making Barrier In The Treatment Of Periodontal Defects
Textured Top Surface Provides Scaffold For Soft Tissue
Smooth Bottom Surface Against Bone For Easy Retrieval
Pore Size Less Than 0.3 Microns

SALVIN®
CytoSurg®

Non-Resorbable PTFE Membrane






1 Membrane Per Box

Size	Product #	Price	10 or More	20 or More
12mm x 24mm	#CYTOSURG-12X24	54.00/Membrane	45.50/Membrane	39.20/Membrane
25mm x 30mm	#CYTOSURG-25X30	93.00/Membrane	80.90/Membrane	70.20/Membrane






NEEDLE SELECTION GUIDE


	 Reverse Cutting Triangular point with cutting edge on the outer curvature.
	 Premium Reverse Cutting Same as Reverse Cutting, but precision-honed for added sharpness.
	 Taper Point Round shaft, straight or curved, tapered point with no cutting edge. Special tapered point needle for connective tissue grafts. Ideal for delicate tissue.
	 Taper Cutting Four (4) cutting edges with a taper shaft. Useful in microsurgery for excellent penetration through tough skin.



The Suture Book By Dr. Lee Silverstein



Beautifully Illustrated Step-By-Step Instructions



**The Suture Book: The Definitive Guide To Dental Suturing And Surgical Flap Closure
By Dr. Lee Silverstein**

- Detailed, Illustrated Step-By-Step Instructions For Dental Suturing
- Characteristics Of Suture Materials & Needle Designs
- Suture Selection Recommendations For Various Procedures

#SUTURE-BOOK **99.00 EACH**
6 OR MORE **90.00 EACH**



**Complete Selection of Sutures
On Pages 36-39**



**Needle Holders
On Pages 158-161**

SALVIN® CytoSurg®



Non-Resorbable PTFE Sutures

- *Excellent Choice For Bone Grafting & Implant Procedures*
- *Smooth & Flexible – Not Fragile Or Brittle – So Knots Stay Tied*
- *Proprietary Monofilament Is Biologically Inert & Non-Cytotoxic*
- *Longer Length – All PTFE Sutures Now 24" Long*



Black Needle For Higher Visibility & Less Glare


Super Sharp Needle Requires Less Force – Easier Suturing With Less Trauma

*"Racetrack" Rather Than
Zig-Zag Package*




*Little-To-No Package Memory
For Excellent Handling
And Secure Knotting*


Ultra-Soft, Ultra-Smooth Suture Designed With Low Surface Porosity



SALVIN **CytoSurg**



Competitor G
(2000x Magnification)



Competitor O

	Salvin Item#	Size	Needle	Length	Per Box	10 Or More
	#PTFE-20FS2	2-0	19mm 3/8 Circle - Reverse Cutting	24" (60cm)	120.60 Box 12	106.90 Box 12
MOST POPULAR	#PTFE-30FS2	3-0	19mm 3/8 Circle - Reverse Cutting	24" (60cm)	120.60 Box 12	106.90 Box 12
	#PTFE-30FS3	3-0	16mm 3/8 Circle - Reverse Cutting	24" (60cm)	120.60 Box 12	106.90 Box 12
	#PTFE-40TF	4-0	13mm 1/2 Circle - <u>Taper Point</u>	24" (60cm)	120.60 Box 12	106.90 Box 12
MOST POPULAR	#PTFE-40FS3	4-0	16mm 3/8 Circle - Reverse Cutting	24" (60cm)	120.60 Box 12	106.90 Box 12
	#PTFE-40PS2	4-0	19mm 3/8 Circle - <u>Premium</u> Reverse Cutting	24" (60cm)	142.20 Box 12	126.50 Box 12
	#PTFE-40P3	4-0	13mm 3/8 Circle - <u>Premium</u> Reverse Cutting	24" (60cm)	142.20 Box 12	126.50 Box 12
	#PTFE-50P3	5-0	13mm 3/8 Circle - <u>Premium</u> Reverse Cutting	24" (60cm)	142.20 Box 12	126.50 Box 12
	#PTFE-60P3	6-0	13mm 3/8 Circle - <u>Premium</u> Reverse Cutting	24" (60cm)	142.20 Box 12	126.50 Box 12



SALVIN[®] SurgiPoint[®] PGA Polyglycolic Acid Synthetic Absorbable-Braided



- Synthetic Absorbable, Braided Coated Polyglycolic Acid – Similar To Vicryl[™] But Greater Tensile Strength
- Unique PolyGlide Coating For Smooth No-Drag Suture Passage Through Tissue

Resorbs In 21-28 Days

	Salvin Item#	Size	Needle	Length	Color	Vicryl [™] Comparison	Per Box	10 Or More
	#LV-423B	3-0	3/8 Circle - Reverse Cutting	30" (76cm)	Undyed	J423G	91.50 Box 12	84.70 Box 12
	#LV-422B	4-0	3/8 Circle - Reverse Cutting	30" (76cm)	Undyed	J422G	91.50 Box 12	84.70 Box 12
	#LV-MC397NC	4-0	3/8 Circle - Reverse Cutting	30" (76cm)	Violet		134.60 Box 12	125.10 Box 12
	#LV-421B	5-0	3/8 Circle - Reverse Cutting	30" (76cm)	Undyed	J421G	91.50 Box 12	84.70 Box 12
	#LV-500B	5-0	3/8 Circle - Reverse Cutting	18" (45cm)	White	J500G	87.90 Box 12	82.10 Box 12
	#LV-493B	5-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed	J493G	157.50 Box 12	147.00 Box 12
	#LV-492B	6-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed	J492G	157.50 Box 12	147.00 Box 12
	#LV-497B	3-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed	J497G	157.50 Box 12	147.00 Box 12
	#LV-496B	4-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed	J496G	157.50 Box 12	147.00 Box 12
	#LV-495B	5-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed	J495G	157.50 Box 12	147.00 Box 12
	#LV-311B	3-0	1/2 Circle - <u>Taper Point</u>	30" (76cm)	Violet	J311G	82.20 Box 12	74.90 Box 12
	#LV-310B	4-0	1/2 Circle - <u>Taper Point</u>	30" (76cm)	Violet	J310G	82.20 Box 12	74.90 Box 12

SALVIN[®] SurgiPoint[®] PGA Polyglycolic Acid Synthetic Absorbable-Braided Fast Absorbing



- Good Alternative To Plain Gut Sutures
- Synthetic Absorbable, Braided Coated Polyglycolic Acid – Similar To Vicryl Rapide[™] But Greater Tensile Strength
- Unique PolyGlide Coating For Smooth No-Drag Suture Passage Through Tissue

Resorbs In 7-14 Days

	Salvin Item#	Size	Needle	Length	Color	Vicryl [™] Comparison	Per Box	10 Or More
	#LVFA-493N	3-0	3/8 Circle - Reverse Cutting	18" (45cm)	Undyed	VR497	102.10 Box 12	95.10 Box 12
	#LVFA-423N	3-0	3/8 Circle - Reverse Cutting	30" (76cm)	Undyed		102.10 Box 12	95.10 Box 12
	#LVFA-492N	4-0	3/8 Circle - Reverse Cutting	18" (45cm)	Undyed	VC422H	102.10 Box 12	95.10 Box 12
	#LVFA-422N	4-0	3/8 Circle - Reverse Cutting	30" (76cm)	Undyed	VR426	102.10 Box 12	95.10 Box 12
	#LVFA-442N	3-0	3/8 Circle - Reverse Cutting	30" (76cm)	Undyed	VR935	102.10 Box 12	95.10 Box 12
	#LVFA-497N	3-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed		102.10 Box 12	95.10 Box 12
	#LVFA-496N	4-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed		102.10 Box 12	95.10 Box 12
	#LVFA-490C	5-0	3/8 Circle - <u>Prem</u> Reverse Cut	18" (45cm)	Undyed		132.60 Box 12	121.40 Box 12

More Sutures On Next 2 Pages →

SALVIN® SurgiPoint® Chromic Gut *Resorbs In 7-14 Days*

- Natural Resorbable Gut – Excellent Suture Except Where There Is Significant Tension Across The Wound

	Salvin Item#	Size	Needle	Length	Ethicon™ Comparison	Per Box	10 Or More
 FS-2 (C6) <i>MOST POPULAR</i>	#L-560B	3-0	3/8 Circle - Reverse Cutting	27" (68cm)	636G	61.90 Box 12	55.90 Box 12
	#L-559B	4-0	3/8 Circle - Reverse Cutting	27" (68cm)	635G	61.90 Box 12	55.90 Box 12
	#L-525B	3-0	3/8 Circle - Reverse Cutting	18" (45cm)		61.90 Box 12	55.90 Box 12
	#L-558B	4-0	3/8 Circle - Reverse Cutting	18" (45cm)		61.90 Box 12	55.90 Box 12
	#L-557B	5-0	3/8 Circle - Reverse Cutting	27" (68cm)		61.90 Box 12	55.90 Box 12
 P-3 (C3) <i>MOST POPULAR</i>	#L-595B	4-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	1654G	117.60 Box 12	109.60 Box 12
	#L-1248B	5-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	687G	117.60 Box 12	109.60 Box 12
 P-1 (C1)	#1242B-CP	6-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)		191.30 Box 12	178.20 Box 12
 SH-1 (T4)	#L-589B	4-0	1/2 Circle - <u>Taper Point</u> For Connective Tissue Grafts	14" (35cm)		75.40 Box 12	68.50 Box 12
 RB-1 (T28)	#L-590B	5-0	1/2 Circle - <u>Taper Point</u> For Connective Tissue Grafts	14" (35cm)		58.90 Box 12	53.90 Box 12

SALVIN® SurgiPoint® Gut-Plain *Resorbs In 5-10 Days*

- Natural Resorbable Gut – Excellent Suture Except Where There Is Significant Tension Across The Wound

	Salvin Item#	Size	Needle	Length	Ethicon™ Comparison	Per Box	10 Or More
 FS-2 (C6) <i>MOST POPULAR</i>	#L-554B	3-0	3/8 Circle - Reverse Cutting	27" (68cm)	H822G	61.90 Box 12	55.90 Box 12
	#L-551B	4-0	3/8 Circle - Reverse Cutting	27" (68cm)	H822G	61.90 Box 12	55.90 Box 12
	#L-591B	4-0	3/8 Circle - Reverse Cutting	18" (45cm)		61.90 Box 12	55.90 Box 12
	#L-592B	3-0	3/8 Circle - Reverse Cutting	18" (45cm)		61.90 Box 12	55.90 Box 12
 P-3 (C3)	#L-596B	5-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	686G	117.60 Box 12	109.60 Box 12
 X-1 (C31)	#L-594B	3-0	1/2 Circle - Reverse Cutting	18" (45cm)	612G	61.90 Box 12	55.90 Box 12
 PC-1 (CP1)	#L-1415B	6-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	1916G	196.30 Box 12	183.10 Box 12

SALVIN® SurgiPoint® Monoderm™ *Resorbs In 21-28 Days*

- Good Alternative To Chromic Gut Sutures
- Synthetic Resorbable Monofilament Glycolide & E-Caprolactone – Similar to Monocryl™

	Salvin Item#	Size	Needle	Length	Ethicon™ Comparison	Per Box	10 Or More
 FS-2 (C6)	#L-388	3-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)		96.60 Box 12	86.50 Box 12
	#L-392C	4-0	3/8 Circle - Black Premium Reverse	18" (45cm)		96.60 Box 12	86.50 Box 12
	#L-Q919	5-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)		96.60 Box 12	86.50 Box 12
 FS-3 (C22)	#L-395	4-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)		96.60 Box 12	86.50 Box 12
 PS-2 (DGL19) <i>MOST POPULAR</i>	L-YD513N	4-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	Y513G	155.70 Box 12	145.70 Box 12
 P-3 (C3)	#L-464B	4-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)		162.60 Box 12	152.40 Box 12
	#L-463B	5-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)		162.60 Box 12	152.40 Box 12
	#L-463C	5-0	3/8 Circle - Black Premium Reverse	18" (45cm)	Y463G	162.60 Box 12	152.40 Box 12
	#L-492B	6-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)		162.60 Box 12	152.40 Box 12
	#L-401B	6-0	3/8 Circle - Violet Premium Reverse	18" (45cm)		162.60 Box 12	152.40 Box 12

Suture Names Are Trademarks Of Their Respective Companies.

SALVIN® SurgiPoint® Silk **Non-Resorbable**



- Non-Resorbable, Braided, Coated Silk – Excellent Handling & Tying Characteristics

	Salvin Item#	Size	Needle	Length	Ethicon™ Comparison	Per Box	10 Or More
FS (C8)	#L-770B	2-0	3/8 Circle - Reverse Cutting	18" (45cm)	685-G	59.90 Box 12	54.70 Box 12
FS-2 (C6) MOST POPULAR	#L-784B	3-0	3/8 Circle - Reverse Cutting	18" (45cm)		59.90 Box 12	54.70 Box 12
	#L-781B	4-0	3/8 Circle - Reverse Cutting	18" (45cm)	683G	59.90 Box 12	54.70 Box 12
	#L-754B	5-0	3/8 Circle - Reverse Cutting	18" (45cm)	682G	59.90 Box 12	54.70 Box 12
FS-1 (C7)	#L-786B	3-0	3/8 Circle - Reverse Cutting	18" (45cm)	684G	59.90 Box 12	54.70 Box 12
C-2 (C17)	#L-787B	4-0	3/8 Circle - Reverse Cutting	18" (45cm)	736G	59.90 Box 12	54.70 Box 12
X-1 (C31)	#L-785B	3-0	1/2 Circle - Reverse Cutting	18" (45cm)	632G	59.90 Box 12	54.70 Box 12

SALVIN® SurgiPoint® Polypropylene **Non-Resorbable**



- Non-Resorbable Blue Monofilament – Inert, Easy To Handle, Passes Easily Through Tissue

	Salvin Item#	Size	Needle	Length	Ethicon™ Comparison	Per Box	10 Or More
P-3 (C3) MOST POPULAR	#L-1085B	5-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	8698G	174.90 Box 12	161.60 Box 12
	#L-1082B	6-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	8695G	174.90 Box 12	161.60 Box 12
P-1 (C1)	#L-1079B	6-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	8697G	174.90 Box 12	161.60 Box 12
	#L-2468N	7-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	8696G	174.90 Box 12	161.60 Box 12
C-2 (C17)	#L-1011B	6-0	3/8 Circle - Reverse Cutting	10" (25cm)	8889H	174.90 Box 12	161.60 Box 12
FS-2 (C6)	#L-8661B	5-0	3/8 Circle - Reverse Cutting	18" (45cm)	8661G	89.90 Box 12	83.90 Box 12
FS-3 (C22)	#L-8680B	5-0	3/8 Circle - Reverse Cutting	18" (45cm)	8680G	89.90 Box 12	83.90 Box 12

SALVIN® SurgiPoint® Nylon Monofilament **Non-Resorbable**



- Non-Resorbable Black Monofilament – Passes Easily Through Tissue – Similar To Resolon™

	Salvin Item#	Size	Needle	Length	Ethicon™ Comparison	Per Box	10 Or More
FS-2 (C6) MOST POPULAR	#L-922B	4-0	3/8 Circle - Reverse Cutting	18" (45cm)	662G	70.90 Box 12	64.60 Box 12
	#L-918B	5-0	3/8 Circle - Reverse Cutting	18" (45cm)	661G	70.90 Box 12	64.60 Box 12
P-3 (C3)	#L-698BKC	5-0	3/8 Circle - Black <u>Premium</u> Reverse	18" (45cm)	698G	125.20 Box 12	117.70 Box 12
	#L-1698B	6-0	3/8 Circle - <u>Premium</u> Reverse Cutting	18" (45cm)	1698G	125.20 Box 12	117.70 Box 12

SALVIN® SurgiPoint® Polyviolene-Green Braided Polyester **Non-Resorbable**



- Green Braided Uncoated Polyester – Smooth, Easy To Handle With Good Knot Security

	Salvin Item#	Size	Needle	Length	Ethicon™ Comparison	Per Box	10 Or More
FS-2 (C6) MOST POPULAR	#L-1180B	4-0	3/8 Circle - Reverse Cutting	18" (45cm)	R683G	59.90 Box 12	54.70 Box 12

PeriAcryl® by GluStitch

Purified Cyanoacrylate Dental Adhesive

- Will Adhere To Moist Living Tissue Without Toxic Or Foreign Body Reaction
- Deep Violet Color Provides Visibility For Controlled Placement
- One Year Guaranteed Shelf Life



Single-Use Applicators Kit

Includes: 12 0.2mL Pre-Loaded Applicators

#PERIACRYL-SINGLE 134.40 BOX OF 12
2 OR MORE BOXES 126.50 BOX OF 12



MOST POPULAR

5mL Bottle Multi-Use Kit

Includes: 5mL Bottle, 50 Pipettes & 1 Reusable Administration Tray

#PERIACRYL 130.90 EACH
2 OR MORE 119.10 EACH

High Viscosity PeriAcryl® by GluStitch

- *High Viscosity Helps Control Flow For Precision Application*
- *Nine Times Thicker Than Regular PeriAcryl For Easier Handling*



Single-Use Applicators Kit

Includes: 12 0.2mL Pre-Loaded Applicators

#PERIACRYL-HV-SINGLE 140.20 BOX OF 12
2 OR MORE BOXES 133.40 BOX OF 12




5mL Bottle Multi-Use Kit

Includes: 5mL Bottle, 50 Pipettes & 1 Reusable Administration Tray

#PERIACRYL-HV 139.20 EACH
2 OR MORE 130.50 EACH

Micropoint Suture Removal Scissors

Takes The "Wince" Out Of Suture Removal, Even In Sensitive Tissue



*Opposing Spring Tips Maintain Force
Between Blades For Incredible
Cutting Without Pinching Tissue*



MOST POPULAR

Featherlite 14cm Curved

Delicate "Castroviejo Style" Scissors

#F-LITE5CURV 345.00 EACH
3 OR MORE 326.30 EACH

Complete Selection Of Suture Scissors On Pages 155-157

HemCon® Dental DressingPro



- *Easy To Use – Adhesive Properties Help Secure Dressing To Mucosal Tissue*
- *Safe – Protects The Wound Surface With No Known Contra-Indications*



"I have used HemCon® for many years to rapidly and predictably control bleeding during surgery. I highly recommend it to my students at the Misch Resnik Implant Institute."

–Dr. Randy Resnik
Prosthodontist, Pittsburgh, PA

**Clinical Studies On File*

*Package Contains 12 Dressings
Dressings Shown Actual Size*



*10mm x 12mm
Best For Extraction Sockets*



*1" x 3" (25mm x 75mm)
Best For Palatal Wound Sites*

HemCon® Dental DressingPro 10mm x 12mm

#HEMCON-10X12MM 196.00 BOX OF 12
2 OR MORE BOXES 184.90 BOX OF 12

HemCon® Dental DressingPro 1" x 3" (25mm x75mm)

#HEMCON-1X3INCH 287.00 BOX OF 12
2 OR MORE BOXES 271.70 BOX OF 12

Portable Electrocautery System



Complete Set Includes:

- Handle & Batteries
- 2 Fine Tips & 2 Loop Tips
- 10 Sterile Handle Sheaths

#CAUT-KIT-DEL-1 135.00 SET



Replacement "Fine" Tips

Box Of 10 Single Use "Fine" Tips
With Sterile Handle Sheaths

#CAUT-FINE-H101 93.00 BOX OF 10



Replacement "Loop" Tips

Box Of 10 Single Use "Loop" Tips
With Sterile Handle Sheaths

#CAUT-LOOP-H103 93.00 BOX OF 10

OraFOAM™ Absorbable Collagen Hemostatic Sponges



Package Contains 18 Sponges

Absorbable Collagen Hemostatic Sponges

- 0.5" x 1" (1.27cm x 2.54cm)
- Controls Bleeding In 2 To 5 Minutes
- Coherent Sponge Structure Allows Easily Controlled Application

#ORAFOAM-PKG-18 340.00 BOX OF 18
 2 OR MORE BOXES 321.30 BOX OF 18

OraFOAM™ Hemostatic Sponge Is Indicated In Surgical Procedures (Other Than Ophthalmological & Urological Surgery) As An Adjunct To Hemostasis When Control Of Bleeding By Ligature Or Conventional Procedures Is Ineffective Or Impractical. Please See Product IFU For Precautions & Adverse Reactions. IFU Available Online At www.salvin.com/Downloads/Guides-Studies.html

BloodSTOP® Absorbent Hemostatic Gauze



BloodSTOP® Absorbent Hemostatic Gauze

- Made Of 100% Natural Cellulose
- Accelerates Coagulation
- Reduces Healing Time & Recurring Bleeding

BloodSTOP® 0.5" x 2"

#BLOODSTOP-IX-0.5X2 203.40 BOX OF 24
 2 OR MORE BOXES 192.70 BOX OF 24

BloodSTOP® 2" x 2"

#BLOODSTOP-IX-2X2 359.00 BOX OF 12
 2 OR MORE BOXES 342.30 BOX OF 12

Hemostats

Mosquito Hemostat
12.7cm

#MOSQ-CURV-5
 73.90 EACH

MOST POPULAR
Kelly Hemostat
Curved - 14cm

#KELLY-CURV-HEMO
 73.90 EACH

Bone Wax



- Packaged Sterile In 2.5 Gram Pouch For Easy Dispensing
- Made From Beeswax & Other Natural Materials

#BONE-WAX
 BOX OF 12 TUBES 164.70 BOX
 2 OR MORE 156.40 BOX

ORASOOTHE®

Oral Hydrogel Wound Dressing – Used To Promote Optimal Healing And Pain Management Before & After Invasive Dental Procedures



OraSoothe "SockIt" Gel

- For Extraction Sites, Bone Grafting, Implants, Scaling, Laser Treatments, Ulcers, Orthodontic Irritations And All Injuries Of The Oral Mucosa
- Unique Syringe Design Allows For Localized And Direct Placement
- Made From Only All-Natural, Nontoxic, Food-Based Ingredients
- Pleasant Fragrance And Taste, Safe To Swallow

#ORASOOTHE-GEL-PKG-5. 102.99 PKG OF 5
#ORASOOTHE-GEL-PKG-25. 299.99 PKG OF 25

OraSoothe Professional Oral Rinse - 3.4oz

- For Mucositis, Full-Mouth Extractions, Periodontal Surgeries, And Pre / Post Hygiene Procedures
- Adheres To Oral Tissue And Forms A Protective Barrier Around Oral Wounds
- Made From Only All-Natural, Nontoxic, Food-Based Ingredients
- Pleasant Fragrance And Taste, Safe To Swallow

#ORASOOTHE-RINSE-3.4OZ 11.00 EACH

YOUNG® PerioDT® WITH DESICCATION TECHNOLOGY



- *Desiccation Of Dental Plaque, Calculus And Biofilm Present In Periodontal Pockets And Surrounding Tooth Surfaces*
- *Assists With Debridement And Removal Of Plaque, Calculus And Biofilm From Periodontal Pockets And Surrounding Tooth Surfaces*
- *Use As An Adjunctive Rinse On Exposed Dentition Or Dental Prosthetics*

PerioDT Gel Biofilm Debridement Sulfonated Phenolics

Features:

- Safe On Dentin, Enamel And Cementum
- Safe On Porcelain, Titanium And Other Metals
- Easy To Use Delivery System
- One Syringe Can Complete A Full Arch
- Red Color To Show Where It Has Been Placed
- Fast Acting – Desiccation Effect Begins Immediately Upon Contact
- Non-Systemic, Non-Antibiotic

1mL – 2 Pack
#100277 62.29 PKG OF 2

0.5mL – 10 Pack
#100278 152.99 PKG OF 10



**BUYER'S GUIDE -
FOR MONITORS ON THIS PAGE**

▶	Pulse Rate
▶	O ₂ Saturation
▶	Auto Blood Pressure

Mindray VS8
MOST POPULAR
Latexfree
PRESSURE CUFF

- 8.4" Touchscreen Color LCD Display
- Nellcor OxiMax SpO₂
- Built-In Printer
- Three Year Warranty

#MINDRAY-VS8
Call For Lowest Price


Built-In Printer

Roll Stand For Mindray VS8

Quick Release Mount & Cable Management

#MINDRAY-ROLL-STAND-VS-EPM

 **CALL FOR PRICE**

Criticare VitalCare™ 506N3

- **ComfortCuff™** Measures On Inflation For Fast, Comfortable Readings
- Excellent Motion Artifact Rejection
- Internal Back-Up Battery

#CR-506-VITALCARE
Call For Lowest Price
Latexfree
PRESSURE CUFF

Roll Stand For eQuality™ 506 Monitor

Adjusts 37" To 47" High, 20" Diameter Base

#CR-ROLL-506 **CALL FOR PRICE**


▶	Pulse Rate
▶	O ₂ Saturation
▶	Auto Blood Pressure
▶	Electrocardiogram
▶	Expired CO ₂



Criticare nGenuity™ 8100

- **ComfortCuff™** Measures On Inflation For Fast, Comfortable Readings
- Three Lead ECG With Freeze & Cascade
- Built-In Printer

WITH CO₂
#CR-NGENUITY-CO2

WITHOUT CO₂
#CR-NGENUITY-8100

Call For Lowest Price



Roll Stand For nGenuity™ 8100 Monitor

Adjusts 37" To 47" High, 20" Diameter Base

#CR-ROLL-8100 CALL FOR PRICE



Precordial Amplified Sedation Stethoscope

- Chest Piece With Built-In Microphone Amplifies Lung Sounds While Suppressing Room Noise
- Diagnostic Quality Sound, Less Ambient Noise
- Set Includes:
 - Precordial Amplifier
 - Volume Controlled Earbuds
 - Precordial Pickup Cable With Piezo Chest Piece
 - 2 AAA Batteries
 - Manual & Case



Wireless Bluetooth® Integrated Technology

#PRE-STETH-WIRELESS 785.00 SET

▶	Pulse Rate
▶	O ₂ Saturation
▶	Auto Blood Pressure
▶	Electrocardiogram
▶	Expired CO ₂

Mindray ePM 10A **MOST POPULAR**

- 10" Multi-Gesture Touchscreen Display
- Extensive Data Storage Capability For Trend Data, Alarms, And Events
- Integrated Sidestream CO₂
- Built-In Printer
- Three Year Warranty



WITH CO₂
#MINDRAY-EPM-10A-C02



WITHOUT CO₂
#MINDRAY-EPM-10A



Roll Stand For Mindray ePM 10A Monitor

21" Diameter Base

#MINDRAY-ROLL-STAND-VS-EPM. . . . CALL FOR PRICE



Call For Lowest Price

CO₂ Monitor Accessories



CO₂ Monitor Single Cannula

- Single Use
- For CO₂ Use Only

#CANNULA-SINGLE

..... 229.70 BOX OF 50



CO₂ Monitor Divided Cannula

- Single Use
- Works With Any O₂ System

#CANNULA-DIVIDED

..... 101.90 BOX OF 25



CO₂ Monitor Sample Line

- Replace Every 20 Uses
- Connects Cannula & Water Trap

#SAMPLE-LINE

..... 117.30 BOX OF 50



Zoe Water Trap

330 Working Hours

#ZOE-WATER-TRAP 51.40 EACH



Criticare nGenuity™ Water Trap

- Replace Every 20 Uses
- Water Trap & Filter

#CR-938-FNC 577.00 BOX OF 30



Mindray Water Trap

- For Use Up To One Month

#MINDRAY-EPM10A-CO2-WATER-TRAP

..... 71.00 BOX OF 2

Zoll AED PLUS™ Automated External Defibrillator



- One-Piece, Pre-Connected Electrode With CPR Feedback
- Sensor Tells You If You Are Pressing At The Correct Rate
- Pictures, Text Displays & Voice Prompts Ensure Proper Sequence
- Auto Adjusts Shock Based On Detected Body Mass
- Works With Both Adults & Children
- Includes: Carry Case, Adult Electrodes, Batteries, Prep Kit & Video
- Memory Allows Download Of Rescue Data With Included Software
- Automatically Self-Tests Defibrillator, Battery & Pad Daily
- Uses 5-Year Lithium-123 Duracell Camera Batteries

#ZOLL-AED-PLUS

Call For Lowest Price



Zoll Replacement Electrodes

5 Year Adult Electrodes

#ZOLL-AED-ELECTRODE244.00 SET

2 Year Pado Electrodes

#ZOLL-AED-PEDO-ELECTR.163.00 SET



Nasal Cannula & Oxygen Adapter



Latexfree

**MOST POPULAR**

Soft-Touch Adult Nasal Cannula

- Curved Tips For Excellent Patient Comfort
- Crush Resistant Tubing
- Individually Packaged - 7' (2.1m) Long

#NASAL-CANNULA

BOX OF 5086.00 BOX

2 OR MORE BOXES79.20 BOX

Nitrous System Adapter For Nasal Cannula Tubing

FOR OXYGEN USE BY ITSELF ONLY

USE YOUR SCAVENGER CIRCUIT WHEN USING NITROUS

#ACC-CANNULA-ADAPT81.00 EACH

Anesthesia Instrumentation

Magill Endotracheal Forcep

12cm Shaft Length

#MAGILL 148.60 EACH

3 OR MORE 138.10 EACH



Young Forcep

17.1cm

#YOUNG-FORCEP 234.50 EACH

3 OR MORE 221.10 EACH



Replacement Tips for Young Forcep

Set Of 2

#YOUNG-FORCEP-TIPS 18.20 SET OF 2

3 OR MORE 16.20 SET OF 2

**MOST POPULAR**

Yankauer Suction Tip

29.8cm

#YANKAUER. 82.50 EACH

3 OR MORE 76.50 EACH



Disposable Sterile Yankauer Tips

- Blue Translucent Plastic Minimizes Visibility Of Fluids
- Flange Tip & Vent

#YANKAUER-DISPOSABLE

BOX OF 50 141.00 BOX

2 OR MORE 131.50 BOX



To Use These Surgical Tips With The High Volume Suction In Your Office, Use Our Tubing Adapter On Page 96

Emergency Kits

**MOST POPULAR**

Manual Single-Use Resuscitator

- Attached Ventilation Bag
- Connects To Oxygen Outlet For Nasal Cannula

#RESUS-SINGLE50.90 EACH

2 OR MORE46.80 EACH

**MOST POPULAR**

Single-Use Emergency Suction

- Use Metal Aspirator Or Plastic Yankauer
- 1000cc Fluid Capacity

#EMERGENCY-SUCTION 86.90 EACH

2 OR MORE 79.50 EACH

Scan For

Emergency Suction
Technique Guide

CO₂ Capnography Nasal Hoods

- Nasal Hood With Quick Connect To Capnography Monitor Tubing
- For Delivery & Scavenging Of Nitrous Oxide & End-Tidal CO₂ Capnography Monitoring
- Secure Female Luer Lock Sample Line Connection


CO₂ Capnography Nasal Hoods

Unscented Gray - Single Use

#ACC-CO2-HOOD-GRAY 145.00 BOX OF 12



ECG Patches & Wrist Clips



ECG Patches - Wet Gel Foam

- Snap Connectors For ECG Leads
- Wet Gel Foam Electrodes
- Recommended For Use With Zoe ECG Monitors

#CR-LEADLOCK-50

1 BOX 31.90 BOX OF 50
2 OR MORE 29.70 BOX OF 50



ECG Wrist Clips

- Snap Connectors For ECG Leads
- Left Wrist Clip Has 2 Leads & Right Wrist Clip Has 1 Lead Of A 3 Lead ECG
- Use Conductive (#SIGNA-SPRAY) Solution For Best Results

#CR 1166 139.50 PAIR



ECG Patches

- Snap Connectors For ECG Leads
- Packaged In Sets Of 3

#CR-LEADLOCK-20

1 BOX 31.90 BOX OF 20 SETS
2 OR MORE .. 29.70 BOX OF 20 SETS



Conductive Solution For ECG Wrist Clips

- Use With ECG Snap Leads
- 250mL Bottle

#SIGNA-SPRAY..... 17.70 EACH

Pocket CPR Mask & Emergency Airways



Single Use Laryngeal Mask Airways (LMA)

- Set Of 3 Allows Rapid Airway Access
- Less Invasive Than An Endotracheal Tube
- Includes #3 (Pedo To Small Adult), #4 (Adult) & #5 (Large Adult)

#LMA-SET-OF-3 77.50 SET OF 3
2 OR MORE 71.50 SET OF 3



Guedel Airways

- Airway, Bite Block & Tongue Depressor
- Coated Plastic Prevents Crimping & Protects Teeth
- 100mm, 90mm, 80mm

#GUEDEL-SET..... 15.60 SET OF 3
3 OR MORE 13.50 SET OF 3

UNIVERSAL COLOR CODE SYSTEM – IV CATHETER SIZES



Terumo SURFLO® IV Catheter

- Clear Flashback Chamber To Confirm Catheter In Vein
- Opaque Catheter



Color	Size	Product #	Price
Green	18g x 1.25"	#IV-CATH-SURFLO-SROX1832CA	97.90 BOX 50
Pink	20g x 1"	#IV-CATH-SURFLO-SROX2025CA	94.90 BOX 50
Pink	20g x 1.25"	#IV-CATH-SURFLO-SROX2032CA	97.90 BOX 50
Blue	22g x 1"	#IV-CATH-SURFLO-SROX2225CA	100.90 BOX 50
Yellow	24g x 0.75"	#IV-CATH-SURFLO-SROX2419CA	94.90 BOX 50

Terumo SurFlash® IV Catheter

- Clear Flashback Chamber To Confirm Catheter In Vein
- Clear Catheter To Show Blood Drawn



Color	Size	Product #	Price
Pink	20g x 1"	#IV-CATH-SURFLASH-SR-FF2025	102.90 BOX 50
Blue	22g x 1"	#IV-CATH-SURFLASH-SR-FF2225	102.90 BOX 50
Yellow	24g x 0.75"	#IV-CATH-SURFLASH-SR-FF2419	102.90 BOX 50

EXEL Safelet IV Catheter

- Clear Flashback Chamber To Confirm Catheter In Vein
- Clear Catheter To Show Blood Drawn



Color	Size	Product #	Price
Green	18g x 1.25"	#IV-CATH-EXEL-26735	75.90 BOX 50
Pink	20g x 1"	#IV-CATH-EXEL-26741	75.90 BOX 50
Pink	20g x 1.25"	#IV-CATH-EXEL-26742	75.90 BOX 50
Blue	22g x 1"	#IV-CATH-EXEL-26746	75.90 BOX 50
Yellow	24g x 0.75"	#IV-CATH-EXEL-26751	75.90 BOX 50


EXEL IV Extension Sets

- 30"

#IV-SET-EXEL-1730..... 169.50 CASE OF 100



BD Insyte™ Autoguard™ Shielded IV Catheter



Spring Loaded Push Button Mechanism Retracts Needle Into Safety Chamber




Color	Size	Product #	Price
Yellow	24g x 0.75 "	#IV-CATH-BD-INSYTE-381412	218.79 BOX 50
Blue	22g x 1 "	#IV-CATH-BD-INSYTE-381423	218.79 BOX 50
Pink	20g x 1 "	#IV-CATH-BD-INSYTE-381433	218.79 BOX 50
Winged - Yellow	24g x 0.75 "	#IV-CATH-BD-INSYTE-WINGED-381512	218.79 BOX 50
Winged - Blue	22g x 1 "	#IV-CATH-BD-INSYTE-WINGED-381523	218.79 BOX 50
Winged - Pink	20g x 1 "	#IV-CATH-BD-INSYTE-WINGED-381533	218.79 BOX 50

B. Braun Introcath Safety® IV Catheter



Passive Safety Shield – Automatically Covers Needle Tip Upon Needle Withdrawal

Helps Prevent Needlestick Injuries



Color	Size	Product #	Price
Straight – Yellow	24g x 0.75 "	#IV-CATH-BB-INTROCAN-4252500-02	165.00 BOX 50
Straight – Blue	22g x 1 "	#IV-CATH-BB-INTROCAN-4252519-02	165.00 BOX 50
Straight – Pink	20g x 1 "	#IV-CATH-BB-INTROCAN-4252543-02	165.00 BOX 50


B. Braun ULTRASITE® Needle Free Luer Access Device




- Quick, effortless Luer access connection
- High flow promotes rapid fluid delivery
- Microbial barrier design
- Easy to clean and disinfect
- Reduces accidental touch contamination

#IV-VALVE-BB-ULTRASITE-415110.466.40 CASE OF 100


EXEL IV Administration Sets



Description	Product #	Price
15 Drop/mL Luer Slip – Y Injection Site – 78" Tubing	#IV-SET-EXEL-29081	97.90 CASE 50
60 Drop/mL Luer Slip – Y Injection Site – 78" Tubing	#IV-SET-EXEL-29083	97.90 CASE 50



Description	Product #	Price
15 Drop/mL Luer Lock – Y Injection Site – 78" Tubing	#IV-SET-EXEL-29082	97.90 CASE 50
60 Drop/mL Luer Lock – Y Injection Site – 78" Tubing	#IV-SET-EXEL-29085	97.90 CASE 50



Description	Product #	Price
15 Drop/mL Luer Lock – 78" Tubing	#IV-SET-EXEL-27086	181.90 CASE 50
15 Drop/mL Luer Lock/Luer Slip – 78" Tubing	#IV-SET-EXEL-29086	97.90 CASE 50



Description	Product #	Price
15 Drop/mL Luer Lock/Luer Slip – 2 Y Injection Site – 105" Tubing	#IV-SET-EXEL-29088	126.90 CASE 50

Baxter IV Administration Sets



Description	Product #	Price
10 Drops/mL – 2 Injections Sites Duo Vent – 104" Tubing	#IV-SET-ACT5439	279.90 CASE 48

B. Braun SafeDay™ IV Administration Sets



Description	Product #	Price
15 Drops/mL – 106" Tubing	#IV-SET-SAFEDAY-352636	327.00 CASE 50

Syringes & Needle Capper



MOST POPULAR

Aspirating Syringe CW Type

Works With Standard 1.8cc Cartridges

#SYRINGE CW **52.00 EACH**
3 OR MORE **47.90 EACH**



Petite Aspirating Syringe

More Comfort & Control For Smaller Hands

Works With Standard 1.8cc Cartridges
#SYRINGE-PETITE. **54.40 EACH**
3 OR MORE **49.90 EACH**



Aspirating Syringe A Type

Works With Standard 1.8cc Cartridges

#SYRINGE A **48.70 EACH**
3 OR MORE **44.60 EACH**



Needle Capper

#NEEDLE CAPPER. **32.20 EACH**
3 OR MORE **29.70 EACH**



Intraligamentary Syringe

- Positive Pressure Design Creates Extra Leverage Without Fatigue
- One Complete Trigger Squeeze Delivers 0.2cc
- Uses Standard Carpules

#NTRALIG. **521.00 EACH**



Palatal Anesthesia Device

- Uses Pressure Anesthesia & Gate Control Theory To Reduce Or Eliminate Palatal Injection Pain
- Double-Ended For Right & Left Palatal Injections

#PAD. **82.70 EACH**
3 OR MORE **78.60 EACH**

Durable Steel Salvin Implant / Surgical Cart***Accommodates Implant Motors, Piezotomes & Vital Sign Monitors****All Metal Construction**41.5" Working Height For Optimal Visibility Of Equipment Screens***Salvin Implant / Surgical Cart**

- Pull Out Side Extension
- Removable Antimicrobial Plastic Top For Easy Cleaning
- No Slam, Soft Closing Drawers
- Removable Ergonomic Side Handles
- 5" Easy Rolling Wheels With Locks
- Power Strip With Six Outlets Attached For Easy & Neat Wire Management
- Foot Pedal Fits In Bottom Space For Mobility & Storage
- 20.5" Wide x 18.5" Deep x 41.5" Tall (Height Includes Wheels)

#SURG-CART-STEEL 1850.00*Spacious Drawers Accommodate Supplies & Wrapped Cassettes**Pull Out Extension Supports Instruments or Equipment**Surfaces & Instruments Shown Undraped For Clarity*

Mayo & IV Stands



MOST POPULAR

Heavy Duty Mayo Stand

- Tray Size Is 19 1/8" x 12 5/8" (48.6cm x 32cm)
- Adjusts From 33" To 53" (83cm to 135cm)
- Includes Stainless Steel Tray

#MAYO STAND 344.10 EACH



Center Post Mayo & Implant Motor Stand

- Tray Size Is 19 1/8" x 12 5/8" (48.6cm x 32cm)
- Adjusts From 29 1/2" To 50" (75cm to 127cm)
- Includes Stainless Steel Tray

#MAYO-STAND-CP..... 353.30 EACH



Heavy Duty IV Stand

- Adjusts From 51 1/4" To 94" (130cm to 239cm)
- Chrome Plated Heavy Bar Steel

#IV STAND 186.00 EACH

Sterile Mayo Stand Cover

- Individually Packaged
- Fits Regular Mayo Stands
- 24" x 24" (61 x 61cm)



#MAYO-COVER
CASE OF 24..... 324.00 CASE
2 OR MORE CASES 311.80 CASE



Latexfree

Sterile Surgical Table Cover

- Individually Packaged
- Fits Surgical Tables
- 53" x 77"
(134 x 195cm)



#SURG-TABLE-COVER
CASE OF 20..... 114.00 CASE
2 OR MORE CASES 106.40 CASE



Latexfree

Over The Patient Tables

**MOST POPULAR****"Over The Patient" Surgical Table**

- Putty Color Work Top Is 24" x 39" (61cm x 99cm)
- Height Adjusts From 31" To 46" Tall (78.75cm To 116.85cm)
- Dimension Between Legs Is 34" (86.4cm)

#SURG-TABLE-24X39 2090.00 EACH

**Stainless Steel "Over The Patient" Surgical Table**

- Extra Heavy-Duty Hospital Style Stainless Steel
- Work Top Is 24" x 48" (61cm x 122cm)
- Height Easily Adjusts With Hand Crank From 36" to 56" (91.5cm to 142cm)
- Dimension Between Legs Is 40" (101.6cm)

#SURG-TABLE-STAINLESS 3805.00 EACH

Memory Foam Chair Pad For "Long Appointment" Patient Comfort



*Memory Foam For
Ultra Comfortable
Pressure Distribution*

MOST POPULAR

- Fits All Dental Chairs - 61" Long x 18" Wide x 3" Thick
- Non-Skid Bottom Won't Slip Or Move On Your Chair
- Sealed Seams & Fluid Proof Top Cover For Infection Control
- Anti-Microbial Foam & Cover

#CHAIR-PAD-SET 769.00 CUSHION & HEADREST
2 OR MORE 715.00 CUSHION & HEADREST

SALVIN ENLIGHTEN-PRO3™ LED SYSTEM



- Continuously Variable Adjustable Intensity
- 100,000 Hour LED Life
- Dual Lithium Polymer Battery Pack
 - 6.5 Hour Capacity On High
 - 14 Hour Capacity On Medium
 - 50 Hour Capacity On Low
- Spot Size Adjusts To 3.5" At 18" Working Distance
- Total Weight: 11oz (311g)
- Optimized Weight Balance And Comfort
- 2 Year Warranty

MOST POPULAR



Includes Batteries & Charger



Includes Anti-Curing Filter



Use With Your Loupes



Includes Carrying Case

SPECIAL

Enlighten-Pro3™ LED Headlamp System

Includes:

- Enlighten-PRO3™ LED Headlamp
- Secure Carrying Case
- Anti-Curing Filter
- 2 Lithium Batteries, Charger & Cable

#HL-PRO3-PKG..... 2059.00

- **Patented Technology**
Reduces LED Temperature
- **No Cables Or Belt Packs**
- **100,000 LUX At 10" (25cm)**
- **15,000 LUX At 18" (46cm)**
- Preferred Working Distance For Most Doctors
(Important To Measure At Your Working Length
As Luminance Decreases With Distance)

Enlighten-Pro3™ Accessories



Extra Enlighten-Pro3™ Lithium Batteries - Set Of 2

#HL-PRO3-BAT (SET OF 2)..... 315.00



Optional Enlighten-Pro3™ Four-Bay Charging Station

- Charge Up To Four Batteries At Once
- Charge Status Indicators For Each Battery
- Batteries Not Included

#HL-PRO-DOC 572.00

#1090 Salvin Split-Sheet Implant Pack

"Eyes Visible"
Patient Drape



MOST POPULAR

#1090 Salvin Split-Sheet Pack

#1090 SOLD IN CASES OF 6

1 CASE 473.70 CASE (78.95 EA)

2 CASES 435.30 CASE (72.55 EA)



- Cooler
- Lighter

*This Implant Pack Contains
2 Aurora™ Gowns*

Scan For

**60^{Second}
Product
Video™**



Scan For

**60^{Second}
Product
Video™**




BUYER'S GUIDE - SALVIN SPLIT-SHEET PACK

▶ "Eyes Visible" Patient Drape

▶ One Pack With Two Separate Inner Wraps For Easier Set Up

Doctor/Assistant Inner Wrap



- 2 Large Operating Room Surgical Gowns
- 2 17" x 20" (38 x 51cm) Absorbent Towels
- 2 Bouffant Caps (Doctor & Assistant)
- 2 Face Masks With Ties
- Outer Wrapper Becomes 24" x 24" (61 x 61cm) Table Cover


Patient / Room Setup Inner Wrap



- 1 Implant "Split Sheet" Patient Drape 77" x 108" (195.5 x 274cm)
- 1 Bouffant Cap (For The Patient)
- 1 Cotton O/R Towel
- 1 Mayo Stand Or Back Table Cover
- 1 Oral Surgery Head Drape
- 2 18" x 24" (45.7 x 61cm) Utility Drapes With Tape
- 1 Machine or X-Ray Cover
- 4 4" x 6" (10 x 15.25cm) Light Handle/Chair Button Covers
- 4 Rapid Roll® Covers
- 1 Chairside Bag
- 2 6" (15.25cm) Cotton Tip Applicators
- 1 Surgical Aspirator Tip
- 1 Suction Tip With Out Slots
- 1 Saliva Ejector
- 25 Lint Free 2" x 2" (5 x 5cm) Sponges
- 8 Lint Free 4" x 4" (10 x 10cm) Sponges
- Outer Wrapper Becomes 45" x 45" (114 x 114cm) Table Cover

Latexfree

#1033 Salvin "Eyes Covered" Implant Pack




MOST POPULAR

#1033 Basic Implant Pack

#1033 SOLD IN CASES OF 6

1 CASE 468.00 CASE (78.00 EA)

2 CASES 431.10 CASE (71.85 EA)




Sterile Surgical Gloves
On Page 68

BUYER'S GUIDE - SALVIN "EYES COVERED" PACK

▶ "Eyes Covered" Patient Drape


▶ **One Pack With Two Separate Inner Wraps For Easier Set Up**

Doctor/Assistant Inner Wrap



- 2 Large Operating Room Surgical Gowns
- 2 17" x 20" (38 x 51cm) Absorbent Towels
- 2 Bouffant Caps (Doctor & Assistant)
- 2 Face Masks With Ties
- Outer Wrapper Becomes 24" x 24" (61 x 61cm) Table Cover

Patient / Room Setup Inner Wrap



- 1 Implant Fenestrated Patient Drape 38 1/2" x 55" (97.5 x 139.5cm)
- 1 Bouffant Cap (For The Patient)
- 1 Cotton O/R Towel
- 1 Mayo Stand Or Back Table Cover
- 3 18" x 24" (45.7 x 61cm) Utility Drapes With Tape
- 1 Machine or X-Ray Cover
- 4 4" x 6" (10 x 15.25cm) Light Handle/Chair Button Covers
- 4 Rapid Roll® Covers
- 1 Chairside Bag
- 2 6" (15.25cm) Cotton Tipped Applicators
- 1 Surgical Aspirator Tip
- 1 Suction Tip With Out Slots
- 1 Saliva Ejector
- 25 Lint Free 2" x 2" (5 x 5cm) Sponges
- 8 Lint Free 4" x 4" (10 x 10cm) Sponges
- Outer Wrapper Becomes 45" x 45" (114 x 114cm) Table Cover

Latexfree

Resnik Institute Implant Surgery Pack

"Eyes Covered"
Patient Drape




Resnik Institute Implant Surgery Pack

#RESNIK-PACK SOLD IN CASES OF 10
1 CASE 424.00 CASE (42.40 EA)
2 CASES 390.00 CASE (39.00 EA)

BUYER'S GUIDE - RESNIK INSTITUTE PACK

▶ "Eyes Covered" Patient Drape

Resnik Institute Implant Surgery Pack



- 1 Implant Fenestrated Patient Drape
38 1/2" x 55" (97.5 x 139.5cm)
- 5 Cotton O/R Towels
- 1 Chair Cover 41" x 85" (104 x 216cm)
- 1 Bouffant Cap
- 2 Rapid Roll® Covers
- 1 20" x 20" (50.8 x 50.8cm) CSR Wrap
- 2 30" x 30" (76.2 x 76.2cm) CSR Wraps
- 4 4" x 6" (10 x 15.25cm) Light Handle/
Chair Button Covers
- 2 Surgical Aspirator Tips
- 1 Suction Tip Without Slots
- 12 Lint Free 4" x 4" (10 x 10cm) Sponges
- 12 Lint Free 2" x 2" (5 x 5cm) Sponges
- Outer Wrapper Becomes
48" x 48" (121.9 x 121.9cm) Table Cover

Latexfree

Salvin E.N.T. Head & Neck Pack

"Eyes Visible"
Patient Drape




Salvin E.N.T. Head & Neck Pack

#2000 SOLD IN CASES OF 8
1 CASE 468.40 CASE (58.55 EA)
2 CASES 425.60 CASE (53.20 EA)

BUYER'S GUIDE - SALVIN E.N.T HEAD & NECK PACK

▶ "Eyes Visible" Patient Drape

Salvin E.N.T. Head & Neck Pack



- 2 Large Operating Room Surgical Gowns
- 2 17" x 20" (38 x 51cm) Absorbent Towels
- 1 Implant "Split Sheet" Patient Drape
77" x 108" (195.5 x 274cm)
- 1 Oral Surgery Head Drape
- 1 Mayo Stand Or Back Table Cover
- 1 18" x 24" (45.7 x 61cm) Utility Drape
With Tape
- 1 Chairside Bag
- Outer Wrapper Becomes
54" x 54" (137 x 137cm) Table Cover

Latexfree



- Cooler
- Lighter

This Implant Pack Contains
2 Aurora™ Gowns

Salvin Perio/Minor Surgery Pack

Salvin Perio/Minor Surgery Pack

#1034 SOLD IN CASES OF 10

1 CASE 364.00 CASE (36.40 EA)

2 CASES 332.00 CASE (33.20 EA)




Sterile Surgical Gloves
On Page 68

BUYER'S GUIDE - SALVIN PERIO/MINOR SURGERY PACK

► This Pack Includes A Sterile Patient Bib

Salvin Perio/Minor Surgery Pack



- 2 Yellow Isolation Gowns
- 2 17" x 20" (38 x 51cm) Absorbent Towels
- 1 18" x 36" (46 x 91cm) Patient Bib
- 2 Face Masks W/Ties
- 1 Cotton O/R Towel
- 1 Chairside Bag
- 2 18" x 24" (45.7 x 61cm) Utility Drapes With Tape
- 2 4" x 6" (10 x 15.25cm) Light Handle/Chair Button Covers
- 2 6" (15.25cm) Cotton Tipped Applicators
- 1 Suction Tip Without Slots
- 1 Saliva Ejector
- 25 Lint Free 2" x 2" (5 x 5cm) Sponges
- Outer Wrapper Becomes 36" x 36" (91 x 91cm) Table Cover

Latexfree

Individual Sterile Patient Drapes

MOST POPULAR

Split Sheet "Eyes Visible" Sterile Patient Drape

Latexfree

- Individually Packaged
- 77" x 108" (195.5 x 274cm)

#DRAPE-SPLIT SHEET

CASE OF 24 504.00 CASE
2 OR MORE CASES 470.70 CASE

Patient Drape Detail



#DRAPE-SPLIT SHEET

Drape 33 "Eyes Covered" Sterile Patient Drape

Latexfree

- Individually Packaged
- 38 1/2" x 55" (97.5 x 139.5cm)

#DRAPE33

CASE OF 24 437.00 CASE
2 OR MORE CASES 414.30 CASE

Patient Drape Detail



#DRAPE 33

Rapid Roll® Covers & Sterile Drapes

SALVIN
Rapid Roll®

MOST POPULAR

- Individually Packaged
- 65" (165cm) Long

#RAPID-ROLL

BOX OF 25 121.50 BOX
2 OR MORE BOXES 112.10 BOX



Latexfree

Scan For

60 Second Product Video™



Scan For

60 Second Product Video™



Sterile Drapes

- Individually Packaged Poly Lined Paper Drape
- 18" x 26" (45.7 x 66cm)
- Use As Patient Bib, Table Cover, Etc.

#TOWEL-DRAPE-50

BOX OF 50 60.20 BOX
2 OR MORE BOXES 57.60 BOX





Latexfree

Salvin Sterile Gown/Towel Packs



Latexfree



Salvin Sterile Gown / Towel Pack - UNISIZE Fit Gown

- 1 Sterile O/R Gown (Large-XLarge Fit)
- 1 Absorbent Scrub Towel 15" x 30" (38 x 76cm)

#1001

CASE OF 24 499.00 CASE

2 OR MORE CASES 483.20 CASE

Salvin Sterile Gown / Towel Pack - SMALL Size Gown

- 1 Sterile O/R Gown (Small-Medium Fit)
- 1 Absorbent Scrub Towel 15" x 30" (38 x 76cm)

#GOWN-SMALL

CASE OF 30 420.00 CASE

2 OR MORE CASES 404.20 CASE

Sterile Absorbent Scrub Towels

- 80 Individually Packaged Absorbent Scrub Towels 15" x 30" (38 x 76cm)

#TOWEL-SCRUB-80

CASE OF 80 166.00 CASE

2 OR MORE CASES 156.70 CASE

Backhaus Towel Clamp

- Sharp Points
- Use For Cloth Drapes

#BACKHAUS 23.60 EACH

3 OR MORE 21.50 EACH



Aurora™ Gown

- More Protection
- Cooler
- Lighter

O/R & Isolation Gowns – Non-Sterile



Latexfree



Aurora™ Gown

- More Protection
- Cooler
- Lighter

Disposable Blue O/R Gown (Non-Sterile)

- The Feel Of Cotton Without The Weight
- Breathable Gown With Knitted Cuffs
- Tie Back & Velcro Neck Closure
- Generously Sized Unisize Fit

#OR-GOWN

CASE OF 50 460.00 CASE

2 OR MORE CASES 434.40 CASE



Latexfree

MOST POPULAR

Isolation Gown-Yellow (Non-Sterile)

- Elastic Cuff
- Unisize Fit

#ISO-GOWN-50

CASE OF 50 108.00 CASE

2 OR MORE CASES 100.50 CASE

Surgical N95 Respirator

Surgical N95 Respirator

- Approved For Use By The FDA And NIOSH
- Sold In Boxes Of 30

#N95-MASK-5076.60 BOX OF 30



Towels, Caps & Shoe Covers



Surgeon's Caps

#SURG CAP 100
BOX OF 100 41.00 BOX
2 OR MORE BOXES37.70 BOX



MOST POPULAR

Bouffant Caps

#BOUF CAP 100
BOX OF 10035.70 BOX
2 OR MORE BOXES32.70 BOX



Blue 100% Cotton O/R Towels

16" x 27" (40.6 x 68.6cm)
Single Use
#OR TOWEL 50
CASE OF 50135.00 CASE
2 OR MORE CASES 128.70 CASE



Shoe Covers-Unisize

- Fits Over Regular Shoes
- No-Slip Grip On Bottom For Safety In The Office

#SHOE 100
BOX OF 10044.00 BOX
2 OR MORE BOXES41.00 BOX



Surgical Masks W/Ties

Soft Polypropylene Facing
For Comfort & Protection
#SURG-MASK-TIE
BOX OF 5026.50 BOX
6 OR MORE BOXES23.70 BOX



MOST POPULAR

Shoe Covers-Sport

- Fits Over Athletic Shoes
- No-Slip Grip On Bottom For Safety In The Office

#SHOE-SPORT-100
BOX OF 10049.50 BOX
2 OR MORE BOXES46.60 BOX

Patient Skin Prep

ChloraPrep® Surgical Skin Prep Individual Swabs

- Chlorhexidine Gluconate/IPA
- Effective Skin Antiseptic
- Will Not Stain Or Irritate Skin
- Rapid Action Against Microorganisms
- Available In US Only

#CHLORA-PREP
BOX OF 2063.70 BOX
2 BOXES60.10 BOX



MOST POPULAR

Patient Prep Individual Swab Sticks

Povidone Iodine (Iodophor)
#SWABSTICK 50
BOX OF 50 23.40 BOX
2 BOXES 21.20 BOX



Cold Therapy



Plastic Outer Covering

Latexfree**MOST POPULAR**

Instant Cold Pack

- Squeeze To Activate
- 6" x 8.75" (15.2 x 22.2cm)

#INSTANT-COLD

BOX OF 2453.80 BOX
2 OR MORE BOXES50.10 BOX



"Sweatless" Covering For Comfort

Latexfree

Deluxe Instant Cold Pack

- Squeeze To Activate
- 6.75" x 9" (17.1cm x 22.86cm)

#INSTANT-COLD-DELUXE

BOX OF 2476.50 BOX
2 OR MORE BOXES70.10 BOX

**Latexfree**

Bilateral Ice Pack

- No Wetness Or Leakage
- Refillable Ice Pack With Velcro Strap

#BILAT-24

CASE OF 24 231.00 CASE
2 OR MORE CASES 221.00 CASE

Sterillium Rub® Surgical Hand Scrub

- Exceeds FDA Requirements For Surgical Hand Antisepsis
- 80% Ethanol For Superior Germ Killing
- Chlorhexidine Gluconate (CHG) Free To Prevent Hand Irritation



Touchless Activation Dispenser

- Mount With Screws (Included) Or Double-Sided Tape
- Uses 4 "D" Batteries (Not Included)

#SCRUB-STERI-DISPENSE182.00 EACH



Sterillium® Rub Surgical Hand Scrub Refills

- Approximately 850 Uses Per Bottle
- Fits Touchless Dispenser On This Page

#SCRUB-STERI-1000. 249.80 PKG OF 2

SPECIAL

Sterillium® Starter Kit

Includes:

- Touchless Activation Dispenser
- 2 1000mL Bottles

#SCRUB-STERI-KIT
370.00 COMPLETE KIT2 OR MORE STARTER KITS
340.00 COMPLETE KIT

Scrub Brushes

**MOST POPULAR**

"Clens" Type Surgical Scrub Brush / Sponge

PCMX "Clens" Solution Is Kinder To Sensitive Hands

#SCRUB-CLEN-30

BOX OF 30 76.40 BOX
2 BOXES 71.10 BOX

Sterile Equipment Covers

Sterile "Self Stick" **MOST POPULAR** Light Handle Covers

- 20cm x 20cm
- Packaged Sterile In Sets Of 2 (50 Sets Per Case)

#LIGHT-HANDLE-COVER

1 CASE 199.90 CASE
2 OR MORE CASES 189.30 CASE



Sterile Mayo Stand Cover

- Individually Packaged
- Fits Regular Mayo Stands
- 24" x 24" (61 x 61cm)

#MAYO-COVER

CASE OF 24 324.00 CASE
2 OR MORE CASES 311.80 CASE



Latexfree

Sterile Machine Or X-Ray Cover

- Individual Sterile Machine Cover For X-Ray Cassette, Arm Rest, Implant Motor, Etc.

- 24" x 40" (61 x 101.6cm)

#MACHINE-COVER-50

CASE OF 50 219.80 CASE
2 OR MORE CASES 209.60 CASE



Latexfree

Sterile Surgical Table Cover

- Individually Packaged
- Fits Surgical Tables
- 53" x 77" (134 x 195cm)

#SURG-TABLE-COVER

CASE OF 20 114.00 CASE
2 OR MORE CASES 106.40 CASE



Latexfree

Sterilization Wrap For Surgical Kits & Steam Autoclave Tape



18" x 18" Sterilization Wrap

45.7cm x 45.7cm

#INST-WRAP-18 248.30 CASE OF 1000

Latexfree

MOST POPULAR

24" x 24" Sterilization Wrap

61cm x 61cm

#INST-WRAP-24 194.00 CASE OF 500

Latexfree



Steam Autoclave Tape - 1" x 60 YDS

- Changes From Beige To Black To Indicate Proper Sterilization Cycle

#CLAVE-TAPE-1X60 12.00 ROLL
5 OR MORE 10.20 ROLL

Latexfree

Instrument Cleaners



EZ-Zyme Enzyme Instrument Cleaner

- Concentrated - 1 Packet Makes A Gallon
- Dissolves Blood & Tissue
- Instrument Pre-Soak
- Ultrasonic & Evac System Cleaner
- Biodegradable & Neutral PH

#EZ-ZYME

BOX OF 32 PACKETS 102.80 BOX
2 OR MORE 96.70 BOX



Instrument Stain Remover

- Removes Stains & Tarnish
- Will Not Harm Instrument Surfaces

#STAIN 19.10 EACH
2 OR MORE 17.00 EACH



MOST POPULAR

Surgical Instrument Spray Lubricant - 8 oz.

"Instrument Milk" In A Spray

- For Hinged Instruments
- Keeps Scissors Sharp
- Preserves Instrument Life
- Non-Greasy, Anti-Corrosive

#SPRAY LUBE 19.10 EACH
2 OR MORE 17.00 EACH

Stainless Steel Instrument Cassettes



MOST POPULAR

8 Instrument Stainless Steel Cassette

- Silicone Cushions & Locking Crossbar
- 7.3" x 5.2" x 0.9" (18.5 x 13.3 x 2.3cm)

#CASS-SS-8 **191.00 EACH**
 3 OR MORE **181.80 EACH**

MOST POPULAR

16 Instrument Stainless Steel Cassette

- Silicone Cushions & Locking Crossbar
- 10.7" x 7.3" x 0.9" (27.3 x 18.5 x 2.3cm)

#CASS-SS-16 **226.00 EACH**
 3 OR MORE **214.10 EACH**

Plastic Instrument Cassettes



10" x 6" Cassette

- "Finger Tip" Mat
- Autoclave, Dry Heat, Chemical Sterilization
- 25cm x 15.3cm - 0.75" (1.9cm) Height

#CASSETTE-10X6 **141.90 EACH**
 3 OR MORE **134.10 EACH**



Ideal For Large Instruments

MOST POPULAR

10" x 6" Deep Cassette

- "Finger Tip" Mats
- Autoclave, Dry Heat, Chemical Sterilization
- 25cm x 15.3cm - 1.5" (3.7cm) Height

#CASSETTE-DEEP-BASE **169.90 EACH**
 3 OR MORE **158.80 EACH**

Sterilization Pouches For Instrument Cassettes



12" x 15" Self-Sealing Instrument Cassette Sterilization Pouches

- Fits All Salvin Instrument Cassettes
- Fits Other Brand Cassettes Up To 8" x 12" (20.3cm x 30.5cm)
- External & Internal Indicators With See-Through Packaging For Visibility

#POUCH-12X15-STER **122.00 BOX OF 200**
 2 OR MORE **114.40 BOX OF 200**

Latexfree

Salvin "Where Is It?" Small Parts Cassettes

Perfect For Busy Practices – No More Looking For Small Parts!



For Larger Parts



"Where Is It?" Stainless Steel Cassette - Large

- 7.1cm Wide x 5.3cm Deep x 2.3cm Tall

#WHERE-IS-IT-LARGE 49.40 EACH
3 OR MORE 45.60 EACH



For Smaller Parts




MOST POPULAR

"Where Is It?" Stainless Steel Cassette - Small

- 7.1cm Wide x 3.7cm Deep x 1.4cm Tall

#WHERE-IS-IT-SMALL 41.60 EACH
3 OR MORE 38.80 EACH



Includes Removable Silicone Mat

MOST POPULAR

"Where Is It?" - Plastic Cassette

- Removable "Finger Tip" Mat
- 7.3cm Wide x 4cm Deep x 1.9cm Tall

#WHERE-IS-IT-PLASTIC 31.90 EACH
3 OR MORE 29.40 EACH



SPECIAL

Set Of 3 Salvin "Where Is It?" Cassettes



#WHERE-IS-IT-SET-3 90.30 SET

Small Parts Organizing System

- One 10"x 6" Deep Cassette & Nine "Where-Is-It" Plastic Cassettes
- Autoclave, Dry Heat, Chemical Sterilization

#CASS-ORG-KIT 397.00 SET



Bur Blocks On Page 89

Signature™ Latex Essential Surgical Gloves

**Synthetic Polymer Coating For Exceptional Damp Hand Donning**

- Beaded Cuff For Easy Donning Without Roll-Back
- Proprietary Coating For Easy Damp/Dry Donning

Glove Guide	
▶	Powder Free
▶	Anatomically Shaped
▶	Sterile
▶	Latex

Size	Product Number	
6	#GL-MSG-5960	
6 1/2	#GL-MSG-5965	
7	#GL-MSG-5970	
7 1/2	#GL-MSG-5975	
8	#GL-MSG-5980	BOX OF 50 PAIRS 134.00 BOX
8 1/2	#GL-MSG-5985	2 OR MORE BOXES 125.60 BOX

SensiCare® - Synthetic Surgical Gloves

MOST POPULAR**Fit & Feel Of Latex But Without The Latex Proteins & Allergens**

Pure Aloe Vera Coating Is Applied To Interior Glove Surface. Body Heat Activates The Aloe, Coating & Protecting Your Hands.

- Beaded Cuff For Easy Donning Without Roll-Back
- Unique Styrene-Isoprene-Styrene (SIS) Synthetic
- Proprietary Coating For Easy Damp/Dry Donning
- Bisque Finish For A Better Grip

Glove Guide	
▶	Powder Free
▶	Anatomically Shaped
▶	Sterile
▶	Latex Free Synthetic
▶	Aloe Vera

Size	Product Number	
6	#GL-MDS-377060	
6 1/2	#GL-MDS-377065	
7	#GL-MDS-377070	
7 1/2	#GL-MDS-377075	
8	#GL-MDS-377080	BOX OF 25 PAIRS 136.00 BOX
8 1/2	#GL-MDS-377085	2 OR MORE BOXES 127.60 BOX

SensiCare® PI Micro - Synthetic Surgical Gloves

**15% Thinner Than Standard Gloves To Enhance Tactile Sensitivity**

- Beaded Cuff For Easy Donning Without Roll-Back
- Proprietary Coating For Easy Damp/Dry Donning

Glove Guide	
▶	Powder Free
▶	Anatomically Shaped
▶	Sterile
▶	Latex Free Synthetic

Size	Product Number	
6	#GL-MSG-9660	
6 1/2	#GL-MSG-9665	
7	#GL-MSG-9670	
7 1/2	#GL-MSG-9675	
8	#GL-MSG-9680	BOX OF 50 PAIRS 236.50 BOX
8 1/2	#GL-MSG-9685	2 OR MORE BOXES 224.90 BOX

See For Yourself! Please Call For Samples Of Any Of Our Gloves

Blood Collection



Basic Blood Collection Setup

Vacutainer® Needle 21 Ga x 1" (2.54cm)	
#VACUTAIN-NEEDLE	76.00 BOX OF 100
Vacutainer® Needle Holder	
#VACUTAIN-HOLDER.	52.00 BOX OF 50
IV Adapter	
#BIO-PRF-IV-CONNECT	206.00 BOX OF 100
Tourniquet-Latex Free	
#VACUTAIN-TOURNIQUET	38.30 BOX OF 25
Vacutainer® Tube-Red - 10mL - No Additives	
#VACUTAIN-RED-TUBE	152.00 PKG OF 100
Vacutainer® Tube-Yellow - 10mL (8.5mL Draw) Citrate Additive	
#VACU-TUBE-YELLOW	175.90 PKG OF 100
10cc Sterile Syringe	
#SYRINGE-10CC.	73.90 PKG OF 100
3cc Sterile Syringe	
#SYRINGE-STERILE-3CC.	68.10 PKG OF 100
Blunt Catheter Sterile 16 Ga x 3.5" (9cm)	
#BLUNT-CATH-16G-9CM	410.00 BOX OF 100
Blunt Catheter Sterile 18 Ga x 1" (2.54cm)	
#BLUNT-CATH-18G-2.5CM.	129.00 BOX OF 100
Blood Collection Tube Rack	
#VACU-TUBE-RACK	39.80 EACH
3cc Sterile Syringe With 21 Ga x 1.5" (3.81cm) Catheter	
#SYRINGE-3CC-21G	64.50 BOX OF 100
Butterfly To Vacutainer® Adapter 21 Ga x .75" (1.9cm)	
#BUTTERFLY-21G	109.00 BOX OF 50

Vacutainer® is a registered trademark of Becton, Dickinson and Company

BIO-PRF Centrifuge



BIO-PRF Horizontal Centrifuge

- Up To 4X More Platelets And Leukocytes
- Higher Growth Factor Release When Compared To Fixed-Angle Centrifugation
- Better Horizontal Layer Separation Between PRF-Clot And Red Blood Cells
- More Even Distribution Of Regenerative Cells Throughout The PRF-Clot
- Less Cell Damage / Death On The Tube Walls Produced On Fixed-Angle Centrifugation
- Less RBC Accumulation

Includes:

- BIO-PRF Centrifuge
- BIO-PRF Box
- Set Of 6 BIO-PRF Instruments
- BIO-PRF Bowl - For Sticky Bone
- BIO-PRF Tray - For Membranes
- BIO-PRF Mixing Bowl
- BIO-PRF 24 Tube Holders
- BIO-PRF Latex Free Tourniquet
- BIO-PRF Cassette
- Box Of 100 BIO-PRF Red Tubes
- Box Of 100 BIO-PRF Blue Tubes
- Box Of 24 BIO-PRF Butterfly Needles

#BIO-PRF-COMP-KIT 3450.00



BIO-PRF Red Tubes

#BIO-PRF-RED-TUBES
.....137.00 BOX OF 100

BIO-PRF Blue Tubes

#BIO-PRF-BLUE-TUBES
.....117.00 BOX OF 100

BIO-HEAT System



BIO-HEAT System To Extend PRF Resorption

- Extend PRF Resorption From 2-3 Weeks To 4-6 Months

Indicated For:

- Barrier Membrane For GBR Procedures
- Covering The Lateral Window During Sinus Lifts
- Root Coverage
- Mixing With Bone Grafting Materials
- Insertion Into Perio Pockets Non-Surgically For Stability

#BIO-HEAT-KIT 2050.00

Centrifuge



Variable Speed Centrifuge

- Now With Digital Display
- Adjustable Speed Range Of 500rpm - 4,000rpm
- Compact Size 10.2" x 11.8" x 8" (26cm x 30cm x 20cm)
- Very Quiet Sealed Motor
- 1 Year Parts & Labor Warranty
- 110 Volts

#CENTRIFUGE-VARISPEED..... 1520.00



"Rutkowski" Blood Collection Starter Package

Disposables Package For Approximately 50 Patients
Includes:

- Box Of 50 Vacutainer® Needle Holders
- Box Of 25 Latex Free Tourniquets
- 1 Blood Collection Tube Rack
- Pkg Of 100 Red Top Vacutainer® Tubes
- Pkg Of 100 Yellow Top Vacutainer® Tubes
- Pkg Of 100 3cc Sterile Syringes
- Box Of 50 Butterfly To Vacutainer® Adapters
- PRF Box

#BLOOD-COLL-RUTKOWSKI

..... 890.00/PKG FOR APPROX 50 PATIENTS

Disposable Cost Approximately 17.00 Per Patient

PRF Box



PRF Box

- Form Fibrin Membrane Or Plug Shape From Clot
- Weighted Top Creates 1mm Thick Membrane - No Force Needed

#PRF-BOX.....355.00



Form Fibrin Membrane Shape From Clot



Form Fibrin Plug Shape From Clot

Piezotome Cube

Safely Cuts Bone While Preserving Soft Tissue

*Power Adjusts Dynamically For Effortless Cuts
Will Not Damage Soft Tissue*


Automatically
Reduces Power
In Softer Bone

Automatically
Increases Power
In Denser Bone



BUYER'S GUIDE	
▶	Console With Irrigation
▶	LED Piezo Handpiece
▶	Touch Screen Display

PROCEDURE GUIDE	
▶	Atraumatic Extractions
▶	Block Grafting
▶	Osteoplasty
▶	Lateral Sinus Lift
▶	Crestal Sinus Lift
▶	Ridge Splitting
▶	Osteotomies
▶	Crown Lengthening



Call For Package Specials



"There is nothing like the Piezotome for lifting the sinus membrane, harvesting autogenous bone & removing lingual tori while protecting the soft tissue in the floor of the mouth."

—Dr. John Russo
Periodontist, Sarasota, FL



"The use of the Piezotome has been the most revolutionary development in extractions in my career. Most extractions are completed in a matter of minutes with NO loss of alveolar bone."

—Dr. Michael Sonick
Periodontist, Fairfield, CT

Implant Center 2



Piezotome + Implant Motor Combination

Scan For
60th Second Product Video™



Motor Fits Your Existing Implant Handpieces
Cuts Bone While Preserving Soft Tissue
LED Handpiece For Excellent Visibility
60 Watts Of Power For Faster Cuts

BUYER'S GUIDE

▶ Console With Irrigation	▶ LED Implant Handpiece
▶ LED Piezo Handpiece	▶ 100 to 40,000 RPM
▶ Touch Screen Display	

PROCEDURE GUIDE

▶ Implant Placement
▶ Atraumatic Extractions
▶ Block Grafting
▶ Osteoplasty
▶ Lateral Sinus Lift
▶ Crestal Sinus Lift
▶ Ridge Splitting
▶ Osteotomies
▶ Crown Lengthening



Call For Package Specials

Piezotome Accessories



Autoclavable Tubing With Pump

Includes 10 Disposable Spikes

#PIEZO-IRR-AC 222.00 EACH

Disposable Irrigation Spikes



Pkg Of 10 Spikes

#PIEZO-IRR-SPIKE 65.00 PKG OF 10



***Salvin Surgical
Cart On Page 54***




Tips Shown By Procedure On Next Two Pages ➔

Atraumatic Extraction Kit **MOST POPULAR**

For Piezotome 2 or Cube

#PIEZO-2-EXTRACTION-KIT

WITH CASSETTE. 1833.00 KIT



LC1 Tip



LC2 Tip



LC2L Tip



LC2R Tip



Ninja Tip




Piezotome 2 or Cube	#PIEZO2-TIP-LC1	#PIEZO2-TIP-LC2	#PIEZO2-TIP-LC2L	#PIEZO2-TIP-LC2R	#PIEZO2-TIP-NINJA
	333.00 EACH	333.00 EACH	333.00 EACH	333.00 EACH	333.00 EACH

Sinus Lift Kit **MOST POPULAR**

For Piezotome 2 or Cube

#PIEZO-2-SL-TIP-KIT

WITH CASSETTE. 1500.00 KIT



SL1 Diamond Tip



SL2 Diamond Tip



SL3 Tip



SL4 Tip



SL5 Tip



Piezotome 2 or Cube	#PIEZO2-TIP-SL1	#PIEZO2-TIP-SL2	#PIEZO2-TIP-SL3	#PIEZO2-TIP-SL4	#PIEZO2-TIP-SL5
	333.00 EACH	333.00 EACH	333.00 EACH	333.00 EACH	333.00 EACH

Bone Surgery Kit

For Piezotome 2 or Cube

#PIEZO-2-BS-TIP-KIT

WITH CASSETTE. 1833.00 KIT



BS1-Slim Tip



BS2L Tip



BS2R Tip



BS4 Tip



BS5 Tip



BS6 Tip



Piezotome 2 or Cube	#PIEZO2-TIP-BS1-SLIM	#PIEZO2-TIP-BS2L	#PIEZO2-TIP-BS2R	#PIEZO2-TIP-BS4	#PIEZO2-TIP-BS5	#PIEZO2-TIP-BS6
	333.00 EACH	333.00 EACH	333.00 EACH	333.00 EACH	333.00 EACH	333.00 EACH

Essential Kit

For Piezotome 2 or Cube

#PIEZO2-ESSENTIAL-KIT

WITH CASSETTE. 1833.00 KIT



						
	BS1-Slim Tip	BS4 Tip	SL1 Diamond Tip	SL2 Diamond Tip	SL3 Tip	LC2 Tip
Piezotome 2 or Cube	#PIEZO2-TIP-BS1-SLIM	#PIEZO2-TIP-BS4	#PIEZO2-TIP-SL1	#PIEZO2-TIP-SL2	#PIEZO2-TIP-SL3	#PIEZO2-TIP-LC2
	333.00 EACH	333.00 EACH	333.00 EACH	333.00 EACH	333.00 EACH	333.00 EACH





Crest Splitting Kit **MOST POPULAR**

For Piezotome 2 or Cube

#PIEZO-2-CS-TIP-KIT

WITH CASSETTE. 1833.00 KIT



						
	CS1 Tip 0.55mm	CS2 Tip 0.85mm	CS3 Tip 0.5mm - For Side Cuts	CS4 Tip 1.8mm	CS5 Tip 2.75mm	CS6 Tip 3.75mm
Piezotome 2 or Cube	#PIEZO2-TIP-CS1	#PIEZO2-TIP-CS2	#PIEZO2-TIP-CS3	#PIEZO2-TIP-CS4	#PIEZO2-TIP-CS5	#PIEZO2-TIP-CS6
	333.00 EACH	333.00 EACH	333.00 EACH	333.00 EACH	333.00 EACH	333.00 EACH







Osteotomy Kit

For Piezotome 2 or Cube

#PIEZO-2-INTRALIFT-KIT

WITH CASSETTE. 1500.00 KIT



					
	TKW1 Tip	TKW2 Tip	TKW3 Tip	TKW4 Tip	TKW5 Tip
Piezotome 2 or Cube	#PIEZO2-TIP-TKW1	#PIEZO2-TIP-TKW2	#PIEZO2-TIP-TKW3	#PIEZO2-TIP-TKW4	#PIEZO2-TIP-TKW5
	333.00 EACH	333.00 EACH	333.00 EACH	333.00 EACH	333.00 EACH






Crown Lengthening Kit

For Piezotome 2 or Cube

#PIEZO-2-CE-TIP-KIT

WITH CASSETTE. 1225.00 KIT



				
	CE1 Tip	CE2 Tip	CE3 Tip	BS6 Tip
Piezotome 2 or Cube	#PIEZO2-TIP-CE1	#PIEZO2-TIP-CE2	#PIEZO2-TIP-CE3	#PIEZO2-TIP-BS6
	333.00 EACH	333.00 EACH	333.00 EACH	333.00 EACH

For Accelerated Orthodontic Bone Cuts



Rounded Piezo Tip

#PIEZO2-TIP-BS1-ROUND

333.00 EACH

Newtron® P5 XS BLED Ultrasonic Scaler

- *Perform Faster, More Comfortable Treatments With 80+ Unique Tips*
- *Auto-Responsive Handpiece & Precision Tips Ensure Optimal Protection For Teeth & Soft Tissue*
- *Smooth, Linear Vibrations Reduce Gingival Impact & Preserve Cementum*

SPECIAL

Newtron® P5 XS BLED Package With Intro Scaling Set

Includes:

- Newtron® P5 XS BLED Unit
- One LED Handpiece
- Newtron® Intro Scaling Kit With 3 Tips
- Autoclavable Cassette

#NEWTRON-F62117
Call For Lowest Price


Newtron® Intro Scaling Kit

- Includes: Tips – H3, H4L, H4R With 3 Blue Wrenches
- Autoclavable Cassette

#NEWTRON-INTRO-KIT
464.00


Implant Protect Tip Kit

- Includes: Tips – IP1, IP2L, IP2R, IP3L, IP3R With Universal Wrench
- Autoclavable Cassette

#PIEZO-F02120 774.00

Air-N-Go® Air Polishing



Air-N-Go® Easy

- Versatile Nozzles & Tips Provide Effective Plaque, Stain, And Implant Cleaning In Both Supragingival & Subgingival Applications
- Gentle, Comfortable Treatments: Supra And Perio Powders Offer Pain-Free Cleaning
- Easy Setup & Use: Quick Start-Up – Connects Directly To High-Speed Handpieces

Call For Lowest Price


Air-N-Go® Classic Air Polishing Powder Variety Pack

- Box Of Four 250gm Bottles
- Variety Pack Contains One Bottle of Each Flavor: Cola, Raspberry, Peppermint, Lemon
- Also Available Per Flavor In Packs Of Four 250gm Bottles

#AIR-N-GO-F10255 147.00 BOX OF 4

SALVIN® **RotoBrush-Titanium®** For Exposed Implant Threads

- Cleans Infected Implant Surface Without Grinding The Threads
- Maintains Structure & Shape Of The Implant
- Gentler On Implant Surface Than Diamond Burs & Ultrasonics
- Start At 300rpm Working Up To 1,200rpm As Needed



MOST POPULAR

RotoBrush-Titanium® For Vertical Access

Single Use

#ROTOBRUSH-POCKET 42.90 EACH
3 OR MORE 37.90 EACH



RotoBrush-Titanium® For Posterior Access

Single Use

#ROTOBRUSH-SHORT-POCKET 42.90 EACH
3 OR MORE 37.90 EACH



RotoBrush-Titanium® For Wide Area

Single Use

#ROTOBRUSH-OMEGA 42.90 EACH
3 OR MORE 37.90 EACH



RotoBrush-Titanium® Short For Wide Area

Single Use

#ROTOBRUSH-DELTA-S 42.90 EACH
3 OR MORE 37.90 EACH



MOST POPULAR

RotoBrush-Titanium® For Non-Flap Case

Single Use

#ROTOBRUSH-NANO 42.90 EACH
3 OR MORE 37.90 EACH



MOST POPULAR

RotoBrush-Titanium® For Flap Open Case

Single Use

#ROTOBRUSH-2D-LONG 42.90 EACH
3 OR MORE 37.90 EACH



*Fits Latch Type
Implant Handpiece*



*Gently Cleans
Implant Surface*

SPECIAL

Salvin RotoBrush-Titanium® Kit



MOST POPULAR



Includes:

- 2 RotoBrush-Titanium® Pocket
- 2 RotoBrush-Titanium® Nano
- 2 RotoBrush-Titanium® 2D-Long
- **FREE** Autoclavable Bur Block

#ROTOBRUSH-KIT
215.00 KIT WITH BUR BLOCK

Peri-Implantitis Tungsten Carbide Bur Kit

SPECIAL

Peri-Implantitis Bur Kit

- Tungsten Carbide Burs Clean & Smooth Titanium Surfaces
- 32mm Long FG Shanks To Reach Difficult To Access Implants
- Different Shapes & Sizes For Variety Of Individual Cases
- Includes: 4 Egg Shaped Burs & 4 Flame Shaped Burs
- Fits Friction Grip Handpieces

#PERIIMPLANTITIS-BURS **148.00 KIT W/BLOCK**



PERI-IMPLANTITIS NON-SURGICAL TREATMENT – IMPLA-SAFE PROTOCOL

1 Probe The Sulcus & Assess For Calculus Using Impla-Safe Probe To Avoid Tissue Seal Rupture



If Calculus Present



Remove Calculus Using Impla-Safe Soft Curette
(Only On Abutments – Never On Implant Threads)

2 Use Impla-Safe Cement Explorer For Detection Of Residual Cement



If Residual Cement Is Present



Use Impla-Safe Cement Remover
To Chip Off Cement

3 Subgingival Irrigation With Saline

4 Shoe-Shine Implant Surface – Use Impla-Safe Ti Micropliers To Safely Push Gauze Soaked In Saline In The Sulcus



5 The NEW Impla-Safe Rotary Brush minimizes iatrogenic damage during implant cleaning



"Impla-Safe is the first fully implant-specific instrument set for peri-implant maintenance, and for the prevention and non-surgical management of peri-implant inflammation. From the flexible probe to the soft titanium alloy instruments for implant-specific factors, such as excess cement, this set addresses the iatrogenic damage caused to implants by dental instruments and scaling using an evidence-based workflow. Should be integral in the maintenance of patients with dental implants."

Dr. George Kotsakis
Periodontist, Hoboken, NJ



"These exceptionally well crafted implant site specific instruments for detecting and removing excess cement are a brilliant innovation and must have for all implant surgically based clinicians."

Dr. Michael A. Pikos
Oral Surgeon, Palm Harbor, FL



"The Impla-Safe set is a game-changing, revolutionary implant-specific solution!"

Dr. Jennifer Hirsch Doobrow
Periodontist, Cullman, AL



"I'm excited to see an implant specific solution with instruments that adapt to the shape of implant-supported restorations and can help detect excess cement during and after cementation."

Dr. Konstantinos Chochlidakis
Prosthodontist, Rochester, NY



"The Impla-Safe Protocol and instrument set addresses what has been missing in peri-implantitis treatment: an evidence-based, implant-specific solution for our patients."

Dr. Eduardo R. Lorenzana
Periodontist, San Antonio, TX

SPECIAL


Impla-Safe™ Peri-Implant Management Set

INCLUDES:

- Impla-Safe™ Probe
- Impla-Safe™ Mirror Handle
- Box of 6 Impla-Safe™ Mirrors
- Impla-Safe™ Titanium Curette
- Impla-Safe™ Cement Explorer™
- Impla-Safe™ Cement Remover
- Impla-Safe™ Titanium Micropliers
- Impla-Safe™ Rotary Brush
- Instrument Cassette

#IMPLA-SAFE-INSTR-SET

.....715.00 SET





Peri-Implant Management Set

Designed By Dr. George Kotsakis

*For Prevention And Non-Surgical Management
Of Peri-Implantitis*



Depth Marked
At 1-15mm

**Impla-Safe™
Perio Probe**

Titanium
Flexible
Ultra Thin
#IMPLA-SAFE-PROBE
.....103.00 EACH




**Impla-Safe™
Mirror #4**

Double Sided
#IMPLA-SAFE-
MIRROR
.... 99.00 BOX OF 6

**Impla-Safe™
Mirror Handle**

#MIRROR-HNDL-BLK
.....69.50 EACH



**Impla-Safe™
Soft Tissue
Curette 3/4**

Titanium
For Supragingival
Calculus
#IMPLA-SAFE-
SOFT-TI-CU
.....103.00 EACH



**Impla-Safe™
Cement
Explorer™**

Titanium
High Modulus Of
Elasticity For Optimal
Tactile Sensitivity
#IMPLA-SAFE-
EXPLORER
.....103.00 EACH



**Impla-Safe™
Cement Remover**

Titanium
Micro-Serrated For
Efficient Cement
Removal
#IMPLA-SAFE-
CEM-REMOVER
.....103.00 EACH



**Impla-Safe™
Titanium
Micropliers -
Curved**

Titanium
#IMPLA-SAFE-
MICROPLIERS
.....134.00 EACH



**Impla-Safe™
Rotary Brush**

Package of 10
#IMPLA-SAFE-BRUSH
.....49.00 PKG

Salvin Motor Systems For Implantology & Oral Surgery

**Single Motor System****BUYER'S GUIDE - IMPLANT MOTORS**

▶	One Motor Unit
▶	300 to 40,000 RPM Brushless Motor
▶	Pre-Selectable Torque Levels 10 -- 70Ncm
▶	Memory Backup For Program Changes
▶	Autoclavable Motor & Cord
▶	Internal & External Irrigation Pump

SPECIAL**Single Motor – 1 Handpiece Package****Includes:**

Single Motor Console & 1 Motor, 1 Handpiece

MOST POPULAR

With 20:1 Handpiece (25 To 2,000 RPM)

#MOTOR-SINGLE-20-1..... **5490.00**

With 32:1 Handpiece (15 To 1,250 RPM)

#MOTOR-SINGLE-32-1..... **5490.00**

With 16:1 Handpiece (30 to 2,500 RPM)

#MOTOR-SINGLE-16-1..... **5490.00****Dual Motor System****BUYER'S GUIDE - IMPLANT MOTORS**

▶	Two Motor Unit
▶	300 to 50,000 RPM Brushless Motor
▶	Pre-Selectable Torque Levels 10 -- 70Ncm
▶	Memory Backup For Program Changes
▶	Autoclavable Motor & Cord
▶	Internal & External Irrigation Pump

SPECIAL**Dual Motor – 1 Handpiece Package****Includes:**

Dual Motor Console & 2 Motor, 1 Handpiece

With 20:1 Handpiece (25 To 2,500 RPM)

#MOTOR-DUAL-20-1..... **6790.00**

With 32:1 Handpiece (15 To 1,550 RPM)

#MOTOR-DUAL-32-1..... **6790.00**

With 16:1 Handpiece (30 to 3,125 RPM)

#MOTOR-DUAL-16-1..... **6790.00****Irrigation Tubing On Page 82**

Implant & Oral Surgery Handpieces



32:1 Reduction Push Button Latch Contra Angle

**FOR DRILLING, TAPPING & PLACEMENT
RECOMMENDED SPEED FOR MOST IMPLANT SYSTEMS**

- 10 To 1,000 RPM For Drilling, Tapping & Placement
- Internal & External Irrigation
- Autoclave Or Chemiclave Up To 275° F (135° C)

#NOU-32-1 1259.00 EACH



16:1 Reduction Push Button Latch Contra Angle

**ORDER THIS HANDPIECE FOR IMPLANT SYSTEMS
REQUIRING DRILLING SPEEDS EXCEEDING 1,250 RPM**

- 20 To 2,400 RPM For Drilling & Placement
- Internal & External Irrigation
- Autoclave Or Chemiclave Up To 275° F (135° C)

#NOU-16-1 1259.00 EACH



MOST POPULAR

20:1 Reduction Push Button Latch Contra Angle

- 15 To 1,700 RPM
- Internal & External Irrigation
- Autoclave Or Chemiclave Up To 275° F (135° C)
- Works With X-Guide

#NOU-20-1 1259.00 EACH



MOST POPULAR

1:1 40,000 RPM Straight Surgical

For Bone & Jaw Surgery

- 300 To 40,000 RPM
- External Irrigation
- Autoclave Or Chemiclave Up To 275° F (135° C)

#NOU-1-1-STRAIGHT 1209.00 EACH




1:1 40,000 RPM Angled Surgical

For Bone & Jaw Surgery

- 300 To 40,000 RPM
- External Irrigation
- Autoclave Or Chemiclave Up To 275° F (135° C)

#NOU-1-1-ANGLED 1335.00 EACH




1:2 Straight Surgical (Multiplies RPM X 2)

For Bone & Jaw Surgery

- 1,000 To 80,000 RPM
- External Irrigation
- Autoclave Or Chemiclave Up To 275° F (135° C)

#NSK-1-2-HANDPIECE 1590.00 EACH

**Handpieces On This Page Fit All
Standard "E" Type Implant Motors**



Disposable Irrigation Tubing Fits Nouvag & Salvin Implant Motor Systems

- Includes Tubing From Solution Bag To Handpiece
- Works With Internal, External Irrigation, Or Both
- Packaged Sterile For One-Time Use

#NOU-IRR-TUBE-DISP209.00 PKG OF 10
2 OR MORE193.40 PKG OF 10



Nouvag & Salvin Motor Tubing Irrigation Needle

#NOU-IRR-NEEDLE.....60.70 EACH

Fits NOUVAG & SALVIN Single & Dual Motors



Disposable Irrigation Tubing Fits W&H Implant Motor Systems

- Includes Tubing From Solution Bag To Handpiece
- Works With Internal, External Irrigation, Or Both
- Packaged Sterile For One-Time Use


#WH-IRR-TUBE-DISP209.00 PKG OF 10
2 OR MORE193.40 PKG OF 10





Fits W&H Motors Except Elcomed 100 Or 200 System



Motor & Handpiece Manufacturer's Names Are Trademarks Of Their Respective Companies.

Salvin Tri-Spade Intermediate Implant Drills Save You Money



	2.00mm(D) x 38mm(L)	#IMPL-DRILL-2.00	102.60 EACH
	2.75mm(D) x 38mm(L)	#IMPL-DRILL-2.75	102.60 EACH
	3.30mm(D) x 38mm(L)	#IMPL-DRILL-3.30	102.60 EACH
	3.50mm(D) x 38mm(L)	#IMPL-DRILL-3.50	102.60 EACH

SPECIAL

4 Or More Intermediate Drills (YOUR ASSORTMENT).....**92.80 EACH**

Salvin **POWER POINT®** Pilot Drills

Ultra Sharp For Precise Re-Angulation In Any Bone Density

- Easy To Read Depth Markings
- Chatter-Free




MOST POPULAR

Salvin POWER POINT® Ultra Sharp Implant Pilot Drills

1.3mm Diameter - Overall Length 33mm

#PILOT-DRILL 159.80 SET OF 2
2 OR MORE 148.70 SET OF 2




MOST POPULAR

Salvin POWER POINT® Ultra Sharp Implant Pilot Drills-Long

1.3mm Diameter - Overall Length 43mm

#PILOT-DRILL-LONG 159.80 SET OF 2
2 OR MORE 148.70 SET OF 2



Drill Extender & Pilot Drills



Internally Irrigated

Implant Drill Extender - Extends 16mm

- Fits Latch Type Implant Drills
- #DRILL-EXTENDER . . 80.60 EACH



2.0mm Inside
(3.0mm Outside)

2mm Trephine Pilot Drill

- Harvest Bone As You Drill
Primary Osteotomy
- #TREP-2 119.00 EACH



MOST POPULAR

Lance Pilot Drill

- 27mm Length
- #IMPL-DRILL-LANCE . . 62.80 EACH
2 OR MORE 57.20 EACH



Implant Pilot Drill - Round

- 22mm Length #6 Round
- #IMPL-DRILL-ROUND . . 20.20 EACH
5 OR MORE 17.70 EACH

SPECIAL

"Hole In One" Pilot Drill Kit

MOST POPULAR

Includes:

- 2 Salvin **POWER POINT®** Pilot Drills - 43mm
- 2 Salvin **POWER POINT®** Pilot Drills - 33mm
- 1 Lindemann Side Cutting Bur - 27mm
- 1 Round Implant Pilot Drill
- 1 Stainless Steel Bur Block

#HOLE-IN-ONE-KIT
. 415.00 KIT W/BLOCK



This Kit Fits Latch-Type Handpieces

Meisinger Diamond Saw Discs



7mm Diameter Diamond Saw Disc

- 0.3mm Disc Thickness
- Fits Latch Type Handpieces

#SAW-DISC-231DC 82.30 EACH
3 OR MORE 79.20 EACH



10mm Diameter Diamond Saw Disc

- 0.3mm Disc Thickness
- Fits Latch Type Handpieces

#SAW-DISC-231DC/100 85.70 EACH
3 OR MORE 82.30 EACH



13mm Diameter Diamond Saw Disc

- 0.3mm Disc Thickness
- Fits Latch Type Handpieces

#SAW-DISC-231DC/130 99.50 EACH
3 OR MORE 93.60 EACH



Saw Disc Kit For Latch Type Handpieces

- 7mm Saw Disc
- 10mm Saw Discs (2)
- 13mm Saw Disc

MOST POPULAR

#SAW-DISC-KIT 314.50 KIT

Oral Surgery Burs

NEW

**Lindemann Side Cutting Bur-HP**

- Fits Straight Handpieces
- 44mm Long

#LINDEMANN-HP

INDIVIDUALLY85.30 EACH

2 OR MORE79.50 EACH

**Lindemann Side Cutting Bur Non-Irrigated****MOST POPULAR**

- 27mm Long

#LINDEMANN

INDIVIDUALLY85.30 EACH

2 OR MORE79.50 EACH

- 32mm Long

#LINDEMANN-LONG

INDIVIDUALLY85.30 EACH

2 OR MORE79.50 EACH

**Debridement Bur - 25mm**

- FGXL #2 Round - 1.0mm Head

#DEBRIDE-BUR-2

INDIVIDUALLY5.90 EACH

2 OR MORE5.40 EACH

- FGXL #4 Round - 1.4mm Head

#DEBRIDE-BUR-4

INDIVIDUALLY5.90 EACH

2 OR MORE5.40 EACH

- FGXL #6 Round - 1.8mm Head

#DEBRIDE-BUR-6

INDIVIDUALLY5.90 EACH

2 OR MORE5.40 EACH

**Lindemann Side Cutting Bur-SU**

- Fits Friction Grip Handpieces
- 11mm Head Length

#LINDEMANN-SU

INDIVIDUALLY85.30 EACH

2 OR MORE79.50 EACH

**700 XXXL Implant Carbides**

- Fits Friction Grip Handpiece
- 30mm Long

#700 XXXL

SET OF 5 92.60 SET OF 5

2 OR MORE SETS... 85.80 SET OF 5

**700 XXL Implant Carbides****MOST POPULAR**

- Fits Friction Grip Handpiece
- 25mm Long

#700 XXL

INDIVIDUALLY17.10 EACH

12 OR MORE14.90 EACH

**Pikos Scoring Bur****MOST POPULAR**

- #8 Round Scoring Bur
- Depth Marked At 15mm
- Fits Straight Handpieces
- 44mm Long

#PIKOS-SCORINGBUR

INDIVIDUALLY7.80 EACH

10 OR MORE6.10 EACH

**Pikos Osteotomy Bur****MOST POPULAR**

- 702-L Fits Straight Handpieces
- 44mm Long

#PIKOS-OSTEO-BUR-S

INDIVIDUALLY7.80 EACH

10 OR MORE6.10 EACH

- 703 Fits Straight Handpieces

- 44mm Long

#PIKOS-OSTEO-BUR-L

INDIVIDUALLY7.80 EACH

10 OR MORE6.10 EACH

**Zygomatic Trough Bur**

- 2.3mm Pear Shape
- Fits Straight Handpieces
- 44mm Long

#MEIS-HM77MF-023-HP

INDIVIDUALLY 31.00 EACH

10 OR MORE 29.35 EACH

SPECIAL**Pikos Block Graft Kit****MOST POPULAR****Includes:**

- 2 Pikos 702-L Osteotomy Burs
- 2 Pikos #8 Round Scoring Burs
- 1 Pikos Large Round Block Contouring Bur
- 1 0.8mm Diameter Perforation Drill
- 1 Stainless Steel Bur Block

#PIKOS-BUR-KIT..... 153.20 KIT W/BLOCK



Special Purpose Drills & Burs



Pikos Block Contouring Bur

MOST POPULAR

- Fits Latch Type Handpieces
 - 25mm Long
- #PIKOS-BLOCK-CONTOUR-RA**
INDIVIDUALLY 51.40 EACH
2 OR MORE 46.10 EACH
- Fits Straight Handpieces
 - 44mm Long
- #PIKOS-BLOCK-CONTOUR**
INDIVIDUALLY 51.40 EACH
2 OR MORE 46.10 EACH



0.8mm Bone Perforation Decorticating Drill

- Fits Latch Type Handpieces
 - 27mm Long
- #PERF-DRILL-0.8-RA**
INDIVIDUALLY 91.00 EACH
3 OR MORE 80.90 EACH
- Fits Straight Handpieces
 - 43mm Long
- #PERF-DRILL-0.8-HP**
INDIVIDUALLY 91.00 EACH
3 OR MORE 80.90 EACH



Sinus Scoring #8 Diamond

MOST POPULAR

- Fits Latch Type Handpieces
 - 22mm Long
- #SINUS-SCORE-DIAM-RA**
INDIVIDUALLY 19.40 EACH
10 OR MORE 16.90 EACH
- Fits Straight Handpieces
 - 44mm Long
- #SINUS-SCORE-DIAM-HP**
INDIVIDUALLY 19.40 EACH
10 OR MORE 16.90 EACH



Sectioning Bur

- 703L Surgical Length Friction Grip Bur
 - 25mm Long
- #SECTION-BUR-703FG**
INDIVIDUALLY 7.60 EACH
2 OR MORE 6.50 EACH



Ridge Contouring Bone Bur

- Contour & Reshape Bone Anatomy
 - Fits Latch Type Handpiece
 - 25mm Long
- #RIDGE-CONTOUR**
INDIVIDUALLY 144.10 EACH
2 OR MORE 136.50 EACH



1.2mm Step Drill Urban Decorticating Bur

- Fits Latch Type Handpieces
 - 21mm Long
- #MEIS-2035-012-RA**
PACK OF 2 107.00 PKG/2
2 OR MORE PACKS . . . 97.60 PKG/2



Neumeyer Root Planing Bur-FG

- Fits Friction Grip Handpieces
 - 20mm Long
- #NEUMEYER-FG**
INDIVIDUALLY 18.20 EACH
10 OR MORE 15.90 EACH



Wire Passing Bur

- Fits Straight Handpieces
 - 1.5mm Diameter
 - 5cm (2") Long
- #WIRE-PASS-1.5**
INDIVIDUALLY 108.70 EACH
2 OR MORE 100.90 EACH



MetalKut™ Carbides

- Our MetalKut™ Carbides Cut Through Metal Like Butter!
 - Fits Friction Grip Handpiece
 - 20mm Long
- #MKUT-PLUS**
INDIVIDUALLY 13.70 EACH
12 OR MORE 11.70 EACH

SPECIAL

Rasmussen Biologic Shaping Bur Kit

Includes:

- 2 1.6mm Coarse Flat End Taper Friction Grip Burs
- 2 1.6mm Fine Flat End Taper Friction Grip Burs
- 1 1.0mm Round End Needle Tip Diamond
- 1 2.3mm Round Friction Grip Diamond
- 1 2.3mm Fine Football Shaped Friction Grip Bur
- 1 Stainless Steel Bur Block

#RASMUSSEN-BUR-KIT 112.00 KIT W/BLOCK



SPECIAL**Peri-Implantitis Bur Kit**

- Tungsten Carbide Burs Clean & Smooth Titanium Surfaces
- 32mm Long FG Shanks To Reach Difficult To Access Implants
- Different Shapes & Sizes For Variety Of Individual Cases
- Fits Friction Grip Handpieces

Includes:

- 4 Egg Shaped Burs
- 4 Flame Shaped Burs

#PERIIMPLANTITIS-BURS **148.00 KIT W/BLOCK**

**SPECIAL****P.E.T. Partial Extraction Therapy Bur Kit - FG****FG Shanks**

Fits Friction Grip Handpieces

Includes:

- Short 1.2mm Diameter Bone Cutter
- Long 1.2mm Diameter Bone Cutter
- Coarse 1.0mm Round Diamond Bur
- 1.8mm Round Prep Bur
- 2.3mm Round Prep Bur
- 2.9mm Round Prep Bur
- Stainless Steel Bur Block

#PET-BUR-KIT **148.00 KIT W/BLOCK**

**SPECIAL****P.E.T. Partial Extraction Therapy Oral Surgery Bur Kit - HP****HP Shanks**

Fits Straight Handpieces

Includes:

- 701L Bur HP
- 1703 Bur HP
- 1.8mm Round Diamond HP
- 3.1mm Round Bur HP
- Lindemann Bur HP
- 2.3mm Pointed Football Diamond HP
- Stainless Steel Bur Block

#PET-OS-BUR-KIT **148.00 KIT W/BLOCK**



SPECIAL

Degranulation Bur Kit

- Multiuse Burs (Typically 100 Procedures)
- Designed For Quick & Effective Bone Debridement
- Removes Remnants Of Tissue Attached To Bone
- Fits Latch Type Implant Handpieces

Includes:

- 1.0mm Degranulation Bur
- 1.8mm Degranulation Bur
- 2.5mm Degranulation Bur
- 3.0mm Degranulation Bur
- 3.5mm Degranulation Bur

#DEGRANULATION-KIT **242.00 KIT W/BLOCK**



SPECIAL

Alveoplasty Bur Kit

- Smooth & Re-Contour Bone Following Tooth Extraction To Provide Ideal Platform For Denture Preparation Or Full Arch Restoration
- Prepare Bone For Implant Placement

Includes:

- 2 Lindemann Bone Cutters (Fits Friction Grip Handpieces)
- 2 Initial Starter Drills (Fits Friction Grip Handpieces)
- 1 Round Diamond (Fits Straight Handpieces)
- 5 Carbide Alveoplasty Burs With Round Ends To Protect Tissue (Fits Straight Handpieces)

#MEIS-ALV18 **298.00 KIT W/BLOCK**



SPECIAL

Crown Removal Bur Kit

- For Cutting Crowns & Bridges
- Includes Burs For Every Type Of Restorative Material
- Recommended Cutting Angle Of 45°
- Optimum Speed Of 160,000rpm
- Fits Friction Grip Handpieces

#MEIS-BC04 **149.40 KIT W/BLOCK**



SALVIN®
RotoBrush-Titanium®
For Exposed Implant Threads
On Page 77




Salvin Long Trephine Drills

For Failed Implant Retrieval & Bone Harvesting



Internally Irrigated

-  2.0mm Inside Diameter x 22mm Barrel
#TREPH-2 **119.00 EACH**
-  3.5mm Inside Diameter x 22mm Barrel
#TREPH-3.5 **119.00 EACH**
-  4.25mm Inside Diameter x 22mm Barrel
#TREPH-4 **119.00 EACH**
-  5.0mm Inside Diameter x 22mm Barrel
#TREPH-5 **119.00 EACH**
-  6.0mm Inside Diameter x 22mm Barrel
#TREPH-6 **119.00 EACH**
-  8.0mm Inside Diameter x 22mm Barrel
#TREPH-8 **119.00 EACH**

SPECIAL

Long Trephine Kit

MOST POPULAR



Internally Irrigated

Set Of 6 Long Trephines With **FREE** Bur Block

- Fits Latch Type Handpieces
- Marked At 8, 10, 13 & 15mm

#TREPH-KIT

..... **597.00 KIT W/BLOCK**

Salvin Short Trephine Drills

For Bone Harvesting



-  2.0mm Inside Diameter x 10mm Barrel
#TREPH-SHORT-2 **119.00 EACH**
-  4.0mm Inside Diameter x 10mm Barrel
#TREPH-SHORT-4 **119.00 EACH**
-  6.0mm Inside Diameter x 10mm Barrel
#TREPH-SHORT-6 **119.00 EACH**
-  8.0mm Inside Diameter x 10mm Barrel
#TREPH-SHORT-8 **119.00 EACH**
-  10.0mm Inside Diameter x 10mm Barrel
#TREPH-SHORT-10 **119.00 EACH**
-  12.0mm Inside Diameter x 10mm Barrel
#TREPH-SHORT-12 **119.00 EACH**

SPECIAL

Short Trephine Kit



Set Of 6 Short Trephines With **FREE** Bur Block

- Fits Latch Type Handpieces
- Marked At 2, 4, 6, 8 & 10mm

#TREPH-SHORT-KIT

..... **597.00 KIT W/BLOCK**

Autoclavable Bur Blocks & Baskets



MOST POPULAR

6 Bur Block - Latch & Friction Grip

#BUR-BLOCK-6-RA.	60.40 EACH
3 OR MORE	56.20 EACH



MOST POPULAR

6 Bur Block - Straight Handpiece

BUR-BLOCK-6-HP.	60.40 EACH
3 OR MORE	56.20 EACH



8 Bur Block - Latch & Friction Grip

#BUR-BLOCK-8-RA.	60.40 EACH
3 OR MORE	56.20 EACH



12 Bur Block - Latch & Friction Grip

#BUR-BLOCK-12-RA.	60.40 EACH
3 OR MORE	56.20 EACH



Autoclavable Bur/Parts Basket

46mm Diameter


#BUR-BASKET	75.00 EACH
3 OR MORE	71.40 EACH



12 Compartment Autoclavable Parts Cassette

21cm x 6cm x 2cm

#CASS-SS-PARTS-12.	103.90 EACH
3 OR MORE	97.70 EACH




Instrument Cassettes On Page 66

SALVIN AutoTome®-Osteotome


Automatic Osteotome Kit



- *Automatic Osteotome For Crestal Sinus Lifts & Maxillary Bone Spreading*
- *Ensures Precision And Saves Time*
- *Enhanced Patient Safety And Comfort*
- *One-Hand Operation For Visibility & Ease Of Use*
- *Creates An Osteotomy For Most Major Implant Systems*



*Atraumatic, Controlled
Bone Preparation*



*Gentle Auto Tapping
Easily Advances
Osteotomy While
Spreading Bone*




*Fits Your
Existing
Implant Motor*



*Set Your Motor To
12,000-15,000 RPM
Translates To 120-150
Gentle Taps Per Minute
(Max 15,000 RPM)*

SPECIAL

AutoTome® Osteotome Starter Kit With Offset Tips




Includes:

- 1 AutoTome® Handpiece
- 1 Manual Handpiece For Use With Mallet When Desired
- 8 Osteotome Tips (4 Offset Concave, 4 Offset Convex)
(2.0-2.8mm, 2.5-3.3mm, 3.0-3.9mm, 3.5-4.5mm Diameters)
- 1 Sterilization Cassette

#AUTOTOME-OST-STARTER-KIT **2295.00 KIT**

SPECIAL

AutoTome® Combination Osteotome & Extraction Kit



Includes:

- 1 AutoTome® Handpiece
- 1 Manual Handpiece
- 16 Osteotome Tips
- Sterilization Cassette
- 6 Periotome Tips
- 1 Autoclavable Tip Block


#AUTOTOME-COMP-KIT **3750.00 KIT**




Conventional Osteotomes On Pages 136-137 & 142-145

SALVIN
AutoTome®-Extract
Automatic Extraction Kit

- Automatic Periotome For Root Elevation & Extraction
- Ensures Precision And Saves Time
- Enhanced Patient Safety And Comfort
- One-Hand Operation For Visibility & Ease Of Use
- Gentle Tapping Facilitates Atraumatic Controlled Extractions



Bone Preserving
Atraumatic Extraction



Gentle Auto Tapping
Easily Facilitates
Controlled Extractions



Fits Your
Existing
Implant Motor

SPECIAL



AutoTome® Extraction Kit With Handpiece

Includes:

- 1 AutoTome® Handpiece
- 6 Periotome Tips
- 1 Autoclavable Tip Block

#AUTOTOME-EXTRACT..... **1795.00 KIT**

SPECIAL



AutoTome® Extraction Tip Kit

(If You Already Have An AutoTome® Handpiece)

Includes:

- 6 Periotome Tips
- 1 Autoclavable Tip Block

#AUTOTOME-EXTRAC-TIPS..... **855.00 KIT**

Replacement AutoTome® Extraction Tips

Straight Narrow - 2.1mm	#AUTOTOME-EXT-TIP-1	145.00 EACH
Straight Wide - 2.6mm	#AUTOTOME-EXT-TIP-2	145.00 EACH
Angled Vertical - 2.1mm	#AUTOTOME-EXT-TIP-3	145.00 EACH
Angled Horizontal - 2.1mm	#AUTOTOME-EXT-TIP-4	145.00 EACH
Straight Narrow Concave - 2mm	#AUTOTOME-EXT-TIP-5	145.00 EACH
Straight Medium Concave - 3mm	#AUTOTOME-EXT-TIP-6	145.00 EACH



Conventional Periotomes On Pages 100-103

"Rear-Exhaust" 45° High Speed Surgical Handpiece

Creates Maximum Access & Visibility On 3rd Molars

- No Air, Just Water In Surgical Field - Air Exits At Rear Of Handpiece For Safety
- High Speed: 380,000 - 450,000rpm
- High Torque: 21W Of Cutting Power



Push-Button Bur Release Uses
Standard Friction Grip Burs



45° High Speed Handpiece With Swivel 4 Hole Connector

New Version With Higher Torque

#SURG-45-PKG. 1210.00 EACH

Handpiece Without Connector

New Version With Higher Torque

#SURG-45. 995.00 EACH



💡 Fiber Optic 45° High Speed Handpiece With Swivel 5 Hole Fiber Optic Connector

New Version With Higher Torque

#SURG-45-FO-PKG. 1425.00 EACH

Handpiece Without Connector

New Version With Higher Torque

#SURG-45-FO. 1237.00 EACH



💡 Fiber Optic Fits Kavo® MULTiflex® Connector

#SURG-45-FO-K 1237.00 EACH

Lubricant For 45° High Speed Handpiece



Pana Spray Lubricant - 10.5 oz.

Available In US Only

#PANA-SPRAY 79.20 EACH

**Remove Failed Implants Or
Broken Abutment Screws**



SPECIAL

Implant Rescue Kit

#IMPLANT-RESCUE-KIT **1840.00 KIT**




SPECIAL

Narrow Implant Broken Screw Kit

For Broken Screw Removal On Small Diameter &
Bone Level Implants

#IMPLANT-RESCUE-NARROW **860.00 KIT**



Broken Abutment Screw Removal

- Helps To Remove Broken Screws Without Damaging The Implant
- Includes Centering Guides To Fit Most Implant Systems



Place Guide On
Implant & Hold
With Stabilizing
Handle.

Set Motor To
REVERSE At
1,000-1,250RPM &
50-70Ncm Torque.



Drill In **REVERSE**
Using "Up & Down"
Motion To Prepare
1-2mm Deep
Dimple Into Top Of
Broken Screw.

*Use Irrigation To
Minimize Heat.*



Suction To Remove
Metal Shavings.
Then Insert Tap Into
Implant Handpiece

Set Motor To
REVERSE At
70-80RPM &
50-70Ncm Torque.



Insert Tap Into
The 1-2mm Dimple
In The Top Of The
Broken Screw.

Use The Tap
In **REVERSE**
To Remove
Broken Screw.



Broken Screw
Removed.

Failed Implant Removal

- Removes Failed Implants Without Compromising Surrounding Bone Or Damaging Neighboring Structures
- Works With Any Threaded Implant System
- Kit Includes Small, Medium And Large Fixture Removers



Insert Fixture
Remover Screw
Into Hand Adapter.



Turn Counter-
Clockwise Into
The Implant.



Use Wrench
Counter-Clockwise
Until The Implant
Becomes Mobile.



Remove Failed
Implant.



For Difficult Cases:
Trephine Top 3-4mm
Of Bone.

Irrigating Syringes & Catheters



Sterile Irrigating Syringes

- 60cc Luer Lock Syringes
- Sterile - Individually Packaged

#SYRINGE-STERILE-60CC

BOX OF 30 76.90 BOX
2 OR MORE BOXES 72.20 BOX



Sterile Irrigating Catheters

- Sterile Blunt Catheters
- 18 Ga x 1" (2.54cm)
- Individually Packaged

#BLUNT-CATH-18G-2.5CM

BOX OF 100 129.00 BOX
2 OR MORE BOXES 122.50 BOX



Monoject #412 Irrigating Syringes

- 12cc Disposable Syringes
- Curved Tip
- Packaged Non-Sterile




#MONOJECT-412

BOX OF 50 61.50 BOX
2 OR MORE BOXES 57.50 BOX

IrriGUIDE ADAPTER

Developed by Dr. Thomas Jackson

- *Compatible With Any Guided Dental Implant System*
- *Ensures Continuous Direct Water Flow To The Osteotomy Site*
- *Minimizes The Risk Of Thermal Damage To Surrounding Tissues*
- *Improves Visibility During Implant Procedures By Eliminating Obstructions Within The Surgical Field*
- *Eliminates The Need To Stop And Start The Drilling Process For Irrigation*



SPECIAL

IrriGUIDE Adapter Starter Set

Includes:

- IrriGUIDE Adapters - 5 pack
- IrriGUIDE Seating Bur
- IrriGUIDE Trimming Bur
- Cyanoacrylate glue
- Assembly Instructions

#IRRIGUIDE-STARTER-KIT
250.00



External Irrigation Wand Systems

- Available With Or Without Suction
- Uses Pump Up Pressure Cuff
- Includes Autoclavable Silicone Irrigation Tubing & IV Bag Striker Insert

Scan For

Suction Irrigator
Technique Guide



Fingertip Irrigation Control



Included Pressure Cuff
Uses 500mL Saline Bag



MOST POPULAR

Wand Irrigator - Without Suction

Complete Kit

#IRR-SET 216.00 SET



Wand Irrigator - With Suction

Complete Kit Including Built In Suction Aspirator

Order Suction Tubing Separately

#SUCT-IRR-SET 321.00 SET



Replacement Irrigation Tubing

#SILICONE-IRR-TUBING

. 22.90 FOR 7FT



Wand Irrigator With Suction Connects
To 1/4" Disposable Suction Tubing Or
Salvin Autoclavable Suction Tubing
On Pages 96-97

Stainless Steel Bowls



Solution Bowl

- Heavy Duty Stainless Steel (20 oz)
- 5" x 2 3/8" (12.7 x 6.1cm)

#SOLN BOWL 5 36.30 EACH
2 OR MORE 33.70 EACH



MOST POPULAR

Prep Cup

- Heavy Duty Stainless Steel (6 oz)
- 3 1/2" x 2" (8.9 x 5.1cm)

#PREP CUP 25.50 EACH
2 OR MORE 23.20 EACH

Autoclavable Tubing Adapter

Use Autoclavable Surgical Tubing With Your Regular High Volume Suction



Latexfree

Connects To Your
Regular Dental Suction



MOST POPULAR

Autoclavable Tubing Adapter

- Autoclavable Hundreds Of Times
- Non-Crimping "U" Adapter Connects To Regular Dental Suction
- 5'(152.4cm) x 1/4"(0.6cm) ID Silicone Surgical Tubing
- Fits All Aspirator Tips

#TUBING ADAPT 50.90 EACH
2 OR MORE 48.80 EACH

● 1/4"(0.6cm) Inside Diameter

Autoclavable Silicone Surgical Suction Tubing

Remains Flexible...Autoclavable Hundreds Of Times



Latexfree

MOST POPULAR

Silicone Tubing By The Foot


1/4"(0.6cm) Internal Diameter

#TUBING..... 4.80 FOOT
50 FT OR MORE 4.20 FOOT



End To End Connector

#TUBE CONNECT. . . . 17.50 EACH
3 OR MORE 15.90 EACH



MOST POPULAR

Autoclavable Cut Off Clamp

#CUT..... 11.40 EACH
3 OR MORE 9.90 EACH



Connects High
Volume Surgical
Aspirators To
Surgical Tubing

High Volume Adapter

Connects High Volume Surgical
Aspirators To Surgical Tubing
#HVE NECK ADAP . . . 50.10 EACH

Autoclavable Suction Cut Off Valves

- *Smooth On/Off Lever Action*
- *Easy To Remove & Clean*



MOST POPULAR

Flip Saliva Ejector Valve

- Fits 7/32" Saliva Ejector Tubing
- Autoclavable Silicone Head
- Includes Cleaning Brush

#FLIP-SE 71.20 EACH



Flip-D Cut Off Valve

- Fits Standard 1/2" Tubing
- Includes Cleaning Brush

#FLIP-D 71.20 EACH



Flip-A Cut Off Valve

- Fits Standard 1/2" Tubing
- Includes Cleaning Brush

#FLIP-A 71.20 EACH

SALVIN
Rapid Roll®



MOST POPULAR

- Easy To Use Cover Ring Slides Over Long Hoses In Seconds
- Individually Packaged - 65" (165cm) Long

#RAPID-ROLL

BOX OF 25 **121.50 BOX**
2 OR MORE BOXES **112.10 BOX**



Disposable Sterile Suction Tubing & Yankauer Suction Tips



MOST POPULAR

6' Sterile Suction Tubing

- 6' (182.9cm) x 1/4" (0.6cm) Inside Diameter
- Individually Packaged

#SUCTION-TUBE-6

BOX OF 50 **192.00 BOX (2.99 EACH)**
2 OR MORE **181.30 BOX (2.81 EACH)**



Disposable Sterile Yankauer Tips

- Blue Translucent Plastic Minimizes Patient Visibility Of Fluids
- Flange Tip & Vent

#YANKAUER-DISPOSABLE

BOX OF 50 **141.00 BOX**
2 OR MORE **131.50 BOX**

10' Sterile Suction Tubing

- 10' (304.8cm) x 1/4" (0.6cm) Inside Diameter
- Individually Packaged

#SUCTION-TUBE-10

BOX OF 50 **203.00 BOX (3.19 EACH)**
2 OR MORE **191.30 BOX (2.99 EACH)**




MOST POPULAR

Surgical Tubing To Dental Suction Adapter

Adapts Surgical Tubing To Standard Perio & General Dentistry High Volume Suction


#TUBETIP **29.90 EACH**

Surgical Aspirators



	1.5mm Surgical Aspirator
#15P1A.....	42.60 EACH
	2.5mm Surgical Aspirator
#15P2A.....	42.60 EACH
	MOST POPULAR
	3mm Surgical Aspirator
#15P3A.....	42.60 EACH
	MOST POPULAR
	4mm Surgical Aspirator
#16RSA.....	42.60 EACH
	<i>Recommended Cleaning Brushes:</i>
	<i>Medium & Small</i>
	2mm Modified Frazier
#FR2A.....	55.30 EACH
	3mm Modified Frazier
#FR3A.....	55.30 EACH
	MOST POPULAR
	4mm Modified Frazier
#FR4A.....	55.30 EACH
	<i>Recommended Cleaning Brushes:</i>
	<i>Small & Medium</i>
	<i>Aspirators On This Page Fit</i>
	<i>Surgical Tubing</i>
	<i>Or</i>
	<i>Standard High Volume Suction</i>

Autoclavable Cleaning Brushes



#BRUSH-MEDIUM.....	41.20 DOZEN
#BRUSH-SMALL.....	41.20 DOZEN
	MOST POPULAR
#BRUSH-X-SMALL.....	41.20 DOZEN

Silicone Tip Aspirator Kit

- Autoclavable Silicone Tips Prevent Metal Contact With Implant
- Tips Are Easily Modified With Scissors Or Knife



MOST POPULAR

Silicone Tip Aspirator Kit

Stainless Steel Aspirator &
1 Trimmable Silicone Tip

#46YST 42.60 EACH


Autoclavable Silicone Tips

Set Of 3 Trimmable Silicone Tips

#YST 34.00 SET OF 3

Recommended Cleaning Brushes:
Medium

No-Clog Surgical Aspirators



No-Clog Aspirator

Built-In Clog Ejector
Sliding Button Forward Ejects
The Clog


(2.5MM TIP)
#B2 NO CLOG 139.00 EACH

MOST POPULAR

(3.0MM TIP)
#B3 NO CLOG 139.00 EACH

Recommended Cleaning Brushes:
Small & Medium

Micro Aspirators For Endo/Root Tip Suctioning



1mm Micro/Endo Aspirator

Bendable Tubing
Includes Cleaning Stylet

#02BA. 99.00 EACH

MOST POPULAR

2mm Micro/Endo Aspirator

Bendable Tubing
Includes Cleaning Stylet

#03BA. 99.00 EACH

Recommended Cleaning Brushes:
X-Small

Micro/Endo Suction Handle



MOST POPULAR

Micro/Endo Suction Handle

Connects To Saliva Ejector Suction
Includes 10 Non-Sterile
Blunt Catheters

#SUCT-LL 55.60 EACH

Sterile Luer Lock Catheters

#BLUNT-CATH-18G-2.5CM
. 129.00 BOX OF 100



Non-Sterile Luer Lock Catheters

#SUCT-LL-TIP-14GA 42.60 BOX OF 100
#SUCT-LL-TIP-18GA 42.60 BOX OF 100
#SUCT-LL-TIP-20GA 42.60 BOX OF 100



SALVIN

PDL-Evator®

Precise Tip Luxating Elevator

Combines The Benefit Of A Periotome & Elevator In One Instrument**MOST POPULAR** PDL-Evator® Straight - 2.5mm Tip

2.5mm Narrow Tip Is Thin Enough To Insert In Tight Spaces
Between The Root And The Socket Wall

#PDL-STRAIGHT 155.20 EACH

**MOST POPULAR** PDL-Evator® Curved

#PDL-CURVED 155.20 EACH

**MOST POPULAR** PDL-Evator® Bayonet

#PDL-BAYONET 155.20 EACH



PDL-Evator® Angled "In" Small

#PDL-IN-SMALL 155.20 EACH



PDL-Evator® Angled "Out" Small

#PDL-OUT-SMALL 155.20 EACH



PDL-Evator® Angled "In" Large

#PDL-IN-LARGE 155.20 EACH



PDL-Evator® Angled "Out" Large

#PDL-OUT-LARGE 155.20 EACH

SPECIAL

SALVIN

PDL-Evator®

Precise Tip Luxating Elevator

MOST POPULAR


Includes:

- Set Of 7 Precise Tip Luxating Elevators
- **FREE** Stainless Steel Cassette

#PDL-PKG

1067.00 SET OF 7 W/CASSETTE

"Interchangeable Tip" Periotomes



MOST POPULAR

**Periotome Handle With
5 Interchangeable Tips**

#PERIOTOME-5-TIPS295.00 SET



*Use BUR-BLOCK-6-HP On Page 89
To Easily Store Periotome Tips*

Thin Straight Replacement Tip

#PERIOTOME-TIP-1 132.00 PKG OF 3



Wide Straight Replacement Tip

#PERIOTOME-TIP-2 132.00 PKG OF 3



Thin Angled Replacement Tip

#PERIOTOME-TIP-3 132.00 PKG OF 3



Flat Angled Replacement Tip

#PERIOTOME-TIP-4 132.00 PKG OF 3



Flat Wide Angled Replacement Tip

#PERIOTOME-TIP-5 132.00 PKG OF 3



Third Molar Elevators



Third Molar Elevator - Left

#ELEV-THIRD-MOLAR-L. 112.00 EACH

Third Molar Elevator - Right

#ELEV-THIRD-MOLAR-R. 112.00 EACH

Periotomes

**MOST POPULAR**

Straight Periotome With Fixed Straight Tip


#PERIOTOME-STR 109.10 EACH
3 OR MORE 103.40 EACH



Straight Periotome With Fixed Rounded Tip

#PERIOTOME-STR-ROUND 109.10 EACH
3 OR MORE 103.40 EACH

Serrated Tip



Straight Periotome With Serrated Tip


#PERIOTOME-SERR-STR 129.00 EACH
3 OR MORE 121.60 EACH



Lane Elevator

Designed By Dr. Bill Lane

#LANE-ELEVATOR 107.00 EACH
3 OR MORE 100.10 EACH




Periosteal-Tome™

Periosteal Elevator/
Periotome Combo

#PERIOSTEAL-TOME 111.60 EACH
3 OR MORE 105.20 EACH

Double-Ended Periotomes



MOST POPULAR

Double Ended Periotome #1

#PERIOTOME-DE-1	107.50 EACH
3 OR MORE	99.80 EACH



Double Ended Periotome #2

#PERIOTOME-DE-2	107.50 EACH
3 OR MORE	99.80 EACH



Serrated Tips

Serrated Tip Periotome #1

#PERIOTOME-SERR-1	129.00 EACH
3 OR MORE	121.60 EACH



Serrated Tips

Serrated Tip Periotome #2

#PERIOTOME-SERR-2	129.00 EACH
3 OR MORE	121.60 EACH



Serrated Tips



Serrated Tips

SALVIN®
AtrauLux®
 Ergonomic Elevators

*Handle Fits Comfortably In Palm
 For Increased Control*



MOST POPULAR

AtrauLux® Elevator
 2mm Straight

#ATRAU-ELEV-2-STR
98.00 EACH

SALVIN
Ergonomic
 Handle



MOST POPULAR

AtrauLux® Elevator
 3mm Straight

#ATRAU-ELEV-3-STR
98.00 EACH

SALVIN
Ergonomic
 Handle



AtrauLux® Elevator
 4mm Straight

#ATRAU-ELEV-4-STR
98.00 EACH

SALVIN
Ergonomic
 Handle



AtrauLux® Elevator
 2mm Curved

#ATRAU-ELEV-2-CUR
98.00 EACH

SALVIN
Ergonomic
 Handle



AtrauLux® Elevator
 3mm Curved

#ATRAU-ELEV-3-CUR
98.00 EACH

SALVIN
Ergonomic
 Handle



AtrauLux® Elevator
 4mm Curved

#ATRAU-ELEV-4-CUR
98.00 EACH

SALVIN
Ergonomic
 Handle

SPECIAL

Set Of 6
 Salvin AtrauLux®
 Elevators

MOST POPULAR



#ATRAU-ELEV-SET-6
 514.00 SET OF 6

Serrated Elevators



**Serrated
Elevator 2.2mm**

#DEL-SERR-2
..... 106.00 EACH



MOST POPULAR
**Serrated
Elevator 3mm**

#DEL-SERR-3
..... 106.00 EACH



**Serrated
Elevator 4mm**

#DEL-SERR-4
..... 106.00 EACH



**Serrated
Elevator 4mm
- Left**

#DEL-SERR-4-L
..... 106.00 EACH



**Serrated
Elevator 4mm
- Right**

#DEL-SERR-4-R
..... 106.00 EACH

SPECIAL

MOST POPULAR

Set Of 5 Bayonet Elevators #DEL-SERR-SETOF5 **442.00 SET OF 5**

Atraumatic Bayonet Elevators



MOST POPULAR
**Bayonet
Elevator
Straight**

#BAYONET-ELEV-1
.... 106.00 EACH



**Bayonet
Elevator
Offset "In"
- Small**

#BAYONET-ELEV-2
.... 106.00 EACH



**Bayonet
Elevator
Offset "Out"
- Small**

#BAYONET-ELEV-3
.... 106.00 EACH



**Bernard Bayonet
Elevator Offset
- Left**

#BERNARD-L
..... 106.00 EACH



**Bernard Bayonet
Elevator Offset
- Right**

#BERNARD-R
..... 106.00 EACH

SPECIAL

MOST POPULAR

Set Of 5 Bayonet Elevators #BAYONET-SETOF5 **442.00 SET OF 5**

Bone Preserving Atraumatic Luxating Elevators



3mm Luxating Elevator Curved

#LUX-3-CURVED
..... 99.80 EACH



3mm Luxating Elevator Straight

#LUX-3-STRAIGHT
..... 99.80 EACH



5mm Luxating Elevator Curved

#LUX-5-CURVED
..... 99.80 EACH



5mm Luxating Elevator Straight

#LUX-5-STRAIGHT
..... 99.80 EACH

SPECIAL

Set Of 4 Luxating Elevators



MOST POPULAR

#LUX-SETOF4
346.00 SET OF 4

"Anatomically Grooved" Luxating Elevators



MOST POPULAR
2mm Grooved Luxating Elevator

#LUX-GROOVE-2
..... 99.80 EACH



3mm Grooved Luxating Elevator

#LUX-GROOVE-3
..... 99.80 EACH



4mm Grooved Luxating Elevator

#LUX-GROOVE-4
..... 99.80 EACH



5mm Grooved Luxating Elevator

#LUX-GROOVE-5
..... 99.80 EACH

SPECIAL

Set Of 4 Grooved Elevators



#LUX-GROOVE-SETOF4
346.00 SET OF 4

Luxating Elevators With Periotome Tip



MOST POPULAR
3mm Luxating Elevator/Periotome
#LUX-TOME-3
..... 99.80 EACH



4mm Luxating Elevator/Periotome
#LUX-TOME-4
..... 99.80 EACH



5mm Luxating Elevator/Periotome
#LUX-TOME-5
..... 99.80 EACH

SPECIAL

Set Of 3 Luxating Elevators

Three small images of the 3mm, 4mm, and 5mm Luxating Elevators/Periotomes.

#LUX-TOME-SETOF3
259.50 SET OF 3

Thin Handle Luxating Elevators



MOST POPULAR
2.5mm Thin Handle Luxating Elevator
#LUX-THIN-2.5
..... 99.80 EACH



3.5mm Thin Handle Luxating Elevator
#LUX-THIN-3.5
..... 99.80 EACH

SPECIAL

Set Of 2 Luxating Elevators

Two small images of the 2.5mm and 3.5mm Thin Handle Luxating Elevators.

#LUX-THIN-SETOF2
167.00 SET OF 2

Salvin Twist Elevators

- Twist Design Provides Improved Vertical Movement Of The Tooth Over Conventional Elevators
- Sharp Tips Engage Tooth & Provide Deeper Subgingival Access



Twist Elevator Straight Left

#TWIST-ELEV-STR-L
..... 117.90 EACH



Twist Elevator Straight Right

#TWIST-ELEV-STR-R
..... 117.90 EACH



Twist Elevator Angled Left

#TWIST-ELEV-ANG-L
..... 117.90 EACH



Twist Elevator Angled Right

#TWIST-ELEV-ANG-R
..... 117.90 EACH



#TWIST-ELEV-SET-4
396.00 SET OF 4

Elevators



Potts Elevator-Left

#POTTS-L
.....96.00 EACH



Potts Elevator-Right

#POTTS-R
.....96.00 EACH



Potts "T-Handle" Elevator-Left

#POTTS-T-L
..... 102.00 EACH



Potts "T-Handle" Elevator-Right

#POTTS-T-R
..... 102.00 EACH



Coupland Gouge Elevator

Useful For Removal Of Teeth
With Amputated Crowns
#COUPLAND-1C . . .102.30 EACH

Elevators



Your Choice Of Elevators

INDIVIDUALLY . . . 97.00 EACH
3 OR MORE 91.50 EACH


SALVIN®
AtrauLux®
 Extraction Forceps

- *Tapered Tips Adapt Exactly To The Tooth With Excellent Subgingival Contact*
- *Lightweight, Ergonomic Design Fits Both Right & Left Handed Practitioners*

SPECIAL**MOST POPULAR**

Set Of 9 Salvin AtrauLux® Forceps

#ATRAU-FORCEP-SET-9 . . . **2215.00 SET OF 9**



AtrauLux® Forcep Cow Horn

Serrated Beak
 Lower Atraumatic Extraction

#ATRAU-FORCEP-COWHORN . . **294.00 EACH**

**MOST POPULAR**

AtrauLux® Forcep #1

Serrated Beak
 Upper Anterior Atraumatic Extraction

#ATRAU-FORCEP-1 **294.00 EACH**

**MOST POPULAR**

AtrauLux® Forcep #2

Serrated Beak
 Lower Anterior Atraumatic Extraction

#ATRAU-FORCEP-2 **294.00 EACH**

**MOST POPULAR**

AtrauLux® Forcep #3

Serrated Beak
 Upper Premolar Atraumatic Extraction

#ATRAU-FORCEP-3 **294.00 EACH**



AtrauLux® Forcep #4

Serrated Beak
 Lower Premolar Atraumatic Extraction

#ATRAU-FORCEP-4 **294.00 EACH**



AtrauLux® Forcep #5

Lower Root Atraumatic Extraction

#ATRAU-FORCEP-5 **294.00 EACH**

**MOST POPULAR**

AtrauLux® Forcep #6

Upper Root Atraumatic Extraction


#ATRAU-FORCEP-6 **294.00 EACH**



AtrauLux® Forcep #7

Serrated Beak
 Lower Molar Atraumatic Extraction

#ATRAU-FORCEP-7 **294.00 EACH**



AtrauLux® Forcep #8

Serrated Beak
 Lower Atraumatic Extraction

#ATRAU-FORCEP-8 **294.00 EACH**

DEEP GRIP Extended Tip Length Serrated Forceps

BUYER'S GUIDE

▶

Superior Grip & Tactile Feel Even With Limited Contact

▶

Less Time & Force Required, Even For Difficult Extractions

Tips Marked At
1mm - 8mm Depths



SPECIAL

MOST POPULAR

Set Of 4 DEEP GRIP Forceps

#DGA-SET-4 922.00 SET OF 4

DEEP GRIP
Extended Tip
Serrated Forcep

DGA-34

Serrated Beak
Upper Atraumatic Extraction

#DGA-34..... 294.00 EACH

DEEP GRIP
Extended Tip
Serrated Forcep

MOST POPULAR

DGA-35

Serrated Beak
Upper Atraumatic Extraction

#DGA-35..... 294.00 EACH

DEEP GRIP
Extended Tip
Serrated Forcep

DGA-36

Serrated Beak
Lower Atraumatic Extraction

#DGA-36..... 294.00 EACH

DEEP GRIP
Extended Tip
Serrated Forcep

MOST POPULAR

DGA-37

Serrated Beak
Lower Atraumatic Extraction

#DGA-37..... 294.00 EACH

Bayonet Atraumatic Extraction Forceps




MOST POPULAR

Bernard Forcep 153

Upper Atraumatic Extraction

#BERNARD-153 230.00 EACH




MOST POPULAR

Bernard Forcep 146

Lower Atraumatic Extraction

#BERNARD-146 230.00 EACH



German Stainless Steel Extraction Forceps



DEF 150AS

Serrated Beak
Upper Anteriors, Split Beaks


#DEF150AS. 219.80 EACH



DEF 151AS

Serrated Beak
Lower Anteriors, Split Beaks

#DEF151AS. 219.80 EACH



DEF 222

Lower Molars
Excellent Access To 3rd Molars

#DEF-222 219.80 EACH



MOST POPULAR

DEF 23-Lower

Lower Cow Horn
1st & 2nd Lower Molars


#DEF 23 219.80 EACH



DEF 23-Upper

Upper Cow Horn
1st & 2nd Upper Molars

#DEF-23-UPPER 219.80 EACH



Wrap Around Extraction Forcep

Right Angle "Wrap Around" Design
Provides Superior Access

Serrated Beak

#DEF-68-WRAP 219.80 EACH

Tungsten Carbide Tip Extraction Forceps – Better Grip & Tactile Feel – Less Trauma


BUYER'S GUIDE

▶ **Superior Grip & Tactile Feel Even With Limited Contact**

▶ **Less Time & Force Required, Even For Difficult Extractions**



TUNGSTEN CARBIDE



MOST POPULAR

DEF 150 TC

Carbide Serrated Beak
Upper Incisors, Cuspids,
Bicuspid & Roots, Universal

#DEF 150 TC 303.00 EACH



TUNGSTEN CARBIDE



MOST POPULAR

DEF 151 TC

Carbide Serrated Beak
Lower Incisors, Cuspids,
Bicuspid & Roots, Universal

#DEF 151 TC 303.00 EACH



TUNGSTEN CARBIDE



DEF 15TD TC

Carbide Serrated Beak
Upper Incisors, Cuspids

#DEF 15TD TC . . . 303.00 EACH



TUNGSTEN CARBIDE



DEF 17 TC

Carbide Serrated Beak
1st & 2nd Lower Molars,
Universal

#DEF 17 TC 303.00 EACH



TUNGSTEN CARBIDE




DEF 53L TC

Carbide Serrated Beak
1st & 2nd Upper Molars, Left

#DEF 53L TC 303.00 EACH




TUNGSTEN CARBIDE



DEF 53R TC

Carbide Serrated Beak
1st & 2nd Upper Molars, Right

#DEF 53R TC 303.00 EACH



TUNGSTEN CARBIDE




DEF 65 TC

Carbide Serrated Beak
Upper Incisors,
Roots & Fragments

#DEF 65 TC 303.00 EACH



TUNGSTEN CARBIDE




MOST POPULAR

DEF 69 TC

Carbide Serrated Beak
Upper & Lower
Fragments & Roots

#DEF 69 TC 303.00 EACH



TUNGSTEN CARBIDE



DEF 150S TC

Carbide Serrated Beak
Pedo, Upper Primary
Teeth & Roots, Universal

#DEF 150S TC . . . 303.00 EACH



TUNGSTEN CARBIDE




DEF 151S TC

Carbide Serrated Beak
Pedo, Lower Primary
Teeth & Roots, Universal


#DEF 151S TC . . . 303.00 EACH

Tungsten Carbide Tipped Root Tip Forceps

**MOST POPULAR****Root Tip-Lower**

Carbide Serrated Beak
Lower Root Tips & Molars

#ROOT CU 303.00 EACH

**MOST POPULAR****Root Tip-Upper**

Carbide Serrated Beak
Upper Root Tips & Molars

#ROOT TC. 303.00 EACH

"Fish-Hook" Atraumatic Extraction Root Forceps


45° Angled Root Forcep
"Fish-Hook" Serrations

Upper Atraumatic Extraction

#FHA-45 288.00 EACH


90° Angled Root Forcep
"Fish-Hook" Serrations

Lower Atraumatic Extraction

#FHA-90 288.00 EACH

SPECIAL

#FISH-HOOK-SET-2 529.00 SET OF 2

Forceps & Root Tip Picks



Tungsten Carbide Tip
20mm Tip Length
For Deep Access

TUNGSTEN CARBIDE



MOST POPULAR

TC Coated Tip Root Fragment Forcep

Extended Tips For Excellent
Access & Gripping

#FRAGMENT-FORCEP-TC 251.40 EACH
2 OR MORE 237.90 EACH



Warwick Root Tip Elevator - Left

#WARWICK-L
. 83.00 EACH



Warwick Root Tip Elevator - Right

#WARWICK-R
. 83.00 EACH

SPECIAL

**Set Of
Warwick Root Tip
Elevators**

**#WARWICK-SET
138.60 SET OF 2**



Heidbrink H1

Single-Ended

#HEID 1
. 76.50 EACH



MOST POPULAR

Heidbrink H2

Single-Ended

#HEID 2
. 76.50 EACH



MOST POPULAR

Heidbrink H3

Single-Ended

#HEID 3
. 76.50 EACH



Heidbrink H2-H3

Double-Ended

#HEID 2/3
. 76.50 EACH



Meisinger Benex Pro Atraumatic Extraction Kit


- *Screw Action Lifting Makes Thin Buccal Plate Extractions Easy!*
- *Excellent For Use In The Aesthetic Zone*




"I use the Benex Extractor on a daily basis in my office for difficult extractions to save time and preserve bone for implant placement."

—Dr. John Russo


Periodontist, Sarasota, FL



Drill Hole In Root



Insert Extraction Screw



Turn Knob To Gently & Safely Extract Root

SPECIAL


Salvin Benex Atraumatic Kit

MOST POPULAR

Includes:

- Benex Extraction Handle
- Drills, Extraction Screws, Pull Wires
- Organizer Holds Both Parts & Handle
- Support Tray To Stabilize Handle If Required

#BENEX-PRO2520.00 COMPLETE KIT



Replacement Parts For Benex Pro Atraumatic Extraction Kit



Narrow Extraction Screw

1.3mm x 10mm

#MEIS-BE001-060

..... 122.90 EACH

1.3mm x 16mm

#MEIS-BE001-230-070

..... 122.90 EACH



Wide Extraction Screw

2.1mm x 10mm

#MEIS-BE001-230-080

..... 122.90 EACH



Traction String

#MEIS-BE001-230020

..... 166.80 EACH



**Salvin
POWER POINT®
Pilot Drills
On Page 83**

Mallets & Chisels



Cottle Curved Chisels

Depths Marked At
6, 9, 12, 14 & 17mm

MOST POPULAR

#COTTLE-4 (4MM WIDE)

#COTTLE-6 (6MM WIDE)

INDIVIDUALLY
..... 139.00 EACH
3 OR MORE
..... 131.40 EACH



Sheehan Osteotomes

Depths Marked At
6, 9, 12, 14 & 17mm


#SHEE 2 (2MM WIDE)

MOST POPULAR

#SHEE 4 (4MM WIDE)

#SHEE 6 (6MM WIDE)

INDIVIDUALLY
..... 139.00 EACH
3 OR MORE
..... 131.40 EACH



Mead Mallet

#MEAD
..... 173.00 EACH



TG Perio Chisel

#PERIO-CHISEL-TG
..... 110.80 EACH
3 OR MORE
..... 105.10 EACH



MOST POPULAR

TG-O Perio Chisel

#PERIO-CHISEL-TGO
..... 110.80 EACH
3 OR MORE
..... 105.10 EACH



MOST POPULAR


Rhodes Chisel

17.2cm
Back Action
#RHODES
..... 109.00 EACH

Serrated Lucas Bone Currettes

Serrated Tips Help Debride Extraction Sockets & Create Bleeding Points For Faster Healing


Designed By Dr. Jin Kim



Serrated Lucas Bone Currettes

INDIVIDUALLY	84.60 EACH
3 OR MORE	78.70 EACH

Non-Serrated Lucas Currettes



Non-Serrated Lucas Bone Currettes

INDIVIDUALLY	79.50 EACH
3 OR MORE	73.80 EACH

ZAGA Zygomatic Surgical Instruments

Designed By Dr. Carlos Aparicio – Creator Of The ZAGA Method



"It is my pleasure to share this great set of instruments designed to sequentially facilitate the placement of zygomatic implants following the ZAGA philosophy."


Dr. Carlos Aparicio
Periodontist, Barcelona, Spain



Salvin ZAGA 1 - Zygoma Dissector

Detaches The Periosteum Beyond The Angle Formed By Temporal & Frontal Anatomy Of The Zygoma To Prepare For The Anchorage Of The Retractor

#ZAGA-1
..... 138.30 EACH



Salvin ZAGA 2 - Zygoma Retractor

- Increases Field Of Vision
- Directs The Drilling
- Prevents Perforation Of The Orbit
- Protects Soft Tissue While Drilling Superior Cortex

#ZAGA-2
..... 217.40 EACH



Salvin ZAGA 3 - Maxilla Dissector

Detaches The Periosteum From Posterior Part Of The Maxillary Wall & Detaches The Floor Of The Nose


#ZAGA-3
..... 138.30 EACH



Salvin ZAGA 4 - Posterior Dissector

Finishes Detachment Of The Periosteum In More Posterior Superior Sectors

#ZAGA-4
..... 138.30 EACH



Salvin ZAGA 5 - Nasal Protector

Protects The Mucosa Of The Nose In Case Of Perforation Of The Nasal Floor

#ZAGA-5
..... 138.30 EACH



Zygomatic Trough Bur

- 2.3mm Pear Shape
- Fits Straight Handpieces
- 44mm Long

#MEIS-HM77MF-023-HP
INDIVIDUALLY . . 29.90 EACH
10 OR MORE . . . 28.35 EACH

SPECIAL

Salvin ZAGA Instrument Set

Includes:

- Zygoma Dissector
- Zygoma Retractor
- Maxilla Dissector
- Posterior Dissector
- Nasal Protector
- Plus **FREE** Salvin 10" x 6" Deep Instrument Cassette

#ZAGA-INSTR-KIT
..... **745.00 SET**



Graft Material Packers - Titanium

MOST POPULAR

Titanium Graft Packer / Placement Site Depth Gauge

- Titanium For Biocompatibility
- Excellent In Difficult Access Sites
- 3.3mm & 4.0mm Diameter Tips
- Marked At 8, 11, 13, 15 & 18mm

#GRAFT-PACKER-TI 138.40 EACH



MOST POPULAR

Serrated Titanium Tip Graft Packer-Small

- Excellent In Difficult Access Sites
- 3mm & 4mm Diameter Tips


#GRAFT-PACK-SERR-S 138.40 EACH



Serrated Titanium Tip Graft Packer-Large

- Excellent In Difficult Access Sites
- 6mm & 8mm Diameter Tips

#GRAFT-PACK-SERR-L 138.40 EACH



Graft Material Packers - Stainless Steel

Graft Packer

- Excellent In Difficult Access Sites
- 6mm & 8mm Diameter Smooth Tips


#GRAFT-PACKER 116.20 EACH



Pikos Graft Packer-Angled

- Concave Tips
- Excellent For Use In Sinus
- 4mm & 6mm Diameter Tips

#PIKOS-PACKER-ANGLED 130.90 EACH



Curved Graft Packers For Pockets Around Implants - Titanium

MOST POPULAR

Curved Titanium Pocket Packer 1.75mm Radius

3mm x 1mm Tip

#POCKET-PACKER-1.75

..... 106.50 EACH



Curved Titanium Pocket Packer 2.5mm Radius

4mm x 1mm Tip

#POCKET-PACKER-2.5

..... 106.50 EACH



Graft Material Transfer Instruments

Graft Material Spoon With 4mm Serrated Graft Packer

Concave Spoon Tip For Larger Grafts

#BONE-SPOON-PACK-4 126.00 EACH



MOST POPULAR

Graft Material Spoon With 3mm Serrated Graft Packer

Concave Spoon Tip For Larger Grafts

#BONE-SPOON-PACKER 126.00 EACH



Sinus Floor Elevator With 3mm Graft Packer

- Marked At 7, 10, 12, 14, 16 & 19mm
- Designed By Dr. Richard Martin

#MARTIN-PLUGGER. 126.00 EACH



Stainless Steel Graft Material Syringes

MOST POPULAR

5mm Diameter Stainless Steel Syringe-Curved

#GRAFTSYRINGE-CURVED 112.50 EACH

3.6mm Diameter Stainless Steel Syringe-Curved

#GRAFT-SYR-S-CURV 112.50 EACH

2.8mm Diameter Stainless Steel Syringe-Curved

#GRAFT-SYR-2.8-CU 112.50 EACH

Load From Tip To
Prevent Clogging



Scan For
60th Second Product Video™




Autoclavable Plastic Syringes For Smaller Grafts & Perio Defects

Perio Defect Graft Syringes

- Disassembles For Easy Cleaning
- Autoclavable, Chemiclavable

#HA SYRINGE 66.90 SET OF 2
2 OR MORE SETS 61.10 SET OF 2

Load From Tip To
Prevent Clogging



Single Use Plastic Bone Graft Delivery Syringes

Bone Graft Delivery Syringes

- Packaged Sterile
- Single Use
- Blunt 5mm Wide Tip For Ease Of Delivery

#BONE-SYRINGE-1CC 51.90 PACK OF 10
2 OR MORE PACKS 49.30 PACK OF 10

Load From Tip To
Prevent Clogging



Blunt 5mm Wide Tip
For Ease Of Delivery



Sold As A Pack Of 10 Syringes

Stainless Steel Graft Material Syringes




3.5mm Diameter Ring Handle Syringe-Curved
#GRAFT-SYRINGE-3.5. 213.00 EACH



Holds Approximately 1.6 cc


MOST POPULAR

4.5mm Diameter Ring Handle Syringe-Curved
#GRAFT-SYRINGE-4.5. 213.00 EACH



Holds Approximately 2.1 cc

6.0mm Diameter Ring Handle Syringe-Curved
#GRAFT-SYRINGE-6.0. 213.00 EACH



Holds Approximately 2.8 cc



Graft Material Mixing Dishes



MOST POPULAR

Graft Material Mixing Cups

- Deep Dish With Rounded Bottom
- Holds 4 To 6 Grams Of Hydrated Material
- Autoclavable, Chemiclavable
- Will Not Crack During Sterilization

#GRAFT MIX. 55.50 SET OF 2
2 OR MORE SETS 51.20 SET OF 2



Stainless Steel Weighted Graft Material Dish

- Deep Dish With Weighted Bottom To Prevent Tipping
- Holds 20 Grams Of Hydrated Material
- Autoclavable, Chemiclavable

#GRAFT-DISH 60.80 EACH
2 OR MORE 56.70 EACH

OsseoScape™ Angled Single Use Bone Harvester**MOST POPULAR**

- Angled Lightweight Ergonomic Handle
- Collection Chamber Holds Approximately 1cc
- Autoclave Before Use

#OSSEO-SCRAPE 185.00 BOX OF 5
 (Hex Driver Tip & Thumb Knob Sold Separately)



.05" (1.25mm) Hex Long - 12mm

#HTLT 60.90 EACH



Thumb Knob

#TK 82.60 EACH

Sterile Single Use Bone Harvester

- Angled Lightweight Ergonomic Handle
- Collection Chamber Holds Approximately 1cc
- Packaged Sterile

#BONE-SCRAPER-STERILE 52.40 EACH
 (Hex Driver Tip & Thumb Knob Sold Separately)



Angled Tip For Easy Access




*Collection Chamber Holds
Approximately 1cc*



*For A New Sharp Edge,
Use A 1.25mm Hex Driver Tip
To Turn The Blade 90°*

Salvin Stainless Steel Bone Scraper



Stainless Steel Bone Scraper

- Clean & Autoclave For Repeated Use
- No Blade To Change

#BONE-SCRAPER 191.00 EACH

Stainless Steel Bone Scraper SS With Replaceable Blade

- Collects And Stores Approximately 0.5cc Of Harvested Bone
- Blade Can Reattach For Unlimited Volume Per Case
- Low Cost – Replace Blade, Not Entire Graftor

#BONE-SCRAPE-SS. 187.20 EACH



Bone Scrape SS Blade

- Easily Replaceable Blade

#BONE-SCRAPE-SS-BLADE 45.50 EACH



*Thin Profile Head
Provides Easy Access*



*Replaceable Blade For
Economy Of Use*




*Collects Approximately
0.5cc Of Bone*

ACM – Autobone Chip Maker

- *Collect 1cc Of Autogenous Bone In 10-20 Seconds*
- *Stopper Prevents Over Drilling & Helps Contain Bone Chips*



Easily Operated By Mounting On A Contra-Angle Handpiece



- Drilling Speed: 300rpm
- Available In Four Diameters:
4.0mm, 4.5mm, 5.0mm, 6.0mm
- Includes:
One Drill, Seven Stoppers

#ACM-4.0	225.00
#ACM-4.5	225.00
#ACM-5.0	225.00
#ACM-6.0	225.00



Osseous Collection Trap

- *Stop Suctioning Usable Bone!*
- *Simplify Bone Harvesting With Our Osseous Trap*



MOST POPULAR



- Uses Any Aspirator That Fits Into A High Volume Suction
- Fits Regular High Volume Suction As Well As Silicone Or Latex Surgical Suction Tubing



3mm Stainless Steel Aspirator Tip

#15P3A	42.60 EACH
--------	------------

Osseous Coagulum Trap With 3 Single Use Filters

- Use "In Line" With Your Suction
- Osseous Material Is Collected Directly Into A Filter
- Can Be Chemically, Autoclaved, Or Dry Heat Sterilized

#OCT	136.00 EACH
------	-------------



Set Of 3 Replacement Single Use Filters

- Supplied Non-Sterile – Autoclave Before Using

#OCT FILTERS	.27.80 PKG OF 3
3 OR MORE SETS	.26.10 PKG OF 3

Urban Master Core Kit

Simple And Safe Extraction Of Bone Cylinders For Bone Augmentation

Developed By Dr. Istvan Urban

Includes:

- Depth Limited Trephines
- Depth Markings For Additional Depth Cont
- Black Coating To Prevent Glare

Master Core Kit

- 6 7mm Trephines Of 3mm-8mm Lengths
- 5 5mm Trephines Of 4mm-8mm Lengths
- 3 3mm Trephines Of 6mm-8mm Lengths
- 1 3mm Trephine Of 2mm Length
- 1 5mm Trephine Of 2mm Length
- 1 7mm Trephine Of 2mm Length
- 1 Drill Extender

#MEIS-BMC00 1257.00 KIT



Master Core Kit - Basic

- 5 7mm Trephines Of 3mm-7mm Lengths
- 5 5mm Trephines Of 3mm-7mm Lengths
- 1 5mm Trephine Of 2mm Length
- 1 7mm Trephine Of 2mm Length
- 1 Drill Extender

#MEIS-BMCBA 1120.00 KIT



Autogenous Bone Grafting Bur Kits



Bone Chip Maker Bur Kit With Stops

4 Bone Harvesting Trephines (3.4mm, 3.9mm, 4.2mm, 5.2mm), 4 Trephine Stops (5mm, 7mm, 10mm, 11.5mm), 1 Bone Remover & 1 Latch To Square Adapter

#BONE-CHIP-MAKER339.00 KIT



Bone Collector Bur Kit

5 Bone Collecting Drills (2.8mm, 3.3mm, 3.8mm, 4.3mm, 4.8mm), 5 Drill Stops (6mm, 8mm, 10mm, 11.5mm, 13mm)

#BONE-COLLECTOR-KIT339.00 KIT

Tissue Rakes & Flap Holding Forceps

For Greater Flap Release & Bone Graft Site Access



*Separate Fibers
For Increased
Flap Release*

Tissue Rake #1

Designed By Dr. Hom-Lay Wang
For Increased Flap Extension

#TISSUE-RAKE-1103.40 EACH



*Separate Fibers
For Increased
Flap Release*

Tissue Rake #2

Designed By Dr. Hom-Lay Wang
For Increased Flap Extension

#TISSUE-RAKE-2103.40 EACH




*Retract Flap For
Easy Access*

Flap Holding Forcep

For Easy Bone Graft Site Access


#FLAP-FORCEP109.50 EACH

Bone Block Clamps



Screw Fits Through Notched Tip

Photos: Dr. Yvan Poitras



Anterior Cortical Block Clamp

- Ratchet Lock Holds Block In Place
- 19.7cm Long

#BONE-BLOCK-CLAMP371.00 EACH



*Tightens To Hold
Block In Place*



Adjustable Straight Cortical Block Clamp

- Adjustable Width To Hold Block In Place
- 10cm Long

#BONE-CLAMP-ADJ311.00 EACH

Hand-Held Bone Morselizer

MOST POPULAR

Bone Morselizer

- Collected Bone Chips Ready To Use For Augmentation
- No Loss Of Bone Material
- Easy To Clean
- Adjustable Screw Allows You To Achieve Preferred Particle Size

#BONE-MORSELIZER 629.00 EACH



Bone Mill For Fast & Easy Crushing Of Cortical Bone

Stainless Steel Bone Mill

- Easy Assembly & Cleaning
- No Loss Of Bone Material
- Bone Collected In Easy To Use Dish
- Collected Bone Chips Ready To Use For Bone Augmentation

#BONE-MILL-PLUS 1995.00 EACH



Titanium Bone Cruncher™

Bone Cruncher™

1. Place Harvested Bone In Well Of Instrument
2. Cover Bone With Top, Serrated Side Down
3. Use Mallet To Reduce To Desired Particle Size

#BONE-CRUNCHER 349.50 EACH

Mead Mallet


#MEAD 173.00 EACH





Salvin Titanium **Tenting** Screws

Photo: Dr. Nick Shumaker

5mm Diameter Head Provides Broad Surface Area
During GBR Tenting Procedures

Dental Latch Friction Fit Screwdriver 
#FIX-DRIVER-LATCH. 110.70 EACH

Drill Bit-Dental Latch (1.1mm) 
#FIX-DRILL-RA 110.70 EACH



Screws Have 4mm Of Thread
For Stabilization In Bone

Available With 3, 4 or 5mm
Unthreaded Surface For
Each Unique Case

**Tenting****1.5mm Diameter Tenting Screws - 5mm Diameter Head**

1.5mm x 7mm Tenting Screw	#TENT-1.5X7	55.30 EACH
1.5mm x 8mm Tenting Screw	#TENT-1.5X8	55.30 EACH
1.5mm x 9mm Tenting Screw	#TENT-1.5X9	55.30 EACH

SPECIAL**Salvin **Tenting** Screw Kit****Includes:**

- 12 Tenting Screws (Your Assortment)
- Bone Fixation Screwdriver
- Latch Screwdriver
- Latch Drill Bit
- Autoclavable Organizer Block

#TENT-KIT-W/DRIVER. **1120.00 KIT**

*Screwdriver Tips Are Proprietary To Manufacturer's Screw System And Are Not Interchangeable
The Salvin Fixation Screwdriver Works With Both Our Fixation Screws And Our Tenting Screws*

Salvin Titanium **Fixation** Screws






Photo: Dr. Yvan Poitras


*Use To Secure Autogenous
or Allograft Bone Blocks In
The Mandible Or Maxilla*

Dental Latch Friction Fit Screwdriver 
#FIX-DRIVER-LATCH. 110.70 EACH

Drill Bit-Dental Latch (1.1mm) 
#FIX-DRILL-RA 110.70 EACH

Lag Drill Bit-Dental Latch (1.5mm) 
Oversize Drill Bit For Drilling Enlarged Hole In Graft
#FIX-DRILL-LAG 110.70 EACH


Drill Bit-Straight Handpiece (1.1mm)
#FIX-DRILL-HP 110.70 EACH


Lag Drill Bit-Straight Handpiece (1.5mm)
Oversize Drill Bit For Drilling Enlarged Hole In Graft
#FIX-LAG-DRILL-HP 110.70 EACH



*Locking Tip Driver Keeps Screw Secure
Until You're Ready To Release It*



Block Fixation

1.5mm Diameter Screws - 3mm Diameter Head

1.5mm x 4mm Bone Screw	#FIX-1.5X4-PT	55.30 EACH
1.5mm x 6mm Bone Screw	#FIX-1.5X6-PT	55.30 EACH
1.5mm x 8mm Bone Screw	#FIX-1.5X8-PT	55.30 EACH
1.5mm x 10mm Bone Screw	#FIX-1.5X10-PT	55.30 EACH
1.5mm x 13mm Bone Screw	#FIX-1.5X13-PT	55.30 EACH
1.5mm x 15mm Bone Screw	#FIX-1.5X15-PT	55.30 EACH



Block Fixation

2.0mm Diameter Emergency Screws - 3mm Diameter Head

2.0mm x 4mm Bone Screw	#FIX-2.0X4-PT	55.30 EACH
2.0mm x 6mm Bone Screw	#FIX-2.0X6-PT	55.30 EACH
2.0mm x 8mm Bone Screw	#FIX-2.0X8-PT	55.30 EACH
2.0mm x 10mm Bone Screw	#FIX-2.0X10-PT	55.30 EACH
2.0mm x 13mm Bone Screw	#FIX-2.0X13-PT	55.30 EACH
2.0mm x 15mm Bone Screw	#FIX-2.0X15-PT	55.30 EACH

SPECIAL

Salvin Block Fixation Screw Kit



or



Includes:

- 18 Fixation Screws (Your Assortment)
- Latch Screwdriver
- Bone Fixation Screwdriver
- Drill Bit (Latch or Straight HP)
- Lag Bit (Latch or Straight HP)
- Autoclavable Organizer Block

#FIX-KIT-SPECIAL 1475.00 KIT

Meisinger Master-Pin-Control Membrane Fixation Kit

- **Fixation Of Resorbable Membranes, Non-Resorbable Membranes & Titanium Mesh**
- **Hybrid Screw & Pin Easily Removed With Hand Screwdriver**

Designed By Dr. Istvan Urban

Includes:

- 2 1.8mm Latch-Type Tri-Spade Pilot Drills
- 2 0.8mm Latch-Type Twist Perforation Drills
- 2 0.6mm Latch-Type Twist Perforation Drills
- 1 Placement Instrument
- 1 Fixation Handle
- 1 Hand Screwdriver



Basic Kit With 10 Pins

#MEIS-BMPBA 1210.00 KIT

Master Kit With 34 Pins

#MEIS-BMP00 1735.00 KIT

Replacement Pins - 3.5mm Length

10 Pins, Non-Sterile

#MEIS-MP10. 235.40 PKG OF 10



Meisinger Master-Pin-X-Control Upgrade

Upgrade Your Existing Master-Pin Kit With The Master-Pin-X

- Enhanced Pin Design: Biocompatible Titanium Alloy (Grade 5) For Secure, Precise Fixation Of Resorbable Membranes, Non-Resorbable Membranes & Titanium Mesh – Even In Hard-To-Reach Areas
- Self Tapping Thread Allows Pin To be Inserted Easily And Precisely Into Cortical Bone
- Torx Interface For Easy And Safe Transfer and Optimum Force During Insertion

Designed By Dr. Istvan Urban

Upgrade Kit With 20 Master-Pin-X Pins

Includes:

- 20 MP30 Master-Pin-X Pins
- Pin Holder To Fit Into Your Existing Kit
- 1 Master-Pin-X Fixation Handle

#MEIS-BMPXU 1182.20 KIT



Replacement Master-Pin-X – 3.2mm Length

5 Master-Pin-X Pins, Non-Sterile

#MEIS-MP30. 185.25 PACK OF 5



SPECIAL



Salvin Membrane Tack Kit With Mallet & Placement Instrument

MOST POPULAR

Includes:

- 25 Tacks
- Autoclavable Tack Block
- Tack Placement Instrument
- Tack Mallet
- Tack Drill For Dense Cortical Bone

#TACK-KIT-TI **720.00 KIT**



SPECIAL



Salvin Membrane Tack Kit With PowerTack™ Handle

Includes:

- 20 Tacks
- Autoclavable Tack Block
- PowerTack™ Delivery Handle
- Tack Forcep

#POWERTACK-KIT **1120.00 KIT**



If You Already Have Our Regular Tack Kit

Includes:

PowerTack™ Delivery Handle, Autoclavable Tack Block

#POWERTACK **973.00 HANDLE & BLOCK**

Meisinger "Pin Control" Membrane Fixation Kit



Use Tacks To Secure Membrane Over Surgical Site



- Available In 3mm & 5mm Lengths
- Packaged 5 Tacks To A Vial, Non-Sterile
- Head Of Tack Is 2.5mm Diameter

3MM LENGTH

#TACK-3-PT **74.10 PKG OF 5**

5MM LENGTH

#TACK-5-PT **74.10 PKG OF 5**

Membrane Placement Instruments



Membrane Forcep With Ring Tip

Non-Serrated For Membranes

#MEMB-RING-FORC

.....126.00 EACH



MOST POPULAR

Membrane Forcep "Paddle Tip"

Non-Serrated For Membranes

#MEMB-PADDLE-FORCEP

.....93.30 EACH



Membrane Forcep - Long "Paddle Tip"

Non-Serrated For Membranes

#MEMB-PADDLE-LONG

.....93.30 EACH



MOST POPULAR

SuperCut Iris Curved 11cm

Ideal For Delicate Membranes

#SC IRIS CU

.....147.40 EACH



MOST POPULAR

Membrane Placement Instrument

Micro-Point Tip For Precise Placement & Manipulation Of Membranes

#MEMB-PLACER.86.60 EACH



Membrane Positioner

Stabilize & Condense

Designed By Dr. Larry Towning

#MEMB-CONDENSER133.40 EACH

Comprehensive Lateral & Crestal Sinus Kit

Minimizes Risk Of Membrane Perforation



Complete Crestal & Lateral Sinus Kit

#SINUS-COMP-KIT1645.00 EACH



6 Crestal Drills



8 Drill Stops



4 Crestal Burs



2 Lateral
Trephined Drills



2 Lateral Dome
Shaped Burs



Irrigation Syringe
With Tubing



8 Irrigation Tips



Syringe Holder



1 Sinus Depth Gauge/Packer



Replacement Irrigation Syringe With Tubing

#SINUS-MAS2-01-099.70 EACH



Lateral Reamer 6.5mm x 2mm

#SINUS-REAM-6.5X2209.00 EACH



Lateral Reamer 6.5mm x 3mm

#SINUS-REAM-6.5X3209.00 EACH



*Diamond Tip
Minimizes Risk
Of Perforation*



*Trephined Bur
Cuts & Detaches
Bony Wall*



*Dome Shape
Minimizes Risk
Of Perforation*



*Irrigation Lifts
Membrane*



*Syringe Holder Helps Control
Irrigation Flow*

Concave Tip Sinus Lift Osteotomes

A conventional osteotomy is made to the floor of the sinus. The Sinus Osteotome is inserted into the osteotomy & gently tapped apically to separate & lift the cortical bone at the sinus floor 2mm-4mm.

*Adjustable
Depth Stops*



*Tips Marked
8-20mm*



— 8mm
— 10mm
— 13mm
— 15mm
— 17mm
— 20mm

*Concave Tips
Allow Introduction
Of Graft Material*



MOST POPULAR

Offset Sinus Lift Set

- 5 Offset Osteotomes With Adjustable Stops
- FREE Stainless Steel Cassette

#OST-SINUS-OFF-PKG 779.00 SET



Straight + Offset Sinus Lift Set

- 5 Straight Osteotomes With Adjustable Stops
- 5 Offset Osteotomes With Adjustable Stops
- FREE Stainless Steel Cassette

#OST-SINUS-BOTH-PKG 1399.00 SET



Straight Sinus Lift Set

- 5 Straight Osteotomes With Adjustable Stops
- FREE Stainless Steel Cassette

#OST-SINUS-STR-PKG 779.00 SET

Extended Length Osteotomes For Surgical Guides

Designed By Dr. Jack Krauser

Adjustable Depth Stops

Length & Markings Accommodate Surgical Guide Tubes

— 8mm
— 10mm
— 13mm
— 16mm
— 18mm
— 20mm
— 23mm
— 26mm

Extended Length To Fit NobelGuide™

*NobelGuide™ is a
Trademark of NobelBiocare*

2.0mm
3.5mm
4.3mm
5.0mm
6.0mm

Extended Length Osteotome Set

- 5 Offset Osteotomes With Adjustable Stops
- FREE Stainless Steel Cassette


#GUIDE-TOME-PKG 784.00 SET

— 2.0mm
— 4.0mm
— 6.5mm
— 8.0mm
— 10.0mm
— 12.0mm
— 14.0mm
— 16.0mm
— 18.0mm

SALVIN
AutoTome®
Automatic Osteotome Kit
On Page 90

Meisinger Crestal Sinus Lift Kit

- **Depth Stops Prevent Puncturing Of The Sinus Membrane**
- **Concave Tipped Drills**
- **Crestal Drills Collect Bone Chips**
- **Can Be Used With Any Implant System**



Pilot Drill

Lift Membrane

Collect Bone

Depth Gauge

Place Graft



Basic Crestal Sinus Lift Kit 2 **MOST POPULAR**


4 Crestal Drills (3.3, 3.6, 3.8, 4.1mm), 7 Drill Stops (3-9mm), 1 Pilot Drill & 1 Bone Graft Condenser

#MEIS-BCLBA2 1027.00 KIT

Complete Crestal Sinus Lift Kit

6 Crestal Drills, 11 Drill Stops, 4 Pilot Drills, 1 Cylinder Bur, 2 Condensers & 1 Depth Gauge

#CRESTAL-LIFT-CONTROL 2056.00 KIT



Depth Gauge For Crestal Sinus Lift Kit

12mm Depth Gauge

#MEIS-CL021 159.40 EACH

Internal Sinus Manipulation Set

Minimally Invasive Sinus Floor Augmentation

Designed By Dr. Jason M. Yamada



Small Rounded Head Bone Graft Spreader & Packer

#ISM-1 136.40 EACH



Large Rounded Head Bone Graft Spreader & Packer

#ISM-2 136.40 EACH



Sinus Lift Instrument -
Keep Sharp Working End Against Bony Floor

#ISM-3 136.40 EACH



Sinus Lift Instrument -
Keep Sharp Working End Against Bony Floor

#ISM-4 136.40 EACH



Membrane Perforation Repair Instrument

#ISM-5 136.40 EACH

SPECIAL

Internal Sinus Manipulation Set




Set Of 5 Instruments With
FREE Stainless Steel Cassette

#ISM-KIT
599.00 SET

Kerrison Rongeurs For Widening The Sinus Window

Low Profile Head
For Visibility
& Safe Access



Kerrison Rongeur - 4mm Bite

10cm Long Beak

#KERRISON 649.00 EACH



MOST POPULAR

Kerrison Rongeur - 3mm Bite

8cm Long Beak

#KERRISON-SMALL 649.00 EACH



Serrated Tip Sinus Window Bone Tampers

To Precisely Separate Sinus Window From Surrounding Bone

MOST POPULAR

4mm Sinus Window / Bone Tamper

#BONE-TAMPER-4MM. 87.00 EACH



7mm Sinus Window / Bone Tamper

#BONE-TAMPER-7MM. 87.00 EACH



SPECIAL



Set Of 2 Bone Tampers

#BONE-TAMPER-SET 154.00 SET OF 2

Sinus Scoring Diamonds



Sinus Scoring Diamond-Latch Type

Fits Latch Type Implant Handpieces

#SINUS-SCORE-DIAM-RA

INDIVIDUALLY 19.40 EACH

10 OR MORE 17.00 EACH



Sinus Scoring Diamond-HP Type

Fits Straight Implant & Oral Surgery Handpieces

#SINUS-SCORE-DIAM-HP


INDIVIDUALLY 19.40 EACH

10 OR MORE 17.00 EACH

Sinus Lift Instruments



Sinus Lift Instruments



SPECIAL

Salvin Basic Sinus Set

MOST POPULAR



Basic Set Includes:

- 4mm Bone Tamper
- Small Serrated Graft Packer
- 4.5mm Graft Syringe
- Salvin Sinus Curette #1
- Salvin Sinus Curette #5
- Salvin Sinus Curette #6
- Salvin Sinus Freer
- Salvin Sinus Freer Curved
- (Or Your Choice Of Any 5 Curettes)
- Plus **FREE** Cassette

#SINUS-KIT-WINDOW
869.00 SET



Double Sided Ridge Splitting / Expanding Chisels - Offset

Depths Marked At
8, 10, 13, 15 & 18mm



#RIDGE-SPLIT-OFF-3MM155.30 EACH

#RIDGE-SPLIT-OFF-4MM155.30 EACH

#RIDGE-SPLIT-OFF-5MM155.30 EACH

#RIDGE-SPLIT-OFF-6MM155.30 EACH


**Ridge Splitting Offset Chisel Set**

- Set Of 4 Ridge Splitting Offset Chisels
- Stainless Steel Cassette

#RIDGE-SPLIT-OFF-SET 593.00 SET

Double Sided Ridge Splitting / Expanding Chisels - Straight

Depths Marked At
8, 10, 13, 15 & 18mm



#RIDGE-SPLIT-3MM155.30 EACH

#RIDGE-SPLIT-4MM155.30 EACH

#RIDGE-SPLIT-5MM155.30 EACH


#RIDGE-SPLIT-6MM155.30 EACH

**MOST POPULAR****Ridge Splitting Chisel Set**

- Set Of 4 Ridge Splitting Chisels
- Stainless Steel Cassette


#RIDGE-SPLIT-SET 593.00 SET

Single Bevel Ridge Splitting / Expanding Chisels



Meisinger Mandibular Ridge Splitting Kit

Horizontal Widening Up To 5mm



MOST POPULAR

Comprehensive Ridge Splitting Kit

3 Spreaders, 5 Burs, 3 Saw Discs, 2 Thumb Drivers & 1 Ratchet

#CREST-CONTROL 1945.00 KIT

Bone Expander Kit



- Alternative To Osteotomes For The Expansion & Condensing Of The Atrophic Mandible & Maxilla
- Use With Ratchet Wrench Or Low Speed Handpiece
- Improve Stability, Maintain Bone Density & Increase Fixation
- Can Be Used With Implants Up To 16mm Long
- Contains 5 Bone Expanders, 1 Ratchet Wrench, 1 Pilot Drill, 1 Ratchet Wrench Connector, 1 Latch-Type Connector, 1 Thumb Knob, 1 Saw Disc

#BONE-EXPANDER-KIT-2760.00 KIT

Splitting / Spreading Combination Kit

Mandibular Ridge Splitting + Maxillary Bone Spreading Kit



- Comprehensive Ridge Splitting System Components
- Comprehensive Bone Spreading System Components
- Large Organizing Cassette
- Plus **FREE** Motor Driver Attachments For Spreaders

TOTAL PRICE OF COMPONENTS IS 2705.00 **SAVE 200.00**

#MEIS-CREST-SPLIT2630.00 KIT



Bone Spreading Osteotomes For Implant Placement In Soft Bone



MOST POPULAR

Offset Bone Spreading Set

- 6 Offset Spreaders With Adjustable Stops
- FREE Stainless Steel Cassette

#OST-SPREAD-OFF-PKG . . . 928.00 SET



Straight + Offset Bone Spreading Set

- 6 Offset Spreaders With Adjustable Stops
- 6 Straight Spreaders With Adjustable Stops
- FREE Stainless Steel Cassette

#OST-SPREAD-BOTH-PKG 1735.00 SET




Straight Bone Spreading Set

- 6 Straight Spreaders With Adjustable Stops
- FREE Stainless Steel Cassette

#OST-SPREAD-STR-PKG . . . 928.00 SET

Rongeurs

**MOST POPULAR****Friedman Rongeur**

#FRIED

..... 309.00 EACH

**Micro Friedman
Right Angle Rongeur**


#MICRO FRIED R

..... 309.00 EACH

**MOST POPULAR****Micro Friedman
Rongeur**

#MICRO FRIED

..... 309.00 EACH

Tip Laser
Marked At
1-10mm**MOST POPULAR****Blumenthal Rongeur**

#BLUM

..... 309.00 EACH

**Blumenthal Rongeur
Right Angle**


#BLUM-RA

..... 309.00 EACH

**MOST POPULAR****Mini Blumenthal
Rongeur**

#MINI BLUM

..... 309.00 EACH

**Mead 1-A Oral Surgery
Rongeur**

#MEAD-1-A RONG

..... 309.00 EACH

Kerrison Rongeurs For Widening The Sinus Window



Kerrison Rongeur - 4mm Bite

10cm Long Beak

#KERRISON 649.00 EACH



MOST POPULAR

Kerrison Rongeur - 3mm Bite

8cm Long Beak

#KERRISON-SMALL 649.00 EACH

Low Profile Head
For Visibility
& Safe Access



Double Action Rongeurs - Increased Leverage & Effortless Cutting



**Double Action Gives You
A 3:2 Mechanical Advantage**



3.5mm Bite

DOUBLE ACTION

MOST POPULAR

**Beyer Rongeur
17.8cm Long**

Double Action For Increased Leverage

#BEYER. 471.00 EACH



5.0mm Bite


DOUBLE ACTION

**Ruskin Rongeur
17.8cm Long**

Double Action For Increased Leverage

#RUSKIN. 471.00 EACH

Retractors


**MOST POPULAR****Implant Retractor**Use As Upper Or Lower
Vestibular Retractor#IMPLANT-RETRACT
..... 198.00 EACH**Implant Retractor
Ergonomic Handle**Use As Upper Or Lower
Vestibular Retractor#IMPLANT-RETRACT-ERGO
..... 255.00 EACH**MOST POPULAR****Plastic Cheek Retractor
Single Handle
Autoclavable**#CHEEK-SINGLE
SET OF 2.....39.40 SET
2 OR MORE SETS.....36.30 SET**Plastic Cheek Retractor
Double Ended
Autoclavable**#CHEEK
SET OF 2.....58.90 SET
2 OR MORE SETS.....55.00 SET**Molt Mouth Gag**

#MOLT GAG..... 255.00 EACH

Molt Mouth Gag - Small#MOLT-GAG-SMALL
..... 255.00 EACH**Replacement Guards
For Molt Mouth Gag**

Autoclavable Silicone Rubber

#MOLT-GAG-GUARDS..... 15.90 SET OF 2

**Orringer Retractor**Spring Steel Flexes
For Comfortable Fit#ORRINGER..... 91.50 EACH
2 OR MORE..... 85.60 EACH**Orringer Retractor-Small**Spring Steel Flexes
For Comfortable Fit#ORRINGER-SMALL... 91.50 EACH
2 OR MORE..... 85.60 EACH

Retractors



Cawood-Minnesota Retractor

#CAWOOD-RETRACT
INDIVIDUALLY44.00 EACH
2 OR MORE40.10 EACH



MOST POPULAR
Minnesota Retractor

#MINN
INDIVIDUALLY44.00 EACH
2 OR MORE40.10 EACH



**Minnesota Retractor
Ergonomic Handle**

#MINN-ERGO
INDIVIDUALLY145.00 EACH
2 OR MORE134.60 EACH



**Dog Bone Sinus Retractor
- 18cm**

Sinus Flap Retractor
#DOG-BONE.68.00 EACH



MOST POPULAR
Columbia Retractor

Excellent For Surgery & Photography
#COL RETR55.60 EACH
2 OR MORE51.60 EACH



Austin Retractor

#AUSTIN-RETRACTOR72.00 EACH



*Easily Bent By Hand
To Desired Shape*

Bendable Retractor Stainless Steel

20cm Long x 18mm Wide
#RETRACTOR-BENDABLE59.80 EACH



Retractors

**MOST POPULAR****Weider Tongue Retractor**

44mm Wide At Tip

#WEIDER 84.50 EACH
 2 OR MORE 78.50 EACH

**Weider Tongue Retractor - Small**

31mm Wide At Tip

#WEIDER-SMALL 84.50 EACH
 2 OR MORE 78.50 EACH

**Salvin Zygoma Retractor**

- Increases Field Of Vision
- Directs The Drilling
- Prevents Perforation Of The Orbit
- Protects Soft Tissue While Drilling Superior Cortex

#ZAGA-2. 217.40 EACH

**Laster 3rd Molar Retractor**

#LASTER-RETRACT
 INDIVIDUALLY 81.50 EACH
 2 OR MORE 77.20 EACH




Photo: Dr. Thomas Balshi

SALVIN
Comfort
Handle

**Balshi Zygomatic Retractor**

Designed By Dr. Thomas Balshi

#BALSHI-ZYGO-RETRACT. 199.80 EACH




Photo: Dr. Michael Pikos

SALVIN
Comfort
Handle

MOST POPULAR**Pikos Ramus Retractor**

Designed By Dr. Michael Pikos

16mm Wide At Tip

#PIKOS-RETRACTOR 199.80 EACH

Retractors



Henahan Retractor Large

- 19.5cm
#HENAHAN-LARGE
..... 103.00 EACH
#3 OR MORE
..... 97.70 EACH



MOST POPULAR

Seldin 23

- #SELDIN
..... 103.00 EACH
#3 OR MORE
..... 97.70 EACH



Mini Seldin

- Round Handle For Comfort
#MINI-SELDIN
..... 103.00 EACH
#3 OR MORE
..... 97.70 EACH



Misch Sinus Retractor

- Designed By Dr. Carl Misch
#MISCH-SINUS-RETRACT
..... 101.80 EACH



MOST POPULAR

Dual View Photo Mirror

- View Buccal, Occlusal & Lingual Surfaces
At Once
#PHOTO-DUALVIEW
..... 110.90 EACH



MOST POPULAR

Autoclavable Silicone Mouth Props

- Autoclavable To 275°F Or 134°C
Dry Heat To 350°F Or 160°C
#MOUTH-PROP-CHILD
SET OF 2..... 42.90 SET
2 OR MORE SETS..... 39.60 SET
#MOUTH-PROP-ADULT
SET OF 2..... 42.90 SET
2 OR MORE SETS..... 39.60 SET

Periosteal Elevators

**MOST POPULAR****Molt 2/4**


#MOLT 2/4
 88.70 EACH
 3 OR MORE
 83.90 EACH

**Molt 2 Angled
Double-Ended**

#MOLT-2-DE
 88.70 EACH
 3 OR MORE
 83.90 EACH

**Molt 4 Angled
Double-Ended**

#MOLT-4-DE
 88.70 EACH
 3 OR MORE
 83.90 EACH

**Periosteal-Tome™**

Periosteal Elevator/
Periotome Combo

#PERIOSTEAL-TOME
 111.60 EACH
 3 OR MORE
 105.20 EACH

**MOST POPULAR****Prichard Elevator**

#PRICHARD
88.70 EACH
 3 OR MORE
83.90 EACH



Tip Marked
at 1 - 10mm
For Spacing

**MOST POPULAR****Seldin 23**

#SELDIN
 103.90 EACH
 3 OR MORE
 97.60 EACH

**Mini Seldin**

Round Handle
For Comfort

#MINI-SELDIN
 103.90 EACH
 3 OR MORE
 97.60 EACH

**Freer**

#FREER
 88.70 EACH
 3 OR MORE
 83.90 EACH

Periosteal Elevators



MOST POPULAR

Buser

- #BUSER88.70 EACH
- 3 OR MORE83.90 EACH




Goldman-Fox Elevator

- 17cm
- #GOLDMAN-FOX-ELEV88.70 EACH
- 3 OR MORE83.90 EACH



Hirschfeld #20

- #HIRSCHFELD-2088.70 EACH
- 3 OR MORE83.90 EACH



MOST POPULAR

Woodson

- #WOODSON88.70 EACH
- 3 OR MORE83.90 EACH



Quinn-Titanium

- Gold Titanium Tips For Use Around Implants
- #Q-ELEV-TI105.90 EACH
- 3 OR MORE98.40 EACH



MOST POPULAR

Quinn

- #Q-ELEV88.70 EACH
- 3 OR MORE83.90 EACH



24G Elevator

- #EP24G88.70 EACH
- 3 OR MORE83.90 EACH



Benque

- #BENQUE88.70 EACH
- 3 OR MORE83.90 EACH

Periosteal Elevators

**Molt 4**

#MOLT-4
..... 98.00 EACH
3 OR MORE
..... 92.90 EACH

**Molt 2**

#MOLT-2
..... 98.00 EACH
3 OR MORE
..... 92.90 EACH

**West**

#WEST
..... 88.50 EACH
3 OR MORE
..... 83.60 EACH

**Henahan Retractor Large**

19.5cm
#HENAHAN-LARGE
..... 103.00 EACH
3 OR MORE
..... 97.70 EACH

**MOST POPULAR****Molt 9**

#MOLT 9
..... 88.80 EACH
3 OR MORE
..... 83.80 EACH

**Molt 9 With Suture Hole**

#MOLT 9-SUTURE
..... 88.80 EACH
3 OR MORE
..... 83.80 EACH


**Molt 9 With 15 Aspirating Holes**

#MOLT 9-ASP
..... 88.80 EACH
3 OR MORE
..... 83.80 EACH

SuperCut Serrated Scissors Cut Effortlessly Through Tissue Or Suture

BUYER'S GUIDE - SUPERCUT SCISSORS

- ▶ Razor Sharp Specially Machined Blades Mesh Perfectly To Provide Superior Shearing Force At The Tips
- ▶ Micro Serrations On Lower Blade Grip Tissue For Slip-Free Cutting



Goldman Fox Scissors
13.2cm Curved

#SC GF SCISS147.40 EACH
3 OR MORE139.50 EACH




LaGrange Scissors
11.5cm Double Curved

#SC-LAGRANGE . . .147.40 EACH
3 OR MORE139.50 EACH



Iris Scissors
11cm Curved

#SC IRIS CU147.40 EACH
3 OR MORE139.50 EACH



Iris Scissors
11cm Straight

#SC IRIS ST147.40 EACH
3 OR MORE139.50 EACH

SALVIN Black Kryptonite™ Ceramic Coated Serrated Surgical Scissors

BUYER'S GUIDE - SALVIN BLACK KRYPTONITE™ SCISSORS

- ▶ Long Lasting, Razor Sharp Precise Cutting
- ▶ Lasts 4-6 Times Longer Than Conventional Scissors
- ▶ Non-Glare Finish For Better Tissue Visibility



Goldman Fox Scissors
12.5cm Curved

#Kryp-GF-SCISS . . . 268.00 EACH
3 OR MORE 255.30 EACH



LaGrange Scissors
11.5cm Double Curved

#Kryp-LAGRANGE . . .268.00 EACH
3 OR MORE 255.30 EACH



Iris Scissors
11.5cm Curved

#Kryp-IRIS-CU . . .268.00 EACH
3 OR MORE255.30 EACH



Castroviejo Micro Scissors
12cm Curved

#Kryp-CASTRO . . . 433.00 EACH
3 OR MORE404.70 EACH

Surgical & Suture Scissors

**Goldman Fox 12.7cm**

Serrated For Slip-Free Cutting

#GF SCISS. 112.60 EACH
3 OR MORE 107.90 EACH

*Double
Curved***MOST POPULAR****LaGrange 11.5cm**

Serrated For Slip-Free Cutting

#LAGRANGE. 112.60 EACH
3 OR MORE 107.90 EACH

**Iris Scissors Straight
11.5cm**

#IRIS-SCISS. 79.40 EACH
3 OR MORE 75.20 EACH

**MOST POPULAR****Iris Scissors Curved
11.5cm**

#IRIS-SCISS-CURVED
. 79.40 EACH
3 OR MORE 75.20EACH

**Long Suture Scissors
15cm**

Hooked End To Lift Suture

#LONG 148.50 EACH
3 OR MORE 140.10 EACH

**Littauer Angled
Suture Scissors 11.5cm**

Hooked End To Lift Suture

#LITTAUER-ANGLED . 119.00 EACH
3 OR MORE 112.60 EACH

**MOST POPULAR****Littauer Suture Scissors
11.5cm**

Hooked End To Lift Suture

#LITTAUER 95.10 EACH
3 OR MORE 91.30 EACH

**Metzenbaum-Mayo
Dissecting Scissors 14.5cm**

Blunt Tip To Dissect Tissue

#METZ-5. 108.00 EACH
3 OR MORE 100.70 EACH

Surgical & Suture Scissors



**Locklin Scissors Angled
15.5cm**

Serrated For Slip-Free Cutting
#LOCKLIN 147.50 EACH
3 OR MORE 140.30 EACH



**Wagner Scissors Angled
11.5cm**

#WAGNER-SCISS-ANGLED
. 123.00 EACH
3 OR MORE 117.30 EACH



Kelly Curved 15.9cm

Serrated For Slip-Free Cutting
#KELLY CURV 112.00 EACH
3 OR MORE 104.40 EACH



MOST POPULAR

Dean 17.2cm

Serrated For Slip-Free Cutting
#DEAN 147.50 EACH
3 OR MORE 140.30 EACH

Featherlite Micropoint Scissors


*Takes The "Wince" Out Of Suture Removal,
Even In Sensitive Tissue*



MOST POPULAR

Featherlite 14cm Curved

#F-LITE5CURV 345.00 EACH
3 OR MORE 326.30 EACH



Opposing Spring Tips Maintain Force
Between Blades For Incredible
Cutting Without Pinching Tissue

Suture Scissors / Suture Forcep Combination Instrument

Cut & Remove Suture With One Instrument



#SUT-REM-FORCEP 345.00 EACH
3 OR MORE 326.30 EACH



Engage Suture On Right
Side Of Knot...



...Cut, Squeeze & Pull
To Remove Suture...



Tungsten Carbide Tipped Needle Holders For Larger Needle 3-0 & 4-0 Sutures



Designed For 3-0 To 4-0 Sutures

Accommodates Larger Needles Without
Damaging The Locking Clasp & Instrument Tip



Serrated Tips

TUNGSTEN CARBIDE

FOR 3-0 TO 4-0 SUTURES

MOST POPULAR

Castroviejo Needle Holder Accommodates 3-0 To 4-0 Sutures

Spring Handle For Precise Control

14cm

#CASTRO-3.0-4.0 482.00 EACH



Serrated Tips

TUNGSTEN CARBIDE

FOR 3-0 TO 4-0 SUTURES

Jacobson Needle Holder Accommodates 3-0 To 4-0 Sutures

Spring Handle For Precise Control

17.8cm

#JACOBSON-3.0-4.0. 569.00 EACH

Tungsten Carbide Tipped Needle Holder With Side Cutting Scissors

Scissor On Side For Faster Suturing Without Changing Instruments



TUNGSTEN CARBIDE


FOR 2-0 TO 4-0 SUTURES



Hold Suture Taut In Slot...



...Slide Down To Knot...



*...Squeeze To Cut,
Leaving A 2mm Tag*

Side Cutting Mayo Hegar 6"

Carbide Tipped With Scissors

15.3cm

#SIDE CUT MAYO 352.40 EACH

3 OR MORE 335.90 EACH

Tungsten Carbide Tipped Castroviejo & Micro Needle Holders




BUYER'S GUIDE - CARBIDE TIPPED NEEDLE HOLDERS	
▶	Superior Grip Holds Needle Without Slipping
▶	Less Hand Fatigue & Cramping During Long Surgeries
▶	Last 4-6 Times Longer Than Standard Instruments
▶	Can Be Economically Re-Tipped



Castroviejo Needle Holder

Spring Handle For Precise Control - 14cm
#CASTRO 482.00 EACH
#CASTRO-CURVED 482.00 EACH



Jacobson Needle Holder

Spring Handle For Precise Control - 17.8cm
#JACOBSON 569.00 EACH
#JACOBSON-CURVED 569.00 EACH



ThumbLok Snagless Needle Holder

17.8cm - Straight Tips
#NH-ROUND-7-TCL-R. 569.00 EACH



Round Jacobson Needle Holder

Spring Handle For Precise Control - 17.8cm
#JACOBSON-ROUND 569.00 EACH



Micro Jacobson Needle Holder

17.8cm - Straight Tips
#MICRO-JACOBSON-STR 569.00 EACH
 17.8cm - Curved Tips
#MICRO-JACOBSON-CUR 569.00 EACH



Micro Castroviejo Needle Holder

13.3cm - Straight Tips
#MICRO-CASTRO-STR 569.00 EACH
 13.3cm - Curved Tips
#MICRO-CASTRO-CUR 569.00 EACH

Tungsten Carbide Tipped Needle Holders



Serrated Tips

TUNGSTEN CARBIDE

FOR 3-0 TO 5-0 SUTURES

Crile Wood Needle Holder - Curved

15.3cm

#CRILE-CURVED 209.00 EACH
3 OR MORE 199.70 EACH



Serrated Tips

TUNGSTEN CARBIDE

FOR 3-0 TO 5-0 SUTURES

MOST POPULAR

Crile Wood Needle Holder

15.3cm

#CRILE 209.00 EACH
3 OR MORE 199.70 EACH



Serrated Tips

TUNGSTEN CARBIDE

FOR 2-0 TO 4-0 SUTURES

MOST POPULAR

Mayo Hegar Needle Holder

15.3cm

#MAYO HEGAR 209.00 EACH
3 OR MORE 199.70 EACH



Serrated Tips

TUNGSTEN CARBIDE

FOR 3-0 TO 5-0 SUTURES

Crile Wood Needle Holder - Long

17.8cm

For Suturing Extracted 3rd Molar Sockets

#CRILE-LONG 209.00 EACH
3 OR MORE 199.70 EACH



Serrated Tips

TUNGSTEN CARBIDE

FOR 3-0 TO 5-0 SUTURES

Mathieu Needle Holder

"Squeeze Through" Quick Release

12.7cm

#MATHIEU 286.00 EACH
3 OR MORE 273.70 EACH



Serrated Tips

TUNGSTEN CARBIDE

FOR 5-0 TO 7-0 SUTURES

SALVIN
Comfort
HandleCastroviejo Needle Holder
Straight Handle


12.7cm

#CASTRO-STR-HANDLE 530.00 EACH



Sutures On Pages 35-39

Tungsten Carbide Tipped Needle Holders




MOST POPULAR

Par Perio Needle Holder

12.7cm

Delicate Serrated

#PAR..... 207.80 EACH
3 OR MORE 195.20 EACH



Derf Needle Holder

12.1cm

#DERF..... 207.80 EACH
3 OR MORE 195.20 EACH



Olsen Hegar Scissors/Needle Holder

14cm

#OH 5..... 244.50 EACH
3 OR MORE 234.90 EACH



Olsen Hegar Scissors/Needle Holder

16.5cm

#OH 6..... 244.50 EACH
3 OR MORE 234.90 EACH

Salvin Non-Glare Featherweight Microsurgery Instruments

Custom Designed For Periodontal Microsurgery



▶	Light & Long
▶	Smooth Tips
▶	Rounded Handles
▶	Non-Locking
▶	Black Non-Glare Coating



"The Salvin microsurgical instruments exhibit a semicircular grip to enhance handling and haptic feedback; these instruments are light in weight to prevent muscle fatigue and have been designed without locking mechanisms to prevent undesirable traumatic instrument manipulation. The dark surface coating minimizes glaring when working under coaxial illumination. Once you use these finely crafted instruments it will be hard to work with any other set!"

Dr. Diego Velasquez
Periodontist, Fenton, MI



"The Salvin microsurgery set combines the features needed in high quality surgical instrumentation. The non-reflective black coating effectively eliminates glare during surgeries, allowing for better visibility and precision. The lightweight design enhances our tactile control, and the smooth tips are conducive for manipulation of micro sutures. Salvin listened to what we needed and delivered!"

Dr. Irene Marron-Tarrazzi
Periodontist, Miami, FL



"When performing perio plastic surgery under magnification and illumination, using specialized instrumentation is critical to my success. The Salvin microsurgery instruments check all of the requirements I need: Light, long, non-reflective, rounded handles and manufactured with the highest quality. This is exactly what I use in my practice."

Dr. John Russo
Periodontist, Sarasota, FL



Non-Glare Micro Tunneling - 1.2mm

#MICRO-BLK-TUNNEL-1.2
..... 111.20 EACH

Non-Glare Micro Tunneling - 2.0mm

#MICRO-BLK-TUNNEL-2.0
..... 111.20 EACH

Non-Glare Micro Tunneling - 2.8mm

#MICRO-BLK-TUNNEL-2.8
..... 111.20 EACH

Double-Ended Micro
Tunnel Elevators Help Reflect &
Elevate Tissue From Recipient
Site In Precise Micro Flaps



Non-Glare Lucas 85 Bone Curette

Sharp Tips Efficiently
Debride Sockets

#LUCAS-85-BLK
..... 95.00 EACH



Non-Glare Prichard Elevator

Double-Ended Micro
Tunnel Elevators Help
Reflect & Elevate Tissue
From Recipient Site In
Precise Micro Flaps

#PRICHARD-BLK
..... 90.50 EACH



Non-Glare Perio Probe

#PERIO-PROBE-BLK
..... 106.00 EACH



Non-Glare Papilla Elevator

- Double-Ended Elevator
- Rounded Tip For
Blunt Tissue Dissection

#PAPILLA-ELEV-BLK
..... 106.00 EACH

SPECIAL



Salvin Non-Glare Featherweight Microsurgery Set

Includes:

- Non-Glare Micro Tunneling - 2.0mm
 - Non-Glare Micro Tunneling - 2.8mm
 - Non-Glare Tissue Forcep - Curved OR Non-Glare Tissue Forcep - Straight
 - Non-Glare Castroviejo Micro Needle Holder - Curved
 - Non-Glare Micro Scissor - Curved
 - Non-Glare Lucas 85 Bone Curette
 - Non-Glare Prichard Elevator
 - Non-Glare Perio Probe
 - Non-Glare Papilla Elevator
 - Non-Glare Micro Blade Handle
 - Non-Glare Micro Mouth Mirror - 5mm
 - Non-Glare Mirror Handle
 - Plus **FREE** Instrument Cassette
- #NON-GLARE-MICRO-KIT**
..... **.2470.00 SET**



**Non-Glare Castroviejo
Micro Needle Holder -
Curved**

- Light (<20 grams) And Long (18cm) For Easy Handling Of Micro Sutures
- Smooth Tips Easily Grasp Delicate 6.0-10.0 Sutures, While The Rounded Handles Allow For Precise Control

#MICRO-BLK-CASTRO-CUR
..... **544.00 EACH**



**Non-Glare
Micro Scissor - Curved**

Sharp, Slightly Curved, Long Scissors For Fine Resections Of Soft Tissues & Trimming Of Mucosal Grafts

#MICRO-BLK-SCISS-CUR
..... **489.00 EACH**



**Non-Glare
Mirror Handle**

#MIRROR-HNDL-BLK
..... **69.50 EACH**



**Non-Glare
Micro Mouth Mirror - 5mm**

#MICRO-MIRR-5-BLK
..... **156.00 EACH**



**Non-Glare
Micro Mouth Mirror
- 3x6mm**

#MICRO-MIRR-3X6-BLK
..... **156.00 EACH**




**Non-Glare
Tissue Forcep -
Curved**

#MICRO-BLK-FORCEP-CUR
..... **489.00 EACH**



**Non-Glare
Tissue Forcep -
Straight**

#MICRO-BLK-FORCEP-STR
..... **489.00 EACH**



**Non-Glare
Suture Pliers**

#MICRO-BLK-SUTURE-PLIER
..... **437.00 EACH**



**Non-Glare
Micro Blade Handle**

- Fits Micro "Beaver-Type" Scalpel Blades
 - Collet Locks Blade In Place
- #MICRO-BLK-BLADE-HNDL**
..... **144.20 EACH**

Kern Tunneling Instruments

Designed By Dr. Ricardo Kern

- *Minimally Invasive Periodontal Microsurgery*
- *Avoids Need For Large Flap*
- *Anatomical Silicone Handles With Stainless Steel Tips*



"Precise and delicate, these instruments have been designed to ensure that you have no limitations in performing to the best of your ability. In this kit you'll find the sophistication you need to explore new surgical possibilities and maximize the performance of current mucogingival reconstruction techniques!"

Dr. Ricardo Kern
Periodontist, Brazil



SPECIAL

Kern Tunneling Set

Includes:

- Kern Tunneler #1
- Kern Tunneler #2
- Kern Tunneler #3
- Kern Tunneler #4
- Kern Tunneler #5
- 1 Cassette

#KERN-TUNNELING-KIT
805.00



Tunnelling Blades And Handles



Aluminum Tunnel Blade Scalpel Handle

#TUNNEL-BLADE-HNDL-A
..... 81.00 EACH



Titanium Tunnel Blade Scalpel Handle

#TUNNEL-BLADE-HNDL-TI
..... 110.00 EACH



Stingray Tunnel Blade 1

- "Back Action" Greatly Facilitates The Division Of Muscle Fibers In The Mucosal Or Periosteal Plane.
- Pack Of 5 – Individually Wrapped Sterile

#TUNNEL-BLADE-STINGRAY 70.00 PACK OF 5



Stingray Tunnel Blade 2

- Micro-Blade For Access To Sulcular Tunneling. Small & Delicate, Perfect For Thin & Slender Phenotypes Even In The Narrow Contours Of The Lower Incisors.
- Pack Of 5 – Individually Wrapped Sterile

#TUNNEL-BLADE-TB2 70.00 PACK OF 5

Vanessa Frazão Tunneling Set

Designed By Dr. Vanessa Frazão

- *Minimally Invasive Periodontal Microsurgery*
- *Designed To Provide Maximum Accuracy, Efficiency And Safety*

NEW

SPECIAL

Vanessa Frazão Tunneling Set

Includes:

- VF1 Tunneler
- VF2 Tunneler
- VF3 Tunneler
- VF4 Tunneler
- 1 Cassette

#FRAZAO-TUNNELING-KIT
560.00




Bendable Miniature & Micro Scalpel Blades

Uncompromised Sharpness & Cutting Ability

Blades Are Sterile & Individually Wrapped

Bendable Blades

*Use As Straight Blade Or
Bend With Needle Holder To
Desired Angle (Up To 45°)*



MOST POPULAR

MOST POPULAR

MOST POPULAR



Mini Blade #64

#BLAD-MINI-64 71.20 BOX OF 12



Mini Blade #67

#BLAD-MINI-67 71.20 BOX OF 12



Mini Blade #63 - *Full Radius Edge*

#BLAD-MINI-63 71.20 BOX OF 6



Mini Blade #69 - *Full Radius Edge*

#BLAD-MINI-69 71.20 BOX OF 6

Bendable Micro-Miniature Surgical Blades



Micro Blade #64

#BLAD-MICRO-64 55.10 BOX OF 6



Nordland Micro Blade #6900 - *Full Radius Edge*

#BLAD-NORDLAND-69 105.00 BOX OF 6

Blades Shown Actual Size



**Blade Handle
For Miniature & Micro Surgical Blades**

#BLAD-HANDLE 78.60 EACH



**Ergonomic Blade Handle
For Miniature & Micro Surgical Blades**

#TVA-BLAD-HANDLE 147.60 EACH

Tunneling & Vestibular Approach Instrument Set

**Perio Probe**

#TVA-PERIO-PROBE
..... 115.00 EACH

**Perio Elevator #1**

#TVA-PERIO-ELEV-1
..... 115.00 EACH

**Perio Elevator #2**

#TVA-PERIO-ELEV-2
..... 115.00 EACH

**Perio Elevator #3**

#TVA-PERIO-ELEV-3
..... 115.00 EACH

**Tunnel Knife #1**

#TVA-TUNNEL-KNIFE-1
..... 115.00 EACH

**Tunnel Knife #2**

#TVA-TUNNEL-KNIFE-2
..... 115.00 EACH

**Gingivectomy Knife**

#TVA-GINGI-KNIFE
..... 115.00 EACH

**Flap Prep Instrument #1**

#TVA-FLAP-PREP-1
..... 115.00 EACH

**Flap Prep Instrument #2**

#TVA-FLAP-PREP-2
..... 115.00 EACH

**Perio Chisel**

#TVA-PERIO-CHISEL
..... 115.00 EACH

Tunneling & Vestibular Approach Instrument Set



Tissue Plugger

#TVA-TISSUE-PLUGGER
..... 115.00 EACH



Tunnel Spoon #1

#TVA-TUNNEL-SPOON-1
..... 115.00 EACH



Tunnel Spoon #2

#TVA-TUNNEL-SPOON-2
..... 115.00 EACH



Tunnel Spoon #3

#TVA-TUNNEL-SPOON-3
..... 115.00 EACH




Micro Blade Handle

#TVA-BLAD-HANDLE
..... 147.60 EACH

SPECIAL

Tunneling & Vestibular Approach Instrument Set

MOST POPULAR



- Perio Probe
 - Perio Elevator #1
 - Perio Elevator #2
 - Perio Elevator #3
 - Tunnel Knife #1
 - Tunnel Knife #2
 - Gingivectomy Knife
 - Flap Prep Instrument #1
 - Flap Prep Instrument #2
 - Perio Chisel
 - Tissue Plugger
 - Tunnel Spoon #1
 - Tunnel Spoon #2
 - Tunnel Spoon #3
 - Micro Blade Handle
 - **FREE** Stainless Steel Cassette
- #TVA-TUNNEL-KIT . . . **1765.00 SET**

Tunneling & Soft Tissue Instruments



Tunneling Elevator #1

#TUNNEL-ELEV-1
.....98.70 EACH

**MOST POPULAR**

Tunneling Elevator #2

#TUNNEL-ELEV-2
.....98.70 EACH



Tunneling Elevator #3

#TUNNEL-ELEV-3
.....98.70 EACH

PD Miller
Tunneling Elevator

#PDMILLER-ELEV-1
.....98.70 EACH



Micro Forceps

#FORCEP-MICRO
.....123.60 EACH



Tissue Augmenter

#TISSUE-AUGMENT
.....123.60 EACH

**MOST POPULAR**

Corn Pliers - Angled

Suture Forcep With
"Pass-Through" Hole
15.4cm
#CORN-PLIER
.....97.70 EACH

**MOST POPULAR**

Suture Pliers

"Delicate Procedure"
15.8cm Suture Forcep
1.6mm Ring Tip
#SUTURE-PLIER
.....159.00 EACH

Orban & Hirschfeld For Tunneling



MOST POPULAR

Orban Knife

#ORBAN
.....88.00 EACH
3 OR MORE
.....83.40 EACH

MOST POPULAR

Orban 1/2 Short

#ORBAN-1/2-SHORT
.....88.00 EACH
3 OR MORE
.....83.40 EACH

Orban 1/2 Long

#ORBAN-1/2-LONG
.....88.00 EACH
3 OR MORE
.....83.40 EACH

Hirschfeld 16/17 Tunneling File

Rounded Tip Prevents Puncturing Tissue

#HIRSCHFELD-16/17-FIL . . 92.50 EACH
3 OR MORE 86.80 EACH

SPECIAL

Salvin Tunneling Basic Instrument Set #2

MOST POPULAR



Includes:

- Tunneling Forcep
- Micro Jacobson Curved N/H
- Tunneling Elevator #1
- Tunneling Elevator #2
- Tunneling Elevator #3
- Orban 1/2 Short
- Orban 1/2 Long
- Hirschfeld 16/17 Tunneling File
- PD Miller Tunneling Elevator
- Younger Good 7/8
- Suture Pliers
- **FREE** Stainless Steel Cassette

#TUNNEL-KIT-2 **1779.00 SET**

Micro Surgery Instruments

**MOST POPULAR****Rhodes Chisel**

17.2cm
Back Action
#RHODES
..... 109.00 EACH

**MOST POPULAR****Micro Elevator #1**

17.2cm
#MICRO-ELEV-1
..... 109.00 EACH

**Micro Elevator #2**

17.2cm
#MICRO-ELEV-2
..... 109.00 EACH

**Micro Elevator #3**

17.2cm
#MICRO-ELEV-3
..... 109.00 EACH

**Micro Elevator #4**

17.2cm
#MICRO-ELEV-4
..... 109.00 EACH

**MOST POPULAR****Micro Chisel #1**

17.2cm
#MICRO-CHISEL-1
..... 109.00 EACH

**Micro Chisel #2**

17.2cm
#MICRO-CHISEL-2
..... 109.00 EACH

SPECIAL**Salvin Micro Surgery Set****MOST POPULAR**

- Micro Castro Needle Holder Curved
- Micro Adson Tissue Forcep
- Ross Scissors
- Rhodes Chisel
- Micro Chisel #1
- Micro Elevator #1
- Micro Elevator #3
- **FREE** Stainless Steel Cassette



#MICRO-SURG-KIT-1
1219.00 SET

Pikos / Ordonez Soft Tissue Manipulation Instrument Set

Ideal For Tunneling & Pocket Flap Procedures

Collaboratively Designed By Dr. Mike Pikos & Dr. Alvaro Ordonez



Pikos/Ordonez #1
Tunneling
Instrument
#SOFTIMAK-1
..... 101.80 EACH



Pikos/Ordonez #2
Used For
Pocket Flaps
#SOFTIMAK-2
..... 101.80 EACH



Pikos/Ordonez #3
Used For
Pocket Flaps
#SOFTIMAK-3
..... 101.80 EACH



Pikos/Ordonez #5
Used For
Pocket Flaps
#SOFTIMAK-5
..... 101.80 EACH



Pikos/Ordonez #6
Used For
Pocket Flaps
#SOFTIMAK-6
..... 101.80 EACH



Pikos/Ordonez #8
Used For
Reverse Flap
#SOFTIMAK-8
..... 101.80 EACH



Pikos/Ordonez #7
Thin Modified Orban
#SOFTIMAK-7
..... 101.80 EACH



Pikos/Ordonez #4
Thin Microperiosteal
#SOFTIMAK-4
..... 101.80 EACH




**Pikos/Ordonez
Tissue Forcep**
TC Coated
#SOFTIMAK-TF-1
..... 129.20 EACH

SPECIAL

**Pikos/Ordonez Soft Tissue
Manipulation Set**

MOST POPULAR



Set Of 9 Instruments With
FREE Stainless Steel Cassette
#SOFTIMAK-KIT
867.00 SET


Tunneling Forceps

For Pulling Tissue Graft Through Tunnel

- Micro Thin For Uninterrupted Passage Through Tunnel
- 75° Angle For All Areas & Applications
- Pressure Limiting Spring Prevents Damage To Graft
- Diamond Coated For Slip Resistance



Guide Through Tunnel To
Grasp Tissue Graft



Close Forcep & Pull Graft
Through Tunnel



MOST POPULAR

**Tunnel Forcep
Straight**

#TUNNEL-FORCEP
..... 412.00 EACH



**Tunnel Forcep
Curved Left**

#TUNNEL-FORCEP-CU-L
..... 412.00 EACH



**Tunnel Forcep
Curved Right**

#TUNNEL-FORCEP-CU-R
..... 412.00 EACH

SPECIAL

Set Of 3 Tunnel Forceps

#TUNNEL-FORCEP-SET-3 955.00 SET OF 3

Micro Needle Holders



Smooth Tips

TUNGSTEN CARBIDE

FOR 6-0 TO 10-0 SUTURES

**SALVIN
Comfort
Handle**

Micro Castro Needle Holder

13.3cm - Straight Tips

#MICRO-CASTRO-STR 568.00 EACH

13.3cm - Curved Tips

#MICRO-CASTRO-CUR 568.00 EACH



Serrated Tips

TUNGSTEN CARBIDE

FOR 5-0 TO 7-0 SUTURES

**SALVIN
Comfort
Handle**

Micro Needle Holder

13.3cm - Straight Tips

#MICRO-N-HOLDER-STR 568.00 EACH

13.3cm - Curved Tip

#MICRO-N-HOLDER 568.00 EACH



MOST POPULAR



Smooth Tips

TUNGSTEN CARBIDE

FOR 6-0 TO 10-0 SUTURES

**SALVIN
Comfort
Handle**

Micro Jacobson Needle Holder

17.8cm - Straight Tips

#MICRO-JACOBSON-STR 568.00 EACH

17.8cm - Curved Tips

#MICRO-JACOBSON-CUR 568.00 EACH

Micro Surgery Instruments



**Ross Scissors
Curved**

15cm
#ROSS.....391.00 EACH
3 OR MORE371.90 EACH



MOST POPULAR

**Castroviejo Micro Scissors
Curved**

15cm
#MICRO-CASTRO-SCISS.....524.00 EACH



**Micro Goldman Fox Scissors
Curved**

11.5cm
#MICRO-GF-SCISS...151.20 EACH
3 OR MORE144.70 EACH



*Extra Thin
Micro Tip*

**Extra Thin Micro
Iris Scissors - Curved**

11cm
#IRIS-THIN-CURVED...127.50 EACH
3 OR MORE119.20 EACH



**Micro Adson 1x2
Tissue Forcep**

1x2 Teeth - 0.5mm Tips
12.1cm
#MICRO-ADSON-1X2.....97.70 EACH
3 OR MORE92.90 EACH



MOST POPULAR

**Micro Adson Serrated
Tissue Forcep**

0.5mm Serrated Jaws
12.1cm
#MICRO-ADSON.....97.70 EACH
3 OR MORE92.90 EACH



**Bishop Harmon
Micro Tissue Forcep**

1x2 Teeth - 0.3mm Tips
8.6cm
#BISHOP-HARMON.....151.20 EACH
3 OR MORE144.70 EACH

Periodontal Knives & Surgical Files

**MOST POPULAR****Orban Knife**

#ORBAN
.....88.00 EACH
3 OR MORE
.....83.40 EACH

**MOST POPULAR****Orban 1/2 Short**

17.2cm
#ORBAN-1/2-SHORT
.....88.00 EACH
3 OR MORE
.....83.40 EACH

**Orban 1/2 Long**

Extended Length For
Tunneling Procedures
17.2cm
#ORBAN-1/2-LONG
.....88.00 EACH
3 OR MORE
.....83.40 EACH

**Kirkland Knife**

#KIRKLAND
.....88.00 EACH
3 OR MORE
.....83.40 EACH

**1S/2S Sugarman
Periodontal File**

#SUGARMAN-1S/2S
.....92.50 EACH

**MOST POPULAR****3S/4S Sugarman
Periodontal File**

#SUGARMAN-3S/4S
.....92.50 EACH

**Hirschfeld 16/17
Tunneling File**

Rounded Tip
Prevents Puncturing Tissue
During Tunneling Procedures
#HIRSCHFELD-16/17-FIL
.....92.50 EACH

Perio Curettes



MOST POPULAR

Prichard 1/2 Universal Perio Curette


#PRICHARD-1/2 49.30 EACH
3 OR MORE 45.70 EACH



Kramer #3 Universal Perio Curette

#KRAMER-3 49.30 EACH
3 OR MORE 45.70 EACH

Bone Files



Bone File #10

#FIL10 89.30 EACH



Bone File #12

#FIL12 89.30 EACH



MOST POPULAR

Bone File #12CA

#FIL-12CA 89.30 EACH

Tissue Forceps

**MOST POPULAR****Adson Long 1x2
Tissue Forceps**

1x2 Teeth
15cm

#ADSON-LONG-1X2
..... 78.90 EACH
3 OR MORE
..... 73.60 EACH

**MOST POPULAR****Adson Long Serrated
Tissue Forceps**

Serrated Teeth
15cm

#ADSON-LONG-SERR
..... 78.90 EACH
3 OR MORE
..... 73.60 EACH

**Adson Curved
Serrated
Tissue Forceps**

Serrated Teeth
12cm

#ADSON-SERR-CURV
..... 78.90 EACH
3 OR MORE
..... 73.60 EACH

**Debakey Atraumatic
Tissue Forcep -
Straight**

Serrated Teeth
16cm

#DEBAKEY-STRAIGHT
..... 138.00 EACH
3 OR MORE
..... 131.90 EACH

**Debakey Atraumatic
Tissue Forcep -
Curved**

Serrated Teeth
16cm

#DEBAKEY-CURVED
..... 138.00 EACH
3 OR MORE
..... 131.90 EACH

**Chalazion Mucocele
Forcep**

12mm x 16mm Ring
8.5cm Length

#CHALAZION
..... 190.40 EACH
3 OR MORE
..... 181.70 EACH

Tissue Forceps



MOST POPULAR

**Adson 1x2
Tissue Forceps**

1x2 Teeth
12cm

#ADSON 1X2
.....74.80 EACH
3 OR MORE
.....70.50 EACH



MOST POPULAR

**Adson 1x2 Curved
Tissue Forceps**

1x2 Teeth
12cm

#ADSON 1X2 CUR
.....78.90 EACH
3 OR MORE
.....73.60 EACH



MOST POPULAR

**Micro Adson
Serrated
Tissue Forcep**

0.5mm Serrated
12cm

#MICRO-ADSON
.....97.70 EACH
3 OR MORE
.....92.90 EACH



**Adson 7x7
Tissue Forceps**

7x7 Teeth
12cm

#ADSON-7X7
.....97.50 EACH
3 OR MORE
.....92.80 EACH



**Adson 7x7 Curved
Tissue Forceps**

7x7 Teeth
12cm

#ADSON-7X7-CURVED
.....97.50 EACH
3 OR MORE
.....92.80 EACH



**Dressing/Tissue
Forcep**

Serrated Tip
15.2cm

#DRESS-FORCEP-6
.....63.80 EACH
3 OR MORE
.....59.00 EACH



Cotton Plier

Serrated Tip
15cm

#COTPLI
.....43.10 EACH
3 OR MORE
.....40.50 EACH

Tissue Forceps

**Gerald 1x2**

Delicate Tissue Forceps
1x2 Teeth
18cm

#GERALD 1X2
.....95.40 EACH
3 OR MORE
.....90.30 EACH

**Gerald Serrated**

Delicate Tissue Forceps
Serrated Teeth
18cm

#GERALD
.....95.40 EACH
3 OR MORE
.....90.30 EACH

**MOST POPULAR****Gerald Serrated Straight**


Delicate Tissue Forceps
Serrated Teeth
18cm

#GERALD-STRAIGHT
.....95.40 EACH
3 OR MORE
.....90.30 EACH

**MOST POPULAR****Gerald 1x2 Teeth Straight**

Delicate Tissue Forceps
1x2 Teeth
18cm

#GERALD-1X2-STRAIGHT
.....95.40 EACH
3 OR MORE
.....90.30 EACH

**MOST POPULAR****Russian Long**

Tissue Forceps
Serrated Tip
19.7cm

#RUSS-LONG
.....93.40 EACH
3 OR MORE
.....88.20 EACH

**Russian**

Tissue Forceps
Serrated Tip
15.2cm

#RUSS
.....84.80 EACH
3 OR MORE
.....80.30 EACH

Tissue Forceps



Adson Titanium Curved Tissue Forceps

Serrated Teeth
12cm
#ADSON-TI-CURVED
..... 262.90 EACH
3 OR MORE
..... 250.40 EACH



MOST POPULAR

Adson Carbide Tipped Tissue Forcep

Serrated Carbide Tips
For Atraumatic Grasp
Of Tissue
12cm
#ADSON TC
..... 185.10 EACH
3 OR MORE
..... 178.30 EACH



TC Gripper Carbide Tipped Tissue Forcep

Serrated Carbide Tips
For Atraumatic Grasp
Of Tissue
14.5cm
#TC GRIPPER
..... 194.80 EACH
3 OR MORE
..... 185.50 EACH



TC Gripper Long Carbide Tipped Tissue Forcep

Serrated Carbide Tips
For Atraumatic Grasp
Of Tissue
18cm
#TC-GRIPPER-LONG
..... 194.80 EACH
3 OR MORE
..... 185.50 EACH


Specialty Suturing Forceps



MOST POPULAR

Corn Pliers - Angled

Suture Forcep With
"Pass-Through" Hole
15.4cm
#CORN-PLIER
..... 97.70 EACH
3 OR MORE
..... 92.70 EACH



V Forceps

Serrated Teeth
Designed By
Dr. Richard Martin
#V-FORCEP-BLK
..... 95.40 EACH
3 OR MORE
..... 90.50 EACH



MOST POPULAR

Suture Pliers

Delicate Suture Forcep
With 1.6mm Ring Tip
15.8cm
#SUTURE-PLIER
..... 159.00 EACH
3 OR MORE
..... 151.50 EACH

Super Sharp Disposable Tissue Punches

Sterile Disposable Tissue Punches

- Razor Sharp Cutting Edge
- Individually Packaged Sterile - Box Of 50
- Can Be Used In Single Tooth Sites

MOST POPULAR**MOST POPULAR**

● 3.5mm

● 4.0mm

● 5.0mm

● 6.0mm



3.5MM PUNCHES-BOX OF 50

4.0MM PUNCHES-BOX OF 50

5.0MM PUNCHES-BOX OF 50

6.0MM PUNCHES-BOX OF 50

	1 Box	2 Boxes
#PUNCH-3.5	189.80 BOX	176.80 BOX
#PUNCH-4.0	189.80 BOX	176.80 BOX
#PUNCH-5.0	189.80 BOX	176.80 BOX
#PUNCH-6.0	189.80 BOX	176.80 BOX

Short Sterile Disposable Tissue Punches

- Razor Sharp Cutting Edge
- Individually Packaged Sterile - Box Of 50
- Designed For Edentulous Areas - Not For Single Tooth Sites



● 3.5mm



● 4.0mm



● 5.0mm



● 6.0mm

3.5MM SHORT PUNCHES-BOX OF 50

4.0MM SHORT PUNCHES-BOX OF 50

5.0MM SHORT PUNCHES-BOX OF 50

6.0MM SHORT PUNCHES-BOX OF 50

	1 Box	2 Boxes
#PUNCH-SHORT-3.5	189.80 BOX	176.80 BOX
#PUNCH-SHORT-4.0	189.80 BOX	176.80 BOX
#PUNCH-SHORT-5.0	189.80 BOX	176.80 BOX
#PUNCH-SHORT-6.0	189.80 BOX	176.80 BOX

*Disposable Punches Are One Size Per Box
Boxes May Be Combined For "2 or More Box" Pricing*

Stainless Steel Multi-Use Rotary Tissue Punches



● 3.0mm



● 4.0mm



● 5.0mm



● 6.0mm



● 7.0mm

3.0MM ROTARY PUNCH

4.0MM ROTARY PUNCH

5.0MM ROTARY PUNCH

6.0MM ROTARY PUNCH

7.0MM ROTARY PUNCH

#PUNCH-3-ROTARY	133.00 EACH
#PUNCH-4-ROTARY	133.00 EACH
#PUNCH-5-ROTARY	133.00 EACH
#PUNCH-6-ROTARY	133.00 EACH
#PUNCH-7-ROTARY	133.00 EACH

Use With Latch Type
Implant Handpiece


SPECIAL**Rotary Tissue Punch Kit****MOST POPULAR**

Set Of 5 Punches With Autoclavable Block

#PUNCH-ROTARY-KIT

589.00 SET OF 5 W/BLOCK

Stainless Steel Tissue Punches



Stainless Steel Tissue Punches - Offset

	Individually	3 or More
4.0MM OFFSET PUNCH	#PUNCH-SS-4-RA. 110.40 EACH	105.10 EACH
5.0MM OFFSET PUNCH	#PUNCH-SS-5-RA. 110.40 EACH	105.10 EACH
6.0MM OFFSET PUNCH	#PUNCH-SS-6-RA. 110.40 EACH	105.10 EACH
8.0MM OFFSET PUNCH	#PUNCH-SS-8-RA. 110.40 EACH	105.10 EACH

Stainless Steel Punches May Be Assorted For "3 or More" Pricing

Tissue Flap Punch




Garvar Tissue Flap Punch Implant Exposing Forcep

Designed By Dr. Lanny Garvar

- Use To Contour Tissue Around Implant
- 5mm Diameter Ring

#FLAP-PUNCH 308.00 EACH



Tissue Punch Forcep



Tissue Punch Forcep

To Remove Tissue After Using Punch


#TISSUE-PUNCH-FORC. 83.40 EACH
3 OR MORE 78.90 EACH



Start With Your Choice Of Tissue Punch → Remove Tissue With Forcep

Scalpel Handles

Blade Locks At 90°,
60°, 30°, Straight,
330°, 300°,
& 270° Angles



Tighten / Loosen
To Lock / Unlock



MOST POPULAR

6 Way Adjustable Angled Scalpel Handle

#SCALP-6-WAY-ADJ
..... 263.00 EACH
2 OR MORE
..... 249.80 EACH



MOST POPULAR

Siegel Round Scalpel Handle

For Precise Control When
Rotating Angle Of Incision

#SIEGEL 96.20 EACH
3 OR MORE 88.70 EACH



Siegel Round Angled Scalpel Handle

For Precise Control When
Rotating Angle Of Incision

#SIEGEL-ANGLED. 96.20 EACH
3 OR MORE 88.70 EACH




MOST POPULAR

#3 Scalpel Handle

#SCALP. 30.10 EACH
3 OR MORE .. 27.40 EACH



Push Lever Forward
To Remove Blade



Safety Ejector Scalpel Handle

Safe & Easy Blade Removal

#SCALP-SAFETY 229.50 EACH



Paquette Blade Holder

Palatal, Gingival &
Connective Tissue Grafts
Uses Standard Double-Edged
Razor Blade
15cm Length

#SCALP-PAQUETTE
..... 236.00 EACH



Double Bladed Scalpel Handle

1.5mm Blade Separation
Uses #12 or #15 Blades

#SCALP-DUAL-1.5
..... 190.80 EACH

Bendable Miniature & Micro Scalpel Blades

Uncompromised Sharpness & Cutting Ability

Blades Are Sterile & Individually Wrapped



Bendable Miniature Surgical Blades

	Mini Blade #64
#BLAD-MINI-64	71.20 BOX OF 12
	Mini Blade #67
#BLAD-MINI-67	71.20 BOX OF 12
	Mini Blade #63 - Full Radius Edge
#BLAD-MINI-63	71.20 BOX OF 6
	Mini Blade #69 - Full Radius Edge
#BLAD-MINI-69	71.20 BOX OF 6

MOST POPULAR

Bendable Micro-Miniature Surgical Blades

	Micro Blade #64
#BLAD-MICRO-64	55.10 BOX OF 6
	Nordland Micro Blade #6900 - Full Radius Edge
#BLAD-NORDLAND-69	105.00 BOX OF 6

MOST POPULAR

MOST POPULAR

Blades Shown Actual Size



Blade Handle For Miniature & Micro Surgical Blades

#BLAD-HANDLE	78.60 EACH
--------------	------------



Miniature & Micro Surgical Blade Handle

#TVA-BLAD-HANDLE	147.60 EACH
------------------	-------------

Standard Scalpel Blades



Carbon Steel Blades For Regular Scalpel Handles

Extremely Sharp, Sterile & Individually Wrapped

	1 Box	2 Or More Boxes
#BLAD 12	BOX OF 100 69.90	BOX 65.60
#BLAD 12B	BOX OF 100 69.90	BOX 65.60
#BLAD 15	BOX OF 100 69.90	BOX 65.60
#BLAD 15C	BOX OF 100 99.90	BOX 94.60

Blades May Be Assorted For "2 Or More Box" Pricing


Implant Model Package Specials

SPECIAL**Bar Overdenture Model + Crown & Bridge Model****MOST POPULAR**

MOD D4



MOD-C7

**Implant Model Special #1**

Regular Price For Package

Is 903.00 - **Save 129.00**

#MOD-SPECIAL-D4&C7

774.00 FOR BOTH MODELS**SPECIAL****Locator™ Model + Crown & Bridge Model****MOST POPULAR**

MOD LOC-2



MOD-C7

**Implant Model Special #2**

Regular Price For Package

Is 820.00 - **Save 125.00**

#MOD-SPECIAL-LOC2&C7

695.00 FOR BOTH MODELS**SPECIAL****Bar Overdenture + Locator™ + Crown & Bridge + "On 4" Models****MOST POPULAR**

MOD D4



MOD LOC-2



MOD-C7



MOD-COMplete-4

**Implant Model Special #3**

Regular Price For Package

Is 1706.00 - **Save 251.00**

#MOD-SPECIAL-SETOF4

1455.00 FOR ALL 4 MODELS**SPECIAL****Set Of Both Mandible & Maxilla Bone Loss Models****MOST POPULAR**

MOD MM



MOD MMM

**Set Of Both Mandible & Maxilla Bone Loss Models**

#MOD-SET

259.00 FOR BOTH MODELS

Crown & Bridge Implant Models



3 Unit Bridge On 2 Implants
Single Crown On 1 Implant
#MOD C2 334.00



MOST POPULAR
3 Unit Bridge On 3 Implants
Single Crown On 1 Implant
3 Unit Natural Tooth Bridge
#MOD-C7 465.00



Single Posterior Implant Crown
Single Anterior Implant Crown
3 Unit Natural Tooth Bridge
#MOD-C8 465.00



6 Implant High Water Bridge
#MOD HWB-6 362.00



6 Implant Full Arch Cemented
#MOD C4 542.00



MOST POPULAR
Complete "On 4" Restoration
#MOD-COMPLETE-4 448.00



Complete "On 6" Restoration
#MOD-COMPLETE-6 465.00



*Most of Our Models Can Be Customized
With Your Practice Or Study Club Name
Please Call For Quantity Pricing*

More Models On Next Pages →

Hader Bar Implant Models

**MOST POPULAR**

4 Implant Hader Clip Bar
Lower Overdenture

#MOD D4..... 438.00



6 Implant Hader Clip Bar
Upper Overdenture

#MOD-D6-UPPER..... 490.00

"O" Ring Abutment Implant Models



2 Implant "O" Ring
Lower Overdenture

#MOD D2-O..... 288.00

**MOST POPULAR**

4 Implant "O" Ring
Lower Overdenture

#MOD D4-O..... 319.00



*Most of Our Models Can Be Customized
With Your Practice Or Study Club Name
Please Call For Quantity Pricing*

Locator™ Abutment Implant Models



MOST POPULAR

4 Implant Locator™ Abutment
Upper Overdenture
#MOD-LOC-4-UPPER 447.00



4 Implant Locator™ Abutment
Lower Overdenture
#MOD-LOC-4 447.00



3 Implant Locator™ Abutment
Lower Overdenture
#MOD-LOC-3 417.00



MOST POPULAR

2 Implant Locator™ Abutment
Lower Overdenture
#MOD-LOC-2 355.00

Locator™ Is A Trademark of Zest Anchors, Inc.



Most of Our Models Can Be Customized
With Your Practice Or Study Club Name
Please Call For Quantity Pricing

Bone Loss Models

**MOST POPULAR****Mandible Bone Loss Model Set**

- 4 Mandibles Show Progressive Bone Loss
- Shows Benefits Of Augmentation & Consequences Of No Treatment

#MOD MM 141.00 EACH

Maxilla Bone Loss Model Set

- 4 Maxillas Show Progressive Bone Loss
- Shows Benefits Of Augmentation & Consequences Of No Treatment

#MOD MMM 170.00 EACH

Atrophied Ridge Mandible

Atrophied Ridge Mandible Shows Vertical & Horizontal Defects

#MAND-RIDGE-ATROPH 59.80 EACH

SPECIAL

Set Of Both Mandible & Maxilla Bone Loss Models

MOST POPULAR

#MOD-SET 259.00 FOR BOTH



Teaching Models For Hands-On Courses

**MOST POPULAR****Suturable Implant Drilling Mandibles**

Similar In Density To Bone With Suturable Rubber Surface

#MAND SUTURABLE 71.50 EACH

12 OR MORE 68.70 EACH

24 OR MORE 65.80 EACH

**MOST POPULAR****Implant Drilling Mandibles**

Similar In Density To Bone

#MAND TEACH 24.80 EACH

12 OR MORE 22.90 EACH

50 OR MORE 20.90 EACH

**Drillable Sinus Lift Teaching Model**

Similar In Density To Bone With Rubber Sinus Membrane



#SINUS TEACH 125.00 EACH

12 OR MORE 116.20 EACH

Call For Study Club & Implant Course Volume Pricing

Penguin^{RFA} Implant Stability Monitor (Resonance Frequency Analysis)

- *Easy-To-Use, Reliable & Affordable Way To Measure Implant Stability*
- *Cordless & Lightweight*
- *Provides Accurate ISQ Measurement*



Penguin^{RFA} II

- **Includes:** Stand-Up Wireless Charging Base, MultiTegTM Driver, And ISQ Test Peg
- Rubberized Grip For Easy Handling
- 2 Year Warranty

#PENGUIN-22495.00



MultiTegTM Multi-Use Sensing Devices

- Available For All Major Implant Systems
- Made Of Durable, Tissue Friendly Titanium
- Laser Marked With Type Number
- ISQ Standard CalibratedTM
- Autoclavable Up To 20 Times

CALL TO CONFIRM WHICH PEGS YOU NEED



Penguin Sterile Covers (20 Pack)

- Durable, Disposable Covers Designed For Use With The Penguin II RFA ISQ Monitors

#PENGUIN-SLEEVE-20 45.00 PACK OF 20

Spotter Smart Implant Detector



*Accurately Locates
Covered Dental Implants*



- Easily Finds Buried Titanium Implants
- Saves Time And Minimizes Recovery Period



Spotter Smart Implant Detector

Includes One Sterilizable Detachable Sensor (Recommended For 10 Uses)

#SPOTTER 755.00




Spotter Sensor Replacements

Package of 5

#SPOTTER-SENSORS-5 140.00 PACKAGE OF 5

Surgical Rulers & Markers



MOST POPULAR

Surgical Ruler

Stainless Steel 6" (15.25cm) Surgical Ruler

#SURGRULER

INDIVIDUALLY 16.40 EACH

2 OR MORE 14.30 EACH



Sterile Marking Pens

Box of 10 Individually Packaged Markers

#TISSUE-MARK

BOX OF 10 55.90 BOX

2 OR MORE BOXES 52.20 BOX

X-Ray Distortion Markers




X-Ray Distortion Markers

- Precise 5mm Stainless Balls Calibrate For Distortion On Periapical & Panoramic Films
- Incorporate Markers Into A Plastic Stent When Taking Radiograph, Then Use 5mm Diameter As A Reference When Taking Measurements On Your Film

#RADIO MARK. 34.50 BAG OF 50

Salvin Implant Drill Guide Set

- *Helps Provide Proper Spacing, Angulation & Parallelism*
- *Autoclavable*



Regular Drill Guide Set Includes:



- Set Of 3 Color Coded Drill Guides & Autoclavable Cassette

#IMPL-GUIDE-SETOF3 214.00 SET OF 3 WITH CASSETTE
US Patent #6869283B2



Designed By
Dr. Harold Sussman

"Regular Implant" Guide Protocol



1. Select Color Coded Guide Based On Implant Size. Position Guide Against Adjacent Tooth At A Right Angle To The Jaw



2. Use Perio Probe To Locate Crest Of Ridge, Then Drill Pilot Hole With 2.0mm Drill Parallel To The Guide Pin



3. Place Guide Pin Into Osteotomy & Align Guide Over Ridge



4. Drill Second Pilot Hole Through Metal Ring, Use Paralleling Pin To Confirm Parallel Alignment With Guide Pin

Salvin Paralleling Guide




MOST POPULAR

Easy Alignment Of Implants In Edentulous Mandibles

- 1.9mm Diameter Guide Pin
- 7mm Length Guide Pin

#PARALLEL-GUIDE 138.00 EACH

Salvin Paralleling Guide Pins



Paralleling Guide Pins

1.9mm Diameter For 2mm Pilot Drill

#GUIDE-PIN-PKG-4. . . . 107.00 SET OF 4 PINS

#GUIDE-PIN-PKG-8. . . . 205.00 SET OF 8 PINS
Includes **FREE** Organizing Block

SPECIAL

Salvin Guide Pin

Set Of 8 Guide Pins With **FREE** Sterilization / Organizing Block

#GUIDE-PIN-PKG-8
205.00 KIT W/BLOCK



Calipers & Gauges




*Flexes For Increased
Tactile Feel While Probing*

Flexible Tip Implant Depth Gauge


Marked At 6, 8, 10, 12, 14 & 16mm

#DEPTH-FLEX 115.40 EACH

**MOST POPULAR****"Ball Tip" Implant Depth Gauge - Titanium**


- Ball Tip Probe Marked At 7mm To 20mm
- Flat End Of Handle Is 6mm Wide For Reference Measurements

#DEPTH-BALL-TIP 103.50 EACH

**Electronic Digital Caliper**

- 0-150mm Range
- For Extra-Oral Use Only

#CALIPER-DIGITAL 99.90 EACH




*Thumb Control
Locks Tips In Place*

**MOST POPULAR****Sliding Caliper**

- Facilitates Positioning of Implants
- Measures Both Inside & Outside Dimensions
- Sliding Depth Measurement Probe
- 0-90mm Range

#CALIPER-SLIDING 308.00 EACH

**Boley Gauge Caliper**

0-100mm Range

#BOLEY 131.60 EACH

Calipers



MOST POPULAR

Ridge Mapping Caliper

- For Mapping Bone Width
- Measures 0-25mm Range

#RIDGE-CALIPER390.80 EACH



MOST POPULAR

Castroviejo Caliper Short

- Measures 0-40mm Range
- 9cm Long

#CASTRO-CAL355.00 EACH



Castroviejo Caliper Short - 30° Angled

- Measures 0-40mm Range
- 9cm Long

#CASTRO-CAL-ANGLED-S355.00 EACH



Socket Wall Caliper

- Measures Flap & Socket Wall Thickness
- Measures 0-10mm Range

#CALIPER-SOCKET58.70 EACH



Castroviejo Caliper Long

- Measures 0-40mm Range
- 17cm Long

#CASTRO-CAL-LONG452.00 EACH



Castroviejo Caliper Long - 30° Angled

- Measures 0-40mm Range
- 17cm Long

#CASTRO-CAL-ANGLED452.00 EACH

Salvin Bendable "Magic" Mirror Handle For Glass Photo Mirrors

MOST POPULAR

Handle Securely Holds
Glass Photo Mirrors



Easily Bends & Stays In Place
For Your Desired Angle

**Bendable Silicone Handle For Glass Photo Mirrors**

Securely Holds Photography Mirror

#PHOTO-MIRR-MAGIC-HND 81.10 EACH
2 OR MORE 76.60 EACH

Dual View Photography Mirror**Dual View Photo Mirror**

View Buccal, Occlusal & Lingual Surfaces With One Mirror

#PHOTO-DUALVIEW 110.90 EACH

Photo Contrasters**Set Of 2 Contrasters Provide Better Visualization Of Shade & Detail****MOST POPULAR**

For Effective Communication With Patient & Laboratory
Black Background For Better Visualization Of Color, Detail,
Shade Variations & Translucencies



Without Contraster



With Contraster

Photos: Dr. Ross W. Nash

**Operative Mirrors On Page 196****Large & Small Contraster Set**

#CONTRASTER-STR-2 69.50 SET
2 OR MORE SETS 65.70 SET

Front Surface Float Glass For Reflected Image With No Distortion

- **Mirror Surface On Both Sides**
- **Sterilization Up To 428°F (220°C)**



MOST POPULAR



MOST POPULAR

Pedo Occlusal

60mm Wide
85mm Long

#PHOTO OCC PED
INDIVIDUALLY

..... 83.90 EACH

2 OR MORE

..... 79.40 EACH

Adult Occlusal

68mm Wide
100mm Long

#PHOTO OCC AD
INDIVIDUALLY

..... 83.90 EACH

2 OR MORE

..... 79.40 EACH

Adult Occlusal XL

68mm Wide
135mm Long

#PHOTO OCC XL
INDIVIDUALLY

..... 106.80 EACH

2 OR MORE

..... 99.80 EACH

Buccal

57mm Wide
152mm Long

#PHOTO BUCCAL
INDIVIDUALLY

..... 83.90 EACH

2 OR MORE

..... 79.40 EACH

Buccal Narrow

41mm Wide
152mm Long

#PHOTO-BUCCAL-2
INDIVIDUALLY

..... 83.90 EACH


2 OR MORE

..... 79.40 EACH

Autoclavable Photo Mirror Cassette

- Silicone Cushions & Locking Crossbar
- Holds Set Of 4 Mirrors:
Pedo Occlusal, Adult Occlusal, Buccal Narrow, Buccal
- 29 x 18.3 x 3.7cm (11.5" x 7.2" x 1.5")

#CASS-SS-MIR-4. 176.70 EACH



SPECIAL


Mirror, Cassette & "Magic Handle" Set

MOST POPULAR



Includes:

- Pedo Occlusal Mirror
- Adult Occlusal Mirror
- Buccal Mirror
- Buccal Narrow Mirror
- Bendable Silicone "Magic" Mirror Handle
- Stainless Steel Cassette



#PHOTO-KIT-4-MAGIC
494.00 COMPLETE SET

Diagnostic Probes



**Titanium
UNC Probe**

#PROBE-UNC-TITANIUM
..... 55.60 EACH
3 OR MORE
..... 53.90 EACH



**Titanium
CP12 Probe**

#PROBE-CP12-TITANIUM
..... 55.60 EACH
3 OR MORE
..... 53.90 EACH



**Stainless Steel
UNC Probe**

#PROBE-UNC
..... 38.90 EACH
3 OR MORE
..... 35.40 EACH



**PerioWise Plastic
Perio Probes**

- Set Of 6 Perio Probes
- Autoclavable Plastic
Will Not Scratch Implants

#PERIO-PROBE-6
..... 117.90 SET OF 6
2 OR MORE SETS
..... 111.70 SET OF 6

Operative Mouth Mirrors & Handles



#MIRR-HNDL
..... 21.80 EACH
3 OR MORE
..... 19.40 EACH



#MIRR-HNDL-A
..... 45.20 EACH
3 OR MORE
..... 41.80 EACH



**#4 Cone Socket
Front Surface Mirror**

Rhodium-Coated To
Reduce Ghost Imaging

#MIRR-4-CS
..... 56.00 BOX OF 12
5 OR MORE BOXES
..... 52.30 BOX OF 12



**#5 Cone Socket
Front Surface Mirror**

Rhodium-Coated To
Reduce Ghost Imaging

#MIRR-5-CS
..... 56.00 BOX OF 12
5 OR MORE BOXES
..... 52.30 BOX OF 12

Implant Buddy



Implant Buddy Square-Type Dental Implant Driver Set

- Service 99% Of Implant Cases
- Prosthetic Connections Included:
 - Straight Screwdriver
 - SCS Screwdriver
 - UniGrip Driver
 - ASC Angled Access Driver
 - Hex 1.3mm/.050 in.
 - Hex 1.2mm/.048 in.
 - Hex 0.9mm/.035 in.
- Short & Long Drivers For Each Prosthetic Screw Connection
 - Short Drivers: Total Length 23mm; Shaft Length 16mm
 - Long Drivers: Total Length 39mm; Shaft Length 32mm
- Tapered And Color Coded Drivers
- Smart Thumb Knob
- Autoclavable
- Magnetic Base









#IMPLANT-BUDDY-SET-SQUARE 599.00 EACH



Implant Buddy LT Latch-Type Dental Implant Driver Set

- Prosthetic Connections Included:
 - Straight Screwdriver
 - Multi-Unit Driver
 - SCS Screwdriver
 - ASC Angled Access Driver
 - UniGrip Driver
 - Hex 1.3mm/.050 in.
 - Hex 1.2mm/.048 in.
 - Hex 0.9mm/.035 in.
- Seven Latch Drivers:
 - Total Length 24mm; Shaft Length 10.5mm
 - Multi-Unit Driver: 20mm
- Tapered And Color Coded Drivers
- Autoclavable
- Magnetic Base

#IMPLANT-BUDDY-SET-LATCH. 405.00 EACH

DRIVER / PROSTHETIC CONNECTION		SYSTEM
	Straight Driver	Straight Systems
	SCS Driver	Straumann®
	UniGrip Driver	Nobel BioCare®
	Hex 1.3mm/.050 in.	Lifecore®, BioHorizons®, Astra®, Zimmer®
	Hex 1.2mm/.048 in.	Biomet 3i®, MegaGen®, Hiossen®
	Hex 0.9mm/.035 in.	Friadent®, Ankylos®, Biomet 3i® Cover Screw
	Angled Driver	Dynamic Systems
	Square 1.3mm/.050 in.	Gold Screw/Nobel BioCare®

SALVIN
EASYReach®
RIGHT ANGLE PROSTHETIC WRENCH

- Excellent Ergonomics With Great Access In Difficult To Reach Areas
- Secures Driver Tip For Patient Safety

▶ Right Angle
▶ Latch Type



Finger Tip Rotation Easily Tightens Or Loosens Screws

Fits "Latch" Driver Tips

MOST POPULAR

EasyReach® Right Angle Prosthetic Wrench

Autoclavable

#EASY-REACH 590.00 EACH



SALVIN®
AUTO TORQ®
Automatic Electric Torque Driver

- Adjustable Torque From 10Ncm To 35Ncm
- Adjustable Speed From 30rpm To 50rpm
- Includes Reverse Function
- Excellent Ergonomics With Great Access In Difficult To Reach Areas
- Secures Driver Tip For Patient Safety

▶ Torque Wrench
▶ Adjustable Torque
▶ Right Angle
▶ Latch Type

One Button Operation To Achieve Desired Torque



Fits "Latch" Driver Tips



AutoTorq® Automatic Electric Torque Driver

#AUTOTORQ-35 1545.00 EACH

Now Torques
Up To 35Ncm



SalvinTorq® Adjustable Torque Wrench

▶	Torque Wrench
▶	Adjustable Torque
▶	Square Connection



Nine Torque Settings:
15, 20, 25, 30, 35, 40,
45, 50 & 60 Ncm

MOST POPULAR

Adjustable Torque Wrench

#SALVIN-TORQ 629.00 EACH

Use With Driver Tips On Page 197 & 201



Lubrication Hole
At 15Ncm Setting

SPECIAL

SalvinTorq® Package

SalvinTorq® Wrench, Thumb Knob & Cassette
With Your Choice Of Four Driver Tips

#SALVINTORQ-KIT-4 905.00 KIT

SalvinTorq® Wrench, Thumb Knob & Cassette
With Your Choice Of Six Driver Tips

#SALVINTORQ-KIT-6 951.00 KIT

Please Specify Implant System & We'll Make Sure That You Have The Correct Driver Tips



Salvin Torque Wrenches

▶	Torque Wrench
▶	Square Connection



20 Newton Centimeters
#TORQ 20237.90 EACH



30 Newton Centimeters
#TORQ 30237.90 EACH



35 Newton Centimeters
#TORQ 35237.90 EACH

Use These Wrenches With Driver Tips On Page 197 & 201

SALVIN® ACCESSTORQ®-PLUS

Right Angle Variable Torque Driver

Perfect For All-on-4® Procedures

▶	Torque Wrench
▶	Adjustable Torque
▶	Right Angle
▶	Latch Type

MOST POPULAR

AccessTorq®-PLUS
Right Angle Variable Torque Driver

Autoclavable

Lubricate With Pana-Spray (Page 92) After Each Use
& Before Sterilization

#ACCESSTORQ-PLUS 1099.00 EACH



Fits "Latch"
Driver Tips



All-on-4® is a trademark of NobelBiocare®

"Latch Type" Prosthetic Driver Tips

		.048" (1.22mm) Tip Hex Driver #HN-LATCH-DRIVER 60.50 EACH
		.05" (1.25mm) Tip Hex Driver #HT-LATCH-DRIVER 60.50 EACH
		.05" (1.25mm) Tip Hex Driver - Short #HT-LATCH-DRIVER-S 60.50 EACH
		.05" (1.25mm) Tip Square Driver #CER-LATCH-DRIVER 60.50 EACH
		PROC Tip "Star" Driver #PROC-LATCH-DRIVER 83.80 EACH
		PROC Tip "Star" Driver - Short #PROC-LATCH-DRIVER-S 83.80 EACH
		LOC Triangle Driver #TRIRA 60.50 EACH
		ITI Driver #ITI-LATCH-DRIVER 76.50 EACH
		ITI Driver - Short #ITI-LATCH-DRIVER-S 76.50 EACH




"Latch" Tip To Square Driver Adapter

*Use Your Latch Type Driver Tip
With These Salvin Torque Wrenches*














#LATCH-ADAPT 46.90 EACH



Square Connection Abutment Driver Tips - Fit Salvin Square Connection Torque Wrenches

Specify Implant Systems When Ordering – We'll Make Sure You Get The Right Driver Tips!

 <p>.048" (1.22mm) Hex - 6mm #HNT 60.90 EACH</p> <p>.048" (1.22mm) Hex Long - 12mm #HNTL 60.90 EACH</p> <p>For NobelBiocare, 3i, Keystone, MegaGen, Hiossen, Dentium</p>	 <p>.05" (1.25mm) Hex - 6mm #HTT 60.90 EACH</p> <p>.05" (1.25mm) Hex Long - 12mm #HTLT 60.90 EACH</p> <p>For Zimmer, Astra, BioHorizons, Intra-Lock, MIS, SteriOss, Sybron, Implant Direct, PerioSeal</p>	 <p>PROC Driver - 6mm #PROC-S 79.90 EACH</p> <p>PROC Driver Long - 12mm #PROC-L 79.90 EACH</p> <p>"Star" Driver For NobelBiocare Active & Unigrip Screw, Neoss</p>
 <p>Screw Driver - 6mm #SN 60.90 EACH</p> <p>Screw Driver Long - 12mm #SNL 60.90 EACH</p> <p>Flat Blade Screwdriver</p>	 <p>.035" (0.9mm) Hex - 6mm #HC 60.90 EACH</p> <p>.035" (0.9mm) Hex Long - 12mm #HCL 60.90 EACH</p> <p>For NobelBiocare, 3i, Keystone, BioHorizons, Friadent, PerioSeal</p>	 <p>.05" (1.25mm) Square - 6mm #CER-S 60.90 EACH</p> <p>.05" (1.25mm) Square Long - 12mm #CER-L 60.90 EACH</p> <p>For 3i</p>
 <p>ITI Driver - 6mm #ITIS 79.90 EACH</p> <p>ITI Driver Long - 12mm #ITIL 79.90 EACH</p> <p>"ITI" Driver For Straumann SCS</p>	 <p>LOC Triangle Driver - 6mm #TRIS 59.50 EACH</p> <p>LOC Triangle Driver Long - 12mm #TRIL 59.50 EACH</p> <p>For Locator™ Abutment</p>	 <p>.07" (1.77mm) Hex - 6mm #HTSI 59.50 EACH</p> <p>For Paragon TSI Abutment</p>
 <p>.08" (2.03mm) Female Hex - 6mm #HFD 59.50 EACH</p> <p>For NobelBiocare, 3i</p>	 <p>TG Driver #TG 93.40 EACH</p> <p>For 3i TG Abutment, Straumann SCS</p>	

Abutment Driver Adapters & Handles

Use Torque Wrench For Final Seating

Torque Wrench Adapter For NobelBiocare



#TORQ-ADAP-NOBEL 92.50 EACH
Allows Use Of All Salvin Driver Tips
On This Page With Your Wrench

Torque Wrench Adapter For Straumann



#TORQ-ADAP-ITI 92.50 EACH
Allows Use Of All Salvin Driver Tips
On This Page With Your Wrench



Reversible Ratchet

#RATCHET 169.00 EACH



Thumb Knob

Convenient For Starting Screws

#TK 82.60 EACH

Implant & Abutment Names Are Trademarks Of Their Respective Companies

Titanium Tip Implant Currettes

Langer 1/2 Titanium Curette

Angled Like A Gracey But
With Two Cutting Edges

#LANGER-1/2 78.00 EACH

Titanium Tipped



Langer 3/4 Titanium Curette

Angled Like A Gracey But
With Two Cutting Edges

#LANGER-3/4 78.00 EACH

Titanium Tipped



Langer 5/6 Titanium Curette

Angled Like A Gracey But
With Two Cutting Edges

#LANGER-5/6 78.00 EACH

Titanium Tipped



SPECIAL

Set Of 3 Langers #TI-LANG-SET-3 185.00 SET OF 3

H6/H7 Titanium Curette

#TI-H6/H7 78.00 EACH

Titanium Tipped



204 SD Titanium Curette

#TI-204SD 78.00 EACH

Titanium Tipped



SPECIAL

H6/H7 & 204 SD #TI-H6-204-SET-2 131.40 SET OF 2

SPECIAL

5 Instrument Titanium Curette Set



MOST POPULAR

Includes:

- Titanium Langer 1/2
- Titanium Langer 3/4
- Titanium Langer 5/6
- Titanium H6/H7
- Titanium 204 SD
- Stainless Steel Cassette

#TI-CURETTE-SET-5

442.00 SET OF 5 W/CASSETTE



XP² Technology® – No sharpening required!

- **PRECISION-MADE IN THE USA:** 36-step manufacturing process with a hand-crafted touch
- **SHARP, LONG-LASTING EDGES:** Made with premium materials for maximum possible life
- **LIGHTWEIGHT & ERGONOMIC:** Ultra-light EagleLite® handles minimize hand fatigue and enhance control

PART NUMBER	DESCRIPTION	PRICE
DOUBLE GRACEY		
AEDGAXPXE	DBL GRACEY ANTERIOR XP ² PURPLE	79.00 EACH
AEDGPXPXE	DBL GRACEY POSTERIOR XP ² YELLOW	79.00 EACH
AEDGMAXPXE	DBL GRACEY MINI ANTERIOR XP ² GREEN	79.00 EACH
AEDGMPXPXE	DBL GRACEY MINI POSTERIOR XP ² RED	79.00 EACH
GRACEY - STANDARD		
AEG5-6XPXE	GRACEY 5-6 XP ² YELLOW	79.00 EACH
AEG7-8XPXE	GRACEY 7-8 XP ² GREEN	79.00 EACH
AEG11-12XPXE	GRACEY 11-12 XP ² ORANGE	79.00 EACH
AEG13-14XPXE	GRACEY 13-14 XP ² BLUE	79.00 EACH
AEG15-16XPXE	GRACEY 15-16 XP ² ORANGE	79.00 EACH
AEG17-18XPXE	GRACEY 17-18 XP ² BLUE	79.00 EACH
GRACEY +3 - ACCESS (EQUALS MINI FIVE™ *)		
AEGA00-0XPXE	GRACEY ACCESS 00-0 XP ² PURPLE	79.00 EACH
AEGA1-2TXPXE	GRACEY ACCESS 1-2 PRO THIN XP ² PURPLE	79.00 EACH
AEGA1-2XPXE	GRACEY ACCESS 1-2 XP ² PURPLE	79.00 EACH
AEGA5-6XPXE	GRACEY ACCESS 5-6 XP ² YELLOW	79.00 EACH
AEGA7-8TXPXE	GRACEY ACCESS 7-8 PRO THIN XP ² GREEN	79.00 EACH
AEGA7-8XPXE	GRACEY ACCESS 7-8 XP ² GREEN	79.00 EACH
AEGA11-12XPXE	GRACEY ACCESS 11-12 XP ² ORANGE	79.00 EACH
AEGA13-14TXPXE	GRACEY ACCESS 13-14 PRO THIN XP ² BLUE	79.00 EACH
AEGA13-14XPXE	GRACEY ACCESS 13-14 XP ² BLUE	79.00 EACH
AEGA15-16TXPXE	GRACEY ACCESS 15-16 PRO THIN XP ² ORANGE	79.00 EACH
AEGA15-16XPXE	GRACEY ACCESS 15-16 XP ² ORANGE	79.00 EACH
GRACEY +3 - DEEP POCKET (EQUALS AFTER FIVE™ *)		
AEGDP5-6XPXE	GRACEY DEEP POCKET 5-6 XP ² YELLOW	79.00 EACH
AEGDP7-8XPXE	GRACEY DEEP POCKET 7-8 XP ² GREEN	79.00 EACH
AEGDP11-12XPXE	GRACEY DEEP POCKET 11-12 XP ² ORANGE	79.00 EACH
AEGDP13-14XPXE	GRACEY DEEP POCKET 13-14 XP ² BLUE	79.00 EACH
AEGDP15-16XPXE	GRACEY DEEP POCKET 15-16 XP ² ORANGE	79.00 EACH
CURETTES		
AECB5-6TXPXE	BARNHART 5-6 PRO THIN XP ² RED	79.00 EACH
AECB5-6XPXE	BARNHART 5-6 XP ² RED	79.00 EACH
AECC13-14XPXE	COLUMBIA 13-14 XP ² BLUE	79.00 EACH
AECC2L-2RXPXE	COLUMBIA 2L-2R XP ² ORANGE	79.00 EACH
AECC4L-4RXPXE	COLUMBIA 4L-4R XP ² RED	79.00 EACH
AECY7-8XPXE	YOUNGER-GOOD 7-8 XP ² BLACK	79.00 EACH
SCALERS		
AES204SXPXE	204S XP ² YELLOW	79.00 EACH
AES204SDXPXE	204SD XP ² YELLOW	79.00 EACH
AESBJTXPXE	BLACKJACK PRO THIN XP ² BLACK	79.00 EACH
AESBJXPXE	BLACKJACK XP ² BLACK	79.00 EACH
AESB513XPXE	BOGE 513 XP ² PURPLE	79.00 EACH
AESECPXE	EAGLE CLAW XP ² RED	79.00 EACH
AESH6-7XPXE	H6-7 (U6-7) XP ² PURPLE	79.00 EACH
AESM23XPXE	M23 XP ² BLUE	79.00 EACH
AESM23TXPXE	M23 PRO THIN XP ² BLUE	79.00 EACH
AESN128DEXPXE	NEB 128 (DE) XP ² ORANGE	79.00 EACH
AESN1XPXE	N1 XP ² ORANGE	79.00 EACH
AESN2XPXE	N2 XP ² RED	79.00 EACH

Abutment & Prosthetic Parts Holders



MOST POPULAR

Hebel-Melton Abutment Holder

Self Tensioning Handle Keeps Parts Secure During Pick Up & Transfer

#HEBEL-MELTON-TI 166.10 EACH
3 OR MORE 159.20 EACH



Titanium Parts Forcep

For Slip-Free Grasp Of Prosthetic Parts
13cm

#TI-PARTS-GRIPPER 198.00 EACH
3 OR MORE 189.60 EACH



TC Gripper Forcep - Carbide Tipped

For Slip-Free Grasp Of Prosthetic Parts
14.5cm

#TC GRIPPER 194.80 EACH
3 OR MORE 185.50 EACH

Premier® Autoclavable Graphite Implant Scalers

MOST POPULAR

Universal Implant Scaler - Columbia 4L/4R

- Safe For Use Around Implants
- Resharpen With Clean Sharpening Stone
- Autoclavable Up To 140°C

#IMPL-SCALER-UNIV 215.00 SET OF 5
2 OR MORE 204.90 SET OF 5



Implant-Prophy+™ Plastic Hygiene Instruments

Plastic Implant Scalers With Sharpening Stone

- Special Plastic Will Not Scratch Abutment Surface
- Easily Resharpened With Sharpening Stone
- Autoclavable Up To 190°C (375°F)

Set Includes:

- 1 Gracey 5/6
- 1 Gracey 11/12
- 1 Gracey 13/14
- 1 Columbia 13/14U
- Sharpening Stone

#PROPHY-INSTRSET 145.00 SET OF 4



Crown Grippers



Carbide Tipped Crown Gripper

Slip-Free Gripping Of Crowns,
Temporaries & Prosthetic Parts

#CROWN GRIP252.40 EACH
2 OR MORE242.90 EACH



Carbide Tipped Locking Crown Gripper

Slip-Free Gripping Of Crowns,
Temporaries & Prosthetic Parts

#CROWN-GRIP-LOCKING252.40 EACH
2 OR MORE242.90 EACH



MOST POPULAR

Rubber Tipped Crown Plier



Includes 10g Of Silicon Carbide
Gripping Powder

#CROWN-PLIER-II.335.00 EACH
2 OR MORE326.40 EACH

Replacement Silicon Carbide Gripping Powder

10g Container

#CROWN-PLIER-POWDER.47.90 EACH

Replacement Crown Plier Tips



#CROWN-PLIER-TIPS-II 34.80 SET OF 2
2 OR MORE SETS 32.30 SET OF 2

Crown Splitters

- *Slit The Crown With A Diamond*
- *Easily Pry Crown Apart*
- *Reinforced Tip*



Crown Splitter – Straight

#CROWN-SPLIT-STR. 128.40 EACH



MOST POPULAR

Crown Splitter – Offset

#CROWN-SPLIT 159.00 EACH

Improv® Temporary Implant Cement For Implants & Natural Teeth

- 90 Seconds Working Time
- Final Set In 3.5 Minutes



Temporary Implant Cement

- 90 Seconds Working Time
- Compressive Strength: 57 MPa
- Water Sorption: 1.314mg/cm³
- Final Set In 3.5 Minutes
- Film Thickness: <20 Microns
- Shelf-Life: 1.5 Years
- Contains:
 - 2 - 6mL Dual-Barrel Syringes
 - 12 - Mixing Tips
 - 1 - Mixing Pad

#IMPROV-CEMENT. 123.00 BOX
3 OR MORE 113.80 BOX



Extra Mixing Tips

Package Of 30

#IMPL-CEM-TIPS 52.90 BOX

PRECISION. PARTNERSHIP. POSSIBILITY.

For 45 years, Preat has delivered trusted, innovative implant prosthetic solutions for dental professionals – and we're proud to bring their state-of-the-art products directly to you. With the Salvin + Preat partnership, we make it easier than ever to simplify your implant restoration workflows, improve outcomes and deliver smiles.

Why We Love Working with Preat:

- Industry-leading compatibility across implant systems
- Proven attachment and restorative solutions
- Exceptional technical support and education

A Shared Commitment

As a partner and sister company to Salvin, Preat shares the same dedication to quality, service, and advancing dental innovation. Together, we help you keep patients smiling longer!

Learn more at info.salvin.com/preat



45° Surgical Handpiece	92	Collagen Plugs	27	HemCon® Dental DressingPro	41
Abutment Holders	204	Contrasters, Photography	194	Hemostatic Sponges & Gauze	41-42
AccessTorq®-Plus Torque Wrench	200	Cortical Bone Fibers, OsseoFiber®	11	Hemostats	42
ACM Autobone Chip Maker	126	Covers, Light Handle	65	High Speed Handpieces	92
AED	47	Covers, Machine & Mayo Stand	65	Horizontal Centrifuge	70
Air-N-Go® Air Polishing	76	Covers, Shoe	63	Hygiene, Implant	43, 202-203
Allograft Putty In A Syringe, FiberPutty	9	Crest-Control, Meisinger	143-144	Ice Packs	64
AlloSculpt®-3D Plus Bone Putty	10	Crestal Sinus Lift Kit	135, 138	Impla-Safe Peri-Implant Set	78-79
AlloSpark-GFT™ Growth Factor	8	Criticare Monitors	44-45	Implant Abutment Holders	204
Alveoplasty Burs	87	Crown Grippers	205	Implant Brush, Titanium	77
American Eagle Scalers & Curettes	203	Crown Removal Burs	87	Implant Buddy Implant Driver Set	197
AmnioExcite® Amnion/Chorion Membrane	12-13	Crown Removers	205	Implant Center Piezo + Motor	73
Amplifi® Long Term Membrane	30	Curettes, Perio	175	Implant Detector, Spotter	190
Anesthetic Syringes	53	Curettes & Scalers	202-204	Implant Equipment Cart	54
Aspirators, Surgical	98-99	Cyanoacrylate	40	Implant Cement	206
Atraulux® Elevators	104	CytoSurg® PTFE Sutures	36	Implant Drills	82-83
Atraulux® Extraction Forceps	110	CytoSurg® PTFE Membrane	34	Implant Hygiene	43, 202-203
Atraumatic Extraction	100-116	Debridement Burs	84	Implant Positioning Guides	191
Autoclave Tape	65	Defibrillators	47	Implant Models	184-187
Automatic Membrane Tack System	133	Degranulation Burs	87	Implant Motor Stands	54-55
AutoTome® Extraction Kit	91	Depth Gauges	192	Implant Motors & Handpieces	80-81
AutoTome® Osteotome Kit	90	Dermis Allograft	18-19	Implant Packs	58-62
AutoTorq® Electric Torque Driver	198	Diagnostic Instruments	189-196	Implant Removal Kit	93
Balshi Zygomatic Retractor	150	Distortion Markers, X-Ray	190	Implant Scalers	202-204
Bendable Scalpel Blades	183	Drapes, Patient	58-62	Implant Stability, Penguin	189
Benex Extraction System	116	Drill Extender	83	Improv Implant Cement	206
BIO-PRF Horizontal Centrifuge	70	Drill Guides	191	Infection Control	58-68
Bleeding Control	41-42	Drills & Burs	82-88	Instrument Cassettes	66
Blood Collection	69-71	Driver Tips	197, 200-201, 207	Instrument Cleaner	65
Blood Pressure Monitors	44-46	DynaGraft-D and DynaBlast Bone Putty	29	Instrument Wrap & Pouches	65-66
BloodStop® Hemostatic Gauze	42	EasyReach® Prosthetic Wrench	198	Irrigation Needles & Syringes	94
Bone Chip Maker	126	ECG Monitors	45-46	Irrigation Tubing	82, 96
Bone Chisels	117	ECG Patches & Wrist Clips	49	Irrigation, External	94-95
Bone Collector Burs	126-127	Electric Torque Driver	198	Irrigation Syringes	94
Bone Curettes, Lucas	118	Electrocautery	41	Isolation Gowns	62
Bone Expander Kits	144	Elevators	100-109	IV Stand	55
Bone Fibers, OsseoFiber®	11	Elevators, PDL	100	IV Catheters & Admin Sets	50-52
Bone Files	175	Elevators, Periosteal	152-154	Kern Tunneling Instruments	164
Bone Fixation Screws	131	Emergency Equipment	47-49	Knives, Perio	174
Bone Graft Packers	120-121	Enlighten® LED Headlamp	57	Laryngeal Tubes	49
Bone Graft Syringes	122-123	External Irrigation	94-95	Laster 3rd Molar Retractor	150
Bone Mills & Bone Morselizer	129	Extraction Forceps	110-115	LED Headlight	57
Bone Putty, AlloSculpt®-3D Plus	10	Extraction, Atraumatic	100-116	LifeNet Health® Allograft Bone	14-17
Bone Putty, FiberPutty	9	Featherlite Scissors	40, 157	LifeNet Health® Dermis Allograft	18-19
Bone Putty, OsseoFiber®	11	FiberPutty	9	Light Handle Covers	65
Bone Scrapers	124-125	Files, Bone	175	Lindemann Burs	84
Bone Spreaders	142-145	Files, Periodontal	174	Live Cells Bone Allograft	6-7
Bone Trap	126	Fixation Screws	130-132	Locator Implant Models	187
Bone Wax	42	Flap Holding Forcep	181	Long Term Membrane, Amplifi®	30
Bouffant Caps	63	Forceps, Membrane	134	Lucas Bone Curettes	118
Bovine Xenograft, SalvinOss®-B	24	Forceps, Extraction	110-115	Luxating Elevators	100, 104-107
Broken Screw Remover Kit	93	Forceps, Root Tip	114	Machine Covers	65
Bur Blocks	89	Forceps, Suture	179	Mallet, Mead	117
Burs & Drills	82-88	Forceps, Tissue	176-179	Mandibles, Drilling	188
Calipers	192-193	Fragment Forceps	115	Marking Pens, Sterile	190
Cannulas, Nasal	48	Frazao Tunneling Set	165	Masks, N95	63
Capnography Monitors /Accessories	45-46, 49	Furniture, O/R	54-56	Masks, Surgical	63
Cassettes, Instrument	66	Gloves, Sterile	68	Master Core Kit	127
Castroviejo Needle Holders	158-160	Gowns, Sterile & Non-Sterile	62	Mayo Stand Covers	65
Cautery Systems	41	Graft Mixing Dishes	123	Mayo Stands	55
Cement, Implant	206	Graft Packers	120-121	Mead Mallet	117
Centrifuges	70-71	Graft Syringes	122-123	Meisinger Bone Master Core Kit	127
Chair Pads	56	Growth Factor, AlloSpark-GFT™	8	Meisinger Pin Control	132
Cheek Retractors	148-151	Guide Pins	191	Membranes	12-13, 30-34
Chisels, Bone	117	Handpiece Lubricant	92	Membrane Forceps	134
Clamps, Suction	96	Handpiece, Third Molar	92	Membrane Tacks	132-133
CO ₂ Monitors	45-46	Handpieces, Implant	81	Mesh, Titanium	25
CO ₂ Nasal Hoods	49	Hand Scrubs	64	Micro Scissors	157, 163, 173
Cold Packs	64	Headlamp	57	Micro Surgery	162-173

Mindray Vital Sign Monitors	44, 46	Prosthetic Wrench, Right Angle	198-200	Supercut Scissors	155
Mirror Handles	196	Prosthetic Parts Drivers	197-201	Surgeon's Caps	63
Mirrors	194-196	PTFE Membrane, CytoSurg®	34	Surgical Aspirators	98-99
Models, Implant	184-187	PTFE Sutures	36	Surgical Gloves	68
Molt Mouth Gag	148	Pulse Oximeters	44	Surgical Guide Osteotomes	137
Monitors, Vital Sign	44-46	Punches, Tissue	180-181	Surgical Handpieces	81, 92
Motors, Implant	80	Radiographic Marking Balls	190	Surgical Motors	80
Motor Stands	54-55	Ramus Retractors	150	Surgical Rulers	190
Mouth Props	151	Rapid Roll Covers	97	Surgical Tables	56
N95 Respirator Masks	63	Renovix-Plus Collagen Membrane	31	Sussman Implant Positioning Guides	191
Nasal Cannula	48	Rescue Kit, Implants & Broken Screws	93	Suture Book	35
Nasal Hoods, CO ₂	49	Resnik Implant Pack	60	Suture Forceps & Pliers	179
Needle Holders	158-161	Resonance Frequency Analysis (RFA)	189	Suture Scissors	40, 155-157
Needles, Irrigation	94	Respirator Masks, N95	63	Sutures	35-39
Neumeyer Root Planing Bur	85	Retractors, Surgical	148-151	Synthetic Graft Material, Orastruct®	25
Newtron® Ultrasonic Scaler	76	Ridge Contouring Bur	85	Syringes, Anesthetic	53
No-Clog Aspirators	99	Ridge Mapping Caliper	193	Syringes, Bone Grafting	122-123
Non-Glare Microsurgery Kit	162-163	Right Angle Prosthetic Wrench	199-200	Syringes, Irrigation	94
Nordland Micro Blade	183	Ridge Splitting Instruments	142-144	Table Covers	65
O/R Gowns	62	Rongeurs	146-147	Tables, Surgical	56
O/R Towels	62	Root Planing Bur, Neumeyer	85	Tacks, Membrane	132-133
OraCELL® Dermis Allograft	18-19	Root Tip Forceps	114	Tenting Screws	130
OraFoam™ Hemostatic Sponges	42	Rotary Tissue Punches	180	Third Molar Handpiece	92
OraGEN® Viable Dental Bone Allograft	6-7	RotoBrush-Titanium™	77	Third Molar Retractor	152
OraMEM®+ Collagen Membrane	32	Rulers, Surgical	190	Tissue Forceps	178-181
OraPLUG® Collagen Plugs	27	SalvinMesh-Ti™ Titanium Mesh	25	Tissue Forceps, Micro	177
Orasoothe Wound Dressing	43	SalvinOss® Porcine Xenograft	20-21	Tissue Marking Pens	190
Orastruct® Synthetic Graft Material	25	SalvinOss®-B Bovine Xenograft	24	Tissue Punches	180-181
OraTAPE & OraCOTE Wound Dressings	26	SalvinOss® Collagen Xeno + Bovine	22-23	Tissue Rakes	130
Orban Knife	174	Saw Discs	83	Tissue Scissors	155-157
OsseoFiber® Allograft Bone Putty	11	Scalers & Curettes	202-204	Titanium Implant Brush	77
OsseoScape™ Bone Scraper	124	Scalpel Blades & Handles	182-183	Titanium Mesh	25
Osseous Trap	126	Scissors, Micro	157, 163, 173	Torque Wrenches	197-201
Osteotomes, Bone Spreading	142-145	Scissors, Supercut	155	Towel Clamps	62
Osteotomes, Sinus	136-137	Scissors, Suture	40, 155-157	Towels, O/R	62
Packs, Sterile	58-62	Scissors, Tissue	155-157	Trephines	88
Paralleling Guides	191	Scoring Bur, Sinus	85	Tubing Covers	97
Parts Cassettes	67	Screws, Fixation	130-132	Tubing, Handpiece Irrigation	82, 96
Patient Drapes	58-62	Scrub Brushes	64	Tubing, Suction	96-97
PDL-Evators®	100	Scrub Towels	63	Tunneling Forceps	172
P.E.T Partial Extraction Therapy Burs	86	Sectioning Bur	85	Tunneling Instruments	162-174
Penguin RFA Implant Stability Analysis	189	Sedation Monitors	44-46	Twist Elevators	108
Peri-Implant Management Set	78-79	Serrated Elevators	105	Ultrasonic Scaler, Newtron®	76
Peri-Implantitis Brushes & Burs	77-79	Serrated Periostomes	103	Ultrasonic Surgery, Piezotome	73-75
Periacryl	40	Shoe Covers	63	Urban Master Core Kit	127
PerioDT Gel Biofilm	43	Side Cutting Needle Holder	158	Vacutainer®	69
Perio Curettes	175	Silverstein Suture Book	35	Vital Sign Monitors	44-46
Perio Knives	174	Sinus Crestal Kit	135, 138	"Where Is It?" Parts Cassettes	67
Perio Probes, Plastic & Stainless	196	Sinus Lateral Burs	135, 139	Wire Cutter	25
Periodontal Surgical Files	174	Sinus Lift Instruments	135-141	Wrenches, Torque	197-201
Periosteal Elevators	152-154	Sinus Osteotomes	136-137	Wrist Clips, ECG	49
Periostomes	100-103	Sinus Rongeur	139, 147	X-Ray Distortion Markers	190
Photography Contrasters	194	Sinus Scoring Diamond	85	Xenograft, Bovine, SalvinOss®-B	24
Photography Mirror Handle	194	Socket Caliper	193	Xenograft, Porcine, SalvinOss®	20-21
Photography Mirrors	194-196	Soft Tissue Manipulation Kit	173	Xeno + Bovine, SalvinOss® Collagen	22-23
Piezotome Ultrasonic Surgery	73-75	Solution Bowl	95	Xymphony No Memory Membrane	33
Pikos Burs	84-85	Split-Control, Meisinger	143-144	Yankauer Suction	48, 95
Pikos Ramus Retractor	150	Sponges, Collagen	42	ZAGA Zygomatic Instrument Kit	119
Pilot Drills	83	Spotter Smart Implant Detector	190	Zygomatic Instruments	119
Pin-Control, Meisinger	132	Stands, Implant Motor	54-55		
Plugs	27	Stands, Mayo	55		
Pocket Packers	121	Sterile Packs	58-62		
Porcine Xenograft, SalvinOss®	20-21	Sterilization Cassettes	66		
PowerTack™ Automatic Tack System	133	Sterilization Wrap & Pouches	65-66		
Preat Implant Prosthetics	207	Sterillium Rub Surgical Hand Scrub	64		
Precordial Stethoscope	45	Stethoscope, Precordial	45		
Prep Cup	95	Suction Cut-Off Valves	96		
PRF Box	71	Suction Tubing	96-97		
Probes, Plastic & Stainless	196	Suction, Emergency	48		

Salvin Dental Specialties, LLC

3450 Latrobe Drive
Charlotte, NC 28211 U.S.A.

CHANGE SERVICE REQUESTED



Toll Free 800.535.6566
www.salvin.com

BOUND PRINTED MATTER
US POSTAGE PAID
PROGRESS PRINTING PLUS

**Whether You Are New To Implant Dentistry Or Performing
Advanced Procedures, We've Got The Solutions You Need To
Make Your Life Easier And Your Practice More Productive.**

SALVIN®

**Everything For Your Implant Practice
But The Implants®**